

# Detect the IMPOSSIBLE

## Protein-Protein Interaction for Flow Cytometry

Traditional flow cytometry has been limited in the ability to detect protein interactions and low abundant proteins events — until now. We have combined Duolink® Proximity Ligation Assay (PLA), [SigmaAldrich.com/Duolink](http://SigmaAldrich.com/Duolink), with flow cytometry in a convenient kit, making the analysis of protein-protein interactions with flow cytometry readouts a reality. With the enhanced sensitivity of Duolink®

flowPLA Kits, detection of low-abundant protein targets in cell populations is now possible. This easy-to-follow protocol works with any flow cytometry instrument to achieve publication-ready results.

The Duolink® flowPLA Starter Kits provide enough material to perform ~40 reactions, each with 100,000 cells in 100 µL reaction volume.

### Duolink® flowPLA Starter and Detection Kits

Starter Kits include PLA probes and necessary reagents

Cat. No.	Description	Excitation (nm)	Excitation Laser Line (nm)*	Emission (nm)
DUO94102	Duolink flowPLA Starter Kit Mouse/Rabbit-Green	495	488	527
DUO94104	Duolink flowPLA Starter Kit-Mouse/Rabbit-Far Red	644	638, 640, 642	669
DUO94105	Duolink flowPLA Starter Kit Mouse/Rabbit-Violet	390	405	476
DUO94001	Duolink® flowPLA Detection - Red	594	532, 561	624
DUO94002	Duolink® flowPLA Detection - Green	495	488	527
DUO94003	Duolink® flowPLA Detection - Orange	554	5,32,561	579
DUO94004	Duolink® flowPLA Detection - FarRed	644	638, 640, 642	669
DUO94005	Duolink® flowPLA Detection - Violet	390	405	476

\*The excitation laser line represents commonly used lasers that excite the fluorophore. It does not necessarily reflect the lasers available for each instrument.

It is recommended that Duolink flowPLA Compensation Beads be used when performing a multiplexed experiment to save time, reagent cost and precious sample.


### Duolink flowPLA Compensation Beads

Cat. No.	Description
DUO84011	Red
DUO84012	Green
DUO84013	Orange
DUO84014	Far Red
DUO84015	Violet
DUO84010	Unstained



# Protein-Protein Interaction for Flow Cytometry

**Figure 1.** Increased amplification time during Duolink® PLA can aid in the detection of low-abundant protein targets by flow cytometry. Duolink® PLA was performed to detect the trimethylation of lysine 27 on histone 3 (H3K27me3) mediated by EZH2. A) Few PLA signals (red) in the nuclei (blue) of DU145 cells were detected by fluorescence microscopy after 100 min amplification. FITC-Phalloidin-stained actin (green) was used as a counterstain. B) Extended amplification times enhanced the detection of low-abundant protein events, such as EZH2-H3K27me3 interactions, by conventional flow cytometry. C) Combining Duolink® PLA with imaging flow cytometry allows localization of proteins or protein events (interactions or modifications) in large cell populations.



For more information, visit [SigmaAldrich.com/flowpla](https://www.sigmaaldrich.com/flowpla)

To place an order or receive technical assistance in Europe, please call Customer Service:  
 France: 0825 045 645 Spain: 901 516 645 Option 1  
 Germany: 069 86798021 Switzerland: 0848 645 645  
 Italy: 848 845 645 United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840  
 Or visit: [MerckMillipore.com/offices](https://www.merckmillipore.com/offices)  
 For Technical Service visit: [MerckMillipore.com/techservice](https://www.merckmillipore.com/techservice)

Merck KGaA  
 Frankfurter Strasse 250  
 64293 Darmstadt, Germany

[MerckMillipore.com](https://www.merckmillipore.com)

© 2020 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Sigma-Aldrich and Duolink are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

