



THE RESEARCHER'S GUIDE TO CRISPR SCREENING



1. Ensure the cell line is a good model in terms of relevance, biological process & genotype
2. Do you need a primary, transformed or stem cell platform?
3. Determine if the cell line can be adapted to your workflow
4. Consider the doubling time and ploidy of the cell line

1. Are you interested in the whole genome or a more focused pathway?
2. Lentivirus or ribonucleoprotein?
3. Pooled or arrayed?
 - Pooled: maximize the number of gRNA per gene target
 - Arrayed: optimize the gRNA design
4. Controls: use non-targeting guides and consider controls for enrichment and depletion depending on your screening approach
5. Use optimal designs for gRNA and—if designing your own libraries—spread them to avoid clusters in inaccessible genomic regions

1. Perform a kill curve to determine the concentration of selection antibiotic needed to kill the untransfected or untransduced cells
2. Determine the functional titer in your intended cell line using
 - A colony forming unit assay based on antibiotic resistance, or
 - A vector containing a fluorescence marker like GFP
3. Use a control vector to optimize the multiplicity of infection (MOI). Use the lowest MOI that offers one gRNA per cell

1. Cas9 expressing cell lines: perform clonal isolation or use the mixed population of Cas9 expressing cells for screening.
2. All-in-one vectors: deliver both the Cas9 effector and gRNA by introducing one construct
3. Considerations for the optimal Cas9 source:
 - Ensures constant expression levels in a uniform genetic background
 - Eliminates concerns about co-transduction of gRNAs
 - Supports high-throughput sgRNA applications

STEP 5

Perform Your Screen

Pooled and arrayed screens have similar workflows with some differences:

STEP	POOLED	ARRAYED
Library Preparation	1000s gRNAs per tube	1 gRNA per well
Library Delivery	Lentivirus required	Multiple feasible formats
Screen Duration	Efficient whole genome screening	Time to screen increases with the number of clones
Screen Capability	<i>in vivo</i> screening possible	<i>in vivo</i> screening not possible
Analysis	Deep sequencing/deconvolution required to analyze data/identify hits	NGS is not required to understand results
Readout	Limited options (e.g. cell death or proliferation) but can be coupled with single cell analysis	Multiple options e.g. fluorescence, luminescence, high content, live cell imaging

[SigmaAldrich.com/Screening](https://sigmaaldrich.com/screening)

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. MK_PS8196EN 36160 072021

Sigma-Aldrich®
Lab & Production Materials