Alternative Method of Hybrid SPE Sample Preparation for High Recovery LC-MS Analysis of Pharmaceutical Compounds in Plasma

T409228

Authors: Xiaoning Lu, and Michael Ye Presenter: Carmen T. Santasania Supelco/Sigma/Aldrich



#### Outline:

#### Introduction to HybridSPE-PPT

Alternative method – Methanol-based method, two examples.

**Summary** 



### Introduction of HybridSPE PPT Technique


A sample cleanup method that removes both proteins and phospholipids in a simple platform.

- Protein are removed by precipitation with addition of organic solvents (e.g. acetonitrile, methanol).
- Phospholipids are removed by proprietary zirconia particles.


SIGMA-ALDRICH

The operation is both simple and fast, and is amenable to high throughput.

### How Are Proteins and Phospholipids Removed?



## Interaction of Phospholipids with Zirconia coated Silica Surface



#### **Standard Protocol:**

Protein precipitation agents including organic solvent and additives

1% formic acid acetonitrile


primary recommended procedure.

 exhibits high recovery for broad range of test compounds in terms of recovery and protein PPT efficiency.  Precipitate Proteins by adding 100 μL plasma or serum to the HybridSPE-PPT plate followed by 300 μL 1% formic acid in acetonitrile. Add I.S. as necessary.

 Mix by vortexing/shaking HybridSPE-PPT plate or by aspirating/dispensing with 0.5-1 mL pipette tip (e.g., TOMTEC Quadra liquid handler)

3) Apply vacuum. The packed-bed filter/frit assembly acts as a depth filter for the concurrent physical removal of precipitated proteins and chemical removal phospholipids. Small molecules (e.g., pharma compounds and metabolites) pass through unretained.

4) Resulting filtrate/eluate is free of proteins and phospholipids and ready for immediate LC-MS-MS analysis; or it can be evaporated and reconstituted as necessary prior to analysis



1. Load sample and protein precipitation agent

2. Mix

#### 3. Apply Vacuum

4. Samples ready for analysis

#### SIGMA-ALDRICH®

#### Standard Protocol: Acetonitrile -Based Protein Precipitation

•In some cases, strongly basic compounds can exhibit lower recovery from the HybridSPE-PPT.

•Secondary ion-exchange interactions with Zirconia coated Silica surface and strongly basic compounds can cause lower recovery

•In some cases, limited solubility of analyte in acetonitrile can be cause of lower recovery

SIGMA-ALDRICH


Investigate alternative method for these types of compounds

#### Alternative Method: Methanol-Based Protein Precipitation

- •1% ammonium formate (NH4FA) in methanol
- Milder condition for unstable metabolites when compared to 1% formic acid acetonitrile
- In cases, can exhibit higher recovery for basic compounds.
- Methanol can exhibits increased solubility for some pharmaceutical compounds.
- •Methanol can disrupt secondary hydrogen bonding of analytes with silica surface



### Alternative Method: Methanol-Based Protein Precipitation



### Protein Precipitation Comparison between Acetonitrile and Methanol




acetonitrile 1% formic methanol 1% fo acetonitrile metha

1% formic acid 1 methanol f

1% ammonium formate methanol

#### Comparable or better protein precipitation with Methanol.


### Methanol Method – Complete Phospholipid Removal from 100 μL of Rabbit Plasma




SIGMA-ALDRICH®

# Examples of Improved Recovery of Basic Compounds

Recovery of Standard without matrix plasma					
Analyte	1% formic acetonitrle	methanol	1% ammonium formate methanol	1% ammonium chloride methanol	150mM NaCl methanol
Mirtazapine (266/195)	0.0	13.2	96.0	38.2	99.0
Risperidone (411/191)	0.0	10.4	99.1	111.6	64.0
Olanzapine (313/256)	0.0	13.6	89.4	NO experiment	74.0



Mirtazapine



Risperidone



Olanzapine,

SIGMA-ALDRICH®


Potential Interactions between basic compounds and the Zirconia coated Silica

- Ion-exchange interaction is significant for basic compounds having at least one tertiary amine group.
- The ion-exchange interactions were effectively suppressed by the addition of salts such as ammonium formate, ammonium chloride and NaCl.
- Additions of additive raise effective pH of sample, decreasing ionization of analytes
- Ammonia ions bind with silanol surface thus decreasing surface activity
- The salt additives are readily dissolved in methanol, limited in acetonitrile.

SIGMA-ALDRICH

Improved Recovery and Reproducibility of Neutral Hydrophobic Compounds

> Customer experienced low an irregular recovery using 1% formic acid acetonitrile protocol
> Variation of recovery 15-30%
> Analyte had limited solubility in acetonitrile




Drospirenone, C24H30O3, MW 366.227 (mono), ACD/LogD (pH 7.4): 3.16 MRM: 367.2→97.1 or 91

SIGMA-ALDRICH

## Recovery and Reproducibility using 1% ammonium formate methanol

	Recovery of Drospirenone Spiked in Rat Plasma		
Replicate	10 ng/mL spike	60 ng/mL spike	
1	89.9	91.9	
2	87.0	93.7	
3	91.8	88.4	
4	96.9	93.5	
5	91.9	90.4	
6	91.0	87.4	
7	94.2	87.5	
8	97.0	88.0	
9	86.1	87.5	
10	87.0	87.9	
11	84.0	87.8	
12	86.9	87.7	
Avg	90.3	89.3	
STD	4.3	2.4	
%CV	4.7	2.7	

Why better recovery and reproducibility?



Drospirenone has a limited solubility in acetonitrile causing increase interaction with precipitated proteins. This resulting precipitate caused fluctuation in analyte recovery



#### Summary

•An alternative and improved method using 1% ammonium formate in methanol

•Demonstrate improved recovery of model basic compounds over standard acetonitrile method

•Can improve recovery and reproducibility for analytes with low solubility in acetonitrile.

•The more basic conditions of the 1% ammonium formate in methanol is less aggressive towards unstable metabolites.