

How to Protect Your HPLC Column and Prolong Its Life

Noneluting sample components accumulating in HPLC columns can alter a column's chromatographic properties. Large particles ($>2\mu\text{m}$) in mobile phases or samples accumulate on the frit at the column inlet, disrupting the uniform flow of samples onto the packing bed. Smaller particles enter the column and obstruct spaces in the packing bed, increasing back pressure. This bulletin describes sample and in-line filters, guard columns, and other protective devices an HPLC analyst should use to protect high performance analytical columns.

Key Words:

- guard columns
- HPLC column care
- mobile phase

Only materials that will be eluted by the mobile phase, within a specified time, should (ideally) be introduced onto an analytical HPLC column. Yet, samples, mobile phases, and even normal wear on the system can introduce particles and noneluting compounds onto a column. These materials degrade column performance.

To keep an analytical column performing at its best, protective devices must be incorporated into the HPLC system. This bulletin will first evaluate ways to protect a column from contaminants introduced by the system itself (including the mobile phase), then consider those that originate in the sample.


Routine Operation of an HPLC System Can Introduce Impurities Onto the Column

Particles and Dissolved Air in the Mobile Phase

Particles in the mobile phase can accumulate in the column, increasing back pressure. They can also damage pump seals and check valves. Consequently, only HPLC grade solvents and reagents should be used, and all water and buffers should be filtered through a $0.45\mu\text{m}$ filter (as shown in Apparatus 1 or 2 on page 7) before they are introduced into the system reservoir. If the system can be used to mix mobile phases from components in two or more reservoirs, you should determine in advance that the components will not form emulsions or precipitates on mixing. Furthermore, air might be more soluble in the separate liquids than in the mixture. If so, gas bubbles can form as the components mix, or after the mobile phase has passed through the column and is no longer under pressure. In the detector, bubbles increase noise or cause spikes in the signal. Figure A clearly indicates the importance of filtering and degassing mobile

Figure A. Filtering the Mobile Phase Reduces Baseline Noise

Column: SUPELCO SIL LC-8, 15cm x 4.6mm, $5\mu\text{m}$ particles
 Cat. No.: 58220-U
 Temp.: 30°C
 Mobile Phase: methanol:water (50:50) (v/v)
 Flow Rate: 2.5mL/min
 Det.: 254nm UV, 0.005 AUFS



713-0023

phases. Both air and particles can be conveniently removed from the mobile phase by passing the material through a $0.45\mu\text{m}$ filter into a vacuum flask.

A filtered, degassed mobile phase may reacquire particles from the air, through interaction with the reservoir, or through slow crystallization of its components. Thus, the mobile phase should be refiltered immediately before it passes through the pump. A typical mobile phase filter is a porous stainless steel element attached to the reservoir end of the tube between the mobile phase reservoir and the pump inlet (suction side). Filters designed for this purpose exclude most particles larger than $2\mu\text{m}$ and can be cleaned by periodic backflushing. Examples of these filters appear on page 8.

High pH Mobile Phases

Silica-based HPLC columns generally are most stable in the pH range of 2-7. More strongly acidic mobile phases may remove a bonded phase from its silica base, whereas alkaline mobile phases tend to dissolve silica. Contact between the mobile phase and the silica is minimized if the bonded phase uniformly coats the silica surface. When an alkaline mobile phase must be used, the analytical column may be further protected by saturating the mobile phase with silica (1,2). This is accomplished by inserting a $7.5\text{cm} \times 4.6\text{mm}$ precolumn filled with porous silica particles ($18\mu\text{m}$ diameter, $500\text{m}^2/\text{g}$ surface area) between the pump and the sample injection valve. Silica saturator kits are found on page 12.

Particles from the Injection Valve and Pump Seals

Normal wear of the injection (sampling) valve rotor seal and pump seals generates particles of various sizes. Small particles that pass through the pores in the column inlet frit will eventually clog the packing bed, causing high back pressures. Such material should be trapped in a guard column (detailed later in this bulletin.) Particles larger than the pores accumulate on the frit, causing symptoms that mimic column failure.

SUPELCOSIL™ HPLC columns filled with 3µm and 5µm packings have 0.5µm and 2µm frits, respectively, to confine the packing bed. The column end fittings are designed to distribute the stream over the entire frit. However, in extended experiments we have found that large particles derived from the injection valve rotor seal accumulate on the center of the inlet frit. Back pressure does not always increase when the frit is partially blocked, but the uniform distribution of samples onto the top of the column bed is disrupted. As a result, peak asymmetry increases and column efficiency is reduced.


The effects of partially obstructing the inlet frit are illustrated in Figure B. Figure B1 shows the sharpness and symmetry of a toluene peak eluted from a new SUPELCOSIL LC-18 column. After 3000 simulated sample injections, the shape of the peak had deteriorated (Figure B2), suggesting that column efficiency had been reduced. The column was then removed from the system and reversed. The former column outlet was connected to the injection valve, and the column was briefly flushed with mobile phase at a high flow rate. The column was not connected to the detector during this operation, to prevent accumulated particles on the frit from being flushed into the detector. Subsequently, the column was reinstalled in its original orientation. Injection of a test sample revealed (Figure B3) that much of the peak distortion evident in Figure B2 was due to disruption of the mobile phase flow profile, not to column damage *per se*.

Any well-packed column (i.e., one having ≥70,000 plates/meter for a 5µm packing or ≥120,000 plates/meter for a 3µm packing) should be operable with flow in either direction. We recommend periodically reversing the flow direction to keep the frits clean. This procedure is preferable to changing frits because the top of the column bed can be damaged during frit replacement.

Particles larger than 2µm can be prevented from reaching the frit in the first place by introducing a zero dead volume filter with a replaceable screen between the sample injection valve and the column inlet. As shown in Table 1, a zero dead volume filter does not detract from column efficiency, as measured in terms of N (theoretical plates) or AF (asymmetry at 10% of peak height, calculated for toluene in our standard reversed phase test mixture). Furthermore, it is far easier to change a filter or a screen than it is to backflush a column. Filters for this purpose are described on the products pages of this bulletin.

Figure B. Large Particles Obstructing the Column Inlet Frit Can Imitate Column Deterioration

Column: SUPELCOSIL LC-18, 15cm x 4.6mm, 5µm particles
 Cat. No.: 58230-U
 Temp.: 30°C
 Mobile Phase: methanol:water (66:34) (v/v)
 Flow Rate: 2mL/min
 Sample: 10µL of test mixture (Cat. No. 58278) for Figures B1, B2, and B3; 3000 injections of 500µL mobile phase (60 injections/min. at 2mL/min., at 60°C) between Figures B1 and B2



713-0022

Table 1. A Replaceable Screen Filter Does Not Affect Column Efficiency (N) and Peak Asymmetry (AF) as Functions of Flow Rate[▲]

| Flow Rate | Column Alone | | Column + Filter | |
|-----------|----------------|-----------------|-----------------|------|
| | N [▲] | AF [■] | N | AF |
| 1mL/min | 10340 | 1.00 | 10525 | 0.98 |
| 2mL/min | 8320 | 1.01 | 8310 | 1.01 |
| 3mL/min | 6885 | 0.97 | 6740 | 1.03 |

▲ Mean for three measurements


▲ N = theoretical plates = $5.54 \left(\frac{t_R}{W}\right)^2$, where t_R = retention time and w = peak width at 50% of peak height.

■ Asymmetry factor for toluene (width of right-hand part of peak/width of left-hand part of peak at 10% of peak height).

Pressure Surges


Routine operation of a sample injection valve produces transient interruptions in flow to the column. Pressure falls at the head of the column after flow is interrupted, then surges when the flow is reinstated. When the column is operated at relatively high flow rates and pressures, and with a slow injection valve (0.5-1 sec. switching time), repeated surges can eventually cause disruptions in the packing bed. This potential problem can be avoided by installing a Rheodyne model 7725 or similar valve, which incorporates Rheodyne's patented pressure relief technology. Alternatively, you can use a fast, automatic switching valve (<0.1 seconds).

Figure C. Refillable 5cm Guard Column and Connecting Hardware



713-0024

Figure D. Disposable 2cm Cartridge Guard Column, Cartridge Holder, and Connecting Hardware



713-0025

Table 2. Efficiency of an HPLC Column without and with a 5cm Guard Column Filled with a Pellicular Packing

| Analytical Column | Column Efficiency (plates/15cm) ^o | | |
|-------------------|--|-------------------|------------|
| | Without Guard Column | With Guard Column | Difference |
| SUPELCOSIL LC-8 | 11,020 | 10,580 | -4% |
| SUPELCOSIL LC-18 | 11,530 | 10,850 | -6% |

Column: **SUPELCOSIL LC-8 (58220-U) or LC-18 (58230-U), 15cm x 4.6mm, 5µm particles, with Pelliguard™ LC-8 (59643) or LC-18 (59644), 5cm x 4.6mm, 40µm packing**
 Temp.: ambient
 Mobile Phase: methanol:water (60:40) (for LC-8) or (66:34) (for LC-18)
 Flow Rate: 1mL/min
 Sample: 10µL test mix (Cat. No. 58278)

^oCalculated for toluene (k' approximately 3.0).

The Sample as a Contributor to Column Damage

Particles Larger Than 2µm

Measures to protect the analytical column from particles include carefully cleaning the samples through filtration, solid phase extraction (SPE), or other procedures. Sometimes, organic solvents are used to extract a sample from the original material, after which the sample is subjected to centrifugation to remove particles. Often the sample volume is small, and this complicates particle removal. SPE generally gives cleaner samples than comparable liquid-liquid extractions. The extracts, however, may need to be filtered to remove traces of the SPE packing. A small volume, disposable 2µm filter conveniently traps particles in samples obtained by either method. The filters can be used with samples from less than a milliliter to several milliliters in volume. Disposable 0.2µm or 0.45µm pore disk-type filters are useful for removing fine impurities from samples. This type of filter usually has a void volume of 0.1mL, so large samples are needed.

Late-Eluting and Noneluting Sample Components

To minimize run time, analysts prefer to have sample components elute from the column at a k' of between 1 and 10 (i.e., t_R for the last compound of interest = 2-11 x t_0). If k' for an analyte is less than 1, the peak usually will be too close to the solvent front (t_0 peak), and if k' is greater than 10, the analysis time is too long.

Very rapidly eluted compounds ($k' < 0.5$) that are not important usually have no undesirable effects on the analysis. On the other hand, late-eluting materials from one sample may appear as broad peaks or baseline drift in subsequent analyses. To minimize this problem, you can use gradient analysis or periodically flush the column with a stronger solvent, preferably with the direction of flow reversed. Alternatively, you can use a guard column to isolate undesirable sample fractions, via a column switching procedure described later in this bulletin.

Noneluting materials progressively cover the surface of the packing at the column inlet, decreasing sample retention and capacity. These changes may eventually alter the chromatographic properties of the packing. For purposes of column protection, small (<2µm) particles and noneluting sample components can be considered together. These materials, like late eluters, can be flushed from the column by reversing the direction of flow and rinsing with a strong solvent. It is more convenient to trap all of these materials on a guard column. Either the entire guard column or the packing can be changed periodically.

^o $k' = (t_R - t_0)/t_0$, where t_0 = elution time for an unretained peak and t_R = elution time for the peak of interest.

Guard Columns as Protective Devices

Guard Columns with Pellicular Packings

High pressure equipment and advanced technology are required for filling analytical or guard columns with 3 and 5 μ m packings. Consequently, analysts have had a strong incentive to prepare guard columns from pellicular (38-44 μ m) packings. Such materials can be poured in dry form into a guard column blank, then vibrated and tapped into a densely packed bed. Guard columns filled with pellicular packings are available commercially as refillable 5cm x 4.6mm columns and disposable 2cm x 4.6mm cartridges.

Refillable 5cm guard columns, filled with pellicular packings, have been widely used and are suitable for many purposes. A refillable guard column has a fixed frit at the end facing the analytical column and a removable frit at the end facing the injector. The guard column is connected to the analytical column, as illustrated in Figure C. When the guard column becomes contaminated, it is usually sufficient to remove only the top 2-5mm of packing and add fresh material. By maintaining an inventory of interchangeable guard columns, you can reduce down time when the guard column in the system must be cleaned and refilled. A drawback to using guard columns with pellicular packings is the accompanying small (about 5%) loss in apparent theoretical plates. This loss is brought about by the large size of the particles and by the extra connections and fittings required to add the guard column to the system (Table 2). A test sample should be analyzed before and after any guard column is installed in the system, to ensure that performance remains within acceptable limits.

In many cases disposable cartridge-type guard columns (Figure D) are more convenient and economical to use than refillable guard columns. Disposable 2cm x 4.6mm cartridge-type guard columns filled with a pellicular packing minimize the loss of operating time during replacement. These 2cm columns, like the 5cm refillable columns with pellicular packings, reduce the number of theoretical plates observed by about 5% (Table 3). Peak asymmetry also increases slightly. These effects are independent of flow rate.

Table 3. Effect of a 2cm Pellicular Guard Column on Analytical Column Efficiency and Peak Symmetry

| Flow Rate | Column Alone | | Column + Guard Column | |
|-----------|--------------|------|-----------------------|------|
| | N | AF | N | AF |
| 1mL/min | 10525 | 0.98 | 10035 | 1.11 |
| 2mL/min | 8310 | 1.01 | 7830 | 1.11 |
| 3mL/min | 6740 | 1.03 | 6650 | 1.10 |

Column: **SUPELCO SIL LC-8 (58220-U), 15cm x 4.6mm, 5 μ m particles, with Pelliguard™ LC-8 (59643), 5cm x 4.6mm, 40 μ m packing**

Temp.: ambient

Flow Rate: 1, 2, or 3mL/min

Mobile Phase: methanol:water (60:40)

Sample: 10 μ L test mix (Cat. No. **58278**)

Data represent mean for 3 measurements.

Table 4. Effects of a Pellicular Guard Column on t_0 , Capacity Factors (k'), and Selectivity (α) of an Analytical Column

| Parameter | Analytical column | Analytical Column + Guard Column |
|---------------------------------|-------------------|----------------------------------|
| t_0 Uracil (min.) | 1.67 | 1.85 |
| k' Acetophenone | 0.80 | 0.73 |
| k' Benzene | 1.59 | 1.45 |
| k' Toluene | 2.85 | 2.61 |
| α (Benzene/Acetophenone) | 1.98 | 1.99 |
| α (Toluene/Benzene) | 1.80 | 1.80 |

Column: **SUPELCO SIL LC-8 (58220-U), 15cm x 4.6mm, 5 μ m particles, with Pelliguard™ LC-8 (59643), 5cm x 4.6mm, 40 μ m packing**

Temp.: 30°C

Mobile Phase: methanol:water (60:40)

Flow Rate: 1mL/min.

Sample: 10 μ L test mixture (Cat. No. **58278**)

Because the primary function of a guard column is to retain materials that the mobile phase will not elute from the analytical column, the retention characteristics of the analytical and guard column packings should be similar (i.e., match silica with silica or octadecyl bonded phase with octadecyl bonded phase). However, a 38-44 μ m pellicular packing has only 1/100th of the surface area a porous 5 μ m analytical packing has, making the capacity of the pellicular packing much smaller than that of the analytical packing. Consequently, all sample components transverse the guard column in essentially equal time. Coupling a 2cm pellicular guard column to a 15cm analytical column increases t_0 and the elution times for retained compounds almost equally. As a result, k' for retained compounds decreases, but selectivity or separation factor ($\alpha = k'_2/k'_1$) is unaffected (Table 4 and Figure E). Because a guard column with a pellicular packing has little or no effect on analyte retention, the chromatographic properties of the guard column and analytical column packings need not match perfectly.

Guard Columns with Five Micron Packings

Guard columns filled with 5 μ m packings offer distinct advantages, compared to those filled with pellicular packings. For instance, a perfect match can be obtained with the analytical packing, and the sample-resolving capability of the system is enhanced. When a 2cm disposable, cartridge-type guard column with a 5 μ m packing is used, the theoretical plate loss attributed to guard column connections and fittings is more than offset by the theoretical plates the packing contributes (Table 5 and reference 3).

As is the case with pellicular guard columns, a small increase in peak asymmetry might be observed when adding a 5 μ m packing guard column to the system. Operating pressure will increase in proportion to the increase in total column length, i.e., 17/15 or 27/25, respectively, for a 15 or 25cm analytical column with a 2cm guard column. Since t_0 and t_r for retained compounds increase proportionately (Figure F), there will be no observable change in k' and α values. Poorly packed or improperly connected guard columns, however, reduce column efficiency and increase peak asymmetry (4).

Figure E. A Guard Column Containing a Pellicular Packing Has Little Effect on Analyte Retention

Column: **SUPELCOSIL LC-8 (58220-U)**, 15cm x 4.6mm, 5µm particles (Figure E1), with Pellicguard LC-8 guard column (59643), 2cm x 4.6mm, 40µm packing (Figure E2)
 Temp.: ambient
 Mobile Phase: methanol:water (60:40)
 Flow Rate: 1mL/min
 Chart Speed: 1cm/min for 6 min., then 10cm/min
 Sample: 10µL of test mixture (Cat. No. 58278)





Figure F. A Guard Column Containing a 5µm Packing Increases t_0 and t_R Proportionately

Column: **SUPELCOSIL LC-8 (58220-U)**, 15cm x 4.6mm, 5µm particles (Figure F1), with Supelguard™ LC-8 guard column (59554), 2cm, 5µm packing (Figure F2)
 Temp.: ambient
 Mobile Phase: methanol:water (60:40)
 Flow Rate: 1mL/min
 Chart Speed: 1cm/min for 6 min., then to 10cm/min
 Sample: 10µL of test mixture (Cat. No. 58278)





Table 5. Effects of a 2cm Guard Column on Efficiency and Peak Symmetry of an Analytical Column Filled with the Same 5µm Packing

| Flow Rate | Column Alone | | Column + Guard Column | |
|-----------|--------------|------|-----------------------|------|
| | N | AF | N | AF |
| 1mL/min | 10525 | 0.98 | 10785 | 1.15 |
| 2mL/min | 8310 | 1.01 | 8850 | 1.13 |
| 3mL/min | 6740 | 1.03 | 7245 | 1.16 |

Column: **SUPELCOSIL LC-8 (58220-U)**, 15cm x 4.6mm, 5µm particles, with Supelguard LC-8 guard column (59562), 2cm, 5µm packing
 Temp.: ambient
 Flow Rate: 1, 2, or 3mL/min
 Mobile Phase: methanol:water (60:40)
 Sample: 10µL test mix (Cat. No. 58278)
 Data represent mean for 3 measurements.

Our experience shows that extra-column effects and dead volumes introduced by pellicular packing guard columns, although small, reduce the efficiency of an analytical column with a 3µm packing. For this reason, only 5µm guard columns should be used with 3µm analytical columns. Basic chromatographic principles explain the problem. The height equivalent of a theoretical plate (HETP) varies with particle diameter, whereas operating pressure varies with the inverse of the square of the particle diameter. Compared to a 5µm packing, therefore, a 3µm packing offers a 1.67 fold ($5/3$) increase in plates at a 2.78 fold ($5^2/3^2$) increase in operating pressure. Pressure on 3µm packings is usually kept within acceptable limits by restricting analytical column length to 7.5 or 15cm. The higher efficiency of columns filled with 3µm packings, combined with short column length, produce smaller peak volumes for sample components than are obtained from columns filled with 5µm packings. Thus, factors peripheral to the column, such as those contributed by a guard column, have a greater effect on the efficiency of columns with 3µm packings.

In addition, 0.5µm frits are used to confine 3µm packings, so the guard column must efficiently retain small (<2µm) particles that otherwise will obstruct the inlet frit. Pellicular packings cannot do this. In contrast, only a small reduction in theoretical plates and increase in peak asymmetry are observed when a 2cm x 4mm disposable guard column with a 5µm packing is used in conjunction with a 7.5cm x 4.6mm, 3µm packing analytical column.

Supelco HPLC guard columns are available in internal diameters of 2.1mm, 3.0mm, and 4.0mm. The 4.0mm ID guards are designed to protect analytical columns of 4.0 and 4.6mm ID.

Figure G. A Sampling Valve Two Pump System for Backflushing a Guard Column




Figure H. A Sampling Valve Single-Pump System for Backflushing a Guard Column





Figure I. A Sample Enrichment System Employing a Sampling Valve and a Guard Column



Using Multiport Valves in Conjunction with Guard Columns

Multiport sampling valves are used to make and break connections between two or more columns and a detector. These valves are equally valuable for backflushing a guard column (to clean it) without manually disconnecting it from the analytical column.

A 10-port sampling valve and an auxiliary high pressure pump, incorporated into the HPLC system along with a guard column, conveniently allows you to protect the analytical column and prolong its useful life. The complete system is illustrated in Figure G. A sample is injected into the sample loop (Figure G1), then the valve is switched so that the sample moves through the guard column into the analytical column (Figure G2). Once the analytes of interest have entered the analytical column, the valve is switched again, isolating the guard column from the system (Figure G1). While the analysis continues, the guard column is backflushed by the second pump. In a variation of this procedure, a specific fraction of the sample is directed from the guard column to the analytical column, while the early portion of the eluant is vented to waste. Similarly, a 10-port valve and one pump can be set up to clean the guard column, but in this system material flushed from the guard column is voided through the detector (Figure H).

A system that includes a guard column (5µm porous packing), a 10-port valve, and two pumps can also be used to concentrate trace analytes (5,6). Generally, the sample is flushed from the sample loop onto the guard column in a weak mobile phase (e.g., water, in reversed phase HPLC). Compounds of interest are retained on the guard column, while less retentive materials are eluted to waste. When the valve is switched, the second pump flushes the analytes from the guard column onto the analytical column, in a stronger mobile phase. The system is illustrated in Figure I. In this application, the guard column does not protect the analytical column from particles and noneluting sample components. A second guard column should be used for this purpose.

When very complex samples must be analyzed, a guard column can be used to separate compounds by classes for more discriminate separations on the analytical column. In this situation, a sampling valve is used to sequentially apply the peaks from the guard column to a single analytical column, or to divert them to different columns.

Additional uses for multiport sampling valves are described in references 6-8. In all of these procedures the underlying goal remains the same. So far as is practical, only analytes that will be eluted by the specific mobile phase, in an appropriate timespan, should be allowed to enter the analytical column.

References

1. Wehrli, A., J.C. Hildenbrand, H.P. Keller, R. Stampfli, and R.W. Frei, J. Chromatogr., 149, 199 (1978).
2. Atwood, J.G., G.J. Schmidt, and W. Slavin, J. Chromatogr., 171, 109 (1979).
3. Reed, G.D. and C.R. Lascombe, Chromatographia, 15, 16 (1982).
4. Lundanes, E., J. Dohl, and T. Greibrokk, J. Chromatogr. Sci., 21, 235 (1983).
5. Schauwecker, P., R.W. Frei, and F. Emi, J. Chromatogr., 136, 63 (1977).
6. Harvey, M.C., and S.D. Stearns, Am. Lab., May 1982, p. 68.
7. Roth, W., K. Beschke, R. Jauch, A. Zimmer, and F.W. Koss, J. Chromatogr., 222, 13 (1981).
8. Sonnefeld, W.J., W.H. Zoller, W.E. May, and S.A. Wise, Anal. Chem., 54, 723 (1982).

References not available from Supelco.


Rheodyne® multi-port sampling valves are available from Supelco in manual, electrical, or air-actuated designs.

Protection for Your HPLC Columns

Your HPLC column is valuable! To maintain highest column efficiency you should (refer to Figure J):

- Filter all samples before injection
- Filter and de-gas solutions used to prepare mobile phases
- Use a filter (A) in the mobile phase inlet line
- Use pH 2-7 mobile phases with silica-based packings, OR use a saturator column (B) for higher pHs
- Use a screen or frit filter (C) after the injector
- Use a guard column (D) between the filter and the analytical column
- Reverse the analytical column routinely to remove particles from the inlet frit

Figure J. Complete Protection for Your HPLC Column



713-0031

Ordering Information:

Supelco Solvent Filtration Apparatus

To protect your instrument and columns, remove particles and gases from solvents, buffer solutions, ion pairing reagents, and other mobile phase components. Apparatus 1 can be connected to a 1000mL side-arm flask for vacuum filtration. Apparatus 2 can be connected to a vacuum line and includes a tapered top (34/45) 1000mL receiving flask. Nylon-66 membrane filters (0.45µm) used in either apparatus are compatible with all solvents commonly used in HPLC.

For other options, refer to the current Supelco catalog.

| Description | Cat. No. |
|--|--------------|
| Filtration Apparatus 1 (250mL glass reservoir, funnel base and stopper, clamp, SS holder and screen, Teflon® gasket, 50 Nylon-66 filters) | 58061 |


Replacement Parts:

| | |
|--|--------------|
| Glass Reservoir, 250mL | 58063 |
| Funnel Base and Stopper | 58064 |
| SS Filter Holder and Screen | 58065 |
| Teflon Gaskets (pk. of 10) | 58066 |
| Nylon-66 Membranes, 0.45µm (pk. of 50) | 58067 |

| | |
|---|----------------|
| Filtration Apparatus 2 (250mL glass reservoir, 34/45 tapered funnel base, 34/45 tapered-top 1000mL flask and glass cap, clamp, SS holder and screen, Teflon gasket, 50 Nylon-66 filters) | 58062-U |
|---|----------------|

Replacement Parts:

| | |
|--|----------------|
| Glass Reservoir, 250mL | 58063 |
| Tapered Funnel Base, 34/45 | 58068 |
| Tapered-Top Flask, 34/45, 1000mL | 58070-U |
| Glass Cap for Flask | 58071 |
| SS Filter Holder and Screen | 58065 |
| Teflon Gaskets (pk. of 10) | 58066 |
| Nylon-66 Membranes, 0.45µm (pk. of 50) | 58067 |



910-0549



910-0113

Iso-Disc™ Syringe-Tip Filter Units for Organic or Aqueous Solvents and Gases

- Unique angled thread on luer fitting provides a more reliable connection than filters with tabs
- No seams to fail under pressure
- Sample cannot bypass the membrane
- Solvent-resistant, autoclavable membranes and housings
- Extractable-free construction

Specifications

| | |
|------------------|---|
| Housing: | polypropylene |
| Membrane: | nylon - use with aqueous solvents PTFE - for gases or with organic solvents cellulose acetate - use with aqueous solvents |
| Pressure rating: | 3mm - 60psi 25mm - 75psi |
| Process volume : | 3mm - ≤2mL 25mm - 2mL to 100mL |
| Holdup volume: | 3mm - <10µL (with air purge) 25mm - <60µL (with air purge) |
| Connections: | female luer lock inlet male luer taper outlet |
| Sterilization: | compatible with ethylene oxide, autoclavable |



910-0118

Iso-Disc Syringe-Tip Filter Units

| Iso-Disc Filter Unit | Membrane | Diameter (mm) | Pore Size (µm) | Color Coding | Quantity | Cat. No. |
|----------------------|-------------------|---------------|----------------|--------------|----------|----------|
| N-254 | Nylon | 25 | 0.45 | green | 50 | 59230-U |
| N-252 | Nylon | 25 | 0.2 | purple | 50 | 59231-U |
| N-34 | Nylon | 3 | 0.45 | green | 100 | 59238 |
| N-32 | Nylon | 3 | 0.2 | purple | 100 | 59239 |
| P-255 | PTFE | 25 | 5.0 | gray | 50 | 59244 |
| P-254 | PTFE | 25 | 0.45 | yellow | 50 | 59234-U |
| P-252 | PTFE | 25 | 0.2 | blue | 50 | 59235-U |
| P-34 | PTFE | 3 | 0.45 | yellow | 100 | 59240-U |
| P-32 | PTFE | 3 | 0.2 | blue | 100 | 59241 |
| CA-254 | Cellulose acetate | 25 | 0.45 | red | 50 | 59242 |
| CA-252 | Cellulose acetate | 25 | 0.2 | white | 50 | 59243 |

Mobile Phase Filter




910-0114

Protect your expensive HPLC column and pumping system against particles in your solvents. Attach a low pressure mobile phase filter to the reservoir end of the tube carrying mobile phase to the suction side of the pump. Stainless steel filter with a 2µm porous filter element. Easily cleansed by backflushing.

| Description | Cat. No. |
|---|----------|
| Mobile Phase Filter, 1/8" Swagelok® Tube Connection | 58268 |

Valco Frit and Screen Filters*



713-0561


Efficient, low dead volume in-line filters protect your columns from particles without reducing column performance. The replaceable 1/8" frit has 0.5µm pores to protect 3µm or 5µm column packings. Replaceable screen has 2µm pores. Choose the frit filter for higher filtration capacity (most applications) or the screen filter for less dead volume (e.g., with microbore columns). 1/16" fittings included. Use with 1/16" OD tubing.

| Description | Cat. No. |
|---------------------------------|----------|
| Frit Filter | 58420-U |
| Replacement Frits (pk. of 10) | |
| 0.5µm pores | 59037 |
| 2.0µm pores | 59129 |
| Screen Filter | 58279-U |
| Replacement Screens (pk. of 10) | 58284 |

*Note: Frits and screens should not be interchanged in these filters.

Guard Cartridges for Discovery, SUPELCOSIL, and SUPELCOGEL Columns

Supelguard guard column cartridges contain a Discovery™, SUPELCOSIL, or SUPELCOGEL™ packing, in a 2cm stainless steel body, enclosed by PEEK encapsulated stainless steel frits (2µm porosity). A **Supelguard kit** (one cartridge, a stand-alone holder, tubing, and 2 nuts and ferrules) enables you to use the cartridge with any analytical column. A **direct-connect guard cartridge holder** can directly connect a Supelguard cartridge to a Supelco analytical column. 3.0mm cartridges are available in packs of two only – purchase a stand-alone or modular holder separately. Low volume 3.0mm ID cartridges are a good choice for protecting an analytical column containing 3µm particles.



897-0038

2cm Supelguard Cartridges with 5µm Discovery, SUPELCOSIL, or SUPELCOGEL Packings

| Supelguard Phase | Use To Protect | 4.0mm ID Cartridge ¹ | | 3.0mm ID Cartridge | | 2.1mm ID Cartridge | |
|------------------------|------------------------|---------------------------------|-------------------|--------------------|--------------|--------------------|--|
| | | Kit Cat. No. | Pk. of 2 Cat. No. | Pk. of 2 Cat. No. | Kit Cat. No. | Pk. of 2 Cat. No. | |
| Discovery | Discovery | | | | | | |
| RP-AmideC16 | RP-AmideC16 | 505080 | 505099 | custom | 505102 | 505110 | |
| ABZ ⁺ Plus | ABZ ⁺ Plus | 59534-U | 59535-U | 59535-C30 | 59604 | 59605 | |
| LC-ABZ | LC-ABZ | 59544-U | 59545-U | 59545-C30 | 59610 | 59611 | |
| Suplex™ pKb-100 | Suplex pKb-100 | 59531-U | 59541-U | 59541-C30 | 59608 | 59609 | |
| Discovery C18 | Discovery C18 | 505129 | 505137 | custom | 505161 | 505188 | |
| LC-18 | LC-18, LC-PAH | 59554 | 59564 | 59564-C30 | 59612 | 59613 | |
| LC-318 | LC-318 | 59502 | 59512 | custom | custom | custom | |
| LC-18-DB | LC-18-DB | 59555 | 59565 | 59565-C30 | 59616 | 59617 | |
| LC-18-S | LC-18-S | 59629 | 59630 | 59630-C30 | 59161 | 59162 | |
| LC-18-T | LC-18-T, LC-DABS | 59620 | 59621 | 59621-C30 | custom | custom | |
| Discovery C8 | Discovery C8 | 59589-U | 59590-U | custom | 59587-U | 59588-U | |
| LC-8 | LC-8 | 59552 | 59562 | 59562-C30 | 59614 | 59615 | |
| LC-308 | LC-308 | 59501 | 59511 | custom | custom | custom | |
| LC-8-DB | LC-8-DB | 59553 | 59563 | 59563-C30 | 59618 | 59619 | |
| LC-DP | LC-DP | 59556 | 59566 | 59566-C30 | custom | custom | |
| LC-3DP | LC-3DP | custom | 59513 | custom | custom | custom | |
| LC-F | LC-F | 59520 | 59521 | 59521-C30 | custom | custom | |
| LC-304 | LC-304 | 59591 | 59592 | custom | custom | custom | |
| Discovery Cyano | Discovery Cyano | 59585-U | 59586-U | custom | 59583-U | 59584-U | |
| LC-CN | LC-CN | 59557 | 59567 | 59567-C30 | custom | custom | |
| LC-PCN | LC-PCN | 59504 | 59514 | 59514-C30 | custom | custom | |
| LC-1 | LC-1 | 59551 | 59561 | 59561-C30 | custom | custom | |
| TPR-100 | TPR-100 | 59570-U | 59571 | 59571-C30 | custom | custom | |
| ODP-50 | ODP-50 | 59312 | 59313-C40 | custom | custom | custom | |
| LC-Si | LC-Si, LC-3Si | 59550 | 59560 | 59560-C30 | custom | custom | |
| LC-NH ₂ | LC-NH ₂ | 59558 | 59568 | 59568-C30 | custom | custom | |
| LC-NH ₂ -NP | LC-NH ₂ -NP | 59515 | 59516 | 59516-C30 | custom | custom | |
| LC-Diol | LC-Diol, LC-3Diol | 59559 | 59569 | 59569-C30 | custom | custom | |
| LC-SAX1 | LC-SAX1 | 59536-U | 59537-U | 59537-C30 | custom | custom | |
| LC-SCX | LC-SCX | 59509 | 59519 | 59519-C30 | custom | custom | |
| Hisep™ | Hisep | 59639 | 59640-U | 59640-C30 | custom | custom | |

¹ For 4.0mm ID and 4.6mm ID analytical columns.

Guard Columns for SUPELCOGEL Ion Exclusion Columns (5cm x 4.6mm column, 9µm particles)

| Supelguard Phase | SUPELCOGEL Column Protected | Cat. No. |
|------------------|-----------------------------|----------|
| H | H, C-610H | 59319 |
| Ca | 611, Ca | 59306-U |
| K | K | 59344 |
| Pb | Pb | 59345 |
| Ag | Ag | 59316 |

HELPFUL HINTS

To protect microbore (1mm ID) columns, use an in-line screen filter (Cat. No. 58279-U).

To protect preparative columns (6.2mm–50.8mm ID), use a 5cm column packed with the same particles and bonded phase as in the preparative column.

To protect other brands of columns, see our general catalog.

Guard Column Holders

Stand-Alone Holder for 2cm Guard Column - Holds 2.1mm, 3.0mm, or 4.0mm ID Supelguard cartridge. Includes 1/16" nuts and ferrules and 2"/5cm of 0.01" ID x 1/16" OD tubing.

Direct-Connect Guard Column Holder - Directly connects a 2cm Supelguard cartridge to a Supelco™ analytical cartridge after the column endfitting is removed.

| Description | Cat. No. |
|-------------|----------|
|-------------|----------|

Guard Column Holders

| | |
|---|---------|
| Stand-Alone | 59660-U |
| Direct-Connect for 2.1mm ID analytical column | 504254 |
| Direct-Connect for 3.0/4.0/4.6mm ID analytical column | 55205 |

Pelliguard Guard Columns

Packed Cartridges

2cm x 4.6mm ID guard columns with 40µm pellicular packings.

- Ready-to-use, disposable
- Low cost
- Large particles for negligible pressure drop
- Sufficient capacity for most samples
- Use with 5µm SUPELCOSIL columns, other columns having corresponding phases



997-0192

Each kit contains one cartridge guard column (filled with a 40µm Pelliguard packing), a reusable column holder, and hardware for connecting the holder to 1/16" tubing. Replacement cartridges come in packages of four.

| Chromatography Mode | 5µm Column to be Protected | Recommended Pelliguard Column Kit | Cat. No. |
|---------------------|----------------------------|-----------------------------------|----------|
| Normal Phase | Silica | 2cm LC-Si Kit | 59641 |
| | | Replacement Cartridges (pk. of 4) | 59651 |
| Polar Bonded Phase | Cyano | 2cm LC-CN Kit | 59645-U |
| | | Replacement Cartridges (pk. of 4) | 59655 |
| | Amino | 2cm LC-NH ₂ Kit | 59646 |
| | | Replacement Cartridges (pk. of 4) | 59656 |
| Reversed Phase | C8 | 2cm LC-8 Kit | 59643 |
| | | Replacement Cartridges (pk. of 4) | 59653 |
| | C18 | 2cm LC-18 Kit | 59644 |
| | | Replacement Cartridges (pk. of 4) | 59654 |

Bulk Pellicular Packing Kits

5cm guard column hardware and 40µm pellicular packing.

Dry pack your own guard columns for maximum economy.

- Reusable 5cm x 4.6mm ID column hardware
- Sufficient capacity for most samples

Each column kit contains an empty 5cm x 4.6mm ID column, 10g of Pelliguard packing, 10 frits, and hardware for connecting the column to 1/16" tubing. About 1.3 grams of packing is needed to pack one 5cm x 4.6mm column.

| Chromatography Mode | 10µm Column to be Protected | Recommended Pelliguard Packing Kit | Cat. No. |
|---------------------|-----------------------------|------------------------------------|----------|
| Normal Phase | Silica | 5cm LC-Si Kit | 58202 |
| | | Pelliguard LC-Si Packing, 10g | 58291 |
| Polar Bonded Phase | Cyano | 5cm LC-CN Kit | 58234 |
| | | Pelliguard LC-CN Packing, 10g | 58235 |
| Reversed Phase | C8 | 5cm LC-8 Kit | 58222-U |
| | | Pelliguard LC-8 Packing, 10g | 58293 |
| | C18 | 5cm LC-18 Kit | 58232 |
| | | Pelliguard LC-18 Packing, 10g | 58294 |

Guard Column Hardware Kit

5cm column, endfittings, 2 frits (2µm porosity), and 2" (5cm) of 0.01" ID x 1/16" OD SS tubing

| Description | Cat. No. |
|-------------------------------|----------|
| Guard Column Hardware Kit | 58319 |
| Replacement Frits (pk. of 10) | 58264 |

Funnel and Tubing for Easier Column Filling

Pour Pelliguard packing easily into your guard columns with this convenient plastic funnel. Connects to column with the Tygon tubing included.

| Description | Cat. No. |
|-------------------|----------|
| Funnel and Tubing | 20390-U |



910-0004

Evaluation Test Mixes for HPLC Columns

Well defined test mixes enable you to troubleshoot chromatographic problems, optimize system efficiency, and evaluate columns. We ship these test mixes in amber ampuls to prevent photodegradation, and we include instructions for proper use and interpretation of results.

Choose from column-specific or application-specific mixes. The amino phase test mix (sugars) calls for refractive index detection; all others call for UV detection. For additional information, request *HPLC Troubleshooting Guide* (Bulletin 826).

HPLC Column Test Mixes

1 mL unless otherwise specified.

| Test Mix | Use to Test | Solvent | Components (conc./mL) | Cat. No. |
|----------------------|---|--|---|----------|
| Amino Phase | LC-NH ₂ columns | acetonitrile:water, 25:75 | D-fructose (25mg) α-D-glucose (25mg) sucrose (25mg) maltose (25mg) lactose (25mg) | 58424 |
| Cyano Phase | LC-CN, LC-PCN columns, any weakly hydrophobic phase | acetonitrile:water, 40:60 | uracil (7μg) acetophenone (7μg) benzene (750μg) toluene (775μg) | 58299 |
| Diphenyl Phase | LC-DP columns | methanol:water, 50:50 | uracil (7μg) benzene (750μg) toluene (775μg) | 58421 |
| Normal Phase Mix 1 | LC-Si (silica) columns | methylene chloride: methanol:water, 99:9:1 | benzene (600μg) benzaldehyde (200μg) acetaldehyde (200μg) | 58281 |
| Normal Phase Mix 2 | LC-Si, LC-CN, LC-NH ₂ columns | ethanol:hexane, 5:95 | toluene (1mg) diethyl phthalate (1mg) dimethyl phthalate (1mg) | 47640-U |
| Nucleosides | LC-18-S columns | water, 10mg/mL sodium formate | 12 nucleosides (10-100μg) | 47310-U |
| LC-PAH | LC-PAH columns | methanol: methylene chloride, 50:50 | 16 PAHs (100-2000μg) | 48743 |
| Peptide Standard | reversed phase columns used for peptide separations (e.g., 300Å phases) | none (dried film) | Gly-Tyr (~0.125mg) Val-Tyr-Val (~0.5mg) Met enkephalin (~0.5mg) Leu enkephalin (~0.5mg) angiotensin II (~0.5mg) | H2016 |
| Reversed Phase Mix 1 | hydrophobic RP columns (e.g., LC-8, LC-18, ABZ ⁺ Plus) | methanol:water, 60:40 | uracil (7μg) acetophenone (7μg) benzene (750μg) toluene (775μg) | 58278 |
| Reversed Phase Mix 2 | hydrophobic RP columns (e.g., LC-8, LC-18, ABZ ⁺ Plus) | acetonitrile:water, 58:42 | uracil (5μg) phenol (700μg) N,N-diethyl-m-toluamide (600μg) toluene (4mg) | 47641-U |
| Chiral 1 | chiral phases | hexane:ethyl acetate, 80:20 | toluene (70μg) (+) TFAE* (25μg) (-) TFAE* (25μg) | 48250-U |
| Chiral 2 | chiral phases | chloroform | benzene (500μg) (+) N-PDBA** (50μg) (-) N-PDBA** (50μg) | 48251 |

*2,2,2-trifluoro-1-(9-anthryl)ethanol

**N-(1-phenylethyl)-3,5-dinitrobenzamide

Custom Test Mixes

For information on made-to-order standards and test mixes, call our Technical Service chemists, or request our Custom Chemical Reference brochure (Publication No. 196905).

Saturator Column Kits



713-0580

Use a silica saturator column with alkaline mobile phases, to prevent dissolution of the silica packing base. Use a C18 saturator column to prevent stripping of phase by acidic mobile phases. Each kit contains a 7.5cm x 4.6mm column, 10g of porous silica (18µm spherical), two frits, and two 1/16" nuts and ferrules.

| Description | Cat. No. |
|--------------------------------------|--------------|
| Silica Saturator Column Kit | 58410 |
| Spherical Silica Packing (18µm), 10g | 58411 |
| C18 Saturator Column Kit | 58418 |
| C18 Packing, 10g | 58419 |
| Funnel and Tubing | 20390 |

Trademarks

Discovery, Hisep, Iso-Disc, Pelliguard, Supelco, SUPELCOGEL, SUPELCOSIL, Supelguard, Suplex — Sigma-Aldrich Co.
 Rheodyne — Rheodyne, Inc.
 Swagelok — Crawford Fitting Co.
 Teflon — E.I. du Pont de Nemours & Co., Inc.
 Valco — Valco Instruments Co., Inc.

Free Technical Literature From Supelco

Supelco produces a wide variety of technical literature, all of which is free. You may request this literature several different ways:

- Use the business reply card located at the back of our catalog.
- Call our Ordering and Customer Service department (800-247-6628 or 814-359-3441).
- Use our ChromFax™ service.

How to Use ChromFax

ChromFax is Supelco's automated digital data system, which allows you to instantly obtain technical literature any time of the day or night. Supelco technical bulletins, application notes, lab hints, and product specification sheets can be sent to your Fax machine.

1. To use ChromFax you must dial from a touch-tone phone. If you do not have a touch-tone phone, call our Ordering and Customer Service department and we will Fax the publication to you.
2. Have the phone number of your Fax machine ready, and dial **800-652-8069 or 814-359-5748**. Remember you must have the *complete* Fax number (011 + country/city code + Fax number) when you are calling from outside the US or Canada, otherwise the document cannot be transmitted.
3. When ChromFax answers, follow the simple recorded instructions. Select the **ChromFax number** of the document you want. Be certain you order documents by their ChromFax numbers and NOT their literature numbers.
4. Enter your Fax number and hang up. The requested information will be sent to you within minutes. To report a problem with the ChromFax system contact our Technical Service department.

Technical Literature Available on ChromFax

The listing of documents available on ChromFax is updated every month. We encourage you to request a list of available ChromFax documents each time you use the ChromFax system.

Contact our Technical Service Department (phone 800-359-3041 or 814-359-3041, Fax 800-359-3044 or 814-359-5468) for expert answers to your questions.

BULLETIN 781

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 **AUSTRALIA** · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 **AUSTRIA** · Sigma-Aldrich Handels GmbH · A-1110 Wien
BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem **BRAZIL** · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP **CANADA** · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 **DENMARK** · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand **FINLAND** · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki
FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex **GERMANY** · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen **GREECE** · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens
HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest **INDIA** · Sigma-Aldrich Co. · Bangalore 560 048 **IRELAND** · Sigma-Aldrich Ireland Ltd. · Dublin 24 **ISRAEL** · Sigma Israel Chemicals Ltd. · Rehovot 76100
ITALY · Sigma-Aldrich s.r.l. · 20151 Milano **JAPAN** · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 **KOREA** · Sigma-Aldrich Korea · Seoul **MALAYSIA** · Sigma-Aldrich (M) Sdn. Bhd. · Selangor
MEXICO · Sigma-Aldrich Química S.A. de C.V. · 50200 Toluca **NETHERLANDS** · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht **NORWAY** · Sigma-Aldrich Norway · Torshov · N-0401 Oslo
POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań **PORTUGAL** · Sigma-Aldrich Quimica, S.A. · Sintra 2710 **RUSSIA** · Sigma-Aldrich Russia · Moscow 103062 **SINGAPORE** · Sigma-Aldrich Pte. Ltd.
SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 **SPAIN** · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid **SWEDEN** · Sigma-Aldrich Sweden AB · 135 70 Stockholm
SWITZERLAND · Supelco · CH-9471 Buchs **UNITED KINGDOM** · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH
UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: supelco@sial.com

H

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

BKN