
Preparation, Separation, Filtration & Monitoring Products



Envirocheck® Dip Slides

Overview

Envirocheck[®] Dip Slides provide semi-quantitative detection of microorganisms on surfaces and in liquids. The flexible paddles, coated with relevant nutrient media, allow for sampling on all surface types. After sampling and incubation, two parameters can be determined in one easy step.

Envirocheck[®] Dip Slides are available for the following parameters:

- Disinfection control (DC)
- Enumeration of total viable count (TVC)
- Enumeration of yeasts and molds (YM)
- Enumeration of Enterobacteriaceae (E)
- Enumeration of total coliforms /*E. coli* (C)

Product & Handling


Main benefits

- Flexible: the dip slides are easy to use, with a pliable paddle allowing sampling even in difficult to access areas.
- **Safe:** before and after sampling, the agar paddle is safely protected against contamination in the plastic tube
- **Economical:** testing in-house provides an efficient, low-cost alternative to independent test lab services
- Easy interpretation: following incubation, results can be read or checked against the provided comparison chart

Molds









Density Charts

Bacteria and Yeasts





CFU/ml, for liquids

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.



Specifications:

Product description	Incubation conditions	Composition	Targeted organisms
Envirocheck [®] Contact Slides DC	35-37 °C for 24-48 hours	• TSA • TSA + LTHTh	 Total count without preservatives
			 Total count in the presence of preservatives or disinfectant residues
Envirocheck [®] Contact Slides TVC	35-37 °C for 24-48 hours	• Nutrient Agar + TTC • Nutrient Agar	 Total viable count with red colony dye
			 Total viable count without dye
Envirocheck [®] Contact YM	30 °C: control TSA Agar after 24–48 hours; control Rose Bengal Agar after 5–7 days	• TSA + TTC • Rose Bengal Chloramphenicol	 Total count with red colony dye
			 Yeast and mold count with bacterial inhibition
Envirocheck [®] Contact E	35–37 °C for 24–48 hours	 Plate Count Agar VRBD Agar 	 Total count (typical food formulation PCA)
			 Enterobacteriaceae count (typical food/meat formulation)
Envirocheck [®] Contact C	35–37 °C for 24–48 hours	 Plate Count Agar Chromocult[®] Coliform Agar 	 Total count (typical food formulation PCA)
			 Simultaneous detection of coliforms and <i>E. coli</i> in water and food samples

Ordering information:

Description	Qty/pk	Cat. No.	Shelf Life	Storage conditions
Envirocheck [®] Contact Slides DC	10	1021470001	8 months	3 – 15 °C
Envirocheck [®] Contact Slides TVC	10	1021490001	8 months	3 – 15 °C
Envirocheck [®] Contact YM	10	1021390001	6 months	3 – 15 °C
Envirocheck [®] Contact E	10	1021370001	8 months	3 – 15 °C
Envirocheck [®] Contact C	10	1021360001	6 months	3 – 15 °C

To place an order or receive technical assistance

In Europe, please call Customer Service: France: 0825 045 645 Germany: 069 86798021 Italy: 848 845 645 Spain: 901 516 645 Option 1 Switzerland: 0848 645 645 United Kingdom: 0870 900 4645


For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: MerckMillipore.com/offices

For Technical Service visit: MerckMillipore.com/techservice

MerckMillipore.com

Merck KGaA Frankfurter Strasse 250 64293 Darmstadt, Germany



© 2019 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Chromocult and Envirocheck are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. MK_DS4152EN Ver. 2.0 2019-23504 10/2019