



## New Stevia Extract Reference Material

### Analysis of Steviol Glycosides using HPTLC

Extracts from the leaves of the *Stevia rebaudiana* plant have a long tradition of being used as a sweetener. Stevia extracts are approximately 300 times sweeter than sucrose while only having a negligible effect on blood glucose and has therefore been increasingly used as an alternative to artificial sweeteners.

### Regulations

Stevia is FDA approved as a dietary supplement and one of the active constituents, rebaudioside A, is considered "Generally Recognized As Safe (GRAS)". The European Community has allowed the use of steviol glycosides as food additives since December 2011, and WHO defined the acceptable daily intake of steviol glycosides to be 4 mg per kg body weight.<sup>1</sup>

### Analysis

HPTLC (High-Performance Thin-Layer Chromatography) is a fast and efficient tool to create molecular fingerprints of complex chemical mixtures. It is therefore particularly well suited for the analysis of plants and plant-derived products.


A new Stevia extract reference material is now available which is provided with comprehensive documentation, including a quantitative value for the major component Stevioside as well as qualitative confirmation of various other constituents (Rebaudiosides A, B, C and D, Dulcoside A, Rubusoside, Steviolbioside, Stevioside). In addition to HPLC with assigned peak identities, the certificate of analysis includes HPTLC fingerprint results.

### HPTLC fingerprint method for the analysis of *Stevia rebaudiana* leaves [2, 3]

#### Samples


Extract: 50 mg were suspended in 50 mL of methanol and sonicated for 10 min. The suspension was centrifuged, and the supernatant used.

Leaf: 0.5 g of powdered leaf were suspended in 30 mL of water and boiled for 10 min. The solution is filtered into a 50 mL volumetric flask and the volume is made up with water.



## Chromatogram after derivatization under UV 366 nm

Track 1: rebaudioside D  
 Track 2: rebaudioside A  
 Track 3: rebaudioside C  
 Track 4: stevioside  
 Track 5: rebaudioside B  
 Track 6: dulcoside A  
 Track 7: steviolbioside  
 Track 8: Stevia extract reference material  
 Track 9: *Stevia rebaudiana* leaf 1  
 Track 10: *Stevia rebaudiana* leaf 2



In the fingerprint of the dry *Stevia rebaudiana* leaf extract (track 8), zones corresponding in color and position to those of the standards Rebaudioside A, C, Stevioside, Rebaudioside B and Dulcoside A (which are

co-eluting), and steviolbioside are seen. The fingerprint is similar to those of *S. rebaudiana* leaf (tracks 9 and 10). Rebaudioside D is only seen in the fingerprint of the leaf, particularly under UV 366 nm (very faint zone).

Description	Quantified Components	Qualitatively Confirmed Components	Cat. No.	Package Size
Stevia extract	Stevioside	Rebaudiosides A, B, C and D, Dulcoside A, Rubusoside, Steviolbioside, Stevioside	<b>06295001</b>	500 mg

**Table 1:** Stevia extract reference material

Description	Format / cm	Pack Size/ No. Plates	Cat. No.
HPTLC Silica Gel 60 F <sub>254</sub>	20 x 10	50	<b>1056420001</b>

Visit our website and find our complete TLC portfolio at [SigmaAldrich.com/tlc](http://SigmaAldrich.com/tlc)

Cat. No.	Product	Package Size
<b>90378</b>	Dulcoside A	10 mg
<b>92273</b>	Isosteviol	10 mg
<b>38462</b>	Rebaudioside A	10 mg
<b>49747</b>	Rebaudioside B	10 mg
<b>30987</b>	Rebaudioside C	10 mg
<b>19189</b>	Rebaudioside D	10 mg
<b>19345</b>	Steviol	10 mg
<b>59754</b>	Steviolbioside	10 mg
<b>50956</b>	Stevioside	10 mg

**Table 3:** Analytical Standards for *Stevia rebaudiana* constituents

To place an order or receive technical assistance in Europe, please call Customer Service:  
 France: 0825 045 645 Spain: 901 516 645 Option 1  
 Germany: 069 86798021 Switzerland: 0848 645 645  
 Italy: 848 845 645 United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840  
 Or visit: [MerckMillipore.com/offices](http://MerckMillipore.com/offices)  
 For Technical Service visit: [MerckMillipore.com/techservice](http://MerckMillipore.com/techservice)

Merck KGaA  
 Frankfurter Strasse 250  
 64293 Darmstadt, Germany

[MerckMillipore.com](http://MerckMillipore.com)

