

## 83931 RNase AWAY® (Decontamination Reagent for Rnase)

### Application

RNase decontamination reagent is a precisely formulated, ready-for-use, ribonuclease decontaminant. When used properly, RNase Decontamination Reagent will effectively remove RNA and RNase contamination from plastic, stainless steel and glass lab. surfaces. Do not use RNase Decontamination Reagent on aluminium, soft metal or gaskets and seals. This product is not intended as reactant. Do not add to reaction.

### Guidelines for use

The following are suggested guidelines only and usage will vary depending upon the surface and the level of contamination present.

Apply RNase Decontamination Reagent liberally to the surface to be decontaminated. It is important to ensure that the solution contacts the entire surface. Vigorously rub the wet surface with a RNase-free laboratory wipe, then dry with a fresh RNase-free wipe. It is important that the entire surface is completely dry as residual RNase Decontamination Reagent may degrade an RNA or DNA sample.

Alternatively, you may soak items overnight in RNase Decontamination Reagent, then rinse thoroughly with RNase-free water and dry.

When cleaning pipettors, make sure that the solution does not come in contact with any seals or gaskets.

Note: Gloves should be worn when handling this product. RNase Decontamination Reagent is alkaline in nature and will cause irritation if prolonged contact with the skin is allowed. In case of contact with eyes, immediately flush with water for fifteen minutes and contact a physician. If swallowed, Do not induce vomiting. Give plenty of water and contact a physician immediately.

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The vibrant M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.  
© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

