Sigma-Aldrich®

Instructions for Use

Silver Stain (Modified GMS) Kit

Procedure No. HT100

Intended Use

The Silver Stain (Modified GMS) kit is intended for use in the histologic demonstration of fungi, basement membrane and some opportunistic organisms. The Silver Stain (Modified GMS) Kit is for "In Vitro Diagnostic Use." For professional use only. The data obtained from this manual, qualitative procedure identifies fungi, basement membrane and some opportunistic organisms in tissue samples of human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of fungal infections.

Silver methenamine-borate procedures are well documented.^{1,2} Generally, these require elaborate solution preparation prior to test performance. Solution stability is limited and results vary remarkably due to the capricious nature of metal impregnation and photographic development.

The Silver Stain (Modified GMS) Kit incorporates a stable working silver methenamine salt along with buffer, toning reagent, and developer. This provides the laboratory a unique package for visualizing fungi, basement membrane and opportunistic organisms.³⁻⁵ Included is an application for rapid staining using a microwave oven.^{3,6}

In brief, microbial cell wall and basement membrane polysaccharides are oxidized to aldehydes by treatment with periodic acid. The aldehyde group, at alkaline pH, reduces silver ion to metallic silver. Rinsing with gold salts forms a more stable gold complex and excess silver is removed by a thiosulfate wash.

Reagents

Periodic Acid Solution (Cat. No. HT1001-100ML)

Periodic acid, 1 g/dL, in deionized water

Borax Solution (Cat. No. HT1002-100ML)

Borax, 5 g/dL, in deionized water. Danger. Causes severe skin burns and eye damage. Wear protective gloves/protective clothing/eye protection/face protection.

Silver Methenamine (Cat. No. HT1003-1VL; HT1003-6VL)

Silver methenamine, 110 mg/vial. Danger. Harmful if swallowed. Causes severe skin burns and eye damage. May cause an allergic or asthma symptoms or breathing difficulties if inhaled. Suspected of causing cancer.

Gold Chloride Solution (Cat. No. HT1004-100ML)

Gold chloride, 200 mg/dL, in deionized water.

Sodium Thiosulfate Solution (Cat. No. HT1005-100ML)

Sodium thiosulfate, 2 g/dL, in deionized water.

Special Materials Required but Not Provided

- Plastic forceps or paraffin coated metal forceps.
- Positive control slides, such as Fungi TISSUE-TROL™ (Cat. No. TTR004) should be included in each run.

Counterstains (optional)

(choice depends on specimen and individual preferences)

- Tartrazine Solution, Cat. No. HT302 [HT3024-120ML; HT3028-250ML]
- Harris Hematoxylin Solution, Cat. No. HHS [HHS16-500ML; HHS32-1L; HHS80-2.5L; HHS128-4L]
- Eosin Y Solution, Alcoholic, Cat. No. HT1101 [HT110116-500ML; HT110132-1L; HT110180-2.5L; HT1101128-4L]
- Eosin Y Solution, Aqueous, Cat. No. HT1102 [HT110216-500ML; HT110232-1L; HT110280-2.5L; HT1102128-4L]
- Eosin Y Solution, Alcoholic, With Phloxine, Cat. No. HT1103 [HT110316-500ML; HT110332-1L; HT110380-2.5L; HT1103128-4L]
- Fast Green FCF
- Light Green SF Yellowish

For Standard Procedure Only:

- Water bath, 62°C

For Microwave Procedure Only:

- Plastic Coplin jars with lids
- Microwave Oven

Storage and Stability

Store Periodic Acid Solution, Working Silver Methenamine, Gold Chloride Solution, and Sodium Thiosulfate Solution refrigerated (2–8°C).

Store Borax Solution at room temperature (18–26°C). NOTE: Borax may form precipitate during refrigeration. Bring to room temperature and redissolve prior to use.

Reagent labels bear expiration date.

Preparation

Prepare Working Silver Methenamine Solution by combining 8 mL Borax Solution, 100 mL deionized water and contents of one Silver Methenamine vial. Mix until dissolved. Use once, then discard.

NOTE: Do not allow Working Silver Methenamine Solution to come into contact with metals.

NOTE: To increase tests per vial to two, Silver Methenamine may be dissolved in 100 mL of

deionized water and stored tightly capped. This solution is stable when refrigerated (2–8°C) for 1 month. When ready to use, take 50 mL of one Silver Methenamine Solution and add 4 mL of Borax Solution. Return unused portion of Silver Methenamine Solution to refrigerator.

Prewarm all reagents to room temperature before use.

Precautions

The IVDs included in this kit are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Fungi TISSUE-TROL™ control slides are paraffin embedded animal tissue containing fungi and should be considered potentially infectious.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Fungus and Opportunistic Organisms such as *P. carinii*

Fix in 10% neutral buffered formalin or Bouin's Solution, embed in paraffin and cut 5 micron sections.¹

Basement Membrane

Fix in 10% neutral buffered formalin, Bouin's or Zenker's Solution and cut paraffin sections 1–4 microns.¹

Notes

Staining dishes should be washed in Dishbrite liquid laboratory glassware cleaning concentrate, Cat. No. S4107, then bleach solution, rinsed well in tap water and then in at least 6 changes of deionized water.

If a microwave oven is used, please see the Owner's Manual for instructions.

Standard Procedure

1. Prepare Working Silver Methenamine Solution as described in "Preparation" section and place in 62°C water bath.
2. Deparaffinize tissue sections and rehydrate to deionized water.
3. Place rehydrated sections in Periodic Acid Solution for **5 minutes** for fungi and opportunistic organisms or **11 minutes** for basement membrane.
4. Remove slides and rinse in 6 changes of deionized water.
5. **For Demonstration of Fungi and Opportunistic Organisms**, place slides in prewarmed Working Silver Methenamine Solution. Examine microscopically after 20 minutes. Fungi and opportunistic organisms should be dark brown against a pale yellow background.
For Demonstration of Basement Membrane, place slides in prewarmed Working Silver Methenamine Solution. Examine microscopically after 30 minutes. Basement membrane should be black against a dark yellow-brown background.
NOTE: While incubating sections at 62°C, place a Coplin jar containing deionized water in 62°C water bath. Use this water for rinsing prior to evaluation in Step 5. Examine sections microscopically for adequate staining at stated intervals. Return slides to warm Working Silver Methenamine Solution until desired intensity is achieved.
6. Remove slides from Working Silver Methenamine Solution after adequate structural staining is achieved and rinse 6 times in room temperature deionized water.
7. Tone sections in Gold Chloride Solution for **30 seconds**.
8. Rinse in room temperature deionized water 6 times.
9. Place sections in Sodium Thiosulfate Solution for **2 minutes**.
10. Wash well in running tap water.
11. Counterstain according to personal preference with Tartrazine Solution, Fast Green FCF, Light Green SF Yellowish, Harris Hematoxylin Solution or Eosin Y Solution.
12. Dehydrate in alcohol, clear in xylene, mount and examine microscopically.

Microwave Procedure

1. Deparaffinize slides and hydrate to deionized water.
2. Place slides in **40 mL** of Periodic Acid Solution contained in a plastic Coplin jar. Loosely cover Coplin jar with lid before placing in microwave oven, or use Coplin jars with holes drilled into the lids. Microwave on **800 watts for 10 seconds**. ***Periodic Acid Solution may be reused if used at room temperature for **5 minutes** instead of microwaving. ***For basement membranes Periodic Acid Solution Incubation time is **11 minutes at room temperature**.
3. Remove slides from Periodic Acid Solution and rinse in 6 changes of deionized water.
4. Place slides in **40 mL** of Working Silver Methenamine Solution contained in a plastic Coplin jar.
5. **For Fungal Organisms** microwave on **600 watts for 35 seconds**. Gently mix solution with a beral pipette or applicator stick. Let incubate for **2–3 minutes**. After no longer than **2 to 3 minutes**, place slides in warm deionized water and check development microscopically. If organisms are not sufficiently developed, return slides to Working Silver Methenamine Solution and let sit for **30 seconds to 1 minute** or until desired tone is achieved.

For Basement Membrane microwave on **600 watts for 35 seconds**. Gently agitate solution. Let slides incubate for **5 minutes**. Microwave a second time at **600 watts for 10 seconds**. Let slides develop for **2 minutes**. Check slides microscopically. If slides are not sufficiently developed, repeat last step until desired tone is achieved.

6. Rinse slides in 6 changes of deionized water.
7. Tone slides in Gold Chloride Solution for **30 seconds at room temperature**.
8. Rinse slides in 6 changes of deionized water.
9. Place slides in Sodium Thiosulfate Solution for **2 minutes at room temperature**.
10. Wash well in running tap water.
11. Counterstain according to personal preference with Fast Green FCF, Light Green SF Yellowish, Tartrazine Solution, or Harris Hematoxylin Solution, or Eosin Y Solution.
12. Dehydrate in alcohol, clear in xylene, mount and examine microscopically.

Performance Characteristics

Target Structure	Staining Result
Fungi	Purplish-brown to black
Opportunistic Organisms	Purplish-brown to black
Basement Membrane	Black
Background	Dependent upon counterstain

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT1001	Periodic Acid Solution	Fungal organisms	3 of 3	3 of 3	3 of 3	3 of 3
HT1002	Borax Solution	Fungal organisms	3 of 3	3 of 3	3 of 3	3 of 3
HT1003	Silver Methenamine	Fungal organisms	3 of 3	3 of 3	3 of 3	3 of 3
HT1004	Gold Chloride Solution	Fungal organisms	3 of 3	3 of 3	3 of 3	3 of 3
HT1005	Sodium Thiosulfate Solution	Fungal organisms	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT100A:

- H271: May cause fire or explosion; strong oxidizer.
- H302: Harmful if swallowed.
- H314: Causes severe skin burns and eye damage.
- H317: May cause an allergic skin reaction.
- H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.
- H351: Suspected of causing cancer.
- 360FD: May damage fertility. May damage the unborn child.
- H373: May cause damage to organs (Thyroid) through prolonged or repeated exposure if swallowed.
- H412: Harmful to aquatic life with long lasting effects.
- P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- P273: Avoid release to the environment.
- P280: Wear protective gloves/ protective clothing/ eye protection/ face protection.
- P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.
- P304 + P340 + P310: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/ doctor.
- P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980

2. Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981

3. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986

4. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986

5. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983

6. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986

7. Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990

8. Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 4.0	2016
Rev. 5.0	2022
Rev. 6.0	2022
	Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards section.
Rev. 7.0	2025
	Updated naming convention. Revised Intended Use statement. Removed references to specific organisms. Added CH-REP and importer information.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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The Initial M, TISSUE-TROL, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

1. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
2. Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
3. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
4. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
5. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
6. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
7. Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago, 1990.
8. Kok & Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon, Leiden, 1992.

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Références

- Sheehan D.C., Hrapchak B.B.: Theory and Practice of Histotechnology, 2nd ed., C.V. Mosby Co. St. Louis (MO), 1980
- Koski J.P.: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J. Histotechnol 4:115, 1981
- Churukian C.J., Schenk E.A., Clark G.: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
- Leong A.S.-Y., Milios J.: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J. Pathol 148:183, 1986
- Brinn N.T.: Rapid metallic histologic staining using the microwave oven. J. Histotechnol 6:125, 1983
- Valle S.: Special stains in microwave oven. J. Histotechnol 9:237, 1986
- Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
- Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

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Bibliografia

1. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
2. Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
3. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
4. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
5. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
6. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
7. Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
8. Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

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Referencer

- Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
- Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
- Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
- Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
- Kok and Boon, Microwave Cookbook for Microscopits, Coulomb Press Leydon Leidon 1992

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Referenser

- Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2:a utg, CV Mosby Co. St. Louis (MO), 1980
- Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
- Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986

- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
- Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
- Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

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Referências

1. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
2. Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
3. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
4. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
5. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
6. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
7. Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
8. Kok and Boon, Microwave Cookbook for Microscopits, Coulomb Press Leydon Leidon 1992

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- 1. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
- 2. Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
- 3. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
- 4. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- 5. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- 6. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
- 7. Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
- 8. Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

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Reference

- Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
- Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J. Histotechnol 4:115, 1981
- Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
- Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
- Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

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Referanslar

- Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
- Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981

- Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
- Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
- Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leiden 1992

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Lähteet

1. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
2. Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
3. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
4. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
5. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
6. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
7. Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
8. Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

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Referenties

1. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2e editie, CV Mosby Co. St. Louis (MO), 1980.
2. Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
3. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
4. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
5. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
6. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
7. Carson Frieda, Histotechnologie: A Self Instructional Text, ASCP Press, Chicago 1990.
8. Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leiden 1992

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Reference

- Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
- Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
- Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
- Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
- Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

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Reference

1. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
2. Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
3. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
4. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986.
5. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983.
6. Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
7. Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
8. Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

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Atsauces

- Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
- Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
- Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986

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- Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
- Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
- Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

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Odkazy

- Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd ed., CV Mosby Co. St. Louis (MO), 1980
- Koski JP: Silver methenamine-borate (SMB): Cost reduction with technical improvement in silver nitrate-gold chloride impregnations. J Histotechnol 4:115, 1981
- Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Churukian CJ, Schenk EA, Clark G: Dilute ammoniacal silver as a substitute for methenamine silver to demonstrate Pneumocystis carinii and fungi. Lab Med 17:87, 1986
- Carson Frieda, Histotechnology: A Self Instructional Text, ASCP Press, Chicago 1990
- Kok and Boon, Microwave Cookbook for Microscopists, Coulomb Press Leydon Leidon 1992

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Ak chcete urobiť objednávku, navštívte našu webovú lokalitu na [SigmaAldrich.com](https://sigmaaldrich.com). Informácie o technickom servise nájdete na stránke technického servisu na našej webovej lokalite na [SigmaAldrich.com/techservice](https://sigmaaldrich.com/techservice).

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