


MicropreSure™

In-Line Filtration Sampler



User Guide

**For research use only.
Not for use in clinical applications.
Single use only.**


MSHA WGT 48 (48/pk)
MSHA WGS 48 (48/pk)

MILLIPORE

Introduction

Millipore’s MicropreSure In-Line Filtration Sampler is used for in-line sampling and filtration of liquids within a closed system, minimizing the risk of false positives from adventitious contamination while sampling in a production environment. Large volumes can be filtered at once, or slowly over an extended period of time. For instance, filtration time can be adjusted to the duration of a work shift or a batch. Designed for use with pressurized process lines or reservoirs sampling, the MicropreSure device handles pressures up to 3 bar. Liquids enter the MicropreSure device and pass through the microporous membrane within its filtration chamber. After the desired sample volume has been filtered, the device’s inlets and outlets can be sealed for transportation back to the lab without risk of contamination. The membrane is accessed simply by removing the MicropreSure dome and is transferred to media for incubation and microbial enumeration.

Diagram of MicropreSure Device Components



Usage Guidelines

- Do not use with pressures exceeding the recommended 3 bar limit at room temperature.
- Do not use with temperatures exceeding 90°C.
- Flush the sampling port thoroughly before filtering samples by opening the valve then closing the valve.
- Maintain the device in a horizontal position while sampling to ensure homogeneous colony distribution.
- Filtration time and throughput will vary depending upon sample pressure and fouling index.
- Remove excess liquid remaining in the MicropreSure device after sampling within 2 hours to avoid compromising test results.

Required Equipment

- MicropreSure device
- MSOpener™ manifold
- 6 mm I.D. tubing
- Sampling port with male Luer outlet (preferably horizontal)
- Media plate
- Graduated container (to measure sampled volume)
- Laboratory pen
- Forceps
- Incubator
- One sterile Petri dish

Identifying the Sampling Port Outlet and Choosing the Valve Adapter

- If the sampling port outlet is a male Luer port (such as an ESP or Milliflex™-P valve), use the MicropreSure device without a sampling port adapter.
- If the sampling port has a Tri-Clover® or an NPT connection valve, install a MicropreSure sanitary sampling valve. Order the correct valve as follows:

Connection size	Millipore Cat. No.
TC 1 1/2 "	MSES PTC L1
NPT 1/4"	MSES P14 N1
NPT 1/8"	MSES P18 N1

Either leave the MicropreSure valve in place or remove the valve after each sampling. The MicropreSure valve is autoclavable.

- If the sampling valve outlet is simple tubing, measure the tubing's outer diameter and then order the correct valve adapter as follows:

Outer Diameter	Millipore Cat. No.	Outer Diameter	Millipore Cat. No.
3 mm	MS0S T03 03	20 mm	MS0S T20 03
6 mm	MS0S T06 03	3/16"	MSST 316 03
8 mm	MS0S T08 03	1/4"	MS0S T14 03
12 mm	MS0S T12 03	3/8"	MS0S T38 03

- If the sampling valve outlet is a male threaded connection, check the format and measure the tubing's outer diameter and then order the correct valve adapter as follows:

Outer Diameter	Millipore Cat. No.
1/8" Gas	MSFM 18G 03
1/4" Gas	MSFM 14G 03
1/8" NPT	MSFM 18N 03
1/4" NPT	MSFM 14N 03

- If the sampling valve outlet is a female threaded connection, check the format and measure the tubing's outer diameter and then order the correct valve adapter as follows:

Outer Diameter	Millipore Cat. No.
1/8" Gas	MSLM 18G 03
1/4" Gas	MSLM 14G 03
1/8" NPT	MSLM 18N 03
1/4" NPT	MSLM 14N 03

- If the sampling valve outlet is a Tri-Clover connection, measure the valve outlet's outer diameter and then order the correct adapter as follows:

Outer Diameter	Millipore Cat. No.
TC 3/4" (25 mm)	MS0T C34 03
TC Mini Clamp (34 mm)	MS0T CMC 03
TC 1 1/2" (50.5 mm)	MSTC 120 03

Sampling Procedure


This section outlines the procedures for filtering samples on site, accessing the membrane, and transferring the membrane to agar in the laboratory.

Filtering the Sample

1. Open the sanitary sampling valve and flush thoroughly, then close the valve.
2. Remove the blue stopper from the MicropreSure inlet and attach to the base as shown below.


NOTE: Only one of the 2 dome connectors is an inlet. The device inlet is identified by the blue stopper and by the INLET mark.

3. Remove the yellow stopper from the MicropreSure outlet on the bottom of the device and attach to the base as shown below.



Filtering the Sample, continued

4. Fit the MicropreSure outlet to the receiving graduated container. If the MicropreSure device does not fit directly on the receiving container, connect a length of tubing to the device outlet. Use either tubing with a male Luer connection to fit into the device outlet, or connect tubing with an internal diameter of approx. 6mm (1/4 inch) around the device outlet.
5. Insert the sampling port outlet into the MicropreSure inlet.



6. Hold the MicropreSure device in position on the valve without tightening the connection. Open the sampling valve until liquid covers the entire filter surface (dome half full) and the desired flow is obtained.

NOTE: This step ensures that the surface of the membrane is fully wetted. If the valve is installed horizontally so that the MicropreSure unit is oriented parallel to the floor and the liquid flows evenly over the membrane surface, skip this step.


7. Secure the MicropreSure device in position by firmly pushing it onto the valve (with a slightly rotating movement).
8. Close the sampling valve when the desired sample volume or the desired sampling time has been reached.
9. Disconnect the MicropreSure device from the sampling valve by gently pulling it away from the valve, using a rotating movement. Maintain the device in a horizontal position.
10. Place the yellow stopper over the device outlet and use the blue stopper to seal the device inlet.
11. Record the sample location, time, and volume on the side of the MicropreSure dome.

Accessing the Membrane

The liquid remaining in the MicropreSure device will not compromise the test results if it remains in the device for a short period of time (2 hours). Therefore, the device containing the remaining liquid may be carried to the laboratory and purged just before transferring the membrane onto the agar media.


Use the following procedure to remove the excess liquid in the MicropreSure device and to gain access to the membrane. If company procedures require it, the MSOpener manifold can be sterilized by autoclaving before entering the working area.

1. Connect the MSOpener manifold to a vacuum source (from 15 psi to 25 psi).
2. Remove and discard both stoppers (inlet and outlet) from the MicropreSure device.
3. Place the device on top of the MSOpener.
4. Turn the grooved ring to position **I**.
5. After purging, press down on the dome to separate it from the base.




Accessing the Membrane, *continued*

6. Remove the dome.
7. Lift the membrane using sterilized forceps and transfer it onto an agar plate.



8. Twist the MicropreSure base to remove it from the MSopener manifold. Return the grooved ring to position ● in preparation for next use.



9. Incubate agar plate upside down as with any standard 47 mm membrane.

Specifications

Device height:	45 mm (1.75 in)
Filter diameter:	47 mm
Filter area:	9.10 cm ²
Filter type:	0.45 µm pore size, white, gridded
Pressure limit at 25°C:	3 bar
Pressure limit at 90°C:	1 bar, 1 minute, 1 liter
Temperature range:	5 °C–90 °C

Materials of Construction

Filtration chamber/base:	Polycarbonate
Membrane:	Mixed ester cellulose, 0.45 µm, white, gridded
Filtration support:	Tenax® mesh spacer

Technical Assistance

For more information, contact the Millipore office nearest you. In the U.S., call **1-800-MILLIPORE** (1-800-645-5476). Outside the U.S., see your Millipore laboratory catalogue for the phone number of the office nearest you. You can reach us by e-mail at tech_service@millipore.com or visit our web site (www.millipore.com).

Product Ordering Information

This section lists the catalogue numbers for MicropreSure filtration devices and accessories. See the Technical Assistance section for information about contacting Millipore. You can also buy Millipore products on-line at www.millipore.com/purecommerce.

Description	Quantity	Catalogue No.
MicropreSure device, EO sterilized, W filter transfer onto agar media, 0,45µm white gridded	48/pk	MSHA WGT 48
MicropreSure sterile device, indiv. pkd W filter transfer onto agar plate, sterile, 0,45µm white gridded	48/pk	MSHA WGS 48
MSopener manifold, complete vacuum manifold for removing excess liquid and accessing membrane	1/pk	MS0P ENR 01
MSopener tulip for mounting on 3 or 6 place manifold	1/pk	MS0P NMS 01
Maintenance kit for MSOpener manifold	1/pk	MS0P KIT 01

Standard Warranty

Millipore Corporation ("Millipore") warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of Millipore products appearing in Millipore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Millipore's sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies Millipore promptly of any such breach. If after exercising reasonable efforts, Millipore is unable to repair or replace the product or part, then Millipore shall refund to the customer all monies paid for such applicable product or part. **MILLIPORE SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER INDIRECT DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

MILLIPORE

©2001 Millipore Corporation. All rights reserved. This publication or parts thereof may not be reproduced in any form without the written permission of the publishers.

Millipore is a registered trademark of the Millipore Corporation.
MicropreSure, Milliflex, and MSOpener are trademarks of the Millipore Corporation
Tenax is a registered trademark of TENAX SpA Italy.
Tri-Clover is a registered trademark of Ladish Co. Inc.

P36263 Rev. -, 08/01