English

Millex®-MP (25 mm) Sterile Filter Unit

SLMP0255S (0.22 µm, 50/box)	Male Luer-slip outlet
SLMPL255S (0.22 µm, 50/box)	Male Luer-Lok™ outlet
● Single use only	● Sterile
● Non-pyrogenic	

Indications for Use/Purpose

Intended for use as a syringe filter to sterilize and clarify low volume solutions in direct patient care and pharmacy admixture applications.

Cautions: Federal (USA) law restricts this device to sale by or on the order of a physician or pharmacist.

Introduction

This document provides compatibility information, operating steps, and specifications for the Millex®-MP sterile filter unit. The Millex®-MP unit removes microorganisms, particles, precipitates, and undissolved powders larger than 0.22 micron (µm) from aqueous solutions. This single-use product consists of a modified polyethersulfone membrane filter sealed in a polyvinyl chloride (PVC) housing. It is non-pyrogenic and non-toxic.

Applications

Millex®-MP syringe filters are suitable for use in hospital pharmacy and direct patient care applications. In the hospital pharmacy, Millex®-MP filter units can be used to sterile filter small volumes of protein pharmaceuticals, diagnostic imaging agents, chemotherapeutics, aqueous solutions, or water during admixture preparation. Direct patient care applications include sterilization and particulate removal from epidural and other liquid anesthetics as well as from irrigation solutions used in ophthalmic, etc, and other surgical procedures.

Chemical Compatibility

The Millex®-MP filter unit is compatible with most aqueous solutions. Based on information from technical publications, materials suppliers, and laboratory tests, Merck Millipore Ltd. believes that the agents listed in the following chart are safe to use with Millex®-MP filter units. However, because of the effects of variability in temperature, concentrations, duration of exposure, and other factors outside of our control, Merck Millipore Ltd. does not provide or imply a warranty with respect to this information. For example, the drug list refers to the compatibility between the drug and filter materials. It does not refer to the specific drug or protein binding to the filter, or the potential loss of the active drug component during filtration. You should qualify specific drugs for protein binding or loss of drug component prior to use. Agents that are not listed should be tested prior to use.

Chemicals

Acetic acid (5%)	Clorox® bleach (5% solution)	Water
Aqueous buffers	Sodium hydroxide (10%)	
Cell culture media	Sulfuric acid (20%)	

Active Drug Compounds

Drug	Ref. No.	Merck® Index 11th Edition Page No.	Comments	Drug	Ref. No.	Merck® Index 11th Edition Page No.	Comments
Aminophylline	477	76	water soluble	Hydrocortisone	4713	758	water soluble
Ampicillin	621	93	water soluble	21-glycol sodium succinate			
Aspartame	861	132	water soluble	Immunoglobulins	4837	780	water soluble
Biomeprins	1324	201	water soluble	Insulin	4887	789	water soluble
Caffeine	1635	248	water soluble	Isoproterenol	5105	821	water soluble
Cefazolin	1935	294	water soluble	Lidocaine	5359	863	water soluble
Cefoxitin	1938	297	water soluble	Mannitol	5629	901	water soluble
Cephalexin	1978	305	water soluble	Metronidazole	6079	968	water soluble
Cisplatin	2319	361	water soluble	Mitoxantrone	6131	978	water soluble
Colistin	2475	387	suspension + surfactant	Mitomycin	6133	979	water soluble
Cytarabine	2790	437	water soluble	Mitoxantrone	6135	979	water soluble
Dactinomycin	2804	441	water soluble	Moxalactam	6201	991	water soluble
Doxorubicin	2825	445	water soluble	Nitroglycerin	6528	1045	water soluble
Dexamethasone	2922	463	5% alcohol	Norepinephrine	6612	1058	water soluble
Diazepam	2977	472	40% alcohol	Penicillin G	7041	1123	water soluble
Digoxin	3150	499	50% alcohol	Potassium phenobarbital	7201	1149	water soluble
Dobutamine	3396	535	water soluble	Phenytoin	7293	1162	50% alcohol
Dopamine	3415	538	water soluble	Piperacillin	7430	1184	water soluble
Doxorubicin	3428	540	water soluble	Plicamycin	7510	1198	water soluble
Ergonovine	3600	573	water soluble	Prochlorperazine	7768	1231	water soluble
Etoposide	3842	610	30% alcohol	Protamines	7898	1253	water soluble
Factor III	3873	616	water soluble	Streptokinase	8784	1390	water soluble
Factor IX	3874	616	water soluble	Sulfamethoxazole	8888	1407	50% alcohol
Fentanyl	3944	628	water soluble	Tobramycin	9413	1494	water soluble
Fluorouracil	4109	654	water soluble	Trimethoprim	9624	1528	water soluble
Folic acid	4140	680	water soluble	Urokinase	9799	1555	water soluble
Furosemide	4221	674	water soluble	Vidarabine	9881	1569	water soluble
Gentamicin	4284	686	water soluble	Vinblastine	9887	1570	water soluble
Hemin	4563	733	water soluble	Vincristine	9891	1571	water soluble
Heparin	4571	735	water soluble				

How to Use the Millex® Sterile Filter Unit

This section lists warnings and cautions and provides steps to use the Millex® sterile filter unit.


WARNINGS:

- To ensure sterility, do not use this product if the package is damaged.
- Do not use this product as an in-line filter for intravenous fluid administration; it was not designed for long-term continuous use.
- Do not use with syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter unit and/or personal injury.

CAUTIONS:

- Do not use the Millex® filter unit to filter fluids at temperatures above 45 °C (113 °F).
- Do not use the Millex® filter unit to filter emulsions or suspensions because it was not designed for that purpose.
- Do not use the Millex® filter unit to filter solutions containing 5 milligrams (mg) or less of active drug unless binds binding studies have been performed.
- Do not re-sterilize or reuse the Millex® filter unit, as Merck Millipore Ltd. cannot assure the sterility, integrity, and performance beyond a single use.

Procedure for Using the Millex® Sterile Filter Unit

Specifications for Millex®-MP (25 mm) Filter Unit with Millipore Express® Membrane

Materials	Hydrophilic Millipore Express® polyethersulfone
Membrane	Pore size: 0.22 µm
Housing	Polyvinyl chloride (PVC)
Dimensions	
Inlet to outlet	SLMP0255S: 25 mm (0.98 in.)
	SLMPL255S: 23 mm (0.91 in.)
Diameter	29 mm (1.14 in.)
Filtration area	4 cm² (0.62 in²)
Temperature limit	45 °C (113 °F) maximum
Housing Pressure at 25 °C	5.2 bar (75 psi) inlet and differential maximum
Filtration volume	1 mL to 100 mL
Hold-up volume	≤ 0.1 mL after air purge
Sterilization method	Ethylene oxide gas
Connections	SLMP0255S: Female Luer-Lok™ inlet; male Luer-slip outlet
	SLMPL255S: Female Luer-Lok™ inlet; male Luer-Lok™ outlet
Flow rate at 2.1 bar (30 psi), 25 °C	≥ 175 mL/min

Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Merck Millipore Ltd. ("Millipore") or an affiliate. Neither Merck Millipore Ltd. nor any of its affiliates assumes responsibility for any errors that may appear in this document.

Technical Assistance

For more information, contact the office nearest you. In the U.S., call 1-800-MILLIPORE (1-800-645-5476). Outside the U.S., go to our web site at www.millipore.com/offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at www.millipore.com/techservice.

Standard Product Warranty

The applicable warranty for the products listed in this publication may be found at www.millipore.com/terms (within the "Terms and Conditions of Sale" applicable to your purchase transaction).

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© 2013 EMD Millipore Corporation. All rights reserved.
Millipore, Millex, and Millipore Express are registered trademarks of Merck KGaA, Darmstadt, Germany.
The M logo is a trademark of Merck KGaA.
All trademarks of third parties are the property of their respective owners.