

## Product Information

### Detergent Screening Kit

Product Number **66317**

Store at 4°C

## TECHNICAL BULLETIN

### Application

Similar to the additive kit this kit provides a range of detergents, that may be used for optimization. Detergents often interfere with the crystallization process. They are used to influence hydrophobic sample-sample interactions which may cause non-specific aggregation. On the other hand they can also influence the structure of water. The use of detergents is typically required for crystallization of membrane proteins or other proteins, if non-specific aggregation is suspected.

The detergent kit includes 48 defined detergents, that are pre-formulated for ease of use. The concentration of each detergent solution is adjusted to the 10-fold value of critical micelle concentration. The kit includes 100 µl of each pre-formulated detergent.

### Quality of reagents – a key to success

The reagents within this kit are formulated using the highest purity reagents available. All solutions are sterile filtered using 0.22 micron filters. All solutions are available separately as 1 ml bottles. Larger quantities are available on request.

### Precautions and Disclaimer

This product is for laboratory research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

### Storage/Stability

It is recommended that the reagents of this kit be stored at 4 °C. Storage at –20 °C will not adversely affect the kit reagents and the reagents as supplied are stable at room temperature for short-term storage. Kit reagents should not be set under ultraviolet light to protect them from microorganisms.

### Handling and Procedures

The application method described below is for the sitting drop method with 10 µl drop size (e.g. in Corning CrystalEX™ plates, our product no Z709409). For smaller drop sizes (e.g. 2 or 4 µl in Greiner CrystalQuick™ plates, our product no. Z666130 and Z617644) smaller volumes with same ratios are used. The detergents can be used for hanging drop and microbatch methods in similar manner.

1. Pipet crystallization reagent into the reservoir.
2. Pipet 5 microliters of sample onto a sitting drop post (50% of drop volume).
3. Pipet 1 microliter of Detergent 1 into the sample drop (10% of drop volume).

4. Pipet 4 microliters of the crystallization reagent into the sample/detergent drop. (40% of drop volume) When mixing use caution to avoid bubbles and to minimize spreading of the drop. Some drop spreading will occur due to the nature of the detergent.

5. Seal the reservoir. Repeat for next detergent. The sealing membrane can be pierced by standard needles used in automatic liquid handling equipment. If you need to remove the sealing membrane, first check for drops hanging at the bottom side of the membrane. If necessary, centrifuge the plate or knock it on the table to get those drops down into the reservoir.

### **Observation**

Drops are typically observed by a stereo microscope at 10 to 100X. Record all observation by scanning every droplet on the slides.

Scan the focal plane for small crystals and record observations for all droplets. Scan the first time shortly after the screen is set up. Then for the first 5-10 days, information may be recorded daily and, thereafter, on a weekly basis. Records should follow the same scheme as without use of detergents.

### **References**

1. Michel, H. (Ed.) Crystallization of membrane proteins, CRC Press, 1991.
2. Ducruix, A., Giege, R. (Eds.), Crystallization of nucleic acids and proteins. The Practical Approach Series. Oxford Univ. Press, (1992).
3. Cudney, B. et al, Screening and optimization strategies for macromolecular crystal growth, Acta Cryst, D50, 414-423, (1994).

## 66317 Detergent Screening Kit Observation Sheet

Sample: \_\_\_\_\_

Date: \_\_\_\_\_

Sample Concentration: \_\_\_\_\_

Reservoir Volume: \_\_\_\_\_

Sample Puffer: \_\_\_\_\_

Temperature: \_\_\_\_\_

Reservoir Composition: \_\_\_\_\_

Drop Volume: Total \_\_\_\_\_ ul Sample \_\_\_\_\_ ul Reservoir \_\_\_\_\_ ul Detergent \_\_\_\_\_ ul

1 Clear Drop

4 Birefringent Precipitate or Microcrystals

7 Plates (2D Growth)

2 Phase Separation

5 Posettes or Spherulites

8 Single Crystals (3D Growth < 0.2mm)

3 Regular Granular Precipitate

6 Needles (1D Growth)

9 Single Crystals (3D Growth > 0.2mm)

No.	Product No.	Name	Date:	Date:	Date:	Date:	Date:	Date:
1.	50006	Hexyl-beta-D-glucoside 200 mM						
2.	44615	Heptyl-beta-D-glucoside 150 mM						
3.	18842	Octyl-beta-D-glucoside 100 mM						
4.	30887	Octyl-alpha/beta-glucoside 100 mM						
5.	11448	Nonyl-beta-D-glucoside 50 mM						
6.	04978	Methyl-6-O-(N-heptylcarbonyl)-alpha-D-glucopyranoside (HECAMEG) 200 mM						
7.	42387	Triethyleneglycol monoocyl ether 1 mM						
8.	41957	Triethyleneglycol monodecyl ether 1 mM						
9.	07891	Triethyleneglycol monododecyl ether 1 mM						
10.	13172	Tetraethyleneglycolmonooctylether 1 mM						
11.	07888	Tetraethyleneglycol monodecyl ether 1 mM						
12.	16802	Tetraethyleneglycol monododecyl ether 1 mM						
13.	28000	Pentaethyleneglycol monodecyl ether 1 mM						
14.	04556	Pentaethyleneglycol monododecyl ether 1 mM						
15.	03269	Hexaethyleneglycol monodecyl ether 1 mM						
16.	15931	Hexaethyleneglycol monododecyl ether 1 mM						
17.	39376	Heptaethyleneglycol monodecyl ether 1 mM						
18.	49534	Heptaethyleneglycol monododecyl ether 1 mM						
19.	07511	Octaethyleneglycol monodecyl ether 1 mM						
20.	55377	Octaethyleneglycol monododecyl ether 1 mM						
21.	40913	Nonaethyleneglycol monododecyl ether 1 mM						
22.	07509	Decyl-beta-D-maltoside 10 mM						
23.	04836	Undecyl-beta-D-maltoside 10 mM						
24.	54202	Dodecyl-beta-D-maltoside 3 mM						
25.	08386	Hexadecyl-beta-D-maltoside 0.1 mM						
26.	54280	Cyclohexylmethyl-beta-D-maltoside 500 mM						
27.	19164	2-Cyclohexylethyl-beta-D-maltoside 200 mM						
28.	41583	6-Cyclohexylpentyl-beta-D-maltoside 30 mM						
29.	51793	6-Cyclohexylhexyl-beta-D-maltoside 30 mM						
30.	39489	Heptyl-beta-D-thiogluconate 300 mM						
31.	07818	Nonyl-beta-D-1-thiomaltopyranoside 30 mM						
32.	38498	Decyl-beta-D-1-thiomaltopyranoside 10 mM						
33.	19899	CHAPS 100 mM						
34.	39490	MEGA-8 500 mM						
35.	44882	MEGA-9 200 mM						
36.	43104	N,N-Dimethylhexylamine N-oxide 100 mM						
37.	41261	N,N-Dimethyloctylamine N-oxide 100 mM						
38.	49556	N,N-Dimethylnonylamine N-oxide 100 mM						
39.	18705	N,N-Dimethyldodecylamine N-oxide 100 mM						
40.	57915	N,N-Dimethyldodecylamine N-oxide 100 mM						
41.	54185	3-(N,N-Dimethyloctylammonio)propanesulfonate 500 mM						
42.	54196	3-(N,N-Dimethyldodecylammonio)propanesulfonate 300 mM						
43.	16974	3-(N,N-Dimethyldodecylammonio)propanesulfonate 30 mM						
44.	53502	3-(N,N-Dimethylmyristoylammonio)propanesulfonate 3 mM						
45.	44490	3-(N,N-Dimethylpalmitoylammonio)propanesulfonate 0.3 mM						
46.	38544	O-(n-Octyl)-phosphorylcholine 1 M						
47.	55461	O-(n-Decyl)-phosphorylcholine 100 mM						
48.	16448	O-(n-Tetradecyl)-phosphorylcholine 1.5 M						

**Precautions and Disclaimer:**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

©2015 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at [www.sigmaaldrich.com](http://www.sigmaaldrich.com) and/or on the reverse side of the invoice or packing slip.