

Application Note

HumanKine[®] Thermostable bFGF, an engineered recombinant protein with enhanced biological functionality on human iPS and neural progenitor cells

Introduction

Fibroblast growth factors (FGFs) are secreted glycoproteins that regulate several fundamental developmental pathways and help regulate mesoderm and ectoderm patterning in the early embryonic development. In adults, FGFs participate in several physiological and pathological processes including cell proliferation, cell survival, differentiation, wound healing, and regulation of angiogenesis. The FGF signaling process involves their binding to FGF receptors (FGFR) on the cell surface. Twenty-three different FGFs have been identified that exert their effects through four major types of tyrosine kinase-linked FGFRs that are designated as FGFR1, FGFR2, FGFR3, and FGFR4. Following ligand binding, these receptors undergo dimerization, which results in the activation of their tyrosine kinase domains and intermolecular transphosphorylation.

Basic FGF, also known as bFGF or FGF-2, belongs to a group of structurally homologous FGFs that are characterized by heparin binding and mitogenic activity on a variety of cells of mesodermal and neuroectodermal origin. The use of bFGF in stem cell biology is widespread and is a critical growth factor that enables the expansion of many stem cells in their pluripotent

or multipotent states. Several recent studies suggest an important role for FGF/extracellular signal-regulated kinase (Erk) signaling in promoting the transition from a naïve to a primed differentiation state and in preventing primed cells from reverting back to a naïve state. FGF signaling, in effect, stabilizes the primed cell state of human embryonic stem cells.

Due to its rapid degradation at 37 °C, many stem cell researchers add high levels of bFGF to their media and require cumbersome frequent media exchanges. The HumanKine[®] Thermostable bFGF is an engineered recombinant version of human bFGF. It is generated by selectively altering the amino acid sequence to enhance its thermostability compared to wild-type bFGF. This humanized recombinant protein is produced using a proprietary human cell expression system and is designated as a xeno-free product. We have shown that this unique protein has enhanced stability at 37 °C while possessing increased biological functionality on both human iPS cells and human neural stem cell populations. The use of this product will reduce the variation commonly observed within many stem cell cultures and will reduce the time and cost needed to culture these demanding cells.




Figure 1.
 FGF signaling pathways lead to cell survival, proliferation, invasion and migration

Materials and Methods

Determining stability of bFGF cytokines at elevated temperatures

The stability of normal and HumanKine® Thermostable bFGF recombinant proteins was analyzed at 37 °C over a period of 7 days. Normal bFGF (Merck Millipore Cat. No. GF003) or HumanKine® Thermostable bFGF (Merck Millipore Cat. No. GF446) was diluted in 20% KOSR medium to a final concentration of 150 ng/mL. Samples were allowed to remain at 37 °C over a period of one week with stopping at different time points during this period. The level of remaining bFGF was quantified at each time point using a quantitative bFGF ELISA kit (Biologend Cat. No. 434307).

Culture of human iPS cells

Human iPS cells were generated from neonatal human foreskin fibroblasts (Merck Millipore Cat. No. SCC058) using STEMCCA™ lentivirus reprogramming kit (Merck Millipore Cat. No. SCR545) followed by excision of transgenes using a cell-permeable Cre-Recombinase (Merck Millipore Cat. No. SCR508). Transgene-free iPS cells were maintained in KOSR medium with standard bFGF (10 ng/mL) on primary mouse embryonic fibroblasts (Merck Millipore Cat. No. PMEF-CF) for

20 passages using a manual passaging technique until used for experiments. At passage 21, Human iPS cells were manually passed into either KOSR medium with no FGF, KOSR + Normal bFGF (Merck Millipore Cat. No. GF003, 10 ng/mL) or KOSR medium + Thermostable bFGF (Merck Millipore Cat. No. GF446, 10 ng/mL) conditions in duplicate at approximately 50 colonies/well with feeder cells present. Cells were fed every third day for one passage under above conditions and were fixed and exposed to anti-Oct-4 (Merck Millipore Cat. No. MAB4401) and anti-TRA-1-60 (Merck Millipore Cat. No. MAB4360) for immunostaining. The total number of differentiated colonies were manually counted on day 7.

Expansion of ENStem™-A human neural progenitors

ENStem™-A cells (Merck Millipore Cat. No. SCR005) at passage 20 were thawed onto Matrigel® matrix-coated plates (1:50) in ENStem™-A neural expansion medium (Cat. No. SCM004) with standard bFGF (20 ng/mL) or HumanKine® Thermostable bFGF (20 ng/mL). Cells were passaged using Accutase™ reagent (Cat. No. SCR005) and, when 80% confluency was reached for three passages, they were immunostained with neural stem cell markers anti-nestin (Cat. No. MAB5326) and anti-Sox-2 (Cat. No. MAB4423).

To assess proliferation, 5×10^4 cells were seeded onto Matrigel® matrix-coated (1:50) 6-well plates in duplicate in either ENStem™-A neural expansion medium + standard bFGF (20 ng/mL) or ENStem™-A neural expansion medium + Thermostable bFGF (20 ng/mL). Cells were allowed to grow for 6 days without any media exchanges. Cells were detached and triplicate samples were mixed with Trypan blue and counted using a hemocytometer at days 2, 4, and 6.

Results

Thermostability Experiments

We studied the thermal stability of normal bFGF and modified, thermostable bFGF at 37 °C over a period of seven days (Figure 2). The reduction in normal bFGF levels was more pronounced (nearly 95% reduction on day 1) compared to the thermostable bFGF (only ~60% reduction on day 1). Both bFGFs showed gradual decrease over a period of 7 days; however, at a slower rate of decline with the thermostable bFGF.




Figure 2. Enhanced stability of thermostable bFGF at 37 °C vs. normal bFGF.

Human iPS cell experiments

When we used a limited feeding schedule to culture human iPS cells and used thermostable bFGF, we observed elevated levels of the pluripotency markers Oct-4 and TRA-1-60 after 7 days and significant reduction in spontaneous differentiation of human iPS cells compared to when we used normal bFGF in the same feeding schedule (Figures 3 and 4). In fact, we found no significant reduction in the number of spontaneously differentiated colonies when using normal bFGF compared to KOSR medium alone (Figure 4).

ENStem™-A human neural stem cell experiments

To determine if thermostable bFGF could also facilitate the culture of neural progenitor cells, we first cultured them with thermostable bFGF using a normal feeding schedule (Figure 5). Immunocytochemistry and bright-field microscopy showed that the cells, after three passages, displayed proper morphology and expression of neural stem cell markers Sox-2 and nestin.




Figure 3. Thermostable bFGF supports pluripotent growth of human iPS cells using an every-third-day feeding regimen. Human iPS cells grown in HumanKine® Thermostable bFGF following a limited feeding regimen express elevated levels of the pluripotency markers Oct-4 (Green) and TRA-1-60 (Red) after 7 days. Human iPS cells grown with KOSR medium with normal bFGF following a limited feeding regimen spontaneously differentiate within the center of the colonies and do not express the same uniform levels of Oct-4 and TRA-1-60 expression.




Figure 4. Human iPS cells show lower levels of spontaneous differentiation using an every-third-day feeding regimen when used in combination with HumanKine® Thermostable bFGF.



Figure 5. ENStem™-A neural progenitor cells grown for 3 passages in HumanKine® Thermostable bFGF following a normal feeding regimen display proper morphology and express the neural stem cell markers Sox-2 and Nestin. Nestin-stained cells were counterstained with DAPI (blue) nuclear stain.

Next, we measured proliferation of neural progenitor cells while keeping them on a minimal feeding schedule, in which cells were incubated for six days without any intervening medium changes. After six days, there were slightly more cells per well in the cultures containing thermostable bFGF as compared to normal bFGF (Figure 6).





Figure 6. Higher cell proliferation rates of ENStem™-A human neural progenitor cells using a limited feeding regimen when using thermostable bFGF compared to normal bFGF. (A) Time course showing cell counts on days 2, 4 and 6 of incubation. (B) Comparing cell number per well after six days of incubation without medium exchange. Experiments were performed using three replicates of each time point and culture condition. Error bars reflect standard deviation.

Discussion

Fibroblast growth factor is known to play a key role in maintaining embryonic stem cells and iPS cells in a proliferative, undifferentiated state. Basic FGF is generally added to the culture medium to sustain ERK phosphorylation and Nanog expression¹. bFGF levels required to support pluripotency over a period of time are about 10-fold higher than those required for other cell types². Chen et al have shown that acidic FGF is highly unstable at 37 °C and fails to maintain ERK phosphorylation and pluripotency in stem cells¹. A similar phenomenon was suggested for bFGF. However, their studies have indicated that, at 37 °C, bFGF undergoes aggregation, thereby reducing its biological availability for maintaining pluripotency.

Consistent with these indications, our data suggest that HumanKine® Thermostable bFGF can be used with human ESC, iPS and neural stem cells to maintain them in a pluripotent and multipotent state as a more effective alternative to normal bFGF. Furthermore, because of its increased stability, thermostable bFGF permits a less stringent feeding schedule, increasing the experimental flexibility of these cell models.

References

1. Chen G, Gulbranson DR, Yu P, Hou Z, Thomson JA. Thermal stability of fibroblast growth factor protein is a determinant factor in regulating self-renewal, differentiation, and reprogramming in human pluripotent stem cells. *Stem Cells*. 2012 Apr;30(4):623-30.
2. Ludwig TE, Levenstein ME, Jones JM, Berggren WT, Mitchen ER, Frane JL, Crandall LJ, Daigh CA, Conard KR, Piekarczyk MS, Llanas RA, Thomson JA. Derivation of human embryonic stem cells in defined conditions. *Nat Biotechnol*. 2006 Feb;24(2):185-7.

Related Products

Description	Catalogue No.
HumanKine® Thermostable bFGF, Human Recombinant	GF446-10UG GF446-50UG GF446-MG
Fibroblast Growth Factor basic, human recombinant	GF003
Fibroblast Growth Factor basic, human Animal-Free recombinant	GF003-AF GF003AF-MG
Accutase™ Reagent	SCR005
FibroGRO™ Xeno-Free Human Foreskin Fibroblasts	SCC058
Human STEMCCA™ Cre-Excisable Constitutive Polycistronic (OKSM) Lentivirus Reprogramming Kit	SCR545
TAT-CRE Recombinase	SCR508
EmbryoMax® Primary Mouse Embryo Fibroblasts, Strain CF1, Mitomycin C- treated, passage 3	PMEF-CF
ENStem™ Human Neural Progenitor Expansion Kit	SCR055
ENStem™-A Neural Expansion Medium	SCM004
Anti-Oct-4 Antibody, clone 10H11.2	MAB4401
Anti-Oct-4 Antibody, clone 10H11.2, Alexa Fluor® 488 conjugate	MAB4401A4
Anti-Oct-4 Antibody, clone 10H11.2, Cy3 conjugate	MAB4401C3
Anti-TRA-1-60 Antibody, clone TRA-1-60	MAB4360
Anti-TRA-1-60 Antibody, clone TRA-1-60, Alexa Fluor® 488 Conjugate	MAB4360A4
Anti-TRA-1-60 Antibody, clone TRA-1-60, Cy3 conjugate	MAB4360C3
Anti-Nestin Antibody, clone 10C2	MAB5326
Anti-Nestin Antibody, clone 10C2, Alexa Fluor® 488 conjugate	MAB5326A4
Anti-Nestin Antibody, clone 10C2, Cy3 conjugate	MAB5326C3
Anti-SOX-2 Antibody, clone 10H9.1	MAB4423
Anti-SOX-2 Antibody, clone 10H9.1, Alexa Fluor® 488 conjugate	MAB4423A4
Anti-SOX-2 Antibody, clone 10H9.1, Cy3 conjugate	MAB4423C3

To Place an Order or Receive Technical Assistance

In Europe, please call Customer Service:

France: 0825 045 645

Germany: 01805 045 645

Italy: 848 845 645

Spain: 901 516 645 Option 1

Switzerland: 0848 645 645

United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: www.merckmillipore.com/offices

For Technical Service visit:

www.merckmillipore.com/techservice

Get Connected!

Join Merck Millipore Bioscience on your favorite social media outlet for the latest updates, news, products, innovations, and contests!

 facebook.com/MerckMilliporeBioscience

 twitter.com/Merck4Bio



www.merckmillipore.com/offices