Sigma-Aldrich®

Instructions for Use

Leukocyte Acid Phosphatase (TRAP) Kit

Procedure No. 387

Intended Use

Sigma-Aldrich Acid Phosphatase kits are intended as a general laboratory use histological stain. Acid Phosphatase reagents are for professional "In Vitro Diagnostic Use" only. This manual, qualitative, histochemical procedure demonstrates acid phosphatase activity in leukocytes of human blood, bone marrow, and tissue touch preparations. Histological visualization of leukocyte associated acid phosphatase is a unique technique currently commonly used in medicine.

Use of substituted naphthol AS phosphates in conjunction with diazonium salts for detection of acid phosphatase in human leukocytes was first reported by Goldberg and Barka.¹ The Sigma-Aldrich procedure for demonstrating this enzyme employs naphthol AS-BI phosphate and freshly diazotized fast garnet GBC salt. The latter compound was selected because it couples rapidly at acid pH, forming highly insoluble dye deposits.

Most procedures, including that provided by Sigma-Aldrich, employ stable diazonium salts. These are formed by reacting an arylamine with sodium nitrite in an acid medium.² The resulting diazonium chloride (usually unstable) can then be treated with compounds such as zinc chloride, zinc sulfate or naphthalene-1,6-disulfonate, forming stable salts. These stabilizers may exert marked inhibition upon some enzymatic systems, whereas the diazonium chlorides are less inhibitory.² For this reason, Sigma-Aldrich now provides a stable solution of fast garnet GBC base and sodium nitrite for acid phosphatase cytochemistry. To further simplify this method, a stable aqueous solution of naphthol AS-BI phosphate is also included. The availability of these stable solutions allows users to adjust working reagent volumes according to needs, eliminating waste.

According to the Sigma-Aldrich procedure, blood films are incubated in a solution containing naphthol AS-BI phosphoric acid and freshly diazotized fast garnet GBC. Duplicate films may be treated with a solution containing L(+)-tartrate. Naphthol AS-BI, released by enzymatic hydrolysis, couples immediately with fast garnet GBC forming insoluble maroon dye deposits at sites of activity. Cells containing tartaric acid-sensitive acid phosphatase are devoid of activity. Those mononuclear cells containing tartaric acid-resistant phosphatase are not affected by such treatment.

Reagents

Naphthol AS-BI Phosphoric Acid Solution (Cat. No. 3871-10ML)

Naphthol AS-BI phosphoric acid, 12.5 mg/mL

Fast Garnet GBC Base Solution (Cat. No. 3872-10ML)

Fast garnet GBC base, 7.0 mg/mL, C.I. 11160, in 0.4 mol/L hydrochloric acid with stabilizer.

Danger. Harmful if swallowed. Causes severe skin burns and eye damage. May cause an allergic skin reaction. May cause cancer.

Acetate Solution (Cat. No. 3863-50ML)

Acetate buffer, 2.5 mol/L, pH 5.2 ± 0.1

Tartrate Solution (Cat. No. 3873-10ML)

L(+)-Tartrate buffer, 0.335 mol/L, pH 4.9 ± 0.1. Warning. Causes mild skin irritation.

Sodium Nitrite Solution (Cat. No. 914-10ML)

Sodium nitrite, 0.1 mol/L

Citrate Solution (Cat. No. 915-50ML)

Citric acid, 18 mmol/L, sodium citrate, 9 mmol/L, sodium chloride, 12 mmol/L, and surfactant, pH 3.6 ± 0.1

Hematoxylin Solution, Gill No. 3 (Cat. No. GHS3-50ML)

Hematoxylin, certified, 6.0 g/L, C.I. 75290, sodium iodate, 0.6 g/L, aluminum sulfate, 52.8 g/L and stabilizers. Danger. Harmful if swallowed. Causes serious eye damage. Wear protective gloves/eye protection/face protection.

Special Materials Required but Not Provided

- Acetone, ACS Reagent
- Formaldehyde, 37%

Storage and Stability

Store Hematoxylin Solution at room temperature (18–26°C). Protect from light. Reagent label bears expiration date.

Store Naphthol AS-BI Phosphoric Acid Solution, Fast Garnet GBC Base Solution, Acetate Solution, Tartrate Solution, Sodium Nitrite Solution and Citrate Solution at 2–8°C.

Store Fixative Solution refrigerated (2–8°C). Warm to 18–26°C prior to use. Stable up to 2 months, if stored tightly capped in refrigerator. Discard Fixative Solution if evaporation is noted and prepare fresh.

Deterioration

Naphthol AS-BI Phosphoric Acid Solution, Tartrate Solution should be discarded if turbidity develops.

Citrate Solution is suitable for use in the absence of microbial growth.

Hematoxylin Solution should be discarded if solution turns brown (over-oxidized from air) or purple (loss of acidity).

Preparation

Reagents are provided ready for use.

Prepare Fixative Solution by combining 25 mL Citrate Solution, 65 mL acetone and 8 mL of 37% formaldehyde. Place in glass bottle and cap tightly.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Blood, bone marrow films, tissue touch preparations and cytocentrifuge preparations may be used. Either EDTA or heparin will serve as anticoagulants. Blood or bone marrow films may be stored fixed at room temperature (18–26°C) for several weeks or unfixed for several days without appreciable change in activity.³ Do not ship whole blood for assay at other laboratories. Send fixed or unfixed slides. Slides should be kept cool during transit. Allow films to dry at least 1 hour prior to fixation.

Notes

Test performance should be monitored by including blood film prepared from a healthy individual.

Procedure

- Prewarm sufficient deionized water for a day's use to 37°C. Check temperature before use.
- Bring Fixative Solution to room temperature (18–26°C). Fix slides by immersing in Fixative Solution for 30 seconds. Rinse thoroughly in deionized water: Do not allow slides to dry.
- To each of 2 test tubes add 0.5 mL Fast Garnet GBC Base Solution and 0.5 mL Sodium Nitrite Solution. Mix by gentle inversion for 30 seconds. Let stand 2 minutes.
- Label two 100 mL beakers A and B and add the following while mixing:

	Beaker A	Beaker B
Deionized water prewarmed to 37°C	45 mL	45 mL
Diazotized Fast Garnet GBC Solution (Step 3)	1.0 mL	1.0 mL
Naphthol AS-BI Phosphate Solution	0.5 mL	0.5 mL
Acetate Solution	2.0 mL	2.0 mL
Tartrate Solution	—	1.0 mL

- Label Coplin jars A and B and transfer solutions from beakers to appropriate Coplin jar. Warm solutions in jars to 37°C in water bath. Check that temperature is at 37°C before adding slides.
- Add slides to Coplin jars and incubate 1 hour in 37°C water bath protected from light.
- After 1 hour, rinse slides thoroughly in deionized water, then counterstain 2 minutes in Hematoxylin Solution, Gill No. 3.
- Rinse several minutes in alkaline tap water to blue nuclei.
- Air dry and evaluate microscopically. Coverslipping is not recommended since dye fades with time.

Performance Characteristics

Results

Stained blood films are usually evaluated subjectively for the presence or absence of tartrate-resistant enzyme. In the absence of tartaric acid, most leukocytes demonstrate granular sites of activity. When incubated with tartrate, an occasional granule may be observed in lymphocytes and some specialized macrophages.⁴ In blood smears a positive reaction is denoted by the presence of more than two cells with diffuse and intense activity, i.e., more than 40 granules.⁴ To evaluate Golgi staining, characteristic of thymus derived lymphocytes (T-cells), Sigma-Aldrich Procedure No. 181 is recommended.

Acid Phosphatase activity appears as purplish to dark red granules in the cytoplasm of most leukocytes.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
3871	Naphthol AS-BI Phosphoric Acid Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
3872	Fast Garnet GBC Base Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
3873	Tartrate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
3863	Acetate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
914	Sodium Nitrite Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
915	Citrate Solution	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
GHS3	Hematoxylin Solution, Gill No. 3	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

387A:

- H290 May be corrosive to metals.
- H302 Harmful if swallowed.
- H314 Causes severe skin burns and eye damage.
- H317 May cause an allergic skin reaction.
- H350 May cause cancer.
- H373 May cause damage to organs (Kidney) through prolonged or repeated exposure if swallowed.
- H402 Harmful to aquatic life.
- P201 Obtain special instructions before use.
- P202 Do not handle until all safety precautions have been read and understood.
- P234 Keep only in original container.
- P260 Do not breathe mist or vapors.
- P264 Wash skin thoroughly after handling.
- P270 Do not eat, drink or smoke when using this product.
- P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.
- P280 Wear protective gloves, protective clothing, eye protection and face protection.
- P301 + P 312 + P330 IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell. Rinse mouth.
- P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
- P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower.
- P304 + P340 + P310 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/ doctor.
- P305 + P351 + P338 + P310 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/ doctor.
- P308 + P313 If exposed or concerned: Get medical advice/ attention.
- P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.
- P362 +P364 Take off contaminated clothing and wash it before reuse
- P390 Absorb spillage to prevent material damage.
- P405 Store locked up.
- P406 Store in corrosive resistant container with a resistant inner liner.
- P501 Dispose of contents/ container to an approved waste disposal plant.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- Goldberg AF, Barka T: Acid phosphatase activity in human blood cells. Nature 189:297 1962
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88-113
- Janckila A, Li CY, et al: The Cytochemistry of the Tartrate-Resistant Acid Phosphatase – Technical considerations. Am J Clin Pathol 70:45, 1978
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 28, 119

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 3.0	2014
Rev. 4.0	2016
Rev. 5.0	2023
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Removed reference number 5. Added Warnings and Hazards. Added CH-REP information.	
Rev. 6.0	2025
Updated SDS information. Added CH Importer information.	

	Sigma-Aldrich, Inc., 3050 Spruce Street, St. Louis, MO 63103 USA an affiliate of Merck KGaA, Darmstadt, Germany +1(314) 771-5765		MDSS GmbH Schiffgraben 41 30175 Hannover, Germany		Merck KGaA, 64271 Darmstadt, Germany
	MDSS CH GmbH Laurenzenvorstadt 61 5000 Aarau Switzerland		Sigma - Aldrich Chemie GmbH Industriestraße 25 CH-9470 Buchs		

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Saure Phosphotase, Leukozyten (TRAP) – Kit

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Referenzen

- Goldberg AF, Barka T: Acid phosphatase activity in human blood cells. Nature 189:297 1962
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88–113
- Janckila A, Li CY, et al: The Cytochemistry of the Tartrate-Resistant Acid Phosphatase – Technical considerations. Am J Clin Pathol 70:45, 1978
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 28, 119

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Kit de phosphatase acide leucocytaire (TRAP)

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Références

- Goldberg AF, Barka T: Acid phosphatase activity in human blood cells. Nature 189:297 1962
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88–113
- Janckila A, Li CY, et al: The Cytochemistry of the Tartrate-Resistant Acid Phosphatase – Technical considerations. Am J Clin Pathol 70:45, 1978
- Sun T: Atlas of Cytochemistry and Immunochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 28, 119

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Referencer

- Goldberg AF, Barka T: Acid phosphatase activity in human blood cells. Nature 189:297 1962
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88–113
- Jancikla A, Li CY, et al: The Cytochemistry of the Tartrate-Resistant Acid Phosphatase – Technical considerations. Am J Clin Pathol 70:45, 1978
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 28, 119

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Referenser

- Goldberg AF, Barka T: Acid phosphatase activity in human blood cells. Nature 189:297 1962
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88–113
- Jancik A, Li CY, et al: The Cytochemistry of the Tartrate-Resistant Acid Phosphatase – Technical considerations. Am J Clin Pathol 70:45, 1978
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 28, 119

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Referanser

- Goldberg AF, Barka T: Acid phosphatase activity in human blood cells. Nature 189:297 1962
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88-113
- Jancikila A, Li CY, et al: The Cytochemistry of the Tartrate-Resistant Acid Phosphatase – Technical considerations. Am J Clin Pathol 70:45, 1978
- Sun T: Atlas of Cytochemistry and Immunochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 28, 119

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Referanslar

- Goldberg AF, Barka T: Acid phosphatase activity in human blood cells. Nature 189:297 1962
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88-113
- Janckila A, Li CY, et al: The Cytochemistry of the Tartrate-Resistant Acid Phosphatase – Technical considerations. Am J Clin Pathol 70:45, 1978
- Sun T: Atlas of Cytochemistry and Immunochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 28, 119

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Reference

- Goldberg AF, Barka T: Acid phosphatase activity in human blood cells. Nature 189:297 1962
- Burstone MS: In Enzyme Histochemistry and Its Application in the Study of Neoplasms. Academic Press, New York, 1962, pp 88–113
- Jancikila A, Li CY, et al: The Cytochemistry of the Tartrate-Resistant Acid Phosphatase – Technical considerations. Am J Clin Pathol 70:45, 1978
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, 1985, pp 28, 119

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