


LC-MS Contaminants

Avoid, identify, minimize.



The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Supelco®
Analytical Products

Table of Contents

| | |
|---|-----------|
| Introduction: the consequences of contamination | 4 |
| Sample components that can interfere with LC-MS results | 5 |
| General system care, maintenance and laboratory practice | 6 |
| Sample preparation is crucial for minimizing contamination | 7 |
| Proper mobile phase preparation to minimize contaminants | 17 |
| Effects of column choice on LC-MS performance | 24 |
| Effects of Reference Material choice on LC-MS performance | 28 |
| Useful LC-MS resources | 30 |
| Appendix I. Common mass spectrometry contaminants and their sources | 31 |
| Appendix II. Monoisotopic ion masses of commonly observed repeating units in LC-MS | 37 |

How to identify and avoid contaminants in LC-MS

(Liquid Chromatography-Mass Spectrometry)

You thought your LC-MS analysis would be straightforward, but there are peaks you didn't expect.

In this technical bulletin, you will learn some tips on identifying LC-MS contaminants and avoiding contamination. Our 110 years of separation expertise, combined with our precision-manufactured products, give you the greatest chance of obtaining reproducible, clean data.

Introduction

A successful LC-MS assay can be defined as one that has high sensitivity, reproducibility, and efficiency while also providing an increase in sample throughput. The presence of contaminants can affect these key analytical figures of merit as well as potentially compromise the performance of the analytical U/HPLC column and/or the instrument itself.

Table 1. LC-MS performance parameters.

| Parameters negatively affected by contaminants | What it is | How to minimize contaminants |
|--|--|---|
| Sensitivity | Lowest level of analyte detectable above background; sensitivity is reduced by ion suppression | Wash columns to mitigate column bleeding |
| | | Effective sample preparation |
| | | Use certified mobile phase |
| Reproducibility | Can mean column-to-column reproducibility or run-to-run reproducibility | Prepare sample with devices that do not introduce extractable impurities |
| | | Run sufficient controls to verify run-to-run stability within a batch |
| | | Remove sample components that interfere with separation, ionization and fragmentation |
| | | Use high-quality HPLC columns |
| Column Lifetime | Number of injections on a column without change in selectivity and efficiency | Use a robust column with high matrix tolerance: e.g., monolithic columns have high tolerance |
| | | Eliminate sample contaminants that adsorb strongly or ionize easily. Avoid polymers that, when fragmented, result in multiple peaks of varying m/z. Use guard column and proper sample preparation (such as centrifugation, filtration, and extraction) to remove particles and extend column lifetime. |
| | | Fully elute/clean column after each sample |

Sample components that can interfere with LC-MS results

Some biological matrices, such as plasma, contain high amounts of phospholipids. If not removed prior to chromatography, separating phospholipids from analytes of interest can require long chromatography run times and high concentrations of organic solvents. Furthermore, phospholipids can build up on analytical column, and unexpectedly elute in future runs. Drug formulation agents, such as polysorbittans and polyethylene glycol, can also interfere and cause ionization suppression.

Besides sample-derived contaminants, additional sources of contamination are sampling devices, solvent impurities, containers, sample preparation devices, volatile organics introduced as a result of handling personal care products, and even columns themselves.

Plasticizers from labware can interfere with LC-MS, resulting in the need to lengthen the chromatography run in order to resolve these peaks from analyte peaks.

A list of common contaminants, their molecular weight, and possible sources can be found in **Table 2**.

Sample components that can interfere with LC-MS results include:

- Metabolites
- Detergents
- Salts/Buffer components
- Degradation products
- Counterions
- Matrix

Table 2. List of selected contaminants observed in mass spectra (ESI, positive mode, ion mass ≤ 1000 Da). Refer to Appendix I for a more complete list, or one of the databases listed in the resources section on page 28.

| Mono-isotopic ion mass (singly charged) | Ion type | Formula for M or subunit or sequence | Compound ID or species | Possible origin and other comments |
|---|-----------------------------------|---|------------------------|--|
| 74.06059 | [M+H] ⁺ | C ₃ H ₇ NO | Dimethyl formamide | Solvent |
| 102.12827 | [M+H] ⁺ | C ₆ H ₁₅ N | TEA | Triethylamine, buffer |
| 107.0782 | [A ₂ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 123.09222 | [M+H] ⁺ | C ₇ H ₁₀ N ₂ | DMAP | Dimethylaminopyridine, solvent |
| 153.13917 | [M+H] ⁺ | C ₉ H ₁₆ N ₂ | DBU | 1,8-Diazabicyclo[5.4.0]undec-7-ene |
| 214.09018 | [M+H] ⁺ | C ₁₀ H ₁₅ NO ₂ S | n-BBS | n-butyl benzenesulfonamide, plasticizer |
| 242.28477 | M ⁺ | C ₁₆ H ₃₆ N | TBA | Tetrabutylammonium, buffer |
| 279.15964 | [M+H] ⁺ | C ₁₆ H ₂₂ O ₄ | Dibutylphthalate | Plasticizer, phthalate ester |
| 371.1018 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₅ | Polysiloxane | Polysiloxane, followed by m/z 388 |
| 371.31614 | [M+H] ⁺ | C ₂₂ H ₄₂ O ₄ | DEHA | Bis(2-ethylhexyl) adipate, plasticizer |
| 391.28484 | [M+H] ⁺ | C ₂₄ H ₃₈ O ₄ | Diisooctyl phthalate | Diisooctyl phthalate, plasticizer |
| 445.12060 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₆ | Polysiloxane | Polysiloxane, followed by m/z 462 |
| 447.2934 | [M+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 515.41341 | [M+H] ⁺ | C ₃₀ H ₅₈ O ₄ S | DDTDP | Didodecyl 3,3'-thiodipropionate, antioxidant |
| 519.13940 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₇ | Polysiloxane | Polysiloxane, followed by m/z 536 |
| 593.15820 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₈ | Polysiloxane | Polysiloxane, followed by m/z 610 |

General system care, maintenance and laboratory practice

In addition to good laboratory practices, such as wearing powder-free, nitrile gloves and monitoring laboratory air (which can contain siloxanes and phthalates), follow these tips for minimizing contamination in LC-MS.

- Flush HPLC system with organic eluent (preferably isopropanol or methanol; acetonitrile [ACN] can polymerize and block valves if system is stopped for several weeks) regularly to prevent microbial contamination. The interval of flushing depends on the eluents and buffers used and should be between two and four weeks.
- Pump debris is collected in the pump outlet filter. Some of these components can leach and be detected by MS. Replace the filter every 1–2 months or after changing from ACN to methanol (or vice versa) for lower baseline noise and general system protection.
- Filter frits attached to the inlets of the mobile phase tubing to protect the LC system from particulate matter should be made out of stainless steel. Cleaning of glass frits is time-consuming (buffer residue is hard to remove); in addition, silica and alkali are dissolved from the glass filter and form adducts $[M+X]^+$.

Sample preparation is crucial for minimizing contamination

Without sample preparation, samples contain components that are incompatible with HPLC/UHPLC/MS analyses:

- Undissolved particles/precipitates in a sample clog and reduce the life of the chromatography column.
- Sample matrices may contain many impurities, making chromatograms challenging to interpret; for example, sample matrix contains components that either elute at the same point in the LC-MS chromatogram as the analyte (potentially causing ionization suppression) or affect analyte signal intensity.
- Particles held up on the column can leach contaminants into the mobile phase (in the current sample and subsequent samples), thereby increasing background.




Figure 1. Without sample preparation, the presence of contaminants in the sample results in more peaks in the chromatogram (as represented in the left hand schematic), making analysis challenging.

Select a sample preparation method that brings the sample into a solution that is free of particles. Additional points of consideration include concentrating the analyte and reducing sample complexity. For example, a plasma sample might benefit from solid phase extraction, which removes contaminants (proteins, lipids) and also concentrates the sample, whereas fruit and vegetable juices, with their high particle load, might benefit simply from dilution and filtration.

Depending on the method chosen, sample preparation may be used, for example, to selectively enrich analytes, increase analyte concentration, or remove impurities that cause ionization suppression.

How can you tell if your LC-MS analysis is suffering from ionization suppression?

Consider performing the following steps to test for ionization suppression:

1. First, assess the detector response to a calibration standard under conditions of zero ionization suppression.
2. Spike an identical concentration of this standard into prepared sample matrix. Assess the detector response again to determine the effect of ionization suppression.
3. Assess the detector response when the spiked sample prepared in step 2 is processed using the sample preparation method(s) being considered.
4. Finally, add additional calibration standard to determine if the expected increase in signal is observed.

To mitigate the interfering effects of ionization suppression, consider performing these steps:

- Dilute sample or reduce volume injected.
- Reduce ESI flow rate to the nL/min range—this will generate smaller, highly charged droplets that can resist the effects of nonvolatile salts, in case those have not been removed from the sample. Note that you should never use nonvolatile salts in the mobile phase.
- Choose a sample preparation method that removes contaminants causing ionization suppression. Using solid phase extraction instead of protein precipitation, for example, can reduce ionization suppression by phospholipids. Phospholipid removal is discussed further below.
- Change the strength of the mobile phase or the slope of the gradient, so that the analytes of interest may elute further from the solvent front and from the end of the gradient. In these regions of the chromatogram, ionization suppression is most likely to occur.

Types of sample preparation commonly used for LC-MS

Ten of the most popular sample preparation procedures currently in use (as ranked by percentage of survey respondents who reported using each method):

- | | | |
|-------------------|--|------------------------------------|
| 1. Filtration | 5. pH adjustment | 8. Sonication |
| 2. Centrifugation | 6. Vortexing | 9. Solid phase extraction (SPE) |
| 3. Dilution | 7. Concentration (e.g., by ultrafiltration, precipitation) | 10. Liquid-liquid extraction (LLE) |
| 4. Evaporation | | |

According to a 2013 survey of users performed by LCGC magazine, filtration was the most commonly used sample preparation method. For complex samples containing components that contribute to high background and/or interfere with analyte ionization

and fragmentation, filtration alone cannot provide the sample necessary for analysis, but forms an integral part of an overall sample preparation strategy, which involves other sample preparation techniques, like extraction, centrifugation and depletion.

Sample preparation tips

Choose an appropriate membrane filter to remove particles from your sample.

The presence of particles in a sample can reduce the signal-to-noise ratio, reduce column lifetime, and increase backpressure in the LC system, potentially causing system failure. Filtration through a microporous membrane is a simple and effective method for

removing particles from a sample. However, particle retention ability is different between different membranes and between different suppliers. As **Figure 2** suggests, PTFE filters with polypropylene housing consistently deliver high particle retention.




Figure 2. Particle retention ability differs between different membranes and between different suppliers. In principle, 100% of filtration membranes should retain particles. To test this hypothesis, microporous membranes from two different vendors (pore size 0.2 μm) were tested for latex particle retention following filtration of a suspension of 0.3 μm latex particles in water. PTFE=polytetrafluoroethylene; PP=polypropylene; PVDF=polyvinylidene fluoride; PES=polyethylenesulfone.

Minimize extractables (contaminants) from sample prep device

Extractable impurities can generate interfering peaks in a chromatogram or mass spectrum, making it difficult or impossible to identify or quantify analytes of interest. Therefore, it is important to use a sample preparation device that leaches minimal impurities into the sample.

Though a number of syringe filters are certified as “low-extractable” for use in high performance LC (HPLC), most of those filters are certified using HPLC

coupled to detection of ultraviolet (UV) absorbance. Though this method provides information about the levels of UV-absorbing extractables coming from a filter, this information does not necessarily correlate with data obtained from an MS detector.

Select vendors now validate syringe filters using mass spectrometric analysis of extractables, which provides valuable guidance in choosing an appropriate membrane for filtering your sample.



TIP

In general, hydrophilic PTFE syringe filters provide the cleanest samples (with the lowest levels of extractable impurities). Presence of polymeric extractable impurities (such as from polypropylene syringe filters) complicate analysis of small molecular analytes.

Table 3. Overall mass spectral signal intensity for five different types of HPLC-certified syringe filters when tested using eight different solvents. The range of chemical compatibility with solvents may indicate the general level of extractables leached by a particular membrane. Millex® LCR filters, which contain hydrophilic PTFE and have broad compatibility with solvents, show the lowest level of signal intensity (and therefore background noise). On the other hand, polypropylene syringe filters from vendor A as well as nylon syringe filters from vendors A and B all show very high levels of extractables, impacting background signal.

| | Hydrophilic PTFE | Polypropylene (Vendor A) | Polypropylene (Vendor B) | Nylon (Vendor A) | Nylon (Vendor B) |
|--|------------------|--------------------------|--------------------------|--------------------|--------------------|
| Reproducibility | | | | | |
| Range of Compatibility with Organic Solvents | Broad | Broad | Broad | Moderate | Moderate |
| Extractable Level | Low | High | Medium | High | High |
| Nature of Extractables | MW 100–400 Da | Polymeric | Variable | Polymeric-Variable | Polymeric-Variable |

Consider these parameters for evaluating the suitability of a membrane filter for LC-MS:

Solvent compatibility of device

When selecting a filter, determine if constituents in the liquid being filtered will chemically attack the filter. If the filter undergoes chemical degradation, its performance will be compromised, and it may release foulants into the sample stream.

Some solvents may be incapable of dissolving the filter, but could be absorbed into the polymer matrix, causing it to swell over time, altering the effective pore size of the filter and changing its performance.

Lot-to-lot reproducibility of extractables level

This parameter reflects the consistency with which filters are manufactured. Since there are very few MS-certified filters, this parameter helps select the right filter for MS applications and indicates the degree of variation in levels of extractables when different lots of syringe filters are used.

Intensity of signal contribution from extractables: Total Ion Current (TIC) chromatograms

LC-MS-certified membrane filters should be supplied with a total ion current chromatogram that shows the intensity of all peaks detected under a specified set of experimental conditions, normalized to an internal standard. The TIC chromatograms can enable comparisons of extractable profiles between membranes and different filter vendors.

Type of extractables: low molecular weight, discrete peaks vs. polymeric peaks

Any type of extractables can confound downstream analysis, but the discrete peaks from low molecular weight extractables are typically less problematic than peaks from polymeric extractables, which typically show peaks separated by a common mass difference ranging over a wide m/z range. (See Appendix II for a table of mass differences of repeating units derived from common contaminating extractables.) Polymeric extractables are also difficult to remove from the sample or mass spectrometer, even after extensive cleaning of the mass spectrometer.

Adsorption of analyte to device

Because the internal surface area of polymeric microporous membranes is 100–600 times as great as the frontal surface area, there is a vast internal surface area available for nonspecific binding.

Choosing a membrane filter with low nonspecific analyte binding ensures that the overall molecular composition of the filtrate is minimally altered upon passing through the device.

Common extractable contaminants

- Polyethylene glycol (PEG)
- Metal ions (e.g., lithium, sodium, potassium, copper, platinum, iron)
- Phthalates (present in many plastics)
- Slip agents (amides)




Figure 3. Few extractable impurities from Millex® LCR syringe filter (left) containing 0.45 µm pore hydrophilic PTFE membrane as detected by MS. In contrast, a syringe filter containing 0.45 µm pore polypropylene membrane (Vendor C, right) shows significant extractables. Presence of polymeric extractable impurities (from polypropylene syringe filters) complicate analysis of small molecular analytes. Millex® LCR filters showed a highest peak intensity of about 8×10^5 for extractable masses, whereas Vendor C polypropylene syringe filters showed extractable levels about 30 times higher (2.5×10^6). Such high signal intensity, which can be comparable to the signal from the analyte of interest, can make analyte quantitation very challenging.




Figure 4. Nylon syringe filters are a common source of extractables. In this experiment, the extractable peaks seen after filtration through nylon were greatly reduced when the first milliliter of filtrate was discarded and the second milliliter was analyzed.

Another potential source of contamination is the syringe used to filter and/or inject the sample. Table 4 shows the level of zinc contamination from various types of syringes.

TIP
Prerinse the filter with sample/solvent to reduce the extractables.

Table 4. Level of zinc contamination with respect to syringe used.

| Syringe Used | Time of Contact | Zn Contamination (ppb) |
|---|-----------------|------------------------|
| Plastic with air gap | No contact | < 10 |
| Plastic with black piston seal | 15 min. | 96 |
| Plastic with black piston seal | 30 min. | 171 |
| Glass with metal Luer fitting & PTFE piston | 30 min. | 470 |



TIP

Use a plastic syringe with an air gap between the sample and the piston. Any surface that comes in contact with the sample has the potential to introduce extractables as well as contribute to analyte binding.

Ultrafiltration separates free from protein-bound analytes

Centrifugal ultrafilters, particularly devices with regenerated cellulose membrane that have defined nominal molecular weight cutoffs, are ideal for separating free from protein-bound microsolute in serum, plasma, and other biological samples,

as illustrated in **Figure 5**. This sample preparation method has been cited in LC-MS analyses for deproteinizing samples to reduce complexity or matrix interference. The method has also been used for LC-MS analyses of binding studies in new drug investigations.



TIP

If filtering many samples at a time, increase throughput while maintaining consistency by using 96-well or 384-well filter plates and a microplate-compatible vacuum manifold.




Figure 5. LC-MS analysis of total vs. free warfarin present in blank matrix (top) and human plasma samples (bottom). Unbound (free) warfarin was separated from protein-bound warfarin using centrifugal ultrafiltration devices, such as the Centrifree® device. Adapted from Jensen BP, Chin PK, Begg EJ. Quantification of total and free concentrations of R- and S-warfarin in human plasma by ultrafiltration and LC-MS/MS. *Anal Bioanal Chem.* 2011 Oct;401(7):2187-93.

When filtration isn't enough

For more complex sample matrices, use more specific sample preparation methods, such as solvent evaporation, protein precipitation, liquid/liquid extraction, QuEChERS, and SPE to transform samples into forms suited for LC-MS.

Some useful tools for these procedures include:

- Separatory funnel
- EXTrelut® pre-packed columns for extraction of lipophilic compounds from aqueous solutions – for sorbent-supported LLE workflows
- Solvents, acids, bases, salts – for protein precipitation
- LiChrolut® product range – for SPE
- Supel™-Select Polymeric SPE – HLB and ion-exchange phases for a wide range of applications and pH conditions
- Discovery® SPE and Supelclean™ SPE lines for a comprehensive range of reverse-phase, ion exchange, mixed mode, and normal-phase SPE products

For samples with high salt load (e.g., food, body fluids or tissue) a desalting (sample preparation step) using LiChrolut® cartridges is recommended.

Sample characteristics determine your SPE procedure

Is your sample matrix:

Is your analyte of interest more soluble in:

Is your compound:

Is your compound:

Is your compound a:

Do you want to recover your analyte from the SPE packing?




Figure 6.

Visit [SigmaAldrich.com/SPE](https://www.sigmaaldrich.com/SPE) for all of your solid phase extraction needs

Interference Removal Using Supel Select SCX SPE Cartridges to Analyze Illicit Bath Salts in Urine

The analysis of bath salts from urine samples is demonstrated using polymeric SPE sample preparation, followed by hydrophilic interaction liquid chromatography (HILIC) analysis with TOF-MS detection. Supel™-Select SCX SPE is used for the processing and sample cleanup of the urine samples. The figure below illustrates the monitored bath salt ions in a spiked urine sample after SPE cleanup (yellow),

in a diluted spiked urine sample without cleanup (green) and in a urine blank after Supel™ Select SCX cleanup (blue). Notice the chromatogram containing the bath salts in the spiked urine sample after SPE cleanup contains no interfering peaks. Therefore, the effectiveness of the Supel™ Select SCX cleanup is demonstrated and the analysis is more robust and reliable.

Figure 7. LC-MS Analysis of Cathinones (Bath Salts) on the Ascentis® Express HILIC (Si) Column


| | |
|------------------|--|
| sample/matrix: | 1 mL urine spiked to 100 ng/mL of bath salt mixture |
| SPE tube: | Supel™-Select SCX, 30 mg/1 mL (54240-U) |
| conditioning: | 1 mL 1% formic acid in acetonitrile, then 1 mL water |
| sample addition: | 1 mL spiked urine |
| washing: | 1 mL water, 1 mL 1% formic acid in acetonitrile, 1 mL water |
| elution: | 2 mL 10% ammonium hydroxide in acetonitrile |
| column: | Ascentis® Express HILIC (Si), 10 cm x 2.1 mm I.D., 2.7 µm (53939-U) |
| mobile phase: | (A) 5 mM ammonium formate acetonitrile; (B) 5 mM ammonium formate water; (98:2, A:B) |
| flow rate: | 0.6 mL/min |
| pressure: | 127 bar |
| column temp: | 35 °C |
| detector: | MS, ESI+, 100-1000 m/z |
| injection: | 1 µL |
| sample: | 200 ng/mL in acetonitrile |

Spiked Urine Sample
Monitored Ions After
Supel-Select SCX Cleanup

Diluted Spiked Urine
Monitored Ions Without SPE
Cleanup

Urine Blank Monitored Ions
After Supel-Select SCX
Cleanup

- 3,4-Methylenedioxypropylvalerone (MDPV)
- Buphenedrone
- 3-Fluoromethcathinone
- Butylone
- Ethylone
- 4-Fluoromethcathinone
- Mephedrone
- Methylone
- Methedrone



Protein precipitation, followed by filtration, is often an effective, simple way to reduce the complexity of the sample matrix (Figure 8). For this process, it can be advantageous to use a filter plate, which enables

precipitation and filtration in a single device, eliminating the need for sample transfer and thereby improving analyte recovery.




Figure 8. Analyte recovery after protein precipitation. Three different drugs (chlorpheniramine, imipramine and niflumic acid) were spiked into plasma at various concentrations. Protein precipitation was carried out using 1:4 water:acetonitrile as the precipitating solvent. The samples were filtered through various multiwell filter plates with polypropylene, PTFE, and washed PTFE (washed with solvent). Drug recovery in the filtrate was determined using LC-MS/MS analysis of the filtrate.



TIP

Prepare samples for nano LC-MS using ZipTip® pipette tips (Figure 9). This sample preparation microdevice is a 10 µL pipette tip with a 0.6 or 0.2 µL bed of chromatography media fixed at its end with no dead volume. It is ideal for concentrating and purifying samples for sensitive analyses such as nano LC-MS or MALDI-ToF MS.




Figure 9. ZipTip® pipette tips increase sensitivity of mass spectrometric analysis. Plasma sample from rat dosed with 10 mg/kg prazosin was injected into an LTQ/Orbitrap mass spectrometer by nanoelectrospray (a) before and (b) after preparation using a C18 ZipTip® pipette tip. Adapted from Erve JCL et al, Rapid Commun. Mass Spectrom. 2008; 22: 3015–3026.

Phospholipids: a concern or LC-MS analysis of small molecules in biological matrices

Phospholipids are present as a major component of all cell membranes.

They are therefore present in all biological sample matrices including serum, plasma and whole blood and can be a problem in LC-MS analysis of small molecules because they often co-elute and ionize along with the analytes of interest. This co-elution results in ion suppression (an erroneous decrease) of the MS signal that can cause variability and impact LC-MS result accuracy. Even if phospholipids do not co-elute with the analyte of interest, they can accumulate on your analytical column.

Phospholipid removal techniques:

To overcome the problem of phospholipid-induced ion suppression, some analysts try traditional SPE. Traditional SPE often requires time-consuming and complex method development, but still only removes nominal amounts of phospholipids. A variety of products designed specifically for the removal of both proteins and phospholipids are now commercially available, including HybridSPE® plates and cartridges (**Figure 10**).




Figure 10. Removing phospholipids can improve the signal-to-noise ratio in LC-MS.

Table 5. Guide to Sample Preparation Tools

| Description | Cat. No. |
|--|-----------|
| Millex® Syringe Filters | |
| Millex®-LCR Filter, 0.45 µm, Hydrophilic PTFE, 33 mm, non-sterile, 50/pk | SLCR033NS |
| Millex®-LCR Filter, 0.45 µm, Hydrophilic PTFE, 33 mm, non-sterile, 250/pk | SLCR033NB |
| Millex®-LCR Filter, 0.45 µm, Hydrophilic PTFE, 33 mm, non-sterile, 1000/pk | SLCR033NK |
| Millex®-LCR Filter, 0.45 µm, PTFE, 13 mm, non-sterile, 100/pk | SLCRX13NL |
| Millex®-LCR Filter, 0.45 µm, PTFE, 13 mm, non-sterile, 1000/pk | SLCRX13NK |
| ZipTip® Pipette Tips | |
| ZipTip® with 0.6 µL C 4 resin, 8/pk | ZTC04S008 |
| ZipTip® with 0.6 µL C 4 resin, 96/pk | ZTC04S096 |
| ZipTip® with 0.6 µL C 4 resin, 960/pk | ZTC04S960 |
| ZipTip® with 0.2 µL C 18 resin, 8/pk | ZTC18M008 |
| ZipTip® with 0.2 µL C 18 resin, 96/pk | ZTC18M096 |
| ZipTip® with 0.2 µL C 18 resin, 960/pk | ZTC18M960 |
| ZipTip® with 0.6 µL C 18 resin, 8/pk | ZTC18S008 |
| ZipTip® with 0.6 µL C 18 resin, 96/pk | ZTC18S096 |
| ZipTip® with 0.6 µL C 18 resin, 960/pk | ZTC18S960 |
| ZipTip® with 0.6 µL strong cation resin, 8/pk | ZTSCXS008 |
| ZipTip® with 0.6 µL strong cation resin, 96/pk | ZTSCXS096 |
| Samplicity® Filtration Systems | |
| Millex®-LCR Filters for Samplicity G2, 0.45 µm Hydrophilic PTFE, 250/pk | SAMP2LCRB |
| Millex®-LG Filters for Samplicity G2, 20 µm Hydrophilic PTFE, 250/pk | SAMP2LGNB |
| Millex®-HV Filters for Samplicity G2, 0.45 µm Hydrophilic PVDF, 250/pk | SAMP2HVNB |
| Millex®-GV Filters for Samplicity G2, 22 µm Hydrophilic PVDF, 250/pk | SAMP2GVNB |
| Samplicity® G2 Filtration System, Bold Blue | SAMP2SYSB |

| Description | Cat. No. |
|---|-----------|
| Ultrafree®-MC and -CL Centrifugal Microfiltration Units | |
| Ultrafree®-MC Filter, 0.22 µm Hydrophilic PTFE, 25/pk | UFC30LG25 |
| Ultrafree®-MC Filter, 0.45 µm Hydrophilic PTFE, 25/pk | UFC30LH25 |
| Ultrafree®-CL Filter, 0.22 µm Hydrophilic PTFE, 25/pk | UFC40LG25 |
| Ultrafree®-CL Filter, 0.45 µm Hydrophilic PTFE, 25/pk | UFC40LH25 |
| Centrifree® Ultrafiltration Device with Ultracel® Membrane | 4014 |
| MultiScreen® Filter Plates | |
| MultiScreen® Solvinert 96-well Plate, 0.45 µm Hydrophilic PTFE, 50/pk | MSRLN0450 |
| MultiScreen® Deep Well Solvinert 96-well Plate, 0.45 µm Hydrophilic PTFE, 10/pk | MDRLN0410 |
| Solid Phase Extraction | |
| Extrelut® NT 20 pre-packed columns for extraction of lipophilic compounds from aqueous solutions (20 mL sample) | 115096 |
| LiChrolut® RP-18 E (40 – 63 µm) 500 mg 3 mL standard PP-tubes 50 extraction tubes per package | 119849 |

Table 6. HybridSPE Cartridges and 96-well Plates

| Description | Qty. | Cat. No. |
|---|------|----------|
| Well Plates | | |
| HybridSPE®-PLus 96-well Plate, 50 mg/well | 1 | 575659-U |
| | 20 | 575673-U |
| HybridSPE®-PL, Small Vol. 96-well Plate, 15 mg/well | 1 | 52794-U |
| | 20 | 52798-U |
| HybridSPE®-PLus 96-Well Plate Essentials Kit (contains: 96-well Plate, 50 mg/well, 1 cap mat, sealing film, and collection plate) | 1 | 52818-U |
| SPE Cartridges | | |
| HybridSPE®-PL Ultra Cartridge, 30 mg/1 mL | 100 | 55269-U |
| HybridSPE®-PL Cartridge, 30 mg/1 mL | 100 | 55261-U |
| | 200 | 55276-U |
| HybridSPE®-PL Cartridge, 500 mg/6 mL | 30 | 55267-U |

Proper mobile phase preparation to minimize contaminants

For LC-MS, use the highest quality of pure solvents and reagents and avoid further contamination by careful handling. Any impurity could cause signal suppression and/or adduct formation with target molecules and therefore decrease sensitivity (signal-to-noise ratio) and/or increase complexity of the mass spectrum.




Figure 11. Combined TICs of the blank runs of four different acetonitrile qualities. All solvents were delivered to the MS source via an LC system.

Hypergrade and gradient grade solvents minimize contaminant peaks

Figure 9 illustrates the influence of LiChrosolv® acetonitrile quality on the background noise intensity in mass spectra. Supelco® solvents labeled “hypergrade for LC-MS LiChrosolv®” are dedicated for use with MS systems and deliver minimized contaminant peaks,

ion suppression, adduct formation and background noise and therefore maximize sensitivity. Gradient grade solvent quality (labeled “gradient grade for liquid chromatography LiChrosolv®”) are suitable for LC-UV gradient runs.

Use ultrapure water (bottled or freshly purified)

Ultrapure water for LC-MS applications can be either bottled or freshly delivered from a water purification system. The choice is mainly determined by daily consumption. Demineralized tap water is not recommended for use in combination with LC-MS setups because of possible system contamination. The quality of LiChrosolv® bottled water for chromatography and freshly purified ultrapure water produced from a Milli-Q® lab water purification system is consistently high and generally independent of the regularity of use.

Ultrapure water quality is perfectly suitable for the production of mobile phase, buffers, blanks, standards preparation, sample dilution, glassware rinsing or extraction used in these critical applications.

Careful storage and handling of water are critical to prevent contamination during drawing. **Figure 12** displays total ion currents (TICs) of Milli-Q® ultrapure water drawn at different points in time: Directly on Monday (after system standby over the weekend), on the same day after discarding several liters prior to ultrapure water collection, and after four days of daily use. Generally, it is recommended to flush the system every morning by drawing and discarding a few liters prior to water collection.



TIP

For your most sensitive LC-MS analyses, we recommend replacing the 0.22 µm filter at the point-of-delivery (POD) of your Milli-Q® ultrapure water system with an LC-Pak® polisher. The LC-Pak® polisher contains C18 silica and is optimized for the most sensitive organic analyses.




Figure 12. Exploiting the full potential of Milli-Q® ultrapure water systems via proper handling. Top: MS spectra of two samples of ultrapure water delivered at different points of time; bottom: TICs of the same samples and one additional sample. All analyses were performed via direct injection of the pre-concentrated solvents into the MS operated in positive ESI mode.

Solvent storage

Tips for maintaining your solvent purity:

- Store all eluents (water and organic) in surface-treated amber glass bottles (original packaging of all Supelco® LC-MS grade solvents) or in borosilicate glass (if solvents have to be decanted).
- Select a solvent storage system that is appropriate for usage volume and withdrawal frequency (**Table 7**).
- Do not use standard glass bottles; silica and alkali dissolve and form adducts $[M+X]^+$ with analytes.
- Use Supelco® HPLC bottle caps/adapters with tube connections and membrane filter mounted directly on the original brown glass bottle. This protects both solvents and environment.
- Avoid decanting; it is a possible source of contamination.
- Avoid improvised repairs for fixing solvent tubing; this may cause leakages and/or release of contaminants to the eluents.
- Do not use plastic devices (bottles, funnels, etc.) to handle or store solvents, buffers, etc. Solvents extract additives (anti-static agents, stabilizers, plasticizers) from plastic, a source of contaminant ghost peaks and increased background noise.

Table 7. Select solvent storage options based on volume of usage and frequency of withdrawal to minimize contamination

| Solvent storage systems | Storage container volume |
|---|---|
| Bottle top adapters for directly connecting solvent bottles to LC system | 1 L, 2.5 L, 4 L (for infrequent withdrawal) |
| Stainless steel barrel with adapter for direct withdrawal | 10 L, 30 L (for frequent withdrawal) |
| Stainless steel barrel directly connected to LC system | |
| Central storage of stainless steel barrels with adapters to supply solvent to multiple different laboratories | |

Using water as a mobile phase?

Keep in mind these additional considerations:

- Keep 5% organic solvent in your eluent if chromatographic conditions allow. This avoids microbial contamination of bottle, tubing and LC system.
- Keep 5% of aqueous eluent in the organic mobile phase to avoid buffer precipitation in the system, e.g., in valves, and subsequent tedious cleaning procedures.

Solvent container cleaning: avoid the dishwasher

Dishwashers are standard laboratory equipment, but they are operated using chemicals such as strong bases and surfactants. Strong bases can lead to dissolution of silica and alkali from glassware and cause the formation of adducts $[M+X]^+$ with analytes, while traces of surfactants remain on the glass surface after the cleaning process and decrease MS sensitivity by increasing background noise.

The easiest way to avoid dishwashing is “cleaning” of all equipment via simple evaporation of both solvents and additives. All chemicals dedicated to the application in LC-MS are volatile; therefore, this procedure is straightforward as long as chemicals are highly pure and microbial growth can be eliminated. In case of equipment contamination, flushing with LiChrosolv® solvents or Milli-Q® ultrapure water or organic hypergrade solvents has to be performed to achieve sustainable cleaning.

Buffers and additives

When working with buffer to adjust pH of eluents, keep in mind:

- Use volatile salts (such as ammonium formate, ammonium acetate, or trimethylamine). Nonvolatile salts (e.g., phosphates, borates, sulfates or citrates) precipitate in and block the MS source, requiring tedious cleaning procedures.
- Total ionic strength of the eluent should not exceed 20 mM. Adjust buffer concentration in the aqueous solvent accordingly. Buffers for LC-MS

should be prepared using the purest salt and acid/base quality available. If possible, avoid working with an ammonium bicarbonate buffer. The salt is often highly contaminated — see comparison with ammonium acetate (**Figure 13**).




Figure 13. Comparison of MS spectra (top) and TIC chromatograms (bottom) of the two buffer systems, ammonium bicarbonate and ammonium acetate. Both mixtures were prepared using Milli-Q® ultrapure water and the same acetic acid source and were analyzed via direct injection into the mass spectrometer operated in positive ESI mode. Note the high MS background signals observed when utilizing ammonium bicarbonate as a buffer.

Buffer pH is generally adjusted via a titration with the respective acid or base and monitored with a pH electrode. The unavoidable contamination of the buffer solution with alkali ions from the pH electrode can be decreased by using a miniaturized system available from several suppliers. Unlike standard equipment with a diameter of approximately 10 mm, the diameter of miniaturized electrodes is only 3 mm.

Buffers not only adjust the pH and ionize a target molecule [M], they can also form adducts [M+buffer], e.g., with ammonium, alkali, halogens, formate or acetate. This leads to the detection of additional peaks in the MS spectrum. Even a complete suppression of the analyte signal is possible when the vapor pressure of the resulting adduct (mainly alkali) is decreased significantly. As a result of this phenomenon and in order to keep the ESI source clean, volatile buffers are recommended.



TIP

Avoid using TFA. Trifluoroacetic acid (TFA) is widely used as an ion pairing reagent to improve the liquid chromatographic separation of peptides or proteins when using standard UV for detection. However, TFA can cause strong ion suppression in mass spectrometry (mainly in negative ESI mode) and also contaminates the LC-MS system. Formic acid (0.1%) is commonly used instead as a mobile phase modifier that is compatible with LC-MS.

Table 8. Guide to Mobile Phase Preparation Reagents

| Description | Cat. No. |
|--|-----------|
| Milli-Q® Water Purification Solutions | |
| Milli-Q® IQ 7000 ultrapure water system | ZIQ7000T0 |
| Milli-Q® IQ 7003/05/10/15 pure and ultrapure water system | ZIQ7005T0 |
| LC-Pak® Application-Specific Polisher | LCPAK00A1 |
| LiChrosolv® Solvents | |
| Acetonitrile hypergrade for LC-MS LiChrosolv® | 100029 |
| Methanol hypergrade for LC-MS LiChrosolv® | 106035 |
| Ethanol gradient grade for liquid chromatography LiChrosolv® | |
| 2-Propanol gradient grade for liquid chromatography LiChrosolv® | 101040 |
| Toluene for liquid chromatography LiChrosolv® | 108327 |
| Water for chromatography LiChrosolv® (LC-MS) | 115333 |
| Suprapur® Inorganic Acids and Bases | |
| Acetic acid (glacial) 100% Suprapur® | 100066 |
| Ammonia solution 25% Suprapur® | 105428 |
| Formic acid 98–100% Suprapur® | 111670 |
| Hydrochloric acid 30% Suprapur® | 100318 |
| Other Reagents | |
| Acetic acid (glacial) 100% anhydrous for analysis EMSURE® ACS, ISO, Reag. Ph Eur | 100063 |
| Ammonia solution 28–30% for analysis EMSURE® ACS, Reag. Ph Eur | 105423 |
| Ammonium acetate for analysis EMSURE® ACS, Reag. Ph Eur | 101116 |
| Dichloromethane for organic trace analysis UniSolv® | 106454 |
| Formic acid 98–100% for analysis EMSURE® ACS, Reag. Ph Eur | 100264 |
| n-Hexane for organic trace analysis UniSolv® | 104369 |
| n-Pentane for organic trace analysis UniSolv® | 107288 |
| Petroleum benzene boiling range 40–60°C for organic trace analysis UniSolv® | 116740 |
| 2-Propanol for analysis EMSURE® ACS, ISO, Reag. Ph Eur | 109634 |


*Contact your local representative for detailed ordering information.

For even more LC-MS solvents and reagents, visit: [SigmaAldrich.com/lc-ms](https://www.sigmaaldrich.com/lc-ms)



TIP

Avoid equilibrating columns with more than 10 column volumes of mobile phase (or one blank gradient run with subsequent equilibration). Contaminants in solvents and additives can accumulate on a stationary phase. **Figure 14** shows this effect for plasticizers dissolved in the eluent on a reversed phase column after equilibration for 0, 15 and 60 minutes. While these compounds would become eluted as very broad peaks under isocratic conditions (and cause an increased background noise), they elute as distinct, intensive peaks under gradient conditions and can interfere with analyte signals. Instead, run samples immediately after two or three blank runs to ensure that the system is stable prior to sample analysis.



Chromatographic conditions

| | |
|---------------------|--|
| System | Bruker Esquire 6000plus |
| Detection | Pos. ESI-MS, TICs |
| Flow rate | 0.4 mL/min |
| Mobile phase | A: Water (Cat. No. 115333) / acetonitrile (Cat. No. 100029) 95/5 (v/v) + 0.1 % formic acid (Cat. No. 100264) B: Acetonitrile + 0.1% formic acid |
| Gradient | 0 min 100% A, 3 min to 95% B, 5 min at 95% B, then back to 100% A |
| Temperature | 25°C |
| Sample | Plasticizers (*) were added by the immersion of plastic tubing in aqueous solvent B |

Figure 14. Accumulation of contaminants on an HPLC column for various periods of time and elution via a gradient profile.

Effects of column choice on LC-MS performance

For samples available in small amounts (such as plasma or serum) or where analyte concentrations are low, a setup consisting of both highly sensitive separation and MS detection techniques is necessary for proper identification of the target molecules.

Suboptimal column choice and misuse of the column could decrease the signal-to-noise ratio and increase background noise.

Tips for proper column choice:

- 1. Refer to comprehensive column selection guides** for full guidance on selecting a column with optimal stationary phase and dimensions and to match method specifications. Column selection guides can be found at: [SigmaAldrich.com/hplc](https://www.sigmaaldrich.com/hplc)
- 2. Pick column in accordance with eluent pH.** Using an eluent pH that is too high (e.g., >8) can dissolve the backbone of silica-based HPLC columns. Using a pH that is too low (<2) can strip the stationary phase (C18, etc.). Both options can lead to additional signals in your spectrum, increased background noise and/or signal suppression. Both scenarios may decrease column lifetime.
- 3. Consider using polymeric columns for highly alkaline samples.** Polymeric columns are more stable at high pH than silica columns, where the silica may dissolve. However, polymeric columns possess smaller phase ratios and therefore, lower resolution. In addition, they are prone to swelling in organic solvents, leading to changed chromatographic characteristics. Furthermore, due to micropores in the stationary phase, column performance may be lower as compared to silica-based columns.
- 4. Highly endcapped stationary phases** are another good option for sample analysis at high pH. Endcapping leads to more pH-stable columns.
- 5. Column diameter influences the sensitivity of the analysis.** The sensitivity increases with decreasing column internal diameter (or increasing mass of the injected sample). For example, when changing from a 4.6 mm i.d. column to a 0.1 mm i.d. capillary

column, sensitivity theoretically increases by a factor of approximately 2000 (**Table 9**). Hence, a combination of capillary chromatography coupled to mass spectrometry may be the best combination for high sensitivity analysis. However, it is important to note that extra-column effects may impact signal-to-noise ratio when column diameter is decreased—for example, the system dwell volume, dead volume in the system, the ability of the system pump to deliver accurate gradient and the volume of the detector cell can all result in peak broadening and loss of sensitivity.


Table 9. Effect of decreasing column diameter on flow rate and relative sensitivity.

| Column i.d. (mm) | Typical flow rate (µL/min) | Relative sensitivity |
|------------------|----------------------------|----------------------|
| 4.6 | 1000 – 6000 | 1 |
| 2.0 | 200 – 800 | 5.3 |
| 0.2 | 0.5 – 20 | 530 |
| 0.1 | 0.4 – 3 | 2100 |
| 0.05 | 0.1 – 0.8 | 8500 |

- 6. If the sample and analyte allow for HILIC chromatography,** consider using HILIC instead of reversed phase columns. In HILIC chromatography, analysis is performed under highly organic conditions (e.g., 80% acetonitrile and 20% aqueous buffer) and polar compounds elute later than in reverse phase chromatography. Under these conditions, the eluent is vaporized more easily in the MS source, resulting in better sensitivity (signal-to-noise ratio).


Tips for column usage

1. Wash columns after each use with appropriate strong eluent to remove all adsorbed compounds. (See “Column bleeding” and **Figure 16.**)
2. Use the column at proper operating temperature (refer to each column User Guide) in order to avoid loss of stationary phase or dissolution of column backbone, both of which contribute to the appearance of additional signals on the spectrum that may interfere with analysis.
3. You can use a guard column (usually either 5 mm or 10 mm long) directly in front of the main column to protect the column against contamination (such as particles, sample matrix). Guard columns should be changed frequently in order to keep system backpressure low and in order to maximize the lifetime of the analytical column.



Chromatographic conditions

| | | |
|------------------|---|-----------------------------|
| Column | Chromolith® Performance RP18 endcapped 100-4.6 mm | |
| Mobile Phase | Acetonitrile/Water 40/60 (v/v) | |
| Flow Rate | 4.0 mL/min | |
| Pressuer | 86 - 40 = 46 bar | |
| Detection | UV 254 nm | |
| Temperature | 23 °C | |
| Injection Volume | 1 µL | |
| Sample | 1. Biphenyl-4,4'-diol | 3.1 mg |
| | 2. Biphenyl-2,2'-diol | 20.8 mg |
| | 3. Biphenyl-4-ol | 12.5 mg |
| | 4. Biphenyl-2-ol | 11.1 mg in 100 mL A/W 50/50 |



Chromatographic conditions

| | | |
|------------------|---|-----------------------------|
| Column | Chromolith® Performance RP18 endcapped 100-3 mm | |
| Mobile Phase | Acetonitrile/Water 40/60 (v/v) | |
| Flow Rate | 1.7 mL/min | |
| Pressure | 65 - 15 = 50 bar | |
| Detection | UV 254 nm | |
| Temperature | 23 °C | |
| Injection Volume | 0.4 µL | |
| Sample | 1. Biphenyl-4,4'-diol | 3.1 mg |
| | 2. Biphenyl-2,2'-diol | 20.8 mg |
| | 3. Biphenyl-4-ol | 12.5 mg |
| | 4. Biphenyl-2-ol | 11.1 mg in 100 mL A/W 50/50 |

Figure 15. Decreasing column diameter may improve sensitivity. Typical fast separation of four compounds in less than two minutes using a Chromolith® 4.6 mm i.d. column at a flow rate of 4 mL/min (left). The same separation was achieved on a Chromolith® 3 mm i.d. column (right). Both chromatograms exhibit excellent column efficiency and peak resolution, however the 3 mm i.d. column demonstrates improved sensitivity at just 1.7 mL/min, thus saving 57% of solvents.

Column bleeding

The stationary phase of every HPLC column (except for normal phase systems) is made out of covalently bound organic entities altering its physical properties. Depending on the quality of both phase modification and a subsequent washing step, these entities (e.g., octadecyl, cyano, phenyl) can be stripped off the column during a chromatographic run and cause weak to severe interfering signals.

This unwanted phenomenon is referred to as “column bleeding” and leads to a decreased sensitivity in MS. It can be avoided by flushing the column prior to analysis using isopropanol and 0.1% formic acid as a solvent at half optimum flow for one hour. This process removes unbound or weakly bound organic entities, minimizes column bleeding and hence increases sensitivity by decreasing background noise (**Figure 16**).




Figure 16. Column washing can compensate for column bleeding. Total ion current (TIC) of a competitor column after 20 gradient runs: left- unwashed column; right- column washed with 0.1% formic acid in isopropanol at 0.2 mL/min for 48 min.

Overview of U/HPLC columns for LC-MS

For LC-MS, particulate HPLC columns, preferably based on high purity Type B silica, are widely used.

1. Small particles deliver high separation efficiency/ peak capacity
2. Suitable for cleaner samples and MS, after removal of matrix components

Typical fully porous particulate columns:

Purospher® STAR HPLC and UHPLC columns are available with particle sizes of 2 µm, 3 µm and 5 µm in various column modifications providing high efficiency, extended pH stability (pH 10.5) for RP-18e and RP-8e and stability in aqueous mobile phases.

Ascentis® and Discovery® HPLC columns are available with 3 µm and 5 µm particle sizes and a very broad range of column chemistries providing selectivity for the separation of almost every compound.

Titan® UHPLC columns based on monodisperse particles of 1.9 µm particle size provide very high efficiency due to a more consistent packed bed.

HILIC is superior for the separation of polar hydrophilic molecules, i.e., many of the endogenous molecules.

SeQuant® ZIC®-HILIC/cHILIC/pHILIC bonded zwitterionic stationary phases combine perfectly with ESI-MS detection due to the applied solvents and additives. A significant increase in sensitivity in comparison with reversed phase chromatography can be achieved. Strongly retained polar analytes can be removed from HILIC columns by changing to a more polar eluent.

Superficially porous particulate columns for maximum resolution and speed:

Fused-core® columns feature narrower particle size distribution and shorter diffusion path compared to fully porous particles. The result is increased resolution, added sensitivity and faster runs.

Ascentis® Express HPLC and UHPLC columns provide about 40% more efficiency in comparison to columns with fully porous particles of the same size. This performance enhancement is applicable to all HPLC instruments (in addition to UHPLC systems). Particle sizes of 2 µm, 2.7 µm and 5 µm are available with a very broad range of modifications, all excellently suitable for LC-MS use.

BIOshell™ HPLC and UHPLC columns deliver maximum speed and efficiency for the separation of biomolecules on both UHPLC and HPLC systems. The Fused-Core® superficially porous silica particles with pore sizes from 90 Å up to 1000 Å allows superior separation of glycans as well as very large proteins. In particular, a pore size of 1000 Å shows very clear advantages over common 300 Å pores for the separation of very large proteins in biopharmaceutical drug development such as monoclonal antibodies (mAbs) or proteins with molecular weights greater than 100 kDa.

The advantages of Fused-Core® columns are:

- Maximum speed and efficiency on both UHPLC and HPLC systems (particle sizes: 2 µm, 2.7 µm and 5 µm)
- 40% more efficiency in comparison to Fully Porous Particles (FPP) of same particle size
- UHPLC columns with 2 µm particles (pressure stable 1000 bar)
- Column dimensions from 0.075 mm ID (capillary columns) to 4.6 mm ID (analytical HPLC columns)
- Very broad range of column chemistries

Monolithic silica columns have high matrix tolerance

The analysis of samples with high matrix load requires tedious and time-consuming sample preparation steps. For cost-effective investigations, sample handling has to be kept as short as possible and combined with robust LC columns displaying a high matrix tolerance and long lifetime. The 50-2 mm monolithic silica column is well-suited for fast gradient run liquid chromatography, and the applied low flow rates make it the perfect choice for MS detection. Analysis of matrix-rich samples, such as food or tissue, can be performed on this robust column type without the need for a guard column or tedious and complex sample preparation procedures.

The advantages of Chromolith® monolithic silica columns are:

- Exceptional robustness or lifetime—described as number of injections—enabling cost savings
- High matrix tolerance decreases tedious sample preparation steps, speeds up all processes and allows for fast and simple HPLC analyses
- Very low column back-pressure, fast analytical speed and high reproducibility on standard HPLC systems as well as using UHPLC instruments

**Consult a column
selection guide today!**

[SigmaAldrich.com/hplc](https://www.sigmaaldrich.com/hplc)

Effects of Reference Material choice on LC-MS performance

How to choose the correct reference material quality grade for your needs?

Who uses reference materials?

Reference materials are a critical component of the analytical testing workflow and can help to avoid or identify contaminants. Through calibration of measurement systems, validation of methods, and quality control programs, reference materials ensure accuracy in testing. Proper selection of the right reference material for the laboratory's testing application is important, because results are only as accurate as your reference. Therefore it is important to understand the five major quality grades of reference materials.


ISO 17034 and quality grades of standards, reference materials and certified reference materials

The reference material hierarchy includes four major quality grades, from national metrology and other

primary standards to Certified Reference Materials (CRMs), Reference Materials (RMs) and Analytical Standards. Level of certification and traceability requirements increase for each higher level. Where national governments give standardization to the top level, specific ISO guidelines provide standardization for CRMs and RMs. These ISO requirements include ISO 17034, ISO/IEC 17025 and ISO Guide 31.

Reference material producers must meet these ISO requirements to manufacture CRMs or RMs. For both of these quality grades, Certificates of Analysis must be provided, and the information contained within is defined by the aforementioned ISO guidelines. The quality specifications for the last two levels are defined by each individual producer rather than by a national government or ISO accreditations specific to CRMs and RMs.

Figure 17. The Hierarchy of Reference Materials - What are the Different Types?



What is measured in each grade of reference material?

Purity and Identity of the material are typically included in the Certificate of Analysis for each of the five quality grades. Content and Stability are required for the primary standards or ISO-defined CRM and RM.

Analytical standards and research chemicals may or may not include these two parameters as their inclusion is dependent on the producer. Analytical Standards can also in some cases be Quality Control materials compliant with ISO Guide 80.

Homogeneity is required for the primary standards, CRM, and RM, but this parameter will not be found with the lower quality grades. Uncertainty and Traceability information are limited to just the primary standards and CRM. In the pharmaceutical world, secondary standards can be CRMs or RMs, but here, there are two different types of traceability – to the SI unit of measurement for the ISO-defined CRM as well as traceability to the primary compendial standard, which is a requirement specific to pharmaceutical secondary standards.

Table 10. The Hierarchy of Reference Materials – What’s the Difference?

| Parameter | NMI Standard | Compendial Standard | CRM | RM | Analytical Standard | Research Chemical |
|--------------|---|---------------------|--|-----------------------------|---------------------|-------------------|
| Purity | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Identity | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Content | ✓ | ✓ | ✓ | ✓ | maybe | |
| Stability | ✓ | ✓ | ✓ | ✓ | ✓ | |
| Homogeneity | ✓ | ✓ | ✓ | ✓ | | |
| Uncertainty | ✓ | | ✓ | | | |
| Traceability | ✓ | | ✓ | ✓ | | |
| Type | Primary Measurement Standard or Primary Standard (Pharma) | | Primary or Secondary Standard (Pharma) | Secondary Standard (Pharma) | | |

Choose the correct reference material for your testing purpose

For instrument qualifications and calibrations, establishing and maintaining traceability is critical. The selected reference material should help the laboratory achieve this. In daily routine system suitability applications, it might be important to qualify something that is practical and easy to use, yet reliable and cost effective for everyday use. In method validation, it’s

critical to use highly accurate and precise materials to show that the laboratory method is accurate and precise. For identity and screening purposes, important attributes of reference materials include proven authenticity and identity. For quantitation, assays, or stability assessment, stable and accurate reference materials are needed.

Table 11

| Type of test | Use of Ref. Mat. | Examples | Requirements of the Ref. Mat. |
|--|---|---|--------------------------------|
| Instrument qualification / Calibration | Establish system performance Measurement accuracy | Annual qualifications Routine balance calibrations | Traceable |
| Routine calibration / System suitability | Daily / weekly System / method specific Establish routine performance | Pre-use balance calibrations System performance checks for LC-UV/MS; GC-FID... | Qualify as suitable for use |
| Method validation | Accuracy Precision Specificity & interferences LOD/LOQ & Linearity | Pharma QC; Environmental testing Standards of the analyte(s), interferences, impurities | Accurate Traceable |
| Identity | Comparison of unknown to known | Incoming raw materials in pharma, food etc. Screening tests | Authenticity |
| Content or assay | Quantitation of analytes | Pesticide/toxin limits Pharma QC - API content | Certified content Traceable |
| Stability assessment | Monitor product stability | Pharma QC | Stable, homogenous |
| Internal Quality Control | Method accuracy | Routine quantitation of analytes - pharma/pesticides/diagnostics | Certified content Traceable |

To learn more visit [SigmaAldrich.com/quality-grades-crm](https://www.sigmaaldrich.com/quality-grades-crm)

Visit [SigmaAldrich.com/standards](https://www.sigmaaldrich.com/standards) to learn more about the Supelco® family of reference materials.

Which quality grade is the best fit for purpose?

Fit for purpose decisions in selection of reference materials can depend on several factors, from regulatory requirements, availability, and type of testing application to level of accuracy and sample matrix.

Table 12. Fit for Purpose Guidance in Standard Selection

| Type of Test | NMI Standard | Compendial Standard | CRM | RM | Analytical Standard | Reagent Chemical | Attribute |
|---|--------------|---------------------|-----|----|---------------------|------------------|---|
| Instrument qualification / Calibration | ✓ | ✓ | ✓ | | | | Traceability & Accuracy |
| Routine calibration/ System suitability | ✓ | ✓ | ✓ | ✓ | maybe | | Qualified standard (Primary or secondary) |
| Method validation | ✓ | ✓ | ✓ | ✓ | | | Accuracy, Precision, Bias |
| Identity | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | Authenticity |
| Content or assay | ✓ | ✓ | ✓ | ✓ | maybe | | Qualified standard |
| Stability assessment | ✓ | ✓ | ✓ | ✓ | maybe | | Qualified standard |
| Internal Quality Control | ✓ | ✓ | ✓ | ✓ | maybe | | Qualified standard |
| Regulatory / Accreditation | ✓ | ✓ | ✓ | ✓ | | | Qualified standard |

Useful LC-MS resources

Solutions for Mass Spectrometry

Since the development of an aluminum oxide for adsorption chromatography in 1904, Merck has offered chromatography products and solutions of exceptional, application-oriented quality. Our columns, plates and mobile phases are used in mass spectrometry labs throughout the world.

This compendium of selected articles is intended to help and support users of chromatographic separation methods in mass spectrometric laboratories worldwide, by providing an overview of different actual topics in this field of analytical science.

Download at: [SigmaAldrich.com/LCMS](https://www.sigmaaldrich.com/LCMS)

John Dolan's LC Troubleshooting Bible (2016 edition)

For over 30 years, John Dolan has been writing the LC Troubleshooting column of LC/GC Magazine, answering questions from chromatographers around the world. Browse the columns all in one place (1984 – 2016) at the LC Troubleshooting Bible website.

Questions addressed cover all aspects of chromatography, ranging from "Importance of Sample Matrix" to "Gradient Elution Surprises."

Visit: [lcresources.com/tsbible](https://www.lcresources.com/tsbible)

chromacademy.com

Described as the "world's largest e-Learning website for analytical scientists," CHROMacademy is a well-organized, all-in-one-place resource for basic LC/LC-MS training, troubleshooting, learning about advanced topics and even obtaining technical support from consultants. It includes video guides and quick tips, such as in "My LC-MS Isn't Behaving!"

Many features are available by subscription only, but for analysts doing a lot of LC, it's worth it.

Training modules within CHROMacademy include:

- Sample Prep
- HPLC
- MS

LC-MS Workflow Webpage

We created a webpage specifically for LC-MS users analyzing biomolecules, food and beverage, environmental samples and pharmaceutical products. It features many advanced techniques and technologies, including enantiomer separation, phospholipid removal, solid phase microextraction (SPME) and molecularly imprinted polymers (MIP).

Visit: [SigmaAldrich.com/lc-ms](https://www.sigmaaldrich.com/lc-ms)

Appendix I.

Common mass spectrometry contaminants and their sources

This list of potential interfering or contaminant ions in mass spectrometry (ESI positive mode, mass ≤ 1000 Da) is adapted from an excerpt of "Interferences and contaminants encountered in modern mass spectrometry" Bernd O. Keller, Jie Sui, Alex B. Young and Randy M. Whittal Analytica Chimica Acta 627, Issue 1, 3 October 2008, Pages 71-81. However, we have updated the masses listed in the previous publication by calculating the singly charged monoisotopic ion mass of each listed ion based on its molecular formula.

| Mono-isotopic ion mass (singly charged) | Ion type | Formula for M or subunit or sequence | Compound ID or species | Possible origin and other comments |
|---|--|---|------------------------------------|---|
| 33.0340 | [M+H] ⁺ | CH ₃ OH | Methanol | Acetonitrile, solvent |
| 42.0344 | [M+H] ⁺ | CH ₃ CN | ACN | Acetonitrile, solvent |
| 59.0609 | [M+NH ₄] ⁺ | CH ₃ CN | ACN | Acetonitrile, solvent |
| 63.0446 | [A ₁ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 64.0163 | [M+Na] ⁺ | CH ₃ CN | ACN | Acetonitrile, solvent |
| 65.0603 | [M ₂ +H] ⁺ | CH ₃ OH | Methanol | Methanol, solvent |
| 74.0606 | [M+H] ⁺ | C ₃ H ₇ NO | Dimethyl formamide | solvent |
| 74.0606 | [A ₁ B ₁ +H] ⁺ | (CH ₃ CN) _n (CH ₃ OH) _m | Acetonitrile/Methanol | ESI solvents |
| 77.0603 | [A ₁ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 79.0218 | [M+H] ⁺ | C ₂ H ₆ OS | DMSO | Dimethylsulfoxide, solvent |
| 83.0609 | [M ₂ +H] ⁺ | CH ₃ CN | Acetonitrile | ESI solvents |
| 85.0265 | [A ₁ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 85.0594 | [M+H] ⁺ | C ₂ D ₆ OS | d6-DMSO | d6-Dimethylsulfoxide, solvent |
| 88.0399 | [A ₁ B ₁ +H] ⁺ | (CH ₃ CN) _n (HCOOH) _m | Acetonitrile/Formic Acid | ESI solvents |
| 96.0425 | [A ₁ B ₁ +Na] ⁺ | (CH ₃ CN) _n (CH ₃ OH) _m | Acetonitrile/Methanol | ESI solvents |
| 99.0422 | [A ₁ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 101.0841 | [M+H] ⁺ | C ₅ H ₁₀ NO | NMP | N-methyl 2-pyrrolidone; solvent, floor stripper |
| 101.0005 | [A ₁ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 101.0037 | [M+Na] ⁺ | C ₂ H ₆ OS | DMSO | Dimethylsulfoxide, solvent |
| 101.0814 | [A ₂ B ₂ +H] ⁺ | [MeOH] _n [H ₂ O] _m | Methanol/Water | ESI solvents |
| 102.0555 | [A ₁ B ₁ +H] ⁺ | (CH ₃ CN) _n (CH ₃ COOH) _m | Acetonitrile/Acetic Acid | ESI solvents |
| 102.1283 | [M+H] ⁺ | C ₆ H ₁₅ N | TEA | Triethylamine, buffer |
| 103.9561 | [M+ ₆₃ Cu] ⁺ | C ₂ H ₃ N | ACN | Acetonitrile, solvent |
| 104.9928 | [M+Na] ⁺ | C ₂ H ₃ O ₂ Na | Sodium acetate | ESI solvents |
| 105.0429 | [M ₂ +Na] ⁺ | C ₂ H ₃ N | ACN | Acetonitrile, solvent |
| 41.0265 | [M+ ₆₅ Cu] ⁺ | C ₂ H ₃ N | ACN | Acetonitrile, solvent |
| 107.0708 | [A ₂ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 115.0161 | [A ₁ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 115.0871 | [A ₁ B ₁ +H] ⁺ | (CH ₃ CN) _n (C ₃ H ₇ NO) _m | Acetonitrile/ Dimethylformamide | solvent |
| 120.0483 | [M+CH ₃ CN+H] ⁺ | C ₂ H ₆ OS | DMSO | Dimethylsulfoxide, solvent |
| 122.0817 | [M+H] ⁺ | C ₄ H ₁₁ NO ₃ | TRIS | TRIS, buffer |
| 123.0633 | [A ₂ B ₂ +Na] ⁺ | [CH ₃ OH] _n [H ₂ O] _m | Methanol/Water | ESI solvents |
| 123.0922 | [M+H] ⁺ | C ₇ H ₁₀ N ₂ | DMAP | Dimethylaminopyridine, solvent |
| 124.0374 | [A ₁ B ₁ +Na] ⁺ | (CH ₃ CN) _n (CH ₃ COOH) _m | Acetonitrile/Acetic Acid | ESI solvents |
| 129.0528 | [A ₂ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 130.1596 | [M+H] ⁺ | C ₈ H ₁₉ N | DIPEA | Diisopropylethylamine, solvent |
| 132.9054 | M ⁺ | Cs | Cs-133 | Cesium, from Cesium Iodide used as calibrant |
| 133.1076 | [A ₃ B ₂ +H] ⁺ | [CH ₃ OH] _n [H ₂ O] _m | Methanol/Water | ESI solvents |
| 135.1021 | [A ₂ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 137.0749 | [M+CH ₃ CN+NH ₄] ⁺ | C ₂ H ₆ OS | DMSO | Dimethylsulfoxide, solvent |
| 142.0303 | [M+CH ₃ CN+Na] ⁺ | C ₂ H ₆ OS | DMSO | Dimethylsulfoxide, solvent |
| 144.1752 | [M+H] ⁺ | C ₉ H ₂₁ N | TPA | Tripropylamine, solvent |
| 144.9827 | [M ₂ + ₆₃ Cu] ⁺ | CH ₃ CN | ACN | Acetonitrile, solvent, together with m/z 147 |
| 145.0267 | [A ₂ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 146.0694 | [M ₃ +Na] ⁺ | CH ₃ CN | ACN | Acetonitrile, solvent |
| 146.9809 | [M ₂ + ₆₅ Cu] ⁺ | CH ₃ CN | ACN | Acetonitrile, solvent, together with m/z 145 |
| 147.1134 | [A ₂ B ₂ +H] ⁺ | (CH ₃ CN) _n (CH ₃ OH) _m | Acetonitrile/Methanol | ESI solvents |
| 149.0239 | [F+H] ⁺ | C ₈ H ₄ O ₃ | Phthalic Anhydride | fragment ion originating from phthalate esters |

| Mono-isotopic ion mass (singly charged) | Ion type | Formula for M or subunit or sequence | Compound ID or species | Possible origin and other comments |
|---|---|---|------------------------|---|
| 150.1283 | [M+H] ⁺ | C ₁₀ H ₁₅ N | Phenyldiethylamine | solvent |
| 151.0970 | [A ₃ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 153.1392 | [M+H] ⁺ | C ₉ H ₁₆ N ₂ | DBU | 1,8-Diazabicyclo[5.4.0]undec-7-ene |
| 155.0895 | [A ₃ B ₂ +Na] ⁺ | [CH ₃ OH] _n [H ₂ O] _m | Methanol/Water | ESI solvents |
| 157.0357 | [M ₂ +H] ⁺ | C ₂ H ₆ OS | DMSO | Dimethylsulfoxide, solvent |
| 157.0841 | [A ₂ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 158.9646 | [M+Na] ⁺ | C ₂ F ₃ O ₂ Na | NaTFA | Sodium trifluoroacetate, salt |
| 163.0395 | [M-CH ₃ OH+H] ⁺ | C ₁₀ H ₁₀ O ₄ | Dimethyl phthalate | Phthalate esters, plasticizer |
| 163.1334 | [M+H] ⁺ | C ₈ H ₁₈ O ₃ | DGBE | Diethylene glycol monobutyl ether, cpd. In scintillation cocktail |
| 169.0953 | [A ₂ B ₂ +Na] ⁺ | (CH ₃ CN) _n (CH ₃ OH) _m | Acetonitrile/Methanol | ESI solvents |
| 170.1188 | [M ₂ +H] ⁺ | C ₂ D ₆ OS | d6-DMSO | d6-Dimethylsulfoxide, solvent |
| 171.0058 | [F+Na] ⁺ | C ₆ H ₄ O ₃ | Phthalic anhydride | from phthalate esters, plasticizer |
| 173.0580 | [A ₂ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 173.0790 | [A ₃ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 179.0176 | [M ₂ +Na] ⁺ | C ₂ H ₆ OS | DMSO | Dimethylsulfoxide, solvent |
| 181.1229 | [M+H] ⁺ | C ₁₁ H ₁₆ O ₂ | BHA | Butylated hydroxyanisole, antioxidant additives |
| 183.0810 | [M+H] ⁺ | C ₁₃ H ₁₀ O | DPK | Diphenyl ketone |
| 183.1444 | [A ₄ B ₃ +H] ⁺ | [CH ₃ OH] _n [H ₂ O] _m | Methanol/Water | ESI solvents |
| 185.1154 | [M+Na] ⁺ | C ₆ H ₁₈ O ₃ | GE | glycol ether |
| 186.2222 | [M+H] ⁺ | C ₁₂ H ₂₇ N | TBA | Tributylamine, solvent |
| 189.0529 | [A ₃ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 193.1440 | [A ₃ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 195.0657 | [M+H] ⁺ | C ₁₀ H ₁₀ O ₄ | Dimethyl phthalate | Phthalate esters, plasticizer |
| 195.1232 | [A ₄ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 203.1048 | [M+Na] ⁺ | C ₁₁ H ₁₆ O ₂ | BHA | Butylated hydroxyanisole, antioxidant additives |
| 205.1263 | [A ₄ B ₃ +Na] ⁺ | [CH ₃ OH] _n [H ₂ O] _m | Methanol/Water | ESI solvents |
| 214.0902 | [M+H] ⁺ | C ₁₀ H ₁₅ NO ₂ S | n-BBS | n-butyl benzenesulfonamide, plasticizer |
| 215.1259 | [A ₃ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 217.1052 | [A ₄ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 221.1905 | [M+H] ⁺ | C ₁₅ H ₂₄ O | BTH | Butylated hydroxytoluene, Antioxidant |
| 225.1967 | [M+H] ⁺ | C ₁₃ H ₂₄ N ₂ O | DCU | N,N'-Dicyclohexylurea |
| 231.0999 | [A ₃ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 231.1167 | [M+NH ₄] ⁺ | C ₁₀ H ₁₅ NO ₂ S | n-BBS | n-butyl benzenesulfonamide, plasticizer |
| 233.0791 | [A ₄ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 236.0721 | [M+Na] ⁺ | C ₁₀ H ₁₅ NO ₂ S | n-BBS | n-butyl benzenesulfonamide, plasticizer |
| 239.1495 | [A ₅ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 239.2254 | [(M ₃ .H ₃₅ .Cl) ₂ -Cl] ⁺ | C ₆ H ₁₅ N | TEA.HCl | Triethylamine-hydrochloride, buffer |
| 241.2224 | [(M ₃ .Cl) ₂ -Cl] ⁺ | C ₆ H ₁₅ N | TEA.HCl | Triethylamine-hydrochloride, buffer |
| 242.2848 | M ⁺ | C ₁₆ H ₃₆ N | TBA | Tetrabutylammonium, buffer |
| 243.1174 | M ⁺ | C ₁₉ H ₁₅ | Trityl cation | Trityl cation, [Ph ₃ C] ⁺ |
| 243.1725 | [M+Na] ⁺ | C ₁₅ H ₂₄ O | BTH | Butylated hydroxytoluene, Antioxidant additives |
| 251.1858 | [A ₄ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 251.2011 | [AB ₁ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 257.0316 | [M ₃ +Na] ⁺ | C ₂ H ₆ OS | DMSO | Dimethylsulfoxide, solvent |
| 261.1314 | [A ₅ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 265.2168 | [AB ₁ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 267.1725 | [M+H] ⁺ | C ₁₂ H ₂₇ O ₄ P | TBP | Tributylphosphate |
| 273.1279 | M ⁺ | C ₂₀ H ₁₇ O | MMT | Monomethoxytrityl cation |
| 273.1678 | [A ₄ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 273.1830 | [AB ₁ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 277.1053 | [A ₅ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 279.0939 | [M+H] ⁺ | C ₁₈ H ₁₅ OP | TPO | Triphenylphosphine oxide |
| 279.1596 | [M+H] ⁺ | C ₁₆ H ₂₂ O ₄ | Dibutylphthalate | Plasticiser, phtalate ester |
| 279.2300 | [AB ₁ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 281.0517 | [M+H-CH ₄] ⁺ | [C ₂ H ₆ SiO] ₄ | Polysiloxane | Polysiloxane, (neutral methane loss from m/z 297) |
| 282.2797 | [M+H] ⁺ | C ₁₈ H ₃₅ NO | Oleamide | Slip agent in polyethylene films |
| 283.1757 | [A ₆ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 284.2953 | [M+H] ⁺ | C ₁₈ H ₃₇ NO | Stearamide | Slip agent in polyethylene films |
| 287.1987 | [AB ₁ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 288.2539 | [M+H] ⁺ | C ₁₆ H ₃₃ NO ₃ | n,n-DDA | n,n-bis(2-hydroxyethyl) dodecanamide, anti-static agent |
| 289.1417 | [A ₄ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |

| Mono-isotopic ion mass (singly charged) | Ion type | Formula for M or subunit or sequence | Compound ID or species | Possible origin and other comments |
|---|---|---|-------------------------------|---|
| 293.2457 | [AB ₁ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton [®] , reduced | 101R Detergents |
| 295.2273 | [AB ₂ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton [®] | X-100, X-114, X-405, or X-45 Detergents |
| 297.0830 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₄ | Polysiloxane | Polysiloxane, followed by m/z |
| 301.1416 | [M+Na] ⁺ | C ₁₆ H ₂₂ O ₄ | Dibutylphthalate | Dibutylphthalate, plasticizer |
| 304.2616 | [M+Na] ⁺ | C ₁₈ H ₃₅ NO | Oleamide | Slip agent in polyethylene films |
| 305.1576 | [A ₆ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 306.2773 | [M+Na] ⁺ | C ₁₈ H ₃₇ NO | Stearamide | Slip agent in polyethylene films |
| 309.2277 | [A ₅ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 309.2430 | [AB ₂ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton [®] | 101 Detergents |
| 315.2535 | [M+H] ⁺ | C ₁₈ H ₃₄ O ₄ | DBS | Dibutyl sebacate, plasticizer |
| 317.1155 | [M+K] ⁺ | C ₁₆ H ₂₂ O ₄ | Dibutylphthalate | Dibutylphthalate, plasticizer |
| 317.2093 | [AB ₂ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton [®] | X-100, X-114, X-405, or X-45 Detergents |
| 321.1316 | [A ₆ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 323.2562 | [AB ₂ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton [®] , reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 325.2590 | [M ₂ +H] ⁺ | C ₈ H ₁₈ O ₃ | DGBE | Diethylene glycol monobutyl ether, cpd. In scintillation cocktail |
| 327.0786 | [M+H] ⁺ | C ₁₈ H ₁₅ O ₄ P | TPP | Triphenyl phosphate, flame retardant in plastics |
| 327.2019 | [A ₇ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 331.2097 | [A ₅ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 331.2249 | [AB ₂ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton [®] | 101 Detergents |
| 337.1190 | [M+H] ⁺ ; (₁₂₀ Sn) | C ₁₃ H ₂₈ O ₂ Sn | Tributyl tin formate | Tributyl tin formate, catalyst |
| 337.2719 | [AB ₂ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton [®] , reduced | 101R Detergents |
| 338.3423 | [M+H] ⁺ | C ₂₂ H ₄₃ NO | Erucamide | Erucamide, (Cis-13-docosenoic amide) |
| 339.2535 | [AB ₃ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton [®] | X-100, X-114, X-405, or X-45 Detergents |
| 347.1836 | [A ₅ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 349.1838 | [A ₇ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 353.2692 | [AB ₃ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton [®] | 101 Detergents |
| 355.0705 | [M+H-CH ₄] ⁺ | [C ₂ H ₆ SiO] ₅ | Polysiloxane | Polysiloxane, (neutral methane loss from m/z 371) |
| 355.3688 | [M-Cl] ⁺ | C ₂₂ H ₄₇ N ₂ OCl | PATC | Palmitamidopropyl-trimonium chloride, personal care products additive |
| 360.3242 | [M+Na] ⁺ | C ₂₂ H ₄₃ NO | Erucamide | Erucamide, (Cis-13-docosenoic amide) |
| 361.2355 | [AB ₃ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton [®] | X-100, X-114, X-405, or X-45 Detergents |
| 365.1578 | [A ₇ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 367.2696 | [A ₆ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 367.2824 | [AB ₃ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton [®] , reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 368.4256 | [M-Cl] ⁺ | C ₂₅ H ₅₄ NCl | BTAC-228 | Behentrimonium chloride, personal care product additive |
| 371.1018 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₅ | Polysiloxane | Polysiloxane, followed by m/z 388 |
| 371.2281 | [A ₈ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 371.3161 | [M+H] ⁺ | C ₂₂ H ₄₂ O ₄ | DEHA | Bis(2-ethylhexyl) adipate, plasticizer |
| 371.3161 | [M+H] ⁺ | C ₂₂ H ₄₂ O ₄ | DOA | Diocetyl adipate, plasticizer |
| 375.2511 | [AB ₃ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton [®] | 101 Detergents |
| 381.2981 | [AB ₃ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton [®] , reduced | 101R Detergents |
| 383.2797 | [AB ₄ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton [®] | X-100, X-114, X-405, or X-45 Detergents |
| 388.1284 | [M+NH ₄] ⁺ | [C ₂ H ₆ SiO] ₅ | Polysiloxane | Polysiloxane, (see m/z 371) |
| 389.2515 | [A ₆ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 391.2848 | [M+H] ⁺ | C ₂₄ H ₃₈ O ₄ | Diisooctyl phthalate | Diisooctyl phthalate, plasticiser |
| 393.2101 | [A ₈ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 397.2954 | [AB ₄ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton [®] | 101 Detergents |
| 405.2255 | [A ₆ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 405.2617 | [AB ₄ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton [®] | X-100, X-114, X-405, or X-45 Detergents |
| 409.1840 | [A ₈ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 411.3086 | [AB ₄ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton [®] , reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 413.2668 | [M+Na] ⁺ | C ₂₄ H ₃₈ O ₄ | Diisooctyl phthalate | Diisooctyl phthalate, plasticiser |
| 415.2543 | [A ₉ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 419.2773 | [AB ₄ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton [®] | 101 Detergents |
| 425.3114 | [A ₇ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 425.3243 | [AB ₄ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton [®] , reduced | 101R Detergents |
| 427.3060 | [AB ₅ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton [®] | X-100, X-114, X-405, or X-45 Detergents |
| 429.0893 | [M+H-CH ₄] ⁺ | [C ₂ H ₆ SiO] ₆ | Polysiloxane | Polysiloxane, (neutral methane loss from m/z 445) |
| 429.2407 | [M+K] ⁺ | C ₂₄ H ₃₈ O ₄ | Diisooctyl phthalate | Diisooctyl phthalate, plasticiser |
| 437.2363 | [A ₉ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 441.3216 | [AB ₅ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton [®] | 101 Detergents |
| 445.1206 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₆ | Polysiloxane | Polysiloxane, followed by m/z 462 |

| Mono-isotopic ion mass (singly charged) | Ion type | Formula for M or subunit or sequence | Compound ID or species | Possible origin and other comments |
|---|--|---|--------------------------|--|
| 447.2934 | [A ₇ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 449.2879 | [AB ₅ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 449.3856 | [M ₂ +H] ⁺ | C ₁₃ H ₂₄ N ₂ O | DCU | N,N'-Dicyclohexylurea |
| 453.2102 | [A ₉ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 453.3441 | [M+H] ⁺ | C ₂₄ H ₄₄ N ₄ O ₄ | nylon | Cyclic oligomer of polyamide 66, (adipic acid/hexamethylene diamine condensation) |
| 454.2933 | [M+CH ₃ CN+Na] ⁺ | C ₂₄ H ₃₈ O ₄ | Diisooctyl phthalate | Diisooctyl phthalate, plasticiser |
| 455.3349 | [AB ₅ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 459.2805 | [A ₁₀ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 462.1471 | [M+NH ₄] ⁺ | [C ₂ H ₆ SiO] ₆ | Polysiloxane | Polysiloxane (see m/z 445) |
| 463.2673 | [A ₇ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 463.3036 | [AB ₅ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 469.3505 | [AB ₅ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 471.3322 | [AB ₆ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 481.2625 | [A ₁₀ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 483.3533 | [A ₈ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 485.3478 | [AB ₆ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 493.3141 | [AB ₆ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 494.5665 | [M-Cl] ⁺ | C ₃₄ H ₇₂ NCl | DPDMA | Dipalmityltrimethylammonium chloride, catalyst, personal care products additive |
| 497.2364 | [A ₁₀ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 499.3611 | [AB ₆ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 484.3376 | [M+H-CH ₄] ⁺ | [C ₂ H ₆ SiO] ₇ | Polysiloxane | Polysiloxane, (neutral methane loss from m/z 519) |
| 503.3068 | [A ₁₁ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 505.3353 | [A ₈ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 507.3298 | [AB ₆ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 513.3767 | [AB ₆ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 515.3584 | [AB ₇ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 515.4134 | [M+H] ⁺ | C ₃₀ H ₅₈ O ₄ S | DDTDP | Didodecyl 3,3'-thiodipropionate, antioxidant |
| 519.1394 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₇ | Polysiloxane | Polysiloxane, followed by m/z 536 |
| 521.3092 | [A ₈ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 522.5978 | [M-Cl] ⁺ | C ₃₆ H ₇₆ NCl | SPDMA | Stearyl-palmityltrimethylammonium chloride, catalyst, personal care product additive |
| 525.2887 | [A ₁₁ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 529.3740 | [AB ₇ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 531.4083 | [M+H] ⁺ | C ₃₀ H ₅₈ O ₅ S | DDTDP | Didodecyl 3,3'-thiodipropionate oxidized to sulfoxide, antioxidant |
| 531.4777 | [M+H] ⁺ | C ₃₅ H ₆₂ O ₃ | Irganox | Irganox 1076, antioxidant in synthetic polymers, antioxidant |
| 536.1659 | [M+NH ₄] ⁺ | [C ₂ H ₆ SiO] ₇ | Polysiloxane | Polysiloxane (see m/z 519) |
| 537.3403 | [AB ₇ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 537.8796 | [M ₆ -6H+ ₃ Fe+O] ⁺ | C ₂ H ₄ O ₂ | Acetic acid-Fe-O-complex | during ESI with metal tips and acetic acid complex |
| 541.2626 | [A ₁₁ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 541.3952 | [A ₉ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 543.3873 | [AB ₇ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 547.3330 | [A ₁₂ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 547.4032 | [M+H] ⁺ | C ₃₀ H ₅₈ O ₆ S | DDTDP | Didodecyl 3,3'-thiodipropionate oxidized to sulfone, antioxidant |
| 550.6291 | [M-Cl] ⁺ | C ₃₈ H ₈₀ NCl | DSDMA | Distearyltrimethylammonium chloride, catalyst, personal care products additive |
| 551.3560 | [AB ₇ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 553.3903 | [M+Na] ⁺ | C ₃₀ H ₅₈ O ₅ S | DDTDP | Didodecyl 3,3'-thiodipropionate oxidized to sulfoxide, antioxidant |
| 553.4597 | [M+Na] ⁺ | C ₃₅ H ₆₂ O ₃ | Irganox | Irganox 1076, antioxidant in synthetic polymers, antioxidant |
| 555.8902 | [M ₆ -6H+H ₂ O+3Fe+O] ⁺ | C ₂ H ₄ O ₂ | Acetic acid-Fe-O-complex | during ESI with metal tips and acetic acid complex |
| 557.4029 | [AB ₇ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 559.3846 | [AB ₈ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 563.3771 | [A ₉ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 569.3149 | [A ₁₂ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 573.4003 | [AB ₈ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 577.1269 | [M+H-CH ₄] ⁺ | [C ₂ H ₆ SiO] ₈ | Polysiloxane | Polysiloxane, (neutral methane loss from m/z 593) |
| 579.3511 | [A ₉ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 581.3666 | [AB ₈ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |

| Mono-isotopic ion mass (singly charged) | Ion type | Formula for M or subunit or sequence | Compound ID or species | Possible origin and other comments |
|---|---|---|-----------------------------|---|
| 585.2888 | [A ₁₂ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 587.4135 | [AB ₈ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 591.3592 | [A ₁₃ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 593.1582 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₈ | Polysiloxane | Polysiloxane, followed by m/z 610 |
| 595.3822 | [AB ₈ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 597.9007 | [M ₇ -6H+3Fe+O] ⁺ | C ₂ H ₄ O ₂ | Acetic acid-Fe-O-complex | during ESI with metal tips and acetic acid |
| 599.4370 | [A ₁₀ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 601.4292 | [AB ₈ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 603.4108 | [AB ₉ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 610.1847 | [M+NH ₄] ⁺ | [C ₂ H ₆ SiO] ₈ | Polysiloxane | Polysiloxane (see m/z 593) |
| 613.3411 | [A ₁₃ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 615.4043 | [M+H] ⁺ | C ₃₂ H ₅₈ N ₂ O ₇ S | CHAPS | 3-[(3-Cholamidopropyl)dimethylammonio]-1-propanesulfonate |
| 617.4265 | [AB ₉ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 621.4190 | [A ₁₀ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 621.9735 | [M ₆ -6H+3Fe+O] ⁺ | C ₃ H ₆ O ₂ | Propionic acid Fe-O complex | during ESI with metal tips and acetic acid |
| 625.3928 | [AB ₉ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 629.3151 | [A ₁₃ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 631.4397 | [AB ₉ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 635.3854 | [A ₁₄ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 637.3929 | [A ₁₀ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 639.4084 | [AB ₉ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 645.4554 | [AB ₉ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 647.4370 | [AB ₁₀ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 651.1456 | [M+H-CH ₄] ⁺ | [C ₂ H ₆ SiO] ₉ | Polysiloxane | Polysiloxane, (neutral methane loss from m/z 667) |
| 657.3673 | [A ₁₄ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 657.4789 | [A ₁₁ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 661.4527 | [AB ₁₀ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 667.1769 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₉ | Polysiloxane | Polysiloxane, followed by m/z 684 |
| 669.4190 | [AB ₁₀ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 673.3413 | [A ₁₄ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 675.4659 | [AB ₁₀ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 679.4116 | [A ₁₅ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 679.4608 | [A ₁₁ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 679.5122 | [M+H] ⁺ | C ₃₆ H ₆₆ N ₆ O ₆ | nylon | Cyclic oligomer of polyamide 66, (adipic acid/hexamethylene diamine condensation) |
| 683.4346 | [AB ₁₀ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 684.2035 | [M+NH ₄] ⁺ | [C ₂ H ₆ SiO] ₉ | Polysiloxane | Polysiloxane (see m/z 667) |
| 689.4816 | [AB ₁₀ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton, reduced | 101R Detergents |
| 691.4633 | [AB ₁₁ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 695.4348 | [A ₁₁ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 701.3936 | [A ₁₅ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 705.4789 | [AB ₁₁ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 713.4452 | [AB ₁₁ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 715.5208 | [A ₁₂ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 717.3675 | [A ₁₅ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 719.4921 | [AB ₁₁ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 723.4378 | [A ₁₆ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 725.1644 | [M+H-CH ₄] ⁺ | [C ₂ H ₆ SiO] ₁₀ | Polysiloxane | Polysiloxane, (neutral methane loss from m/z 741) |
| 727.4608 | [AB ₁₁ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 733.5078 | [AB ₁₁ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 735.4895 | [AB ₁₂ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 737.5027 | [A ₁₂ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 741.1957 | [M+H] ⁺ | [C ₂ H ₆ SiO] ₁₀ | Polysiloxane | Polysiloxane, followed by m/z 758 |
| 745.4198 | [A ₁₆ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 749.5051 | [AB ₁₂ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 753.4766 | [A ₁₂ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 757.4714 | [AB ₁₂ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 758.2223 | [M+NH ₄] ⁺ | [C ₂ H ₆ SiO] ₁₀ | Polysiloxane | Polysiloxane (see m/z 741) |
| 761.3937 | [A ₁₆ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 763.5184 | [AB ₁₂ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 767.4640 | [A ₁₇ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |

| Mono-isotopic ion mass (singly charged) | Ion type | Formula for M or subunit or sequence | Compound ID or species | Possible origin and other comments |
|---|---|---|------------------------|---|
| 771.4871 | [AB ₁₂ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 773.5626 | [A ₁₃ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 777.5340 | [AB ₁₂ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 779.5157 | [AB ₁₃ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 789.4460 | [A ₁₇ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 793.5313 | [AB ₁₃ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 795.5446 | [A ₁₃ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 798.5884 | [M ₂ +NH ₄] ⁺ | C ₂₄ H ₃₈ O ₄ | Diisooctyl phthalate | Diisooctyl phthalate, plasticiser |
| 801.4976 | [AB ₁₃ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 803.5438 | [M ₂ +Na] ⁺ | C ₂₄ H ₃₈ O ₄ | Diisooctyl phthalate | Diisooctyl phthalate, plasticiser |
| 805.4199 | [A ₁₇ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 807.5446 | [AB ₁₃ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 809.4875 | [AB ₁₀ +Na] ⁺ | [C ₁₈ H ₃₄ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 20 |
| 811.4903 | [A ₁₈ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 811.5185 | [A ₁₃ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 815.5133 | [AB ₁₃ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 819.5177 | [M ₂ +K] ⁺ | C ₂₄ H ₃₈ O ₄ | Diisooctyl phthalate | Diisooctyl phthalate, plasticiser |
| 821.5602 | [AB ₁₃ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 823.5419 | [AB ₁₄ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 831.6045 | [A ₁₄ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 833.4722 | [A ₁₈ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 837.5575 | [AB ₁₄ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 845.5238 | [AB ₁₄ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 849.4461 | [A ₁₈ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 851.5708 | [AB ₁₄ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 853.5137 | [AB ₁₁ +Na] ⁺ | [C ₁₈ H ₃₄ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 20 |
| 853.5864 | [A ₁₄ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 855.5165 | [A ₁₉ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 859.5395 | [AB ₁₄ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 865.5501 | [AB ₁₀ +Na] ⁺ | [C ₂₂ H ₄₂ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 40 |
| 865.5864 | [AB ₁₄ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 867.5681 | [AB ₁₅ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 869.5604 | [A ₁₄ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 877.4984 | [A ₁₉ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 881.5838 | [AB ₁₅ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 889.5501 | [AB ₁₅ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 889.6464 | [A ₁₅ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 891.5657 | [AB ₁₀ +Na] ⁺ | [C ₂₄ H ₄₄ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 80 |
| 893.4724 | [A ₁₉ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 893.5814 | [AB ₁₀ +Na] ⁺ | [C ₂₄ H ₄₆ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 60 |
| 897.5399 | [AB ₁₂ +Na] ⁺ | [C ₁₈ H ₃₄ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 20 |
| 899.5427 | [A ₂₀ B+H] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 903.5657 | [AB ₁₅ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 905.6803 | [M+H] ⁺ | C ₄₈ H ₈₈ N ₆ O ₈ | nylon | Cyclic oligomer of polyamide 66, (adipic acid/hexamethylene diamine condensation) |
| 909.5763 | [AB ₁₁ +Na] ⁺ | [C ₂₂ H ₄₂ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 40 |
| 909.6127 | [AB ₁₅ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 911.5943 | [AB ₁₆ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 911.6283 | [A ₁₅ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 921.5246 | [A ₂₀ B+Na] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 925.6100 | [AB ₁₆ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 927.6022 | [A ₁₅ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 933.5763 | [AB ₁₆ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 935.5919 | [AB ₁₁ +Na] ⁺ | [C ₂₄ H ₄₄ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 80 |
| 937.4986 | [A ₂₀ B+K] ⁺ | [C ₂ H ₄ O] _n H ₂ O | PEG | Polyethylene glycol, ubiquitous polyether |
| 937.6076 | [AB ₁₁ +Na] ⁺ | [C ₂₄ H ₄₆ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 60 |
| 939.6232 | [AB ₁₆ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton®, reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 941.5661 | [AB ₁₃ +Na] ⁺ | [C ₁₈ H ₃₄ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 20 |
| 947.5919 | [AB ₁₆ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |
| 947.6882 | [A ₁₆ B+H] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 953.6025 | [AB ₁₂ +Na] ⁺ | [C ₂₂ H ₄₂ O ₆][C ₂ H ₄ O] _n | Tween® | Tween® 40 |
| 953.6389 | [AB ₁₆ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton®, reduced | 101R Detergents |
| 955.6205 | [AB ₁₇ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton® | X-100, X-114, X-405, or X-45 Detergents |
| 969.6362 | [AB ₁₇ +H] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton® | 101 Detergents |

| Mono-isotopic ion mass (singly charged) | Ion type | Formula for M or subunit or sequence | Compound ID or species | Possible origin and other comments |
|---|-------------------------------------|---|-------------------------------|---|
| 969.6702 | [A ₁₆ B+Na] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 977.6025 | [AB ₁₇ +Na] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton [®] | X-100, X-114, X-405, or X-45 Detergents |
| 979.6181 | [AB ₁₂ +Na] ⁺ | [C ₂₄ H ₄₄ O ₆][C ₂ H ₄ O] _n | Tween [®] | Tween [®] 80 |
| 981.6338 | [AB ₁₂ +Na] ⁺ | [C ₂₄ H ₄₆ O ₆][C ₂ H ₄ O] _n | Tween [®] | Tween [®] 60 |
| 983.6494 | [AB ₁₇ +Na] ⁺ | [C ₁₄ H ₂₈ O][C ₂ H ₄ O] _n | Triton [®] , reduced | X-100R, X-114R, X-405R, or X-45R Detergents |
| 985.5923 | [AB ₁₄ +Na] ⁺ | [C ₁₈ H ₃₄ O ₆][C ₂ H ₄ O] _n | Tween [®] | Tween [®] 20 |
| 985.6441 | [A ₁₆ B+K] ⁺ | [C ₃ H ₆ O] _n H ₂ O | PPG | Polypropylene glycol, ubiquitous polyether |
| 991.6181 | [AB ₁₇ +Na] ⁺ | [C ₁₅ H ₂₄ O][C ₂ H ₄ O] _n | Triton [®] | 101 Detergents |
| 997.6287 | [AB ₁₃ +Na] ⁺ | [C ₂₂ H ₄₂ O ₆][C ₂ H ₄ O] _n | Tween [®] | Tween [®] 40 |
| 997.6651 | [AB ₁₇ +Na] ⁺ | [C ₁₅ H ₃₀ O][C ₂ H ₄ O] _n | Triton [®] , reduced | 101R Detergents |
| 999.6470 | [AB ₁₈ +H] ⁺ | [C ₁₄ H ₂₂ O][C ₂ H ₄ O] _n | Triton [®] | X-100, X-114, X-405, or X-45 Detergents |

Appendix II.

Monoisotopic ion masses of commonly observed repeating units in LC-MS

Positive ion

| Mass difference | Origin |
|-----------------|--|
| 14.0157 | -[CH ₂]-, alkane chains, waxes, fatty acids, methylation |
| 15.9949 | O, oxidation |
| 18.0106 | H ₂ O, water clusters |
| 28.0313 | -[C ₂ H ₄]-, natural alkane chains such as fatty acids |
| 32.0262 | CH ₃ OH, methanol clusters |
| 41.0266 | CH ₃ CN, acetonitrile clusters |
| 42.0470 | -[C ₃ H ₆]-, propyl repeating units, propylation |
| 44.0262 | -[C ₂ H ₄ O]-; polyethylene glycol, PEG, and related components such as Triton [®] - and Tween [®] -containing buffers |
| 49.9968 | -[CF ₂]-, from perfluoro compounds |
| 53.0032 | NH ₄ Cl salt adducts/clusters |
| 56.0626 | -[C ₄ H ₈]-, butyl repeating units, butylation |
| 57.9586 | NaCl, sodium chloride clusters |
| 58.0419 | -[C ₃ H ₆ O]-; polypropylene glycol and related compounds, PPG, and related compounds |
| 63.0320 | CHOONH ₄ , ammonium formate adducts/clusters |
| 67.9874 | NaHCO ₂ , sodium formate clusters |
| 67.9874 | CHOONa, sodium formate adducts/clusters |
| 72.0395 | -OH replacement with -OSi(CH ₃) ₃ , (= [C ₃ H ₆ Si]), trimethylsiloxane, endcapping reagent |
| 73.9326 | KCl adducts/clusters |
| 74.0188 | -[O-Si(CH ₃) ₂]-, polysiloxane, silicone rubber polymer (typical series at m/z's 355, 429, 503, 593, 667, 741, 815) |
| 78.0139 | C ₂ H ₆ OS, DMSO adducts/clusters, dimethylsulfoxide solvent |
| 82.0031 | NaCH ₃ CO ₂ , sodium acetate clusters |
| 84.0516 | C ₂ D ₆ OS, deuterated DMSO adducts/clusters, NMR solvent |
| 135.9748 | NaCF ₃ CO ₂ , sodium trifluoroacetate clusters |
| 162.0528 | -[C ₆ H ₁₀ O ₅]-, polysaccharides residues |
| 226.1681 | -[C ₁₂ H ₂₂ N ₂ O ₂]-, cyclic oligomers from polyamide 66 (series observed with m/z 453, 679, 905) |
| 259.8099 | CsI, cesium iodide clusters, used as calibration |

Supelco®

Analytical Products

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

SigmaAldrich.com

To place an order or receive technical assistance

Order/Customer Service: [SigmaAldrich.com/order](https://www.sigmaaldrich.com/order)

Technical Service: [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice)

Safety-related Information: [SigmaAldrich.com/safetycenter](https://www.sigmaaldrich.com/safetycenter)

© 2020 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Supelco, Chromolith, Durapore, Suprapur, EMSURE, Millex, Millipore Express, Smplicity, SeQuant, ZipTip, ZIC-HILIC, Titan, LiChrosolv, LiChropur, Ascentis, Discovery, HybridSPE, BioShell, Supel, Milli-Q, LC-Pak, Purospher, UniSolv, Multiscreen, Ultrafree, Centrifree, Ultracel, EXTrelut, and LiChrolut are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK_BR6428EN Ver. 1.0
32480
09/2020