


purification perfection

Quality materials for absorption,
adsorption & filtration



The life science
business of Merck
operates as
MilliporeSigma in
the U.S. and Canada.

Supelco[®]
Analytical Products

purification perfection

Purification is one of the most important applications in analytical laboratories. To ease your daily work, we offer a complete range of absorption and adsorption reagents, as well as filtration and clarification materials – all with excellent take-up properties. Our products are suitable for a wide variety of applications, such as absorbing or binding substances, as well as for decolorization, clarification and filtration. Regardless of the purpose, they deliver quality perfected for your intended use.



Your benefits:



Reliability:

All natural products used are tested for organic impurities, and various anions and cations. The products are specified and offer excellent **batch-to-batch consistency**.



Convenience:

Comprehensive portfolio allows successful implementation of a wide variety of purification methods.



Variety:

Our products are available in **different grades** to meet **individual application requirements**.



Flexibility:

Products are available in various pack sizes to suit all our customers – from small labs to large testing or production facilities.



Sustainability:

Most of our absorption, adsorption and filtration materials are natural reagents which are **not harmful to the environment**.

Overview of methods

Adsorption [from Latin "adsorbere": to add, to attach] describes the take-up of a substance on the surface of a solid without a chemical reaction. Due to adhesive forces, gaseous, liquid and, in rare cases, solid substances accumulate on the surface or interface of the adsorbent. This releases the bond energy in the form of heat. The larger the surface, the greater the adsorption capacity of a substance. The specific surface (in "m²/g") is therefore often given for specification of the adsorbent. Adsorption strength is reduced by warming and increased by cooling.

Absorption [from Latin "absorbere": to devour, swallow up] describes a physical or chemical process of taking-up molecules or atoms by solid or liquid materials. This process differs from adsorption as molecules are taken up by the volume, not by the surface.

Filtration is a mechanical separation process for purification of substances. The filter and clarification materials used for this are usually reagents with absorbent or adsorbent properties, although materials which work in a similar way to frits or which form capillaries as filter cake, can also be used to clarify solutions, and decolorize or remove solid or colloiddally dissolved substances.

FROM A TO Z

Quality materials for absorption, adsorption & filtration at a glance

Absorption tubes	5
Aluminium oxide fibers	6
Calcium oxide	6
Charcoal activated	7
Charcoal wood	7
Glass wool	8
Graphite	8
Molecular sieves	9
Sea sand	11
Sodalime	11

Explore our complete range on:
SigmaAldrich.com/absorption-filtration



Absorption tubes

Application advice

Drying of incoming or outgoing air or general gas absorption

Additional information

Absorption tubes are 15 cm long and have a diameter of 2 cm. Hoses of various dimensions can also be connected: hoses with 3 mm external diameter can be pushed into the tube, while hoses with 8–12 cm internal diameter can be pushed over the opening. The tube must be opened on both sides before use by removing the two stoppers.

Absorption tube for H₂O Application advice

H₂O absorption tubes are used to keep reagents or solvents dry, for example during Karl Fischer titration or when working with solvent from withdrawal systems.

Additional information

In H₂O absorption tubes, the filling also contains an indicator displaying the degree of saturation.

Absorption tube for CO₂ Application advice

CO₂ absorption tubes are used to bond carbon dioxide, for example to keep the titer stable during titration with NaOH and to prevent the formation of Na₂CO₃ through the CO₂ contained in the incoming air. The absorbent, which is placed on an inert carrier, also prevents caking after absorption. The air or gas can still flow through the absorption tube without a change in resistance.

Additional information

CO₂ absorption tubes are filled with sodium hydroxide on an inert carrier material. The reagent in these absorption tubes has a significantly higher capacity than caustic alkali pellets.

Product	CAS-No.	Content	Packaging	Cat. No.
Absorption tubes				
Absorption tube for H ₂ O (molecular sieve 0.3 nm with indicator)	-	3 units	Plastic can	1.06107.0003
Absorption tubes for CO ₂ (sodium hydroxide on support)	-	3 units	Plastic can	1.01562.0003



Aluminium oxide fibers

Application advice

Due to their high melting point of 1700 to 1800 °C, aluminium oxide fibers are ideal for use as filtration and insulation agents.

Information on sustainable protection

In contrast to asbestos fibers, our aluminium oxide fibers are non-carcinogenic and can be easily disposed of.

Product	CAS-No.	Content	Packaging	Cat. No.
Aluminium oxide fibers				
Aluminium oxide fibers for gooch crucibles	142844-00-6	100 g	Plastic bottle	1.15754.0100
		1 kg	Plastic bottle	1.15754.1000



Calcium oxide

Application advice

Calcium oxide is used as a drying agent and absorbent for CO₂ and for the creation of calcium bonds, e.g. in the production of glass, paper, rubber and soda.

Information on sustainable protection

Calcium oxide is obtained from selected natural marble.

Product	CAS-No.	Content	Packaging	Cat. No.
Calcium oxide				
Calcium oxide from small marble lumps ~ 3-20 mm	1305-78-8	1 kg	Plastic bottle	1.02109.1000
		25 kg	Fibre carton	1.02109.9025

Charcoal activated

Application advice

Activated charcoal is used for a wide range of applications, such as adsorption, decolorization and purification of gases and liquids, and as a carrier for catalysts. Thanks to its very large internal surface area (800-1500 m²/g), activated charcoal is able to absorb molecules of diverse sizes. Since chemical substances must diffuse into the pores during adsorption, powdered carbon is easier to 'load'. In contrast, the granulated form offers lower resistance, making it more suitable for flow-through processes with gases.

Additional information

Activated charcoal is used in many different fields, including the chemical industry or for environmental protection. Activated charcoal looks the same as charcoal, but, due to a special 'activation process', it has a porous, sponge-like structure with a larger internal surface area.

Information on sustainable protection

Activated Charcoal is gained from natural raw materials.



Product	CAS-No.	Content	Packaging	Cat. No.
Charcoal activated				
Charcoal activated granular 1,5 mm extra pure	7440-44-0	1 kg	Plastic bag	1.02514.1000
		5 kg	Fiber carton	1.02514.5000
		25 kg	Fiber carton	1.02514.9025
Charcoal activated powder pure	7440-44-0	1 kg	Plastic bag	1.02183.1000
		17.5 kg	Fiber carton	1.02183.9018
Charcoal activated powder extra pure	7440-44-0	1 kg	Plastic bag	1.02184.1000
		5 kg	Fiber carton	1.02184.5000
		20 kg	Fiber carton	1.02184.9020
Charcoal activated for analysis	7440-44-0	250 g	Metal can	1.02186.0250
		1 kg	Metal can	1.02186.1000
		20 kg	Fiber carton	1.02186.9020

Charcoal wood

Application advice

Charcoal wood can be used for adsorption, filtration or clarification for many technical purposes.

Additional information

Coking or carbonization of beech wood (heating up to 400 °C without air) produces charcoal wood with a very large internal surface area.

Information on sustainable protection

Natural production by coking beech wood.

Product	CAS-No.	Content	Packaging	Cat. No.
Charcoal wood				
Charcoal wood powder	7440-44-0	25 kg	Fiber carton	1.02204.9025


Glass wool

Application advice

Due to its chemical resistance even against strong acids (except hydrofluoric acid), glass wool is used as a filter material in analytical processes, and as a supporting substance for drying agents in drying towers, for example to prevent caking in phosphorus pentoxide. It can also be used for insulation up to around 500 °C.

Additional information

Glass wool is made of soda lime silic acid glass fibers.



Product	CAS-No.	Content	Packaging	Cat. No.
Glass wool				
Glass wool	65997-17-3	250 g	Metal can	1.04086.0250
		1 kg	Fiber carton	1.04086.1000

Graphite

Application advice

Graphite is used as an absorbent and lubricant at temperatures of up to around 350 °C.

Information on sustainable protection

Graphite consists mainly of pure carbon.

Product	CAS-No.	Content	Packaging	Cat. No.
Graphite				
Graphite	7782-42-5	2.5 kg	Fiber carton	1.04206.2500
		25 kg	Fiber carton	1.04206.9025

Molecular sieves

Application advice

Molecular sieves can be used in desiccators and drying tubes; for keeping absolute solvents dry; for drying gases or solvents in drying columns; for selective adsorption (e.g. phosgene from chloroform)

Advantages

- Easy to use: practically chemically inert, non-toxic, easy disposal, dried liquids can be decanted
- High adsorption capacity even with low water content in substance to be dried
- High adsorption capacity even at high temperatures
- High adsorption affinity for polar and unsaturated organic molecules (however, H₂O is always preferentially adsorbed)
- Selective adsorption: only molecules that can pass through the pores are adsorbed

Capacity

15–24% at 25 °C

Indicator

Brown gel indicator changes from brown to yellowish at H₂O uptake of approximately 7–10 g/100 g molecular sieve

Regeneration

Maximum regeneration temperature is 450 °C. Molecular sieves can be dried in a drying oven above 250 °C to a water content of 2–3 g/100 g molecular sieve. The remaining water can be removed at 300–350 °C using a vacuum oil pump (10–1–10–3 mbar), whereby a cold trap containing carbon dioxide coolant or liquid air should be connected. Due to their high partial water vapor pressure, water pumps are not suitable for this purpose. For safety reasons, molecular sieves that have been used to dry solvents should be freed from possible solvent by mixing water prior to regeneration. Molecular sieves with a moisture indicator should not be heated above 150 °C.

Chemical and physical properties

Molecular sieves are crystalline, synthetic zeolites. Their crystal gratings are similar to a cage with numerous hollow spaces. The cavities are accessible from all sides by pores of exactly defined dimensions. Depending on the type of molecular sieve, these can be 0.3, 0.4 or 1.0 nm in diameter. If the water on the hollow spaces is removed due to heating, the material becomes an extremely active absorbent. However, this only occurs with molecules adsorbed that are small enough to pass through the pores (sieve effect).



Pore diameter	Type	Composition	Structure
0.3 nm	3A	Potassium sodium aluminium silicate	Zeolite
0.4 nm	4A	Sodium aluminium silicate	Zeolite
1.0 nm	13A/X	Sodium aluminium silicate	Zeolite



Physical properties

The molecular sieve crystallites obtained by hydrothermal manufacture are formed into rods and beads using 1–2% clay as binding agent. Vibration caused by transport may bring about some abrasion which collects in the first fraction during dynamic drying.

Bulk density	0.75 kg/l
Surface (BET)	800 m ² /g
Form supplied	Beads (~2 mm), rods (~1.6 mm)
Effective pore diameter depending on type	0.3, 0.4 or 1.0 nm
Hollow space volume	0.3 cm ³ /g
Specific heat	>0.8 kJ/kg
Heat of absorption per kg adsorbed water	4,200 kJ



Product	CAS-No.	Content	Packaging	Cat. No.
Molecular sieves				
Molecular sieve 0.3 nm beads ~2 mm (suitable for use in Karl Fischer titration)	1318-02-1	250 g	Plastic bottle	1.05704.0250
		1 kg	Plastic bottle	1.05704.1000
		10 kg	Plastic bucket	1.05704.9010
Molecular sieve 0.3 nm beads, with moisture indicator ~2 mm (suitable for use in Karl Fischer titration)	-	250 g	Plastic bottle	1.05734.0250
		1 kg	Plastic bottle	1.05734.1000
Molecular sieve 0.3 nm rods ~1.6 mm (1/16")	1318-02-1	250 g	Plastic bottle	1.05741.0250
		1 kg	Plastic bottle	1.05741.1000
Molecular sieve 0.4 nm beads ~2 mm Reag. Ph Eur	1318-02-1	250 g	Glass bottle	1.05708.0250
		1 kg	Glass bottle	1.05708.1000
		10 kg	Plastic bucket	1.05708.9010
Molecular sieve 0.4 nm beads, with moisture indicator ~2 mm	1318-02-1	250 g	Glass bottle	1.05739.0250
		1 kg	Glass bottle	1.05739.1000
Molecular sieve 0.4 nm rods ~1.6 mm (1/16")	1318-02-1	1 kg	Plastic bottle	1.05743.1000
Molecular sieve 0.5 nm beads ~2 mm	1318-02-1	250 g	Glass bottle	1.05705.0250
		1 kg	Glass bottle	1.05705.1000
Molecular sieve 1.0 nm beads ~2 mm	1318-02-1	1 kg	Glass bottle	1.05703.1000

Sea sand

Application advice

Sea sand is used in laboratories to filter and clarify extremely contaminated water or solutions. It is also used as a grinding aid, and to create a heating bath (sand bath). To remove contaminations which can be dissolved in acid, sea sand is boiled in hydrochloric acid and ignited in order to carbonize organic material and achieve the quality required for laboratory use.

Information on sustainable protection

Sea sand is a natural product which consists predominantly of silicon oxide and small quantities of other metal oxides.



Product	CAS-No.	Content	Packaging	Cat. No.
Sea sand				
Sea sand, extra pure	14808-60-7	1 kg	Plastic bottle	1.07711.1000
		5 kg	Plastic bottle	1.07711.5000
		10 kg	Fiber carton	1.07711.9010
Sea sand, purified by acid and calcined, for analysis	14808-60-7	1 kg	Plastic bottle	1.07712.1000
		5 kg	Plastic bottle	1.07712.5000
		10 kg	Plastic bottle	1.07712.9010
		25 kg	Fiber carton	1.07712.9025

Sodalime

Application advice

Soda lime is used to absorb CO₂ from air and gases, offering a large internal surface area for this purpose. It is suitable for small or large absorption systems in laboratories and technical facilities.

Additional information

Soda lime contains an indicator which changes from beige to violet. The color may change further from violet to blue when in intermittent use or after a period of shutdown. This is due to a back or balance reaction, in which hydroxide groups that have become active again are

available on the surface. However, no CO₂ is released when this occurs. Even if the indicator returns to beige, the soda lime does not have more absorption capacity than at the start, because the CO₂ bonding which has already taken place reduces the remaining capacity available.

Product	CAS-No.	Content	Packaging	Cat. No.
Sodalime				
Sodalime with indicator, granules ca. 1–2.5 mm	-	500 g	Plastic bottle	1.06733.0501
		2.5 kg	Plastic bottle	1.06733.2500
Sodalime pellets with indicator for analysis	-	1 kg	Plastic bottle	1.06839.1001
		5 kg	Plastic bottle	1.06839.5001
		25 kg	Fiber carton	1.06839.9025

Supelco®

Analytical Products

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany
Phone +49 (0)6151 720

SigmaAldrich.com

© 2019 Merck KGaA, Darmstadt, Germany and its affiliates. All Rights Reserved.
Merck, the vibrant M, and Supelco are trademarks of Merck KGaA, Darmstadt,
Germany or its affiliates. All other trademarks are the property of their respective owners.
Detailed information on trademarks is available via publicly accessible resources.

Lit. No. MK_BR5034EN Ver. 1.0 12/2019