


Product Information

94716 Abberior® STAR 470SXP, NHS ester

Absorption & Fluorescence Spectrum



Key Features

- Ideal depletion behaviour in STED microscopy ~750 nm
- 2-color labeling partner with STAR 635P for 2-color STED microscopy
- Tested in the Leica TSC 2-color STED (Ti:Sa) system

Description

STAR470SXP is a new member of an entirely new class of dyes introduced exclusively by Abberior – **first commercially available phosphorylated fluorescent dye**. They turn out to be **extremely suitable for STED applications**. STAR470SXP is the latest development of long-Stokes-shift dyes for STED microscopy. The dye can be excited from 450 to 480nm.

Abberior STAR470SXP can substitute dyes like Chromeo® 494.


Chemical Data : Abberior® Star 470SXP

Chemical Structure:	on request
Molecular Formula:	C ₃₁ H ₂₉ F ₃ N ₃ O ₁₀ P
Molecular Weight:	691.5 g/mol
Exact Mass:	691.15 Da
Solubility:	PBS, pH 7.4; DMF; DMSO; aq. acetonitrile; MeOH
Polarity:	polar (anionic)
Net Charge (at PH 7.4):	-2
Content:	> 90 %

Photophysical Data : Abberior® Star 470SXP

Absorption Maximum, λ_{max} , nm:	467 (PBS, pH 7.4) 468 (water) 472 (MeOH)
Fluorescence Maximum, λ_{fl} , nm:	598 (PBS, pH 7.4; water) 597 (aq. acetonitrile) 587 (MeOH)
Extinction Coefficient, ϵ , $\text{M}^{-1}\text{cm}^{-1}$:	25 000 (PBS, pH 7.4) 21 500 (water) 22 500 (aq. acetonitrile) 23 500 (MeOH)
Correction Factor, $\text{CF}_{260} = \epsilon_{260}/\epsilon_{\text{max}}$:	0.55 (PBS, pH 7.4, water) 0.46 (aq. acetonitrile; MeOH)
Correction Factor, $\text{CF}_{280} = \epsilon_{280}/\epsilon_{\text{max}}$:	0.39 (PBS, pH 7.4, water) 0.37 (aq. acetonitrile; MeOH)
Recommended STED Wavelength, λ_{STED} , nm:	740 – 770
Fluorescence Quantum Yield, η :	0.12 (PBS, pH 7.4)
Fluorescence Lifetime, τ :	0.8 ns (PBS, pH 7.4)

Applications



Two-color STED imaging of the Golgi apparatus using the Leica TCS STED microscope. HeLa cells were fixed and treated with anti-p230, labeled with Abberior STAR635P, and anti-Gpp130, labeled with Abberior STAR440SXP, to stain the trans- and cis- sides of the Golgi respectively. Scale bars 1 μ m. Image courtesy Francesca Bottanelli, Rothman Laboratory, Yale university, New Haven, USA.

Abberior STAR 470SXP is particularly designed and tested for 2-color STED microscopy in combination with our **STAR635P** using a single STED wavelength. **The dye is our recommendation for usage in the Leica TCS STED (Ti:Sa) 2-color system.** For more information see our **2-color dye package** section.

Literature

1. C.A. Wurm et al. "Novel red fluorophores with superior performance in STED microscopy", *Optical Nanoscopy* (2012) 1:7

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.