

Product Information

Phenazine Methosulfate

P9625

Product Description

Molecular Formula

 $C_{14}H_{14}N_2O_4S$

Molecular Weight

306.3 g/mol

CAS Number

299-11-6

Extinction coefficient

E^M for oxidized phenazine methosulfate (PMS) at maximum wavelength 387 nm is 26,300 (pH 2-8). The difference in molar extinction coefficients between oxidized and reduced PMS is 25,000.¹

Melting Point

158-160 °C

PMS is an electron acceptor and carrier in enzyme systems. The oxidized form is yellow, and the reduced form is colorless. Since the reduced PMS is easily oxidized by oxygen, it is used in assays as an electron carrier between enzymes and oxygen, cytochrome c, indophenols, or tetrazolium salts. The reduced PMS is used as an electron donor to reduce cytochrome c or in photosynthetic experiments. PMS is reduced by flavoproteins such as succinic dehydrogenase. PMS is reduced non-enzymatically by NADH and NADPH. PMS is also reduced by dithionite, sodium borohydride, ascorbic acid, reduced ubiquinones, and reduced vitamin K. PMS reacts with sulfhydryl groups including those of enzymes.¹ Usage of PMS for detection of specific dehydrogenases has been reported.^{2,3}

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Storage/Stability

Solutions should be prepared in deionized water, and not in neutral buffers. Solutions in water are stable frozen and protected from light for several months. Decomposition of solutions (with formation of pyocyanine and other products) may take hours when they are exposed to diffuse light, but only 5-10 minutes when exposed to sunlight.

Solutions are stabilized by polyvalent cations. Cations such as Ca^{2+} , Mg^{2+} , Mn^{2+} , Co^{2+} , Cu^{2+} , Ba^{2+} , Zn^{2+} , and Sn^{3+} at 100 mM each, prevented decomposition of PMS solutions when stored for 1 month, but were not protective when solutions were exposed to sunlight.^{1,2}

References

1. Data for Biochemical Research, 3rd ed., Dawson, R. M. C., et al., Oxford University Press (New York, NY: 1986), pp. 356-357.
2. Nachlas, et al., J. Biol. Chem., 235, 499-503 (1960).
3. Babson, A. L., and Babson, S. R., Clin. Chem., 19(7), 766-769 (1973).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.
All other trademarks are the property of their respective owners. Detailed information on
trademarks is available via publicly accessible resources.

© 2002-2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

P9625pis Rev 11/24

The Merck logo, consisting of the word "MERCK" in a bold, red, sans-serif font.