

Data sheet

Aldrich® glass reference electrode

Double junction, reference Ag/AgCl

Z113107

Protocol

1. Wet lower end of electrode in tap water. Carefully remove the lower cap.
NOTE: The cap is used only for shipping and storage. You might want to keep it for later use.
2. Use a wet paper towel and clean the wick area by wiping several times.
3. For first time usage, or after long storage, soak the lower end of the electrode in saturated KCl solution for five minutes.
4. Thoroughly rinse the electrode tip in distilled water, then place with standard glass pH electrode for use.

Storage and Handling

1. For short-term storage (up to one week), the electrode may be placed in a beaker containing approximately 1mL saturated KCl solution.
2. For long-term storage (over one week), fill storage cap with KCl solution and place the cap back on the electrode tip. The electrode may then be placed back in the box and stored.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Aldrich®, Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
Z113107-dat .Rev 11/23

The Merck logo, consisting of the word "MERCK" in a bold, red, sans-serif font.