



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sial.com  
sigma-aldrich.com

## Product Information

### Lambda DNA Mixed Digest

Catalog Number **D2916**  
Storage Temperature  $-20^{\circ}\text{C}$

#### Product Description

The mixed digest contains uncut bacteria phage Lambda plus Lambda completely digested by *Apa* I, *Kpn* I, *Xba* I, and *Xho* I in separate reactions. The fragments were filled in using DNA polymerase I Klenow fragment to prevent reannealing. Suitable for size determination of double-stranded DNA using DNA electrophoresis.

The digest is supplied in a solution of 10 mM Tris-HCl, pH 8.0, with 1 mM EDTA and contains 10 fragments, 1,503–48,502 bp (see Table 1).

**Table 1.**

Fragment Sizes: Base Pairs (bp)

48,502  
38,416  
33,498  
29,946  
24,508  
23,994  
17,053  
15,004  
10,086  
1,053

#### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### Storage/Stability

Store the digest at  $-20^{\circ}\text{C}$ .

#### Procedure

A sample of the marker should be diluted with gel loading buffer to an appropriate loading concentration. The Loading Mixture may be prepared as follows:

- 1  $\mu\text{g}$  Lambda Mixed Digest
- 6  $\mu\text{l}$  Gel Loading Solution (Catalog Number G2526, 0.05% w/v bromophenol blue, 40% w/v sucrose, and 0.1 M EDTA, pH 8.0)
- 2  $\mu\text{l}$  10 $\times$  TAE Buffer (Catalog Number T9650)
- 11  $\mu\text{l}$  Water, Molecular Biology Reagent (Catalog Number W4502)

The Loading Mixture was heated to  $65^{\circ}\text{C}$  for 5 minutes and quick cooled on ice. 0.5  $\mu\text{g}$  was loaded on a 0.4% agarose gel. After electrophoresis, the gel was stained with ethidium bromide. The ethidium bromide background can be reduced by destaining 30–45 minutes in 1 $\times$  TAE buffer.

#### References

1. Daniels, D.L., *et al.*, Appendix II: Complete Annotated Lambda Sequence in Lambda-II, Hendrix, R.W., *et al.*, eds. Cold Spring Harbor Laboratory (Cold Spring Harbor, NY: 1983).

JB,KH,MAM 05/07-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.