

# Workflow for Cannabinoids Analysis in Cannabis using a High Resolution Monolithic Silica HPLC Column Providing Low Backpressure and Extended Column Lifetime

Sunil Badal, Senior Scientist; Benjamin Peters, Lab Head Instrumental Analytics R&D; Uma Sreenivasan, Head of Reference Materials R&D

## Introduction

The legal use of recreational and medical cannabis is expanding globally along with hemp-based products (based on cannabidiol), for health and wellness. Hemp is defined legally in various geographies as cannabis varieties with limits on total tetrahydrocannabinol (THC) content. To ensure consumer safety, cannabis and hemp products need to be tested to determine accurate potency of the active cannabinoid constituents. Cannabis products in the market range from plants to distillates, and edibles to cosmetics. This broad variety of matrices underscores the need for robust columns and high throughput analytical methods.


This work provides a complete HPLC-DAD (high performance liquid chromatography-diode array detection) workflow for cannabinoids analysis using robust Chromolith<sup>®</sup> HighResolution (HR) HPLC columns which are based on monolithic silica. Chromolith<sup>®</sup> HPLC columns enable fast and cost-efficient separations due to low column backpressure and the very high robustness of the column. The low backpressure allows fast separation at high flow rates with the same mobile phase consumption per sample

compared to slower low flow-rate methods. The workflow offers the following:

- Detailed hemp bud sample preparation for HPLC-UV analysis.
- Fast and cost-efficient use of low back pressure Chromolith<sup>®</sup> HPLC columns to determine potency of a hemp bud sample.
- Demonstration of robustness of the Chromolith<sup>®</sup> column.
- Separation of 14 cannabinoids within 10 minutes.
- Calibration curve preparation using Certified Reference Materials (CRMs).

## Chemical Structures of 14 Cannabinoids

There are more than 100 distinct cannabinoids that have been isolated from cannabis. Delta-9-Tetrahydrocannabinol ( $\Delta^9$ -THC) is the primary psychoactive compound and cannabidiol (CBD) is another major non-psychoactive constituent in cannabis. Structures of  $\Delta^9$ -THC, CBD, and some other cannabinoids analyzed by the method are shown in **Figure 1**.



**Figure 1.** Chemical structures of fourteen cannabinoids included in this study.

## Experimental

In this work, hemp bud samples were analyzed to determine their potency. Sample preparation involved ethanol extraction of cannabinoids from plant material. The extract was then analyzed applying an HPLC-UV method and using a Chromolith® HR RP-18e monolithic silica HPLC column. Quantitation was performed utilizing a 6-point calibration curve obtained from HPLC-UV analysis of standard solutions prepared from CRMs. Peaks were identified using the retention times from a chromatogram of a 14 cannabinoids mix. Cannabinoid peaks were also verified by comparing UV absorption spectra of both samples and standards. Furthermore, robustness of the monolithic silica based Chromolith® column was demonstrated via retention time stability and separation efficiency after 1400 injections.

### Preparation of Mobile Phases

For mobile phase A, 0.1%  $\text{H}_3\text{PO}_4$  (aq.) was prepared by adding 1 mL  $\text{H}_3\text{PO}_4$  to 1000 mL of water. Pure methanol was used as mobile phase B.

### Preparation of Standard Solutions

Standard solutions containing six major analytes were prepared using Supelco® CRMs as shown in the **Table 1**.

**Table 1.** Preparation of standard solutions for 6 cannabinoids determined

Step	Instructions			
1	To a 1.5 mL autosampler vial, add the following volumes of CRM solutions to prepare the stock solution:			
#	Cannabinoid	Cat. No.	Conc. (mg/mL)	Volume ( $\mu\text{L}$ )
1	CBD	C-045	1.0	100
2	CBG	C-141	1.0	100
3	CBDA	C-144	1.0	100
4	CBN	C-046	1.0	100
5	$\Delta^9$ -THC	T-005	1.0	100
6	THCA	T-093	1.0	100
2	Add 400 $\mu\text{L}$ of methanol and mix well. Final concentration is 100 $\mu\text{g}/\text{mL}$ for all cannabinoids.			
3	Prepare the solution for 6-point calibration curve following the dilution scheme below using methanol as a diluent:			
Conc. ( $\mu\text{g}/\text{mL}$ )	Solution ID	Dilution Factor	Source Solution	
100	Stock	1	Stock	
25	A	4	Stock	
5	B	5	A	
1	C	5	B	
0.5	D	2	C	
0.25	E	2	D	

## Preparation of Peak Identification Solution

A peak identification solution containing 14 cannabinoids was prepared using CRMs, as shown in **Table 2**.

**Table 2.** Preparation for peak identification solution for 14 compounds

Step	Instructions				
1	To a 1.5 mL autosampler vial, add the following volumes of CRM solutions:				
	#	Cannabinoid	Cat. No.	Conc. (mg/mL)	Volume (µL)
	1	8-mix (neutrals)	C-219	0.5	100
	2	CBDVA	C-152	1.0	50
	3	CBDA	C-144	1.0	50
	4	CBNA	C-153	1.0	50
	5	THCA	T-093	1.0	50
	6	CBLA	C-171	0.5	50
	7	CBCA	C-150	1.0	50
2	Add 600 µL of methanol and mix well. Final concentration is 25 µg/mL for CBLA and 50 µg/mL for all other cannabinoids.				

## Preparation of Hemp Bud Extract

Cannabinoids were extracted from hemp buds using ethanol extraction as explained below:

- Homogenize 1 g hemp bud (particle size <1 mm). (Low temperature homogenization such as frozen ball-milling is the preferred method of homogenization without sample degradation.<sup>1</sup>)
- Transfer the homogenized sample to a 50 mL polypropylene centrifuge tube.
- Dispense 20 mL ethanol and vortex for 30 s.
- Incubate sample on horizontal shaker for 30 min.
- Centrifuge sample at 4000 rpm for 5 min to pellet plant material.
- Transfer the supernatant into amber 100 mL volumetric flask and keep the pellet for second extraction.
- Perform second extraction with 20 mL ethanol and add the supernatant to amber 100 mL volumetric flask containing contents of the first extraction.
- Fill flask to 100 mL mark with ethanol and mix well.
- Perform 1:10 and 1:100 dilution of sample with ethanol.
- Filter samples into HPLC vials with 0.2 µm PTFE membrane. Here, syringeless filter-vials were used.

Subsequent analysis was performed applying a 2 mm I.D. Chromolith® HR RP-18e HPLC column using conditions described in **Table 3**.

**Table 3.** Chromatographic conditions for determination of 14 cannabinoids by HPLC-UV

HPLC Parameters			
<b>Instrument:</b>	Agilent 1290 Infinity II LC System; Quaternary Pump; 0.12 mm ID tubing; 10 mm Max-light Cartridge Cell 1.0 µL		
<b>Column:</b>	Chromolith® HighResolution RP-18e 50x2 mm I.D. (1.52321)		
<b>Mobile phase:</b>	[A] Water with 0.1% H <sub>3</sub> PO <sub>4</sub> ; [B] methanol		
<b>Gradient:</b>	<b>Time (min)</b>	<b>%A</b>	<b>%B</b>
	0.0	38	62
	1.0	38	62
	8.0	10	90
	10.0	10	90
	10.1	38	62
<b>Flow:</b>	0.5 mL/min		
<b>Max Pressure:</b>	150 bar		
<b>Column temp.:</b>	40 °C		
<b>Detection:</b>	DAD, UV 228 nm		
<b>Injection:</b>	1.0 µL		

## Results and Discussion


Hemp bud sample was homogenized at low temperature to prevent analyte degradation using cryo-cup grinder followed by double extraction with ethanol. Resulting solution was diluted, filtered, and subjected to HPLC-DAD analysis. Calibration curves were obtained by analyzing solutions prepared from CRMs. Cannabinoids in hemp bud extract were identified based on retention time match with standards and cross verified with UV absorption spectra.

### System Suitability: Peak Identification Solutions


CRMs as 1.0 mg/mL or 0.5 mg/mL solutions in methanol or acetonitrile, were used to prepare calibration and peak identification solutions. Separation of 14 cannabinoids was demonstrated with good resolution and analyte signal reproducibility (**Table 4**). Separation of 14 cannabinoids was achieved in less than 10 minutes (**Figure 2**).

**Table 4.** Peak resolution and system reproducibility

Analyte	Resolution	%RSD of area count (n=5)	Analyte	Resolution	%RSD of area count (n=5)
CBDV	-	0.20	Δ <sup>9</sup> -THC	4.5	0.18
CBDVA	6.3	0.17	Δ <sup>8</sup> -THC	2.1	0.23
THCV	7.7	0.18	CBC	6.9	0.16
CBD	1.7	0.19	CBNA	2.8	0.20
CBG	2.1	0.19	THCA	4.6	0.18
CBDA	2.7	0.17	CBLA	1.8	0.17
CBN	6.9	0.18	CBCA	2.3	0.19



**Figure 2.** Chromatogram of 14 cannabinoids mixture obtained with a Chromolith® HR RP-18e 50×2 mm I.D. column at 228 nm.




**Figure 3.** Chromatogram of calibration mixture at 100 µg/mL obtained with a Chromolith® HR RP-18e 50×2 mm I.D. column at 228 nm.

## Quantitation

Calibration curves were obtained for six major cannabinoids (CBD, CBG, CBDA, CBN, Δ<sup>9</sup>-THC, and THCA); see chromatogram in **Figure 3**. Linearity of  $R^2 > 0.995$  was obtained for all analytes within the range of 0.25-100 µg/mL (**Figure 4**).

Results showed that the “as is” hemp bud sample contained 7.37% (wt/wt) total CBD and 0.25% (wt/wt) total THC (**Table 5**). Stable retention time for cannabinoids was observed during 1400 injections, demonstrating the robustness of the column towards a complex matrix like hemp bud extract.



**Figure 4.** Calibration curves of six cannabinoid analytes obtained with a Chromolith® HR RP-18e 50×2 mm I.D. column at 228 nm. Calibration curve ranges from 0.25 to 100 µg/mL. Linearity:  $R^2 > 0.995$  for all six analytes.


**Table 5.** Potency of hemp bud sample (Total CBD = CBD + 0.877xCBDA, and Total THC =  $\Delta^9$ -THC + 0.877xTHCA).

Analyte	Conc.* (% by weight)
CBD	0.37
CBG	0.00
CBDA	7.37
$\Delta^9$ -THC	0.04
CBN	0.05
THCA	0.21
Total CBD	7.74
Total THC	0.25


\*Note: The stated cannabinoid concentrations were calculated on “as is” basis and were not adjusted to dry weight. For many countries “official” potency testing the THC content needs to be referred to a dry weight sample base.<sup>2</sup> For example: the USDA definition is that samples should be dried to a consistent weight (typically 5-12% moisture content). Alternatively, Karl Fischer titration can be applied to determine moisture content.

In HPLC-UV analysis, identity of analytes depends on retention times and can be compromised by co-eluting peaks. Therefore, it is necessary to ensure that no co-elution of matrix compounds with the peak of interest is taking place. Here, we checked for the effects of

matrix impurities by comparing the UV absorption spectra of the analytes identified in the sample with those of the standards. As can be seen in **Figure 6**, most analytes in the hemp bud extract display absorption profiles similar to those of the standard. Among them, the spectra of CBN seems to contain an impurity which is visible as an extra peak in the chromatogram in **Figure 5** as well. This additional verification with UV absorption spectra further ensures the identity of detected analytes.



**Figure 5.** Overlay of chromatograms of peak identification solution and hemp bud extract obtained with a Chromolith® HR RP-18e 50 x 2 mm I.D. column at 228 nm.




**Figure 6.** Overlay of UV spectra from hemp bud sample (green) and standard solution at same retention time (blue) for six analytes. Slight differences in CBN spectra between sample and standard might be due to the overlap with another peak as can be seen in **Figure 5**.

## Robustness of Chromolith® Columns

The separation of matrix-rich samples, such as herbs, food, or biological samples tends to reduce the lifetime of particulate columns if insufficient sample preparation/cleanup is performed before HPLC separation. Chromolith® monolithic silica HPLC columns due to their bimodal pore structure allow the separation of matrix-rich samples with extended column lifetime, with no or very reduced sample preparation required. Extended column lifetime and reduced sample preparation significantly reduces the overall cost of operation. **Figure 7** shows the stability of the retention factors for the analysis of cannabinoids in a hemp extract sample with a Chromolith® HR RP-18e 100×2 mm column over 1400 injections. Results demonstrate the robustness of monolithic silica based Chromolith® columns.

In between cannabinoid sample analyses, the separation efficiency of the Chromolith® column was tracked with a performance test using anthracene.

Separation efficiency was tracked by calculating plate count (N). Results were compared with those obtained for two other HPLC columns with small particle size: Ascentis® Express 2 µm (superficially porous particles, SPP) and Purospher® STAR 2 µm (fully porous particles, FPP) (**Table 6**). In addition, retention time stability was compared using CBD. For the CBD retention time comparison, the cannabinoid separation was optimized for the Chromolith® column and was then transferred to Purospher® STAR and Ascentis® Express columns in order to create comparable stability data. No significant change in the CBD retention time was observed after 1400 injections. The efficiency for the monolithic silica based Chromolith® column only slightly decreased (2.7%), while for the fully porous and superficially porous, a reduction of 25.3% and 11.3%, respectively, was determined. These results again show the robustness of both bimodal pore structure and rigid monolithic silica skeleton over an extended analysis period.



**Figure 7.** Retention factor stability performance test/hemp sample on a Chromolith® HighResolution RP-18e 100×2 mm I.D. (Column temp.: 25 °C, mobile phase A: 0.1% H<sub>3</sub>PO<sub>4</sub>, mobile phase B: methanol, flow rate: 0.38 mL/min, injection volume: 0.2 µL, gradient: 72% B for 0.1 min, 72-90% B in 7 min, hold at 90% B for 3 min).

**Table 6.** Separation efficiency for anthracene and retention time of CBD using three different HPLC column types

Column	Separation Efficiency Anthracene			Retention Time CBD		
	Initial (N/m)	After 1400 injections (N/m)	Relative Change	Initial (min)	After 1400 injections (min)	Relative Change
Chromolith® HR RP-18e	165,990	161,500	-2.7%	4.19	4.14	-1.2%
Purospher® STAR RP-18e 2 µm (FPP)	188,450	140,770	-25.3%	6.42	6.34	-1.2%
Ascentis® Express C18 2 µm (SPP)	241,690	214,310	-11.3%	5.28	5.14	-2.7%

## Conclusion

This work demonstrates an HPLC-DAD workflow, using a monolithic silica based Chromolith® HR RP-18e HPLC column, for the determination of cannabinoids in hemp bud samples. Sample homogenization, use of accurate CRMs, separation of 14 cannabinoids with good selectivity, and robustness of Chromolith® columns are important elements of the workflow. Hemp bud samples were homogenized at low temperature to prevent analyte degradation using a cryo-cup grinder, followed by double extraction with ethanol. The resulting solution was diluted, filtered, and subjected to HPLC-DAD analysis. Calibration curves were obtained by analyzing solutions prepared from CRMs. Cannabinoids in hemp bud extract were identified based on retention time match with standards and cross verified with UV absorption spectra. Results showed that the hemp bud samples contained 7.37% (wt/wt) total CBD and 0.25% (wt/wt) total THC on as is basis without determining dry weight data. The robustness of Chromolith® columns was also demonstrated by the analysis of hemp bud and stable retention factors for cannabinoids over 1400 injections, proving once more the suitability of these columns for matrix-rich samples.

## References

1. Vaclavik L, Benes F, Fenclova M, Hricko J, Krmela A, Svobodova V, Hajslova J, Mastovska K. Quantitation of cannabinoids in Cannabis dried plant materials, concentrates, and oils using liquid chromatography–diode array detection technique with optional mass spectrometric detection: Single-laboratory validation study, first action 2018.11. Journal of AOAC International. 2019;102(6):1822–1833. <https://doi.org/10.1093/jaoac/102.6.1822>
2. AOAC SMPR 2019.003. Standard Method Performance Requirements (SMPRs) for quantitation of cannabinoids in Plant Materials of Hemp (Low THC Varieties Cannabis sp.). [https://www.aoac.org/wp-content/uploads/2019/10/SMPR-2019\\_003.pdf](https://www.aoac.org/wp-content/uploads/2019/10/SMPR-2019_003.pdf)

## Featured Products

Description	Cat. No.
<b>HPLC Columns</b>	
Chromolith® HighResolution RP-18e 50x2 mm I.D. HPLC Column	1.52321
Chromolith® HighResolution RP-18e 100x2 mm I.D. HPLC Column	1.52322

Description	Cat. No.
Ascentis® Express C18, 10 cm x 2.1 mm I.D., 2 µm	50813-U
Purospher® STAR RP-18 endcapped (2 µm) Hibar® HR 100-2.1, 10 cm x 2.1 mm I.D.	1.50648
<b>Sample Preparation, Solvents &amp; Reagents</b>	
Ultrapure Water from Milli-Q® IQ 7005 water purification system or bottled water	ZIQ7005T0C or 1.01262
Phosphoric acid, suitable for HPLC, LiChropur™, 85%	49685
Methanol, suitable for HPLC, ≥99.9%	34860
Ethanol 200 proof, HPLC/spectrophotometric grade	459828
Filter Vial, PTFE 0.2 µm, Whatman® Mini-UniPrep® G2 standard septum	WHAGN203NPEORG
<b>Certified Reference Materials and System Suitability Reagents (Cerilliant®)</b>	
Cannabichromene (CBC), 1.0 mg/mL in methanol	C-143
Cannabichromenic Acid (CBCA), 1.0 mg/mL in acetonitrile	C-150
Cannabicyclic Acid (CBLA), 0.5 mg/mL in acetonitrile	C-171
Cannabidiol (CBD), 1.0 mg/mL in methanol	C-045
Cannabidiolic Acid (CBDA), 1.0 mg/mL in acetonitrile	C-144
Cannabidivarin (CBDV), 1.0 mg/mL in methanol	C-140
Cannabidivarinic Acid (CBDVA), 1.0 mg/mL in acetonitrile	C-152
Cannabigerol (CBG), 1.0 mg/mL in methanol	C-141
Cannabinol (CBN), 1.0 mg/mL in methanol	C-046
Cannabinolic Acid (CBNA), 1.0 mg/mL in acetonitrile	C-153
Δ <sup>9</sup> -Tetrahydrocannabinol (Δ <sup>9</sup> -THC), 1.0 mg/mL in methanol	T-005
Δ <sup>8</sup> -Tetrahydrocannabinol (Δ <sup>8</sup> -THC), 1.0 mg/mL in methanol	T-032
Tetrahydrocannabinolic Acid (THCA), 1.0 mg/mL in acetonitrile	T-093
Tetrahydrocannabinavarin (THCV), 1.0 mg/mL in methanol	T-094
8 Cannabinoids Mix (Neutrals)	C-219
<b>Accessories</b>	
Clear Glass Vial, 2 mL, pk of 100, Natural PTFE/silicone septa	29652-U
Amber Glass Vial, 2 mL, pk of 100, Natural PTFE/silicone septa	29653-U

Find more information on cannabis testing at [SigmaAldrich.com/cannabis](https://www.sigmaaldrich.com/cannabis)

Learn more about Chromolith® columns at [SigmaAldrich.com/chromolith](https://www.sigmaaldrich.com/chromolith)

## To place an order or receive technical assistance

Order/Customer Service: [SigmaAldrich.com/order](https://SigmaAldrich.com/order)  
Technical Service: [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice)  
Safety-related Information: [SigmaAldrich.com/safetycenter](https://SigmaAldrich.com/safetycenter)

[SigmaAldrich.com](https://SigmaAldrich.com)

We have built a unique collection of life science brands with unrivalled experience in supporting your scientific advancements.

**Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®**

Merck KGaA  
Frankfurter Strasse 250  
64293 Darmstadt, Germany

