


Overnight Express™ Autoinduction Systems



Can Merck Chemicals make my protein expression process easier?

Yes! With Overnight Express™
Autoinduction Systems, the protein expression
process begins without having to
manually add induction agent.


We've done it for you, so that
you don't have to!!

That's what's in it for you. Merck Chemicals

The **Overnight Express™ Autoinduction Systems** enable regulated protein expression in *E. coli*, without monitoring the culture or adding inducer during cell growth. The simplified protocol offers great convenience, allowing you to focus on your research while it does its job.

Overnight Express™ Autoinduction Systems

With Overnight Express™ Autoinduction Systems, a period of cell growth is followed by spontaneous induction of protein expression – without monitoring cell density and without conventional induction with IPTG. The method is based on media components that are metabolized differentially to promote growth to high density and automatically induce protein expression from lac promoters.



Analysis of eukaryotic target proteins purified from cultures induced with Overnight Express System 1 versus IPTG
 pET recombinants encoding the indicated His•Tag® fusion proteins were transformed into BL21(DE3). Overnight Express System 1 induction was accomplished by inoculating 5 ml medium in 10-ml × 24-well plates with a single colony and incubating overnight (approximately 16 h) at 30°C with shaking at 250 rpm. IPTG induction was accomplished by inoculating 5 ml medium in 10-ml × 24-well plates with a single colony and incubating at 16°C with shaking at 250 rpm to an average OD600 of 1.0 followed by addition of IPTG to 1 mM final concentration and incubating an additional 16 h prior to harvest. Target proteins were extracted and purified using the RoboPop Ni-NTA His•Bind® Purification Kit. The indicated sample volumes were analyzed by SDS-PAGE (10–20% gradient gel) and Coomassie blue staining.




Overnight Express™
 systems...
 for a good
 night's sleep



Overnight Express™ Autoinduction Systems

We've done it for you, so that you don't have to!

The Overnight Express Autoinduction System is extremely convenient for routine expression of proteins in multiple cultures and is ideal for high-throughput parallel analysis of protein expression, solubility, and purification from multiple expression clones



Overnight Express™ Autoinduction Medium Selection Guide

Please reference the selection guide below for the medium for your specific application.

Overnight Express™ Instant LB Medium		6
Medium	Complete autoinduction media in granulated Luria-Bertani formulation	
Application	Most routine recombinant protein expression applications	
Overnight Express™ Instant TB Medium		6
Medium	Complete autoinduction media in granulated Terrific Broth formulation	
Application	Most routine recombinant protein expression applications	
Overnight Express™ Autoinduction System 1		7
Medium	Autoinduction media to be added to glucose-free medium (eg. 2X YT, SOC, LB and TB)	
Application	Most routine recombinant protein expression applications	
Overnight Express™ Autoinduction System 2		9
Medium	Chemically defined autoinduction media	
Application	Compatible with selenomethionyl (Se-Met) labeling of proteins for x-ray crystallography	
Overnight Express™ NMR Medium – Optimization		11
Medium	Chemically defined autoinduction media	
Application	Determine optimal culture conditions for high-level protein expression before isotopic protein labeling. It can also be used for ¹⁵ N protein labeling when user provides ¹⁵ N-ammonium chloride.	
Overnight Express™ NMR Medium – ¹⁵N		11
Medium	Chemically defined autoinduction media	
Application	High level incorporation of ¹⁵ N for initial NMR analysis, to assess suitability for structure determination	
Overnight Express™ NMR Medium – ¹⁵N, ¹³C		11
Medium	Chemically defined autoinduction media	
Application	High level incorporation of ¹⁵ N and ¹³ C for backbone and side-chain assignments and for restraint measurements in structure determination	




Overnight Express™ Instant LB and TB Medium

The complete solution for high level protein expression without the need to monitor cell growth

- Autoinduction media ▶ Simplifies protocol by eliminating the monitoring and induction steps
- Rich culture media ▶ High cell densities and protein expression level
- Granulated format ▶ Rapid and uniform dissolution in water
Minimize the chance of inhaling hazardous powdered ingredients
- Singe use EasyPak packaging ▶ Minimal sample handling, significantly reduce the chance of contamination
Convenient, no need to spend extra time on weighting out different ingredients

Granulated Instant LB Medium and TB Medium formulations are combinations of Overnight Express™ System 1 with either Luria-Bertani Broth (LB) or Terrific Broth (TB). These extremely convenient, rich-culture autoinduction media provide high-level protein production in pET and other IPTG-inducible bacterial expression systems.



Overnight Express™ Instant LB and Instant TB Medium are ideal for routine expression of proteins in multiple cultures and for high-throughput parallel analysis of protein expression, solubility and purification from multiple expression clones.

Two packaging formats are available. The EasyPak aluminum foil pouch contains enough granulated medium to prepare 1 L culture. Just add the EasyPak contents to 1 L sterile water, supplement with 10 ml glycerol, and bring it to a boil in a microwave oven for 2 minutes.

Ordering Information

Product	Size	Cat. No.	Product	Size	Cat. No.
Overnight Express™ Instant LB Medium	1 EasyPak (1 x 45 g)	71757-3	Overnight Express™ Instant TB Medium	1 EasyPak (1 x 60 g)	71491-3
	5 EasyPak (5 x 45 g)	71757-4		5 EasyPak (5 x 60 g)	71491-4
	1 kg	71757-5		1 kg	71491-5
	10 kg	Bulk.		10 kg	Bulk.

Overnight Express™ Autoinduction System 1

Self induced protein expression system compatible to your specific culture medium

- Autoinduction media  Simplifies protocol by eliminating the monitoring and induction steps
Ideal for protein production in pET expression system or other IPTG-inducible bacterial system
- Sterile stock solutions  Minimal sample handling, significantly reduce the chance of contamination
Convenient, no need to spend extra time on weighting out different ingredients
- Compatible to any conventional glucose-free bacterial growth medium  Just add it into your existing culture medium; offers you with simplicity and flexibility

The unique Overnight Express™ Autoinduction Medium provides high-level protein production in pET and other IPTG-inducible bacterial expression systems. The Overnight Express™ Autoinduction System 1 kit contains three concentrated sterile solutions:

OnEx™ Solution 1 A blend of carbon sources optimized for tightly regulated uninduced growth to high cell density, followed by high-level induction

OnEx™ Solution 2 A concentrated buffer and nitrogen-source bend that mediates metabolic acid production and provides additional nitrogen for increased protein synthesis

OnEx™ Solution 3 Provides magnesium for maximal cell density.




Save time!
We've done the prep work
for you, so that
you don't have to!



Overnight Express™ Autoinduction System 1 (continued)

Cell mass and target protein yield are often increased several-fold compared to conventional protocols using induction with IPTG.



Expression and purification of target proteins from cultures induced with Overnight Express System 1 (OE1) versus IPTG. pET recombinants encoding B-gal, NusA, and GST His•Tag® fusion proteins were transformed into BL21(DE3). Protein expression was induced in parallel cultures either by Overnight Express System 1 or 1 mM IPTG. Cells were harvested by centrifugation and protein was extracted with BugBuster® HT Protein Extraction Reagent plus rLysozyme® Solution. After removing a sample, the extract was used for robotic affinity purification using Ni-NTA His•Bind® Resin. Samples of the crude protein extract (7 ml) and of the purified target protein (2 ml) were analyzed by SDS-PAGE (10-20% gel) and Coomassie blue staining.

Ordering Information

Product	Size	Components	Cat. No.
Overnight Express™ Autoinduction System 1	1 Liter Kit	OnEx™ Solution 1: 1 x 20 ml	71300-3
		OnEx™ Solution 2: 1 x 50 ml	
		OnEx™ Solution 3: 1 x 1 ml	
	5 Liter Kit	OnEx™ Solution 1: 5 x 20 ml	71300-4
		OnEx™ Solution 2: 2 x 125 ml	
		OnEx™ Solution 3: 5 x 1 ml	

Overnight Express™ Autoinduction System 2

Complete, chemically defined medium to selenomethionine (Se-Met) label proteins for crystallography

- Chemically defined  Consistent product performance, minimal lot-to-lot variability
Ideal for Se-Met labeling of proteins for crystallography
- Autoinduction media  Simplifies protocol by eliminating the monitoring and induction steps
Ideal for protein production in pET expression system or other IPTG-inducible bacterial system
- Sterile stock solutions  Minimal sample handling, significantly reduce the change of contamination
Convenient, no need to spend extra time on weighting out different ingredients

The Overnight Express™ Autoinduction System 2 features the capability to label proteins with selenomethionine (Se-Met) for downstream crystallization and x-ray diffraction studies. All components is chemically defined, resulting in consistent product performance by minimizing lot-to-lot variability. The Overnight Express™ Autoinduction System 2 kit contains six concentrated sterile solutions:

OnEx™ Solution 1 A blend of carbon sources optimized for tightly regulated uninduced growth to high cell density, followed by high-level induction

OnEx™ Solution 2 A concentrated buffer and nitrogen-source bend that mediates metabolic acid production and provides additional nitrogen for increased protein synthesis

OnEx™ Solution 3 Provides magnesium for maximal cell density.

OnEx™ Solution 4 Provides trace metals to minimize growth limitations associated with mineral deficiencies and satisfies the metal requirements of metal-containing target proteins, even at high expression levels

OnEx™ Solution 5 A mixture of amino acids lacking methionine


OnEx™ Solution 6 A separate methionine solution. Sufficient methionine (Met) is provided to support growth of the Met auxotroph B834 while providing the ability to reduce the level of unlabeled Met for selenomethionine incorporation by Met auxotrphs

Overnight Express™
systems...
for a good
night's sleep




Overnight Express™ Autoinduction System 2 (continued)

The Overnight Express Autoinduction System 2 is the perfect solution for routine expression of proteins in multiple cultures and is ideal for high-throughput parallel analysis of protein expression, solubility screening, and purification from multiple expression clones



Analysis of crude, soluble, and purified proteins from cultures grown in Overnight Express Autoinduction System 2 medium

pET-44b(+) plasmids expressing a His•Tag®/NusA fusion protein were transformed into B834(DE3) cells and grown for 16 h at 37°C in Overnight Express System 2 medium containing Se-Met (125 µg/ml) or Met (OnEx Solution 6). Cultures were centrifuged and cells were resuspended in Ni-NTA Bind Buffer containing AEBSF Hydrochloride, Benzamidine Hydrochloride, and Lysonase™ Bioprocessing Reagent. The suspension was sonicated and centrifuged at 12,000 x g for 10 min. Soluble cell extracts were processed by Ni-NTA His•Bind® chromatography. Two micrograms of purified protein (P), 5 µl of soluble extract (S), and a sample of crude extract (C; standardized to harvest OD₆₀₀) were analyzed by SDS-PAGE (10–20% gradient gel) and Coomassie blue staining. The predicted molecular mass for the pET-44b(+) protein is 67.8 kDa. Lane M: Perfect Protein™ Markers, 10–225 kDa.



Mass spectra

Mass spectroscopy analysis shows Se-Met incorporation into target proteins with Overnight Express System 2 autoinduction as described in figure at left. (Spectra provided by the Mass Spectrometry/Bioanalytical Facility at the University of Wisconsin Biotechnology Center.)

Ordering Information

Product	Size	Components	Cat. No.
Overnight Express™ Autoinduction System 2	1 Liter Kit	OnEx™ Solution 1: 1 x 20 ml	71366-3
		OnEx™ Solution 2: 1 x 50 ml	
		OnEx™ Solution 3: 1 x 1 ml	
		OnEx™ Solution 4: 1 x 1 ml	
		OnEx™ Solution 5: 1 x 20 ml	
		OxEx™ Solution 6: 1 x 20 ml	
	5 Liter Kit	OnEx™ Solution 1: 5 x 20 ml	71366-4
		OnEx™ Solution 2: 2 x 125 ml	
		OnEx™ Solution 3: 5 x 1 ml	
		OnEx™ Solution 4: 5 x 1 ml	
		OnEx™ Solution 5: 1 x 100 ml	
		OnEx™ Solution 6: 1 x 100 ml	

Overnight Express™ Autoinduction NMR Media

Complete, chemically defined medium to ¹⁵N- and ¹³C- label proteins for NMR spectroscopy

- Chemically defined ▶ Consistent product performance, minimal lot-to-lot variability
- Stable isotope-labeled nutrients ▶ Reliable and efficient isotopic labeling of proteins.
Ideal for ¹⁵N- and ¹³C- labeling of proteins for NMR analysis
- Autoinduction media ▶ Simplifies protocol by eliminating the monitoring and induction steps.
Ideal for protein production in pET expression system or other IPTG-inducible bacterial systems
- Sterile stock solutions ▶ Minimal sample handling, significantly reduce the change of contamination.
Convenient, no need to spend extra time on weighing out different ingredients

The Overnight Express™ Autoinduction NMR Media offer convenient, high-yield expression and efficient isotopic labeling of recombinant proteins for NMR analysis. Each kit includes reagents to make non-inducing, chemically defined starter medium. Starter medium can be used for staged growth without expression, prior to labeling, and for stable storage of bacterial strains containing expression plasmids.

Overnight Express™ Autoinduction NMR Medium – Optimization Intended Purpose: To optimize cell growth and protein expression in minimal medium, before preparing isotopically labeled target proteins for NMR spectroscopy

Overnight Express™ Autoinduction NMR Medium – ¹⁵N Intended Purpose: To produce [U-¹⁵N] labeled proteins to assess suitability for further analysis

Overnight Express™ Autoinduction NMR Medium – ¹⁵N, ¹³C Intended Purpose: To produce [U-¹⁵N, U-¹³C] labeled proteins to determine structure when supplemented with ¹³C-labeled glycerol (available separately, Cat. No. CLM 1510-EMD)

Overnight Express Autoinduction NMR Media Components

Solution	Description	NMR Medium–Optimization ml/L	NMR Medium– ¹⁵ N ml/L	NMR Medium– ¹⁵ N, ¹³ C* ml/L	Non-inducing starter medium ml/200 ml
OnEx™ NMR Solution 1	Autoinducing carbon source; single label	20	20		
OnEx™ Dual-NMR Solution 1	Autoinducing carbon source; dual label			20	
OnEx™ Solution 2	Buffer salts with nitrogen				10
OnEx™ NMR Solution 2	Buffer salts without nitrogen	50	50	50	
OnEx™ Solution 3	Magnesium	2	2	2	0.2
OnEx™ Solution 4	Trace metals	0.2	0.2	0.2	0.04
OnEx™ NMR Solution 5	Amino acids				4
OnEx™ NMR Solution 6	Carbon source (non-autoinducing)				5
*1.0 M ammonium chloride (unlabeled)	Nitrogen source for optimization media	50			
*1.0 M ¹⁵ N-ammonium chloride	¹⁵ N source		50	50	
2.5% (w/v) ¹³ C-glucose	¹³ C source*			20	


*For high level incorporation of ¹³C, suitable for NMR structural determination, Overnight Express Autoinduction NMR Medium–¹⁵N, ¹³C must be supplemented with ¹³C3 Glycerol, available separately (Cat. No. CLM 1510-EMD). ¹³C3 Glycerol is manufactured by Cambridge Isotope Laboratories, Inc.

*Solution is prepared using the solid component included in the medium kit.



Overnight Express™ Autoinduction NMR Media (continued)

SDS-PAGE analysis of expression and purification samples of proteins produced using Overnight Express Autoinduction NMR Media



Lane	Sample	Lane	Sample
1	Perfect Protein™ Markers 10-225 kDa (PPM)	10	At1g01470 - (S)
2	At5g22580 - total cellular protein (TCP)	11	At1g01470 - (P) (3 mg)
3	At5g22580 - soluble protein (S)	12	PPM
4	At5g22580 - purified protein (P) (3 mg)	13	At1g01470 - (TCP)
5	PPM	14	At1g01470 - (S)
6	At5g01610 (S)	15	At1g01470 - Ni-NTA His•Bind® elution
7	At5g01610 (P) (3 mg)	16	At1g01470 - after protease cleavage
8	PPM	17	At1g01470 - subtractive Ni-NTA His•Bind flowthrough
9	At1g01470 - (TCP)	18	At1g01470 - gel filtration flowthrough

For additional information on this experiment, including isotopic incorporation and NMR data, please see the *inNovations* 26 article at www.merck4biosciences.com/innovations26OnEx

Ordering Information

Product	Size	Components	Cat. No.
Overnight Express™ Autoinduction NMR Medium – Optimization	1 Liter Kit	OnEx™ NMR Solution 1: 1 x 20 mL OnEx™ Solution 2: 1 x 10 mL OnEx™ NMR Solution 2: 1 x 50 mL OnEx™ Solution 3: 3 x 1 mL OnEx™ Solution 4: 1 x 1 mL OnEx™ NMR Solution 5: 1 x 4 mL OnEx™ NMR Solution 6: 1 x 5 mL Ammonium Chloride: 2.7 g	71760-3
Overnight Express™ Autoinduction NMR Medium – ¹⁵ N	1 Liter Kit	OnEx™ NMR Solution 1: 1 x 20 mL OnEx™ Solution 2: 1 x 10 mL OnEx™ NMR Solution 2: 1 x 50 mL OnEx™ Solution 3: 3 x 1 mL OnEx™ Solution 4: 1 x 1 mL OnEx™ NMR Solution 5: 1 x 4 mL OnEx™ NMR Solution 6: 1 x 5 mL Ammonium Chloride: 2.7 g	71759-3
Overnight Express™ Autoinduction NMR Medium – ¹⁵ N	5 Liter Kit	OnEx™ NMR Solution 1: 1 x 100 mL OnEx™ Solution 2: 1 x 50 mL OnEx™ NMR Solution 2: 2 x 125 mL OnEx™ Solution 3: 1 x 11 mL OnEx™ Solution 4: 2 x 1 mL OnEx™ NMR Solution 5: 1 x 20 mL OnEx™ NMR Solution 6: 1 x 25 mL [¹⁵ N]Ammonium Chloride: 5 x 2.7 g	71759-4
Overnight Express™ Autoinduction NMR Medium – ¹⁵ N, ¹³ C	1 Liter Kit	OnEx™ NMR Solution 1: 1 x 100 mL OnEx™ Solution 2: 1 x 50 mL OnEx™ NMR Solution 2: 2 x 125 mL OnEx™ Solution 3: 1 x 11 mL OnEx™ Solution 4: 2 x 1 mL OnEx™ NMR Solution 5: 1 x 20 mL OnEx™ NMR Solution 6: 1 x 25 mL [¹⁵ N]Ammonium Chloride: 5 x 2.7 g	71789-3
¹³ C3 Glycerol	5 g		CLM 1510-EMD

Additional guidelines

Glycerol stock preparation

When growing cultures to prepare glycerol stocks that will be used to inoculate Overnight Express cultures, we recommend the addition of 0.5% glucose to a glucose-free medium (e.g., TB, LB broth, or 2X YT) to maintain plasmid stability. Grow the cells to an OD₆₀₀ of 0.6–0.8 and add 0.1 vol of sterile 80% glycerol. Mix well and store at –70°C.

Aeration: Efficient growth to saturation and utilization of carbon sources provided in Overnight Express medium requires vigorous agitation and proper aeration. Optimized culture volume:vessel dimension ratio is required for proper aeration.

Temperature and length of incubation

It is important to grow the cells to stationary phase when using the Overnight Express Systems. Stationary phase is usually reached as quickly as 8–10 hours, if the cultures are incubated at 37°C. When lower incubation temperatures are used, saturation may only be reached by incubation for 24 hours or more. Continued incubation for several hours after stationary phase appears to have no deleterious effects.

To export the target protein using the signal sequence leaders present in a number of pET vectors or to improve the yield of soluble protein, growth, and induction at 25°C or 30°C may be optimal.

Bacterial strains: Because lactose is used for induction, expression hosts should produce functional lac permease (encoded by the *lacY* gene) and β-galactosidase (encoded by the *lacZ* gene) for consistent results in both complex and defined media. *lacY* mutant strains will not efficiently transport lactose for induction and *lacZ* mutants will not convert a portion of the transported lactose into the allolactose inducer. Elevated levels of target gene expression in *lacY* and *lacZ* mutant strains may occur as cells approach stationary phase in some complex media, including Overnight Express medium. However, this induction may vary depending upon medium composition, cell growth stage, and nutrient availability, all of which affect pH and the levels of cyclic AMP and acetate (Grossman et al., 1998).

If using a plasmid with a T7lac promoter for expression, a host strain that does not contain a pLysS plasmid is recommended [e.g., BL21(DE3)]. The combination of the T7 lysozyme expressed by the pLysS plasmid and the lac repressor encoded by pET vectors carrying the T7lac promoter results in significantly reduced levels of protein expression when using the Overnight Express Autoinduction Systems. When the “plain” T7 promoter is used, the low level of lysozyme provided by pLysS has little effect on expression of target proteins.

Expression vectors

Overnight Express Autoinduction Systems are compatible with pET bacterial expression vectors and other IPTG-inducible bacterial expression systems.



Related Products

Benzonase® Nuclease

Benzonase® Nuclease is a genetically engineered endonuclease from *Serratia marcescens*. It degrades all forms of DNA and RNA (single-stranded, double-stranded, linear, and circular) while having no proteolytic activity. It is effective over a wide range of conditions and possesses an exceptionally high specific activity. The enzyme completely digests nucleic acids to 5'-monophosphate terminated oligonucleotides 2 to 5 bases in length (below the hybridization limit), which is ideal for removal of nucleic acids from recombinant proteins, enabling compliance with FDA guidelines for nucleic acid contamination. The ability of Benzonase® Nuclease to rapidly hydrolyze nucleic acids makes the enzyme an excellent choice for viscosity reduction, to reduce processing time and to increase yields of protein. For example, the enzyme is compatible with BugBuster® and PopCulture® Protein Extraction Reagents and can be added along with these reagents to eliminate viscosity by removing nucleic acids from *E. coli* extracts.

BugBuster® 10x Protein Extraction Reagent

BugBuster® 10X is a concentrated formulation of the proprietary detergents employed in BugBuster® Protein Extraction Reagent without the addition of salts or buffer components. This concentrated formula provides a flexible alternative to the ready-to-use standard 1X reagent, allowing user-defined dilution and addition of buffer components. BugBuster 10X has all of the bioprocessing benefits of standard BugBuster® Protein Extraction Reagent plus the freedom to control pH, reagent concentration, and buffer additives necessary for maximum extraction and activity of your target protein.

PopCulture® Reagent

PopCulture Reagent is a buffered mixture of concentrated detergents formulated to extract proteins from *E. coli* cells directly in their culture medium. Using this method, cell culture, protein extraction and purification performed in the original culture tube or multiwell plate. PopCulture reagent perforates the *E. coli* cell wall without denaturing soluble protein and protects protein from the pH extremes produced in high density culture media. Recombinant proteins can be assayed directly or purified by adding an affinity matrix, washing the matrix-target protein complex to remove culture medium and cellular contaminants and eluting the purified protein from the matrix. Purification of fusion proteins from total culture extracts has been demonstrated using both IMAC and GST affinity chromatography methods. *This product is covered by EP Patent 1,432,822 owned by EMD Chemicals Inc. or its Affiliates.*

RoboPop™ Ni-NTA His•Bind® Purification Kit

The RoboPop™ Ni-NTA His•Bind® Purification Kit is designed for filtration-based 96-well format purification of soluble His•Tag® fusion proteins directly from *E. coli* cultures without harvesting cells. The kit features PopCulture® Reagent, rLysozyme™ Solution, and Benzonase® Nuclease for centrifuge-free cell lysis and extract preparation in one step. The combination of PopCulture extraction, Ni-NTA His•Bind Resin, and a 2-ml filter plate allows high-throughput processing of up to 5 ml of *E. coli* culture per well. Whereas the magnetic-based His•Mag™ kit purifies up to 125 µg target protein per 1 ml culture, the filtration-based kit purifies up to 1 mg His•Tag fusion protein per 5 ml culture. Bacterial culture, cell lysis, and resin binding steps are carried out in standard 24-well plates (not supplied), which accommodate a maximum volume of 5 ml per well. The reaction slurry is then transferred to a 96-well Filter Plate (included) and the washing and elution steps are carried out on a vacuum manifold. The Filter Plate is compatible with standard filter manifolds for manual sample processing, and the entire purification has been validated for robotic sample processing with the Genesis® Freedom™ Workstation from Tecan and the MultiPROBE® II liquid handling work station from PerkinElmer Life Sciences. A 96-well Collection Plate (1-ml wells) with an air-tight aluminum foil sealer is provided for storage of the purified proteins.

Antibiotics

Antibiotics are natural substances secreted by microorganisms that are toxic to other microorganisms, but are generally non-toxic toward higher organisms. We offer a wide range of antibiotics for your protein production needs.

Ampicillin, Sodium Salt

Ampicillin is a β-lactam antibiotic that kills growing cells by interfering with the terminal reaction in bacterial wall synthesis. Active against Gram-negative bacteria. Inactivated by lactamases. Suitable for use in research that uses ampicillin-resistant plasmids.

Chloramphenicol

Chloramphenicol is a synthetic bacteriostatic antibiotic that inhibits translation of RNA by binding to the 50S ribosomal subunit, blocking the peptidyltransferase reaction of ribosomes. Chloramphenicol resistance is encoded by the resistance gene *cat*, whose product, chloramphenicol acetyltransferase, inactivates the antibiotic via acetylation. Typically used at 30-50 µg/ml for selection of chloramphenicol-resistant *E. coli*.

Kanamycin Sulfate, *Streptomyces kanamyceticus* sp.

Kanamycin Sulfate is an aminoglycoside antibiotic effective against Gram-positive and Gram-negative bacteria. It inhibits protein biosynthesis by acting on the 30S ribosome, causing misreading of the genetic code. In mammals, kanamycin may cause renal damage and is ototoxic.

Rifampicin

Antibiotic that specifically inhibits DNA-dependent RNA polymerase in bacteria by forming an inactive complex. Does not affect mammalian RNA polymerase. Inhibits transcription by preventing the initial transcription complex from entering the elongation mode.

Spectinomycin, Dihydrochloride, Pentahydrate, *Streptomyces* sp.

Spectinomycin is a broad-spectrum aminoglycoside antibiotic containing two glucose moieties. It is effective against Gram-positive and Gram-negative bacteria. It inhibits initiation, elongation, and termination of protein synthesis in prokaryotes and induces misreading of the genetic code. Footprint studies indicate that spectinomycin exerts regional effects on ribosomal structure.

Streptomycin Sulfate, *Streptomyces* sp.

Streptomycin Sulfate is an antibiotic effective against Gram-positive and Gram-negative bacteria. It inhibits initiation, elongation, and termination of protein synthesis in prokaryotes, induces misreading of the genetic code. It is often used in culture media to control growth of microorganisms.

Tetracycline, Hydrochloride

Tetracycline, Hydrochloride blocks protein synthesis by inhibiting aminoacyl tRNA binding to the A-site of ribosomes. It induces a cold shock-response and enhances P450 expression in bacteria.

Ordering Information

Product	Size	Cat. No.
Benzonase® Nuclease, Purity > 99%	10 KU	70664-3
Benzonase® Nuclease HC, Purity > 99%	25 KU	71206-3
BugBuster® Master Mix	100 ml	71456-3
	500 ml	71456-4
BugBuster® Protein Extraction Reagent	100 ml	70584-3
	500 ml	70584-4
BugBuster® 10x Protein Extraction Reagent	10 ml	70921-3
	50 ml	70921-4
	100 ml	70921-5
PopCulture® Reagent	15 ml	71092-3
	75 ml	71092-4
	250 ml	71092-5
RoboPop™ Ni-NTA His•Bind® Purification Kit	1 Kit	71188-3
Ampicillin, Sodium Salt	5 gm	171254-5GM
	25 gm	171254-25GM
Chloramphenicol	25 gm	220551-25GM
	100 gm	220551-100GM
Kanamycin Sulfate, <i>Streptomyces kanamyceticus</i> sp.	5 gm	420311-5GM
	25 gm	420311-25GM
Rifampicin	1 gm	557303-1GM
	5 gm	557303-5GM
Spectinomycin, Dihydrochloride, Pentahydrate, <i>Streptomyces</i> sp.	100 gm	567570-100GM
Streptomycin Sulfate, <i>Streptomyces</i> sp.	100 gm	5711-100GM
Tetracycline, Hydrochloride	10 gm	58346-10GM
	25 gm	58346-25GM
	50 gm	58346-50GM



CUSTOM SERVICES

Need a **larger size** or **special service**?

From package size to documentation to management of your shipments, if you have a custom need, we will partner with you to meet your requirements.

Contact us at:

contact@merck4biosciences.com



Argentina

Merck Quimica Argentina S.A.I.C.
Tel: +54 11 4546 8100
Fax: +54 11 4546 7369
E-mail: mveltri@merck.com.ar
www.merck-chemicals.com.ar

Australia

Toll Free: 1800 335 571
Merck Pty. Limited
Tel: +61 3 9728 7600
Fax: +61 3 9728 1351
E-mail: merck@merck.com.au
www.merck-chemicals.com.au

Brazil

Merck S.A.
Toll Free: 0800 21 9292
Tel: +55 11 3346 8500
Fax: +55 11 3207 5040
E-mail: quimica@merck.com.br
www.merck.com.br

Caribbean

Merck, S.A.
See information under Central America

Central America

Merck, S.A. Guatemala
Tel: +50 2 2410 2300
Fax: +50 2 2250 5193
E-mail: quimicos@merck.com.gt
www.merck.com.gt

Chile

Merck S.A.
Tel: +56 2 3400 000
Fax: +56 2 3400 399
E-mail: mqch@merck.cl
www.merck-chemicals.cl

China

Merck China
Tech Service: 400 820 8872
Tel: +86 21 3222 4788
Fax: +86 21 6247 9680
E-mail: Bioteam@merck-china.com
www.merck-china.com

Colombia

Merck Colombia S.A.
Tel: +57 1 425 4770
Fax: +57 1 425 5407
E-mail: mcsa@merck.com.co
www.merck.com.co

Ecuador

Merck C.A.
Tel: +593 2 298 1677
Fax: +593 2 298 1644
E-mail: sicmerck@merck.com.ec
www.merck.com.ec

Guatemala

Merck S.A.
See information under Central America

Hong Kong

Onwon Trading Limited
Tel: +852 2757 7569
Fax: +852 2757 7211
E-mail: info@onwon.com.hk
www.onwon.com.hk

India

Merck Specialities Private Limited
Tel: +91 22 6663 9884
Tel: +91 22 6663 9800
Fax: +91 22 6663 9893
Fax: +91 22 6663 9894
E-mail: life.science@merck.co.in
www.merck-chemicals.com/in

Indonesia

Pt. Merck Tbk.
Toll Free: 0800 140 1253
Tel: +62 21 2856 5600
Fax: +62 21 2856 5603
E-mail: chemicals@merck.co.id

Israel

Mercury Scientific & Industrial
Products Ltd.
Tel: +972 3 938 7164
Fax: +972 3 938 7174
E-mail: mercury@mercury-ltd.co.il
www.mercury-ltd.co.il

Japan

Merck Ltd.
Tel: +81 0120 189 390
Fax: +81 0120 189 350
E-mail: service@merck.co.jp
www.merck4bio.jp

Korea

Merck Limited
Tel: +82 2 2185 3840, 3842
Fax: +82 2 2185 3870
E-mail: products@merck.co.kr
www.merck.co.kr

Malaysia

Merck Sdn Bhd
Tel: +6 03 7494 3688
Fax: +6 03 7491 0850
Fax: +6 03 7491 0860
E-mail: chemlab@merck-de.com.my
www.merck.com.my

Mexico

Merck S.A. de C.V.
Tel: +52 55 2122 1600 ext. 7502
Fax: +52 55 2122 1703
E-mail: ruben.martinez@merck.com.mx

Mexico City

Tel: +52 55 5208 5197
+52 55 5208 5198
+52 55 5208 8116
Fax: +52 55 5203 6229

New Zealand

Merck Limited
Toll Free: 0800 46 37 25
Tel: +64 06 356 7328
Fax: +64 06 356 7311
E-mail: info@merck.co.nz
www.merck-chemicals.co.nz

Pakistan

Merck (Private) Limited
Tel: +92 21 455 9210
Fax: +92 21 453 5294
E-mail: chem-query@merck.com.pk
www.merck.com.pk

Peru

Merck Peruana S.A.
Tel: +51 1 6187 500
Fax: +51 1 4372 955
E-mail: merck.peruana@merck.com.pe
www.merck.com.pe

Philippines

Merck Inc.
Tel: +63 2 815 4067
Tel: +63 2 815 4208
Fax: +63 2 815 4882
E-mail: tish.aligada@merck.ph
www.merck-chemicals.ph

Singapore

Merck Pte. Ltd.
Customer Hotline: +65 6890 6660
Tel: +65 6890 6660
Fax: +65 6890 6639
E-mail: chem@merck.sg
www.merck-chemicals.sg

South Africa

Merck Laboratory Supplies Division
Tel: +27 11 345 9000
National Contact No. 08600 63725
Fax: +27 11 345 9251
E-mail: labsupply@merck.co.za
www.merck.co.za

Taiwan

Merck Ltd.
Tel: +886 2 2742 2788
Fax: +886 2 2742 2766
E-mail: bioservice@merck.com.tw
www.merckbio.com.tw

Thailand

Merck Ltd.
Tel: +66 2 667 8333
Fax: +66 2 667 8338
E-mail: customercare@merck.co.th
www.merck.co.th

Venezuela

Merck S.A.
Tel: +58 21 2235 1379
Fax: +58 21 2237 7632
E-mail: mven@merck.com.ve
www.merck.com.ve

Vietnam

Merck Representative Office
Tel: +84 8 3842 0100
Fax: +84 8 3842 0130
E-mail: merckvnadmin@hcm.vnn.vn