Millex® (33 mm) Sterile Filter Unit with MF-Millipore™ Membrane

SLG50335S SLHA0335S SLAA0335S
SLG50335B SLHA0335B SLAA0335B

English Deutsch Ελληνικά
Français Español Nederlands
Italiano Português Dansk

PH24062, Rev. A, 01/13

AW-1301-01

Made in Ireland
Merck Millipore Ltd.
Tullagreen,
Carrigrohilly,
Co. Cork, IRL

© 2013 EMD Millipore Corporation. All rights reserved.

Millipore and Millex are registered trademarks of Merck KGaA, Darmstadt, Germany.
The M logo and MF-Millipore are trademarks of Merck KGaA.
All trademarks of third parties are the property of their respective owners.

English		
Millex® 33 mm Sterile Filter Unit with MF-Millipore™ Membrane		
SLG50335S (0.22 µm, 50/boi)	SLHA0335S (0.45 µm, 50/boi)	SLAA0335S (0.8 µm, 50/boi)
SLG50335B (0.22 µm, 250/boi)	SLHA0335B (0.45 µm, 250/boi)	SLAA0335B (0.8 µm, 250/boi)
● Single use only	● Sterile	● Non-pyrogenic
	● Contains no natural latex rubber	

Indications for Use/Purpose

Outside the U.S. and Japan, Millex®-GS/HA/AA filter units (Cat. Nos. SLO50335S, SLO50335B, SLHA0335S, SLHA0335B, SLAA0335S, and SLAA0335B) are intended for use as a syringe filter to sterilize and/or clarify low volume solutions in direct patient care and pharmacy admixture applications.

Within the U.S., the approved Millex® filter units for direct patient care and pharmacy admixture applications are Cat. Nos. SLO5M335S, SLHAM335S, and SLAAM335S.

Introduction

This document provides compatibility information, operating steps, and specifications for the MF-Millipore® (mixed esters of cellulose nitrate and cellulose acetate) family of sterile Millex® filter units. The Millex® filter unit's bidirectional support of the filter membrane enables users to filter aqueous solutions in either direction; forward (from the syringe into the container) or backward (from the container into the syringe). The Millex® filter unit removes microorganisms, particles, precipitates, and undissolved powders larger than the membrane's rated pore size. These single-use filter units consist of a membrane filter sealed in an acrylic housing. They are non-pyrogenic and non-toxic.

Applications

Typical research laboratory applications include the sterile filtration (GS) and/or clarification (GS/HA/AA) of protein solutions, tissue culture media, additives, buffers, and water. In the clinical pharmacy, Millex® filter units with MF-Millipore™ membrane can be used for the sterile filtration (GS) and/or clarification (GS/HA/AA) of small volumes of protein pharmaceuticals, diagnostic imaging agents, chemotherapeutics, aqueous solutions, or water during admixture preparation. Direct patient care applications include sterilization (GS) and/or particulate removal (GS/HA/AA) from epidural and other liquid anesthetics as well as from irrigation solutions used in ophthalmic, etc, and other surgical procedures.

Chemical Compatibility

The Millex® filter unit with MF-Millipore™ membrane is compatible with most aqueous solutions. Based on information from technical publications, materials suppliers, and laboratory tests, Merck Millipore Ltd. believes that the agents listed in the following chart are safe to use with Millex® filter units. However, because of the effects of variability in temperature, concentration, duration of exposure, and other factors outside of our control, Merck Millipore Ltd. does not provide or imply a warranty with respect to this information. For example, the drug list refers to the compatibility between the drug and filter materials. It does not refer to the specific drug or protein binding to the filter, or the potential loss of the active drug component during filtration. You should qualify specific drugs for protein binding or loss of drug component prior to use. Agents that are not listed should be tested prior to use.

Chemicals

Acetic acid (5%)	Hydrogen	Petroleum ether
Alconox® detergent (1%)	HYPQ (aqueous solution)	Phenol (aqueous solution)
Aliphatic ethers	Kerosene	Silicone oils
Ammonium hydroxide (6 N)	Lactic acid (50%)	Sodium carbonate
Benzyl alcohol (1%)	Lubrol® PX emulsifier (aqueous solution)	Sodium dodecyl sulfate
Boric acid (aqueous solution)	Lidocaine	Sodium chloride (2 M)
Formic acid (50%)	Mineral spirits	Sulfuric acid (6 N)
Freon® solvent (TF or PCA)	Nitrogen	Trifluoroacetic acid
Gasoline	Nonidet™-P 40 surfactant (aqueous solution)	Tween® 20 surfactant (aqueous solution)
Glycerine (glycerol)	Ozone	Urea (6 M)
Guanidine hydrochloride (6 M)	Pentane	Perchloroethylene
Guanidine thiocyanate (5 M)	Helium	Petroleum based oils
Helium	Hexane	Water (brine)
Hexane	Hydrochloric acid (1 N)	Water (deionized)

Active Drug Compounds

Drug	Ref. No.	Merck® Index 11th Edition Page No.	Comments	Drug	Ref. No.	Merck® Index 11th Edition Page No.	Comments
Aminophylline	477	76	water soluble	Hydrocortisone	4713	758	water soluble
Ampicillin	621	93	water soluble	21-glycol sodium succinate			
Aspartame	861	132	water soluble	Immunoglobulins	4837	780	water soluble
Bleomycin	1324	201	water soluble	Insulin	4887	789	water soluble
Caffeine	1635	248	water soluble	Isopterenolol	5105	821	water soluble
Cefazolin	1925	294	water soluble	Lidocaine	5359	863	water soluble
Cefoxitin	1938	297	water soluble	Mannitol	5629	901	water soluble
Cephatholin	1978	305	water soluble	Metronidazole	6079	968	water soluble
Cisplatin	2319	361	water soluble	Mitoxantrone	6131	978	water soluble
Colistin	2475	387	suspension	Minomycin	6133	979	water soluble
			+ surfactant	Minoxantone	6135	979	water soluble
Cytarabine	2790	437	water soluble	Moxalactam	6201	991	water soluble
Dactinomycin	2804	441	water soluble	Nitroglycerin	6528	1045	water soluble
Daurorubicin	2825	445	water soluble	Norepinephrine	6612	1058	water soluble
Dexamethasone	2922	463	5% alcohol	Penicillin G potassium	7041	1123	water soluble
Dobutamine	3396	535	water soluble	Phenobarbital	7201	1149	water soluble
Dopamine	3415	538	water soluble	Piperacillin	7430	1184	water soluble
Doxorubicin	3428	540	water soluble	Plicamycin	7510	1198	water soluble
Ergonovine	3600	573	water soluble	Procaine	7768	1231	water soluble
Factor III	3873	616	water soluble	Protamine	7898	1253	water soluble
Factor IX	3874	616	water soluble	Streptokinase	8784	1390	water soluble
Fentanyl	3944	628	water soluble	Tobramycin	9413	1494	water soluble
Fluorouracil	4109	654	water soluble	Trimethoprim	9624	1528	water soluble
Folic Acid	4140	660	water soluble	Urokinase	9799	1555	water soluble
Furosemide	4221	674	water soluble	Vidarabine	9881	1569	water soluble
Gentamicin	4284	686	water soluble	Vinblastine	9887	1570	water soluble
Hemin	4563	733	water soluble	Vincristine	9891	1571	water soluble
Heparin	4571	735	water soluble				

How to Use the Millex® Sterile Filter Unit

This section lists warnings and cautions and provides steps to use the Millex® sterile filter unit.


WARNINGS:

- To ensure sterility, do not use this product if the package is damaged.
- Do not use this product as an in-line filter for intravenous fluid administration; it was not designed for long-term continuous use.
- Do not use with syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter unit and/or personal injury.

CAUTIONS:

- Do not use the Millex® filter unit to filter fluids at temperatures above 45 °C (113 °F).
- Do not use the Millex® filter unit to filter emulsions or suspensions because it was not designed for that purpose.
- Do not use the Millex® filter unit to filter solutions containing 5 milligrams (mg) or less of active drug materials unless binding studies have been performed.
- Do not use the same Millex® filter unit to filter solutions in both directions.
- Do not re-sterilize or reuse the Millex® filter unit, as Merck Millipore Ltd. cannot assure the sterility, integrity, and performance beyond a single use.

Procedure for Using the Millex® Sterile Filter Unit

Specifications for Millex® (33 mm) Filter Unit with MF-Millipore™ Membrane

Materials	
Membrane	Mixed esters of cellulose nitrate and cellulose acetate
	Pore size: Millex®-GS filter unit: 0.22 µm
	Millex®-HA filter unit: 0.45 µm
	Millex®-AA filter unit: 0.8 µm
Housing	Acrylic
Dimensions	
Inlet to outlet	26 mm (1.02 in.)
Diameter	33 mm (1.30 in.)
Filtration area	4.52 cm² (0.70 in²)
Temperature limit	45 °C (113 °F) maximum
Housing Pressure at 25 °C	10.3 bar (150 psi) inlet maximum
Filtration volume	10 mL to 100 mL
Hold-up volume	≤ 0.1 mL after air purge
Sterilization method	Ethylene oxide gas
Connections	Female Luer-Lok® inlet; male Luer-slip outlet
Flow rate at 2.1 bar (30 psi), 25 °C	Millex®-GS filter unit: ≥ 75 mL/min
	Millex®-HA filter unit: ≥ 180 mL/min
	Millex®-AA filter unit: ≥ 360 mL/min

Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Merck Millipore Ltd. ("Millipore") or an affiliate. Neither Merck Millipore Ltd. nor any of its affiliates assumes responsibility for any errors that may appear in this document.

Technical Assistance

For more information, contact the office nearest you. In the U.S., call 1-800-MILLIPORE (1-800-645-5476). Outside the U.S., go to our web site www.millipore.com/offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at www.millipore.com/techservice.

Standard Product Warranty

The applicable warranty for the products listed in this publication may be found at www.millipore.com/terms (within the "Terms and Conditions of Sale" applicable to your purchase transaction).

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