


Millipore®

Preparation, Separation,
Filtration & Testing Products

MOBIUS® SINGLE-USE FINAL FILL SOLUTIONS



The life science business of Merck
operates as MilliporeSigma in the
U.S. and Canada.

MERCK

Mobius® Single-Use Final Fill Solutions

The Mobius® single-use final fill offering allows end users to design robust, easy-to-use assemblies that ensure product quality and patient safety.

We have a comprehensive library of single-use assembly components which includes filters, tubing connectors, needles and pumps. We provide customers with high quality assemblies, technical expertise, and the services and support needed to minimize their risk.





Efficient and Cost-Effective

Whether you are looking to introduce single-use technologies into an existing process or you are building a new single-use facility, Mobius® products and solutions can help meet your evolving manufacturing needs.

Mobius® single-use products can increase your speed and efficiency, while also being cost effective. Our single use specialists will partner with you from initial specification of single-use components to optimization of your single-use operation.

At the heart of our Mobius® solutions are industry leading Durapore® and Millipore Express® filtration technologies.

Flexible Filling Productivity

Mobius® single-use final fill assemblies provide flexibility, allowing you to meet evolving manufacturing needs. They minimize costly cleaning validations and maximize change-over efficiency. In one customer case study, filling campaign setup, installation and production time was reduced from 36 to 12 hours through the use of a Mobius® single-use final fill solution versus a traditional stainless steel system.

Features and Benefits

- Reduced upfront capital investment
- Reduced risk of cross contamination and enhanced operator safety
- Flexibility and increased filling productivity
- Advanced single-use filtration technologies to maximize yields

PROCESS DESIGN CONSIDERATIONS


Sterility Assured

When claiming sterility, regulatory agencies worldwide recommend pre-use integrity testing and require post-use integrity testing of sterilizing-grade filters. Successful filter integrity tests are a critical link between your filter validation, successful processing, and product release.

Sterile filtration is a critical operation in final filling and there are multiple options for system design, either single, dual or redundant filtration. Regulatory requirements, costs, facility fit and risk profile are all important factors when designing your filtration system.

Questions to consider when designing your single-use final filtration assembly:

1. Do you require single, dual or redundant filtration?
2. Do you require the filter in or out of the isolator?




Filter Flushing Design Options: Bag vs. Millipak® Barrier Filter

Performing pre-use and post-use integrity testing (IT) requires careful flushing and manipulation of the assembly to ensure you do not compromise the downstream sterility. The diagrams below show possible filter configurations for in-line integrity testing. Millipak® Barrier filters, containing both hydrophobic and hydrophilic sterilizing-grade membrane, were developed to facilitate in-line integrity testing of sterilized liquid filtration systems.


Sterile Boundary Attribute					
	Sterility of Product Filter Maintained During Pre-Use Test	✓	✓	✓	✓
	Ability to retest	✓	✗	Limited by the size of the flush bag	
	Ability to Dry the Product Filter	✓	✗	✗	✗
	Simple Design	✓	✓	✓	✓

Sterile Filter Placement Considerations: Inside or Outside the Isolator/RABS

Sterile Filter in Grade A/ISO 5/Class 100 Space




Sterile Filter in Grade C/ISO 7/Class 10,000 Space




Complete SUS transferred to Grade A	
Benefits	Considerations
Low risk of fluid path contamination with majority of operations (venting, connections, etc.) occurring inside isolator/RABS	Sterile environment risk if performing open venting Handling challenges <ul style="list-style-type: none"> • Large number of components to transfer from Grade C into Grade A • Pre-use integrity testing challenge

Only filling needles transferred to Grade A	
Benefits	Considerations
Majority of operations (installation, flushing, etc.) occur outside isolator/RABS Potential to replace filter if integrity test fails	More potential points for microbial ingress Demonstration of assembly integrity more critical

CLOSED SYSTEM FILLING TRANSFER SET TO ISOLATOR



No.	Description
1	Sterile air line to attach Integritest® for in-line filter integrity test pre-use or post-use.
2	Closed and sterile NovaSeptum® sampling capabilities with a variety of container types and volumes ranging from 5mL to 1L.
3	Millipak® Final Fill Filter or Opticap® Product Filter
4	Millipak® Barrier filters have both a hydrophobic and hydrophillic membrane, allowing flushing of water and blowdown with air, maintaining sterility during filter integrity testing operations.
5	The Lynx® S2S connector has been subjected to rigorous microbiological validation studies to demonstrate its robust design and functionality. The results indicate that the connector provides sterility assurance during a sterile fluid transfer process. Please refer to tech brief TB1008N00 (Validation of the Lynx® S2S Connector for Sterile Fluid Transfer) for additional information.
6	Hanging header bags are available in a range of volumes from 1L to 50L, and in a range of inlet/outlet port configurations and sizes from 1/4 inch to 1 inch hosebarbs.
7	The header bag outlet uses an overmolded 1.5" TC Multi-tube Manifold or individual bag ports to accomodate 1 to 10 line filling machines or larger.



No.	Description
8	Pump tubings from several suppliers available in various sizes to fit a variety of peristaltic dosing pumps
9	The DPTE BetaBag® is a Rapid Transfer Port (RTP) bag that attaches to an Alpha port of an isolator so the single-use components can be transferred into the Class A space in a sterile manner after the isolator is cleaned and sterilized.
10	Single-use Filling needles are available in a range of sizes and styles to fit a variety of filling machine models and needle holders

FORMULATION

Mobius® Single-Use Mixing System

Key Features

Easy-to-use

Low shear impeller design

No particulate generation

Broad range of sizes available (10 L – 1000 L)

Jacket and load cells available as standard options in sizes 100 L – 1000 L

Standard and custom bag options available





BioReliance® Services

Formulation, Filtration & Final Filling

BioReliance® Services are customized to address your specific needs and backed by an industry leading team of scientists, engineers and validation specialists.

Standard Offerings

- Filter compatibility
- Bacterial retention
- Product specific bubble point
- Standard model solvents extractables
- Mobius® compatibility
- Mobius® extractables
- Mobius® leachables
- Patient toxicological assessment
- Particle shedding

Custom Offerings

- Shipping validation

INDUSTRY LEADING EXPERTISE AND QUALITY TO EASE YOUR IMPLEMENTATION

Our Mobius® team of experts will work with you to determine your specific final fill requirements. The result: a pre-sterilized, single-use final fill system that’s ready to use out of the box, featuring modular, fully-optimized assemblies. The assemblies are connected to create a closed, disposable system that meets your filling needs/requirements. These systems are replicable, scalable, sterile and simplify technology transfer across multiple sites.

As your sterile filtration partner, you will have the peace of mind that comes from working with the industry leader who has over 50 years of expertise in sterile filtration with market leading Durapore® (PVDF) and Millipore Express® (PES) membranes.

QUALITY

Component Qualification

We have a well-defined and documented process for qualifying single-use components into our library, to ensure the most stringent quality requirements are met.

Criteria
Gamma compatibility >40kGy
Functional testing
Regulatory statements (Animal Origin, Latex, BPA, etc.)
USP <88> Class VI
USP <85> Endotoxin
USP <788> Particulates
USP <661> Physiochemical
Shelf life > 2 years
Sterility per ANSI/AAMI/ISO 11737
Bacteriastasis/Fungistasis
Bioburden


Mobius® Assembly Leak Testing

- Gold certified assemblies are 100% tested
- Testing performed using a pressure decay method

Packaging and Transport Validation

Mobius® packaging protects assemblies from the rigors of typical shipping conditions. Representative Mobius® Final Fill assemblies have been proven to be integral after being subjected to a packaging validation per the International Safe Transit Association (ISTA) 2A Procedure.

Particulate Mitigation and Control Strategy





Millipore®

Preparation, Separation,
Filtration & Testing Products

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

MerckMillipore.com

To place an order or receive technical assistance

In Europe, please call Customer Service:

France: 0825 045 645

Spain: 901 516 645 Option 1

Germany: 069 86798021

Switzerland: 0848 645 645

Italy: 848 845 645

United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: [MerckMillipore.com/offices](https://www.merckmillipore.com/offices)

For Technical Service visit: [MerckMillipore.com/techservice](https://www.merckmillipore.com/techservice)

MerckMillipore.com

© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Mobius, Durapore, Millipore Express, Millipak, Millidisk, Lynx, Opticap, NovaSeptum and BioReliance are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

Lit. No. BR1144EN Ver. 1.0
2017 - 02308
10/2018