


Astec Cellulose DMP

Efficient, Rugged, and Economical Columns
for Chiral HPLC & SFC



Chiral selectivity in normal phase, polar organic, and SFC modes

Efficient, rugged, reproducible, and scalable

Low backpressure

Essential component of chiral column screening

Economically priced

Astec Cellulose DMP

Astec Cellulose DMP is a chiral stationary phase (CSP) comprising spherical, high-purity porous silica coated with DMPC (3,5-dimethylphenylcarbamate)-derivatized cellulose, and packed in analytical to preparative size HPLC columns. It separates a wide range of chiral compounds under normal phase, polar organic, and SFC conditions, with high efficiency, high loading capacity, and excellent column lifetime. Performance is comparable to other DMPC-derivatized cellulose CSPs, but the Astec Cellulose DMP columns are offered at a substantially lower price.

Astec Cellulose DMP is complementary to the other Astec CSPs, including CHIROBIOTIC®, CYCLOBOND™, and the P-CAP™ product lines, and a must-have for every chiral HPLC or SFC screening protocol.

Key Features and Application Areas

- Classic DMPC-cellulose chiral selectivity
- Efficient, rugged, reproducible, and scalable
- Low backpressure
- Ideal for chiral analysis in the pharmaceutical industry and for small analytes in chemical and environmental areas
- Routine chiral column method development screening protocols
- Approximately one-half the cost of most DMPC-cellulose columns

What is Cellulose?

The polysaccharide cellulose is a naturally occurring, optically active, linear polymer comprising hundreds to thousands of D-(+)-glucose units joined by β 1,4-glycosidic bonds. The long polysaccharide chains form rope-like bundles held together via multiple hydrogen bonds between proximate hydroxyl groups. In 1973, Hesse and Hagel described the enantioselective properties of microcrystalline cellulose triacetate (1). In the mid-1980's, Okamoto and colleagues published their work that led to the use of derivatized cellulose adsorbed onto silica as chiral HPLC stationary phases (2,3). Since then, the polysaccharides, particularly cellulose and amylose, have become the most commercially successful class of CSPs.





Figure 1. DMPC-Cellulose Structure



Cellulose is a linear polymer of D-(+)-glucose units linked by β 1,4-glycosidic bonds. This figure shows the cellulose tris(3,5-dimethylphenylcarbamate) derivative used in Astec Cellulose DMP.

How do Cellulose-Based CSPs Separate Enantiomers?

Derivatized cellulose-based CSPs, like Astec Cellulose DMP, owe their high enantioselectivity to the large number of chiral centers in the polysaccharide backbone and to its highly-ordered structure. The shape of the pockets formed by the intertwined chains provides chiral discrimination based on molecular shape. Derivatives at the 2, 3, and 6-position hydroxyls confer additional enantioselectivity. The dimethylphenyl carbamate derivative (Figure 1) separates a wide range of enantiomers; the phenyl ring and carbamate groups provide π - π and hydrogen bonding interactions, respectively, with predisposed analytes.

What Makes Astec Cellulose DMP Unique?

Astec Cellulose DMP is unique in offering classic DMPC-cellulose chiral selectivity, but at approximately one-half the cost of most DMPC-cellulose columns. It is an ideal component of any chiral column screening portfolio, and should be investigated as an alternative to higher-priced DMPC-cellulose columns for existing methods. The cost savings are especially dramatic when comparing preparative column dimensions.

Normal Phase Chiral Separations

DMPC-derivatized cellulose is commonly run in normal phase mode. Typical normal phase mobile phases are hexane or heptane with ethanol or IPA as polar modifiers. The performance of Astec Cellulose DMP in terms of selectivity and compatibility under normal phase conditions meets or exceeds competitive phases of similar composition. Figures 2 - 4 show the resolution of three racemates on Astec Cellulose DMPC and two higher-priced competitive phases.

Comparable to Other DMPC-Cellulose Columns

For the compounds we have tested, Astec Cellulose DMP columns provide similar retention and selectivity, but lower backpressure and higher efficiency compared to competitive columns. Astec Cellulose DMP is not a clone of other DMPC-derivatized cellulose columns on the market, but selectivity and retention is similar enough to use it instead of these columns

in your chiral column HPLC and SFC screening protocols. It also should be investigated as a possible replacement for these columns in established methods. The data in **Table 1** compares resolution, pressure, and selectivity on Astec Cellulose DMP and the leading competitive column. Comparison also appears in **Figures 2 - 4** against two competitive columns.

Table 1. Comparison of Astec Cellulose DMP vs. Leading Competitor for Sample Set of Chiral Compounds

Compound	Mobile Phase*	Astec Cellulose DMP				Competitor D			
		Pressure (bar)	Selectivity	N (Peak 1)	Resolution	Pressure (Bar)	Selectivity	N (Peak 1)	Resolution
Alprenolol	A	16	3.0	9,528	12.8	19	1.4	5,963	3.5
Atropine	A	16	1.4	5,768	4.9	19	1.7	3,666	6.1
Benzoin	B	16	1.4	7,690	6.2	19	1.6	6,217	7.1
Diperodon	C	21	3.9	5,915	14.4	26	4.0	4,846	12.8
Etodolac	B	16	2.6	6,323	10.2	19	2.8	5,568	10.1
Hydroxyzine	A	16	1.2	5,477	1.9	19	1.2	4,173	2.1
Ketamine	A	16	1.2	9,506	2.4	19	1.2	8,172	2.3
Metoprolol	A	16	2.3	7,208	7.0	19	2.6	5,343	6.0
Mianserin	A	16	1.2	6,936	1.9	19	1.2	6,078	1.6
Proglumide	B	16	1.8	3,672	5.3	19	2.2	2,963	5.5
<i>trans</i> -Stilbene Oxide (TSO)	D	16	1.9	13,753	9.9	19	1.9	11,871	9.1
Tröger's Base	A	16	1.4	9,398	4.8	19	1.3	7,375	3.0

Columns: 15 cm x 4.6 mm I.D., 5 µm particles, flow rate: 0.5 mL/min., temp.: 25 °C (Note: Separations not optimized on either column)
 * Mobile Phase A: 10:90:0.1, IPA:heptane:DEA
 Mobile Phase B: 10:90:0.1, IPA:heptane:TFA
 Mobile Phase C: 0.1% w/v ammonium formate in methanol
 Mobile Phase D: 10:90, IPA:heptane

Figure 2. Tröger's Base - Competitive Comparison

columns: 15 cm x 4.6 mm I.D., 5 µm particles
 mobile phase: 10:90:0.1, IPA:heptane:DEA
 flow rate: 0.5 mL/min.
 temp.: 25 °C
 det.: UV, 230 nm
 inj.: 2 µL
 sample: Tröger's Base (2 mg/mL)









Figure 3. Metoprolol - Competitive Comparison

Conditions same as Figure 2 except:
sample: metoprolol

**Figure 4. Dipiperidon - Competitive Comparison**


Conditions same as Figure 2 except:
mobile phase: 0.1% ammonium formate in methanol
sample: dipiperidon



Normal Phase Applications

Figure 5. Alprenolol

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 μ m particles (51098AST)
mobile phase: 10:90:0.1, IPA:heptane:DEA
flow rate: 0.5 mL/min.
temp.: 25 $^{\circ}$ C
det.: UV, 254 nm
inj.: 2 μ L
sample: alprenolol (2 mg/mL in mobile phase)

**Figure 6. Atropine**

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 μ m particles (51098AST)
mobile phase: 10:90:0.1, IPA:heptane:DEA
flow rate: 0.5 mL/min.
temp.: 25 $^{\circ}$ C
det.: UV, 254 nm
inj.: 2 μ L
sample: atropine (2 mg/mL in mobile phase)




Figure 7. Benzoin

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 µm particles (51098AST)
mobile phase: 10:90:0.1, IPA:heptane:TFA
flow rate: 0.5 mL/min.
temp.: 25 °C
det.: UV, 254 nm
inj.: 2 µL
sample: benzoin (2 mg/mL in mobile phase)




Figure 10. Etodolac

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 µm particles (51098AST)
mobile phase: 10:90:0.1, IPA:heptane:TFA
flow rate: 0.5 mL/min.
temp.: 25 °C
det.: UV, 254 nm
inj.: 2 µL
sample: etodolac (2 mg/mL in mobile phase)




Figure 8. Hydroxyzine

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 µm particles (51098AST)
mobile phase: 10:90:0.1, IPA:heptane:DEA
flow rate: 0.5 mL/min.
temp.: 25 °C
det.: UV, 254 nm
inj.: 2 µL
sample: hydroxyzine (2 mg/mL in mobile phase)




Figure 11. Ketamine

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 µm particles (51098AST)
mobile phase: 10:90:0.1, IPA:heptane:DEA
flow rate: 0.5 mL/min.
temp.: 25 °C
det.: UV, 254 nm
inj.: 2 µL
sample: ketamine (2 mg/mL in mobile phase)




Figure 9. Proglumide

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 µm particles (51098AST)
mobile phase: 10:90:0.1, IPA:heptane:TFA
flow rate: 0.5 mL/min.
temp.: 25 °C
det.: UV, 254 nm
inj.: 2 µL
sample: proglumide (2 mg/mL in mobile phase)





Figure 12. Warfarin

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 µm particles (51098AST)
mobile phase: 100:0.2:0.1, CH₃OH:acetic acid:TEA
flow rate: 0.5 mL/min.
temp.: 25 °C
det.: UV, 278 nm
inj.: 2 µL
sample: warfarin (2 mg/mL in mobile phase)



Astec Cellulose DMP

Polar Organic Mode Chiral Separations

Polar organic mode (POM) mobile phases comprise methanol or acetonitrile, often with acids, bases, or salts added to control peak shape and retention of certain sensitive analytes. The benefit of POM is realized when dealing with compounds that are poorly soluble in non-polar normal phase mobile phases. For preparative chiral applications, solubility is especially important; analyte concentration per injection influences the throughput. Astec Cellulose DMP operates in POM to permit choice of mobile phase based on analyte solubility. A few of the racemates resolved on Astec Cellulose DMP with different POM mobile phases appear in **Table 2** and **Figures 13-15**.

Not only does Astec Cellulose DMP provide excellent resolution in POM and NP modes, it also is rugged enough to hold up to repeated NP-POM-NP-POM cycles without loss of performance. **Figure 16** demonstrates the same Astec Cellulose DMP column used alternately in normal phase and POM mobile phases. Performance values appear in the figure.

Figure 13. Homatropine

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 μ m particles (51098AST)
 mobile phase: 0.1% v/v DEA in methanol
 flow rate: 0.5 mL/min.
 temp.: 25 $^{\circ}$ C
 det.: UV, 230 nm
 inj.: 2 μ L
 sample: homatropine (2 mg/mL in mobile phase)




Figure 14. Indapamide

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 μ m particles (51098AST)
 mobile phase: 0.1% w/v ammonium formate in methanol
 flow rate: 0.5 mL/min.
 temp.: 25 $^{\circ}$ C
 det.: UV, 230 nm
 inj.: 2 μ L
 sample: indapamide (2 mg/mL in mobile phase)




Figure 15. Ketoconazole

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 μ m particles (51098AST)
 mobile phase: acetonitrile
 flow rate: 0.5 mL/min.
 temp.: 25 $^{\circ}$ C
 det.: UV, 230 nm
 inj.: 2 μ L
 sample: ketoconazole (2 mg/mL in mobile phase)




Table 2. Performance of Astec Cellulose DMP in Polar Organic Mode (POM)

Mode:	POM - Methanol		POM - LC-MS Conditions		POM - Acetonitrile (no additives)	
Mobile Phase:	0.1% v/v DEA in Methanol		0.1% w/v Ammonium Formate in Methanol		100% Acetonitrile	
Compound	t _R Peak 1 (min.)	Selectivity	t _R Peak 1 (min.)	Selectivity	t _R Peak 1 (min.)	Selectivity
Diperodon	6.63	3.94	5.86	3.92	9.76	8.08
Homatropine	4.27	2.02	3.68	2.00		
Indapamide	4.71	2.32	4.66	2.26		
Ketoconazole					18.47	1.19
Mianserin	6.06	1.33	6.15	1.43	5.28	1.26
Tröger's Base	7.37	1.14	7.41	1.17		
Warfarin*	6.16	1.31				

Column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 μ m particles (51098AST), flow rate: 0.5 mL/min., temp.: 25 $^{\circ}$ C


*Warfarin mobile phase: 100:0.2:0.1, methanol:acetic acid:TEA

Figure 16. Stable Performance After Repeated NP-POM Cycles

columns: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 µm particles (51098AST)
 mobile phase (normal phase): 10:90, IPA:heptane
 mobile phase (POM): 0.1% w/v ammonium formate in methanol
 flow rate: 0.5 mL/min.
 temp.: 25 °C
 sample: TSO (normal phase) or mianserin (POM)
 det.: UV, 254 nm
 inj.: 2 µL


Run 1: NP

Pressure: 18 bar
 Selectivity: 1.86
 N (peak 1): 13,753
 Resolution: 9.92




Run 2: POM

Pressure: 23 bar
 Selectivity: 1.33
 N (peak 1): 11,792
 Resolution: 3.30




Run 3: NP

Pressure: 18 bar
 Selectivity: 1.80
 N (peak 1): 13,575
 Resolution: 9.38




Run 4: POM

Pressure: 23 bar
 Selectivity: 1.38
 N (peak 1): 11,898
 Resolution: 3.85




Analyte Structure:

trans-Stilbene Oxide



Mianserin



SFC Chiral Separations

SFC (supercritical fluid chromatography) is gaining in popularity, especially for chiral separations, due to its speed advantages over HPLC. The CO₂ is readily removed from the eluate, which makes it ideal for prep. DMPC-derivatized cellulose is widely used for chiral separations by SFC, both analytical and prep. The Astec Cellulose DMP works well in SFC mode, providing rapid separations with excellent selectivity. For example, **Figure 17** shows the SFC separation of a mixture of six diastereomers of a single compound in less than four minutes on Astec Cellulose DMP. As a testament to its utility as a screening tool, the durability of Astec Cellulose DMP permits rapid (ballistic) gradients of methanol, ethanol or IPA in CO₂ with long column lifetime and low backpressure, but without significant column bleed (**Figure 18**).

Figure 17. Rapid SFC Separations

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 µm particles (51098AST)
 mobile phase: 10% ethanol in CO₂
 flow rate: 3 mL/min.
 temp.: 35 °C
 pressure: 100 bar
 det.: UV, 220 nm
 inj.: 5 µL
 sample: six diastereomers of a single compound (proprietary drug substance)





Figure 18. Stability and Resolution under Ballistic SFC Gradients

Conditions same as Figure 17 except:
 gradient: 5-65% methanol in CO₂; hold 1 min.
 flow rate: 4 mL/min.
 sample: warfarin, 2 mg/mL



SFC data for Figures 17 and 18 kindly provided by Dr. Christina Kraml, Lotus Separations, LLC, Princeton, NJ.



Astec Cellulose DMP

Ideal for Preparative Applications

A characteristic of the polysaccharide-based CSPs that has contributed to their popularity is their utility for preparative applications. Although there are many factors to consider in a preparative scale up, the particle contribution comes primarily from the amount of available stationary phase and to what degree it resolves the enantiomers. When we designed the Astec Cellulose DMP, we sought to achieve the high sample loading and throughput that chiral chromatographers have come to expect. The sample loading (mg per injection) of Astec Cellulose DMP is comparable to Competitor D (Figure 19). An example of the scale-up of an analytical separation (4.6 mm I.D. column) to a preparative scale (21.2 mm I.D.) is shown in Figure 20 for the anti-Alzheimer's drug BAY 73-6691.

Figure 19. Loading Capacity

columns: 15 cm x 4.6 mm I.D., 5 μ m particles
 mobile phase: 10:90, IPA:hexane
 flow rate: 1 mL/min.
 temp.: 28 $^{\circ}$ C
 det.: UV, 210 nm
 inj.: 100 μ L
 sample: *trans*-stilbene oxide (TSO), 0.05 – 50 mg/mL




Figure 20. Scale-Up Example


columns: Astec Cellulose DMP, 5 μ m particles
 mobile phase: 80:10:10, heptane:methyl *tert*-butyl ether (MTBE):ethanol
 temp.: 25 $^{\circ}$ C
 det.: UV, 230 nm
 sample: BAY 73-6691 in mobile phase

Analytical:

dimensions: 15 cm x 4.6 mm I.D. (51098AST)
 flow rate: 0.5 mL/min.
 inj.: 10 μ L (2 mg/mL)

Preparative:

dimensions: 25 cm x 21.2 mm I.D. (51103AST)
 flow rate: 13 mL/min.
 inj.: 5000 μ L (3.3 mg/mL)



Rugged Columns, Stable Phase Chemistry, No Memory Effect

Although Astec Cellulose DMP is a coated phase, the DMPC-cellulose is held securely onto the silica surface. **Figure 21** shows chromatograms after long-term use. The column was exposed to over 13 liters of mobile phase, nearly 2,000 injections and 10 days of continuous operation without significant change in any chromatographic parameter. Also note the discussion about operation in polar organic mode in **Table 2** and **Figure 16**. The Astec Cellulose DMP has little memory effect when switching between mobile phase systems. This stability and lack of memory effect make the Astec Cellulose DMP columns even more of a value; not only are they considerably less expensive than competitive phases, their ruggedness means they maintain their high performance for long-term operation. (Avoid using polar, aprotic, and halogenated solvents.)

Figure 21. Stability Demonstration


column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 μ m particles (51098AST)
 mobile phase: 10:90, IPA:heptane
 flow rate: 1 mL/min.
 temp.: 28 °C
 det.: UV, 210 nm
 inj.: 2 μ L
 sample: TSO, 1 mg/mL

Initial performance

N = 10,800; AF10 = 0.94;
 Rs = 11.1; Press. = 26 bar

After ~5000 column volumes (~13 L of mobile phase)

N = 10,400; AF10 = 0.95;
 Rs = 10.8; Press. = 27 bar





Reproducibility for Reliable Methods

No separation is useful if it is not reproducible. We have designed Astec Cellulose DMP to have the reproducibility you require for method validation. The example in **Figure 22** shows columns from three different production lots of Astec Cellulose DMP in both normal phase and polar organic modes.

Figure 22. Reproducibility Demonstration

column: Astec Cellulose DMP, 15 cm x 4.6 mm I.D., 5 μ m particles (51098AST)
 mobile phase (normal phase): 10:90, IPA:heptane
 mobile phase (POM): 0.1% w/v ammonium formate in methanol
 flow rate: 0.5 mL/min.
 temp.: 25 °C
 sample: TSO (normal phase) or mianserin (POM)
 det.: UV, 230 nm
 inj.: 2 μ L



Astec Cellulose DMP

Complementary to the Other Astec CSPs

Cellulose DMP is the latest addition to the Astec line of high-quality chiral stationary phases manufactured by Sigma-Aldrich. Other Astec HPLC CSPs include CHIROBIOTIC, CYCLOBOND, P-CAP, P-CAP-DP, CLC-L, and CLC-D. Astec CHIRALDEX and Supelco DEX are the market leaders for chiral GC separations. The HPLC CSPs are complementary to each other in terms of selectivity and mobile phase compatibility. This means:

- It is likely that at least one Astec CSP will give the necessary selectivity. Incorporating Astec Cellulose DMP, CHIROBIOTIC and CYCLOBOND in your HPLC or SFC screening protocol will give at least 90% success rate.

- Multiple Astec CSPs may provide the necessary enantioselectivity, but one may operate in a preferred mobile phase system, one that is more compatible with the detection mode, or provides better analyte solubility or shorter retention time, or many other considerations. For example, polar ionic CHIROBIOTIC mobile phases are ideal for LC/ESI-MS.
- Different CSPs may provide reversal of elution order, a useful attribute for prep and for low-level detection of the presence of an unwanted enantiomer in large excess of the opposite enantiomer (trace analysis).

The wide choice of CSPs in the Astec line means they cover many different areas of interest within chiral chromatography. Some of these areas are captured in **Table 3**.

Table 3. Techniques, Applications, and Fields of Use for Astec Chiral Phases

	Astec Cellulose DMP	Astec CHIROBIOTIC®	Astec CYCLOBOND™	Astec P-CAP™	Astec CLC™	Astec CHIRALDEX®	Supelco DEX™
Routine Chiral Column Screening	Highly suitable	Highly suitable	Marginally suitable, or limited to specific applications	Marginally suitable, or limited to specific applications	Not suitable nor recommended	Highly suitable	Highly suitable
Normal Phase HPLC	Highly suitable	Marginally suitable, or limited to specific applications	Marginally suitable, or limited to specific applications	Highly suitable	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended
SFC	Highly suitable	Marginally suitable, or limited to specific applications	Marginally suitable, or limited to specific applications	Marginally suitable, or limited to specific applications	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended
Reversed-Phase HPLC	Not suitable nor recommended	Highly suitable	Highly suitable	Not suitable nor recommended	Marginally suitable, or limited to specific applications	Not suitable nor recommended	Not suitable nor recommended
Hydrophilic Interaction HPLC (HILIC)	Not suitable nor recommended	Highly suitable	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended
Polar Organic Mode (HPLC)	Highly suitable	Highly suitable	Marginally suitable, or limited to specific applications	Highly suitable	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended
Polar Ionic Mode (HPLC)	Not suitable nor recommended	Highly suitable	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended
Ligand Exchange Mode HPLC	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Highly suitable	Not suitable nor recommended	Not suitable nor recommended
Gas Chromatography (GC)	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Highly suitable	Highly suitable
Prep (LC and/or SFC)	Highly suitable	Highly suitable	Marginally suitable, or limited to specific applications	Highly suitable	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended
Polar/Ionic Analytes	Marginally suitable, or limited to specific applications	Highly suitable	Marginally suitable, or limited to specific applications	Marginally suitable, or limited to specific applications	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended
Amino Acids, Peptides	Not suitable nor recommended	Highly suitable	Marginally suitable, or limited to specific applications	Marginally suitable, or limited to specific applications	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended
Non-Aromatic Organic Acids	Not suitable nor recommended	Not suitable nor recommended	Marginally suitable, or limited to specific applications	Not suitable nor recommended	Highly suitable	Not suitable nor recommended	Not suitable nor recommended
Mass Spec (LC/ESI)	Marginally suitable, or limited to specific applications	Highly suitable	Marginally suitable, or limited to specific applications	Marginally suitable, or limited to specific applications	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended
Mass Spec (LC/APCI)	Highly suitable	Highly suitable	Marginally suitable, or limited to specific applications	Marginally suitable, or limited to specific applications	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended
Mass Spec (GC/MS)	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Not suitable nor recommended	Highly suitable	Highly suitable
Bioanalysis (drugs in biological fluids)	Marginally suitable, or limited to specific applications	Highly suitable	Marginally suitable, or limited to specific applications	Marginally suitable, or limited to specific applications	Not suitable nor recommended	Marginally suitable, or limited to specific applications	Marginally suitable, or limited to specific applications

	Highly suitable
	Marginally suitable, or limited to specific applications
	Not suitable nor recommended

Ordering Information

Astec Cellulose DMP Columns

Particle Size (µm)	Length (cm)	I.D. (mm)	Cat. No.
5	15	2.1	51100AST
5	25	2.1	51101AST
5	10	4.6	51097AST
5	15	4.6	51098AST
5	25	4.6	51099AST
5	25	10	51102AST
5	25	21.2	51103AST

Test Mix

Description	Package Size	Cat. No.
Test mix for polysaccharide-based chiral HPLC columns, contains <i>trans</i> -stilbene oxide (TSO) and 1,3,5-tri- <i>tert</i> -butylbenzene (V_0 marker), each 30 µg/mL in hexane	1 mL	40119-U

Guards

Description	Length (cm)	I.D. (mm)	Cat. No.
Guard column kit for 2.1 mm I.D. columns. Contains holder, one cartridge, 2" x 1/16" tubing, two nuts, and two ferrules	2	2.1	51105AST
Replacement guard cartridges for 2.1 mm I.D. columns, Pk of 2 Requires holder 59660-U, sold separately	2	2.1	51104AST
Guard column kit for 4.6 mm I.D. columns. Contains holder, one cartridge, 2" x 1/16" tubing, two nuts, and two ferrules	2	4	51107AST
Replacement guard cartridges for 4.6 mm I.D. columns, Pk of 2. Requires holder 59660-U, sold separately	2	4	51106AST
Replacement guard cartridge for 10 mm I.D. columns. Requires holder 567499-U, sold separately	1	10	51108AST
Replacement guard cartridge for 21.2 mm I.D. columns. Requires holder 581392-U, sold separately	1	21.2	51109AST

Guard Cartridge Holders

Description	Cat. No.
Holder for 2 cm x 2.1 mm & 4 mm I.D. Guard Cartridges	59660-U
Holder for 1 cm x 10 mm I.D. Guard Cartridges	567499-U
Holder for 1 cm x 21.2 mm I.D. Guard Cartridges	581392-U

Please visit sigma-aldrich.com/chiral to see our complete line of products and services for chiral chemistry and chiral separations.

References

- Hesse, G.; Hagel, R. *Chromatographia* **1973**, 6(6), 277 – 280.
- Okamoto, Y.; Kawashima, M.; Yamamoto, K.; Hadata, K. *Chem. Lett. (The Journal of the Chemical Society of Japan)* **1984**, 739 – 742.
- Okamoto, Y.; Kawashima, M.; Hadata, K. *J. Amer. Chem. Soc.* **1984**, 106, 5357 – 5359.

Trademarks

CHIRALDEX, CHIROBIOTIC, CLC, CYCLOBOND, DEX and P-CAP are trademarks of Sigma-Aldrich Biotechnology LP.



Astec Cellulose DMP

Sigma-Aldrich® Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Hungary

Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9063
Fax: (+36) 1 269 6470

India

Telephone

Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8003

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 4087 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8000

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com

*Accelerating Customers'
Success through Innovation and
Leadership in Life Science,
High Technology and Service*

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAFC® (800) 244-1173**
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com