Instructions for Use

Methyl Green-Pyronin Stain

Procedure No. HT70

Intended Use

Methyl Green-Pyronin Stain is for the histological identification of Plasma Cells and RNA in Tissue Sections. Methyl Green-Pyronin Stain is for "In Vitro Diagnostic use." For professional use only. This manual qualitative procedure detects the presence of plasma cells and RNA in tissue sections and cytologic preparations of human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of certain clinical conditions such as plasma cell endometritis or chronic endometritis.

Methyl Green-Pyronin (MGP) is used to demonstrate plasma cells and RNA in tissue sections and cytologic preparations. The procedure is a simplified method for use with formalin fixed tissue or alcohol-ether fixed smears. Pyronin stains the cytoplasm of plasma cells and most nucleoli red. Methyl green stains DNA bluish-green, but the presence of crystal violet in solution makes this solution unsuitable for the reliable demonstration of DNA.

Reagents

Methyl Green-Pyronin (MGP) Solution (Cat. No. HT70116-500ML)

Methyl green (certified), C.I. 42590, 0.012%, pyronin Y (certified), 0.01% methanol, 0.75% in deionized water.

Special Materials Required but Not Provided

- Positive control slides, such as normal lymph node or tonsil should be included in each run
- Microscope
- Microscope slides, coverslips and staining dishes
- Acetone, ACS Reagent
- Xylene or xylene substitute

Storage and Stability

Methyl Green-Pyronin (MGP) Solution should be stored at room temperature (18–26°C). Reagent label bears expiration date.

Preparation

Methyl Green-Pyronin Solution is ready to use.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Formalin fixed paraffin sections cut at 5–6 microns may be used.^{1,2}

Procedure

1. Deparaffinize tissue and hydrate to deionized water.
2. Stain in Methyl Green-Pyronin (MGP) Solution for 5 minutes.
3. Rinse thoroughly in deionized water.
4. Dehydrate and clear using acetone, acetone/xylene (1:1) and xylene.^{4,5} Mount.

Acetone:	2 or 3 quick dips
Acetone/Xylene:	2 changes of 1–2 minutes each
Xylene:	2 changes of 1–2 minutes each

Performance Characteristics

Nuclei	Purple to bluish
Nucleoli	Red
Cytoplasm of plasma cells	Red
RNA	Red

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT70	Methyl Green-Pyronin (MGP) Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm of plasma cells	3 of 3	3 of 3	3 of 3	3 of 3
		RNA	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT70116:

Not a hazardous substance or mixture.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Brown G: A modification of methyl green-pyronin stain for plasma cells and RNA in formalin fixed tissue. J Histotech 2:19, 1979
2. Diagnostic Cytology and Histopathologic Bases, 3rd ed., Vol 2, Edited by LG Koss, Lippincott, Philadelphia, 1979
3. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 151
4. Luna, Lee G., Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, Third Edition, McGraw-Hill Book Company, New York 1968
5. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 106

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 3.0	2014
Rev. 4.0	2016
Rev. 5.0	2023

Transferred to new template with current branding. Specified for professional use in intended use and precautions. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards. Added CH-REP information.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Bibliografia

1. Brown G: A modification of methyl green-pyronin stain for plasma cells and RNA in formalin fixed tissue. J Histotech 2:19, 1979.
2. Diagnostic Cytology and Histopathologic Bases, 3rd ed., Vol 2, Edited by LG Koss, Lippincott, Philadelphia, 1979.
3. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 151.
4. Luna, Lee G., Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, Third Edition, McGraw-Hill Book Company, New York 1968.
5. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 106.

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Referencias

- Brown G: A modification of methyl green-pyronin stain for plasma cells and RNA in formalin fixed tissue. J Histotech 2:19, 1979
- Diagnostic Cytology and Histopathologic Bases, 3rd ed., Vol 2, Edited by LG Koss, Lippincott, Philadelphia, 1979
- Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis, MO, 1980, p 151
- Luna, Lee G., Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, Third Edition, McGraw-Hill Book Company, New York 1968
- Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis, MO, 1980, p 106

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- Brown G: A modification of methyl green-pyronin stain for plasma cells and RNA in formalin fixed tissue. J Histotech 2:19, 1979
- Diagnostic Cytology and Histopathologic Bases, 3rd ed., Vol 2, Edited by LG Koss, Lippincott, Philadelphia, 1979
- Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 151
- Luna, Lee G, Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, Third Edition, McGraw-Hill Book Company, New York 1968
- Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 106

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- Brown G: A modification of methyl green-pyronin stain for plasma cells and RNA in formalin fixed tissue. J Histotech 2:19, 1979
- Diagnostic Cytology and Histopathologic Bases, 3rd ed., Vol 2, Edited by LG Koss, Lippincott, Philadelphia, 1979
- Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 151
- Luna, Lee G., Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, Third Edition, McGraw-Hill Book Company, New York 1968
- Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 106

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Referenties

1. Brown G: A modification of methyl green-pyronin stain for plasma cells and RNA in formalin fixed tissue. J Histotech 2:19, 1979
2. Diagnostic Cytology and Histopathologic Bases, 3rd ed., Vol 2, Edited by LG Koss, Lippincott, Philadelphia, 1979
3. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis, MO, 1980, p 151
4. Luna, Lee G., Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, Third Edition, McGraw-Hill Book Company, New York 1968
5. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis, MO, 1980, p 106

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Reference

1. Brown G: A modification of methyl green-pyronin stain for plasma cells and RNA in formalin fixed tissue. J Histotech 2:19, 1979
2. Diagnostic Cytology and Histopathologic Bases, 3rd ed., Vol 2, Edited by LG Koss, Lippincott, Philadelphia, 1979
3. Theory and Practice of Histochemistry. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 151
4. Luna, Lee G., Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, Third Edition, McGraw-Hill Book Company, New York 1968
5. Theory and Practice of Histochemistry. Edited by DC Sheehan and BB Hrapchak, 2nded. Mosby, St. Louis, MO, 1980, p 106

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MDSS CH GmbH Laurenzenovstadt 61 5000 Aarau Switzerland