Instructions for Use

Mucicarmine Stain

Procedure No. HT30

Intended Use

Sigma-Aldrich Mucicarmine Staining reagents are intended for the histological staining of tissue mucins with epithelial origin. Mucicarmine staining reagents are for "In Vitro Diagnostic Use". For professional use only. The data obtained from this manual, qualitative procedure serves only as an aid to diagnosis of human specimens and should be reviewed in conjunction with other clinical diagnostic tests or information.

Tissue acid-mucopolysaccharides may be demonstrated by a variety of techniques. Mayer's mucicarmine procedure, as modified by Southgate,¹ achieves excellent contrast between red staining mucin and other cellular components. Aluminum, contained in the formulation, is believed to form a chelate complex with mucin, to which carmine attaches by dye-lake formation.² The stain displays a specificity towards mucins of epithelial origin, whereas mucins of fibroblastic origin stain poorly.³ The procedure described by Sigma-Aldrich is similar to that of Southgate. A stable solution of tartrazine has been substituted for metanil yellow. Included is a mucicarmine stain technique for rapid staining in microwave oven.⁴⁻⁶

Reagents

Mucicarmine Stock Solution (Cat. No. HT3018-250ML or HT30116-500ML)

Carmine, certified, 10 g/L, C.I. 75470, aluminum hydroxide, 10 g/L, and aluminum chloride, 5 g/L, in ethanol, 50% v/v.

Special Materials Required but Not Provided

- Positive control slides, such as Mucin TISSUE-TROL™ (Cat. No. TTR008) should be included in each run.
- Tartrazine Solution (Cat. No. HT3024-120ML or HT3028-250ML)
- Reagent Alcohol
- Xylene or xylene substitute

Standard Procedure Only:

- Weigert's Iron Hematoxylin Set (Cat. No. HT1079-1SET)
- Staining Jars (Cat. No. S5641)

Microwave Procedure Only:

- Hematoxylin Solution, Gill No. 3 (Cat. No. GHS316; GHS332; GHS380; GHS3128)
- Scott's Tap Water Substitute Concentrate (Cat. No. S5134-6X100ML)
- Coplin Jar with vented lids
- Microwave Oven

Storage and Stability

Store Mucicarmine Stock Solution in refrigerator (2–8°C). Reagent is stable until expiration date shown on label.

Working Mucicarmine Solution is stable for 2–3 days when stored in refrigerator (2–8°C).

Preparation

Weigert's Iron Hematoxylin solution is prepared by mixing equal parts of Weigert's Iron Hematoxylin Solution, Part A and Weigert's Iron Hematoxylin Solution, Part B.

Working Mucicarmine Solution is prepared by diluting Mucicarmine Stock Solution 1:4 with tap water (e.g., 10 mL Mucicarmine Stock Solution plus 30 mL tap water).

Tartrazine solution and Hematoxylin Solution, Gill No. 3 are ready to use. Solutions are stable until the stated expiration date.

Scott's Tap Water Substitute working solution is prepared by mixing 1 part of Scott's Tap Water Substitute Concentrate with 9 parts deionized water. Solution is stable until the stated expiration date.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Mucicarmine TISSUE-TROL™ control slides are paraffin embedded human tissue containing collagen and muscle tissue and should be considered potentially infectious.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Any well-fixed paraffin sections cut at 5–6 microns may be used. Incorporate appropriate control slides.

Notes

Do not substitute deionized or distilled water for tap water when preparing the mucicarmine working solution.

Masking of mucin may occur if the nuclear stain or counterstains are too intense.

Procedure

Standard Procedure

- Prepare Working Weigert's Iron Hematoxylin Solution.
- Deparaffinize slides to deionized water.
- Stain in Working Weigert's Iron Hematoxylin Solution for **5 minutes**.
- Wash in running tap water for **5 minutes**.
- Stain in Working Mucicarmine Solution for **30 minutes** or longer at **room temperature**. Check control slide under microscope.
NOTE: The time may be decreased if staining is performed at 37°C or 56°C.
- Rinse in deionized water.
- Stain in Tartrazine Solution for **1–5 seconds**.
- Rinse slides, dehydrate through alcohol, clear in xylene and mount.

Microwave Procedure

- Deparaffinize slides and hydrate to deionized water.
- Place slides in **40 mL** of Hematoxylin Solution, Gill No. 3, contained in a plastic Coplin jar. Loosely cover with lid before placing in oven, or use Coplin jar lids with holes drilled into them.
- Microwave on **800 WATTS** for **10 seconds**.
- Rinse in running tap water for **1–2 minutes**.
- Blue in Scott's Tap Water Substitute working solution.
- Place slides in **40 mL** of Mucicarmine Working Solution.
- Microwave on **800 WATTS** for **20 seconds**. Mix gently with beral pipet or applicator stick. Let incubate for **5 minutes**.
- Rinse well in deionized water.
- Counterstain in Tartrazine Solution for **30 seconds** at **room temperature**.
- Rinse quickly in deionized water.
- Dehydrate quickly through alcohols and clear in xylene. Mount slides.

Performance Characteristics

Target Structure	Staining Result
Nuclei	Blue
Mucin	Deep Rose to Red
Capsule of Cryptococci	Deep Rose to Red
Other Tissue Elements	Yellow

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT301	Mucicarmine Stain	Mucin	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT3018, HT30116:

H226: Flammable liquid and vapor.

H315: Causes skin irritation.

H319: Causes serious eye irritation.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P240: Ground and bond container and receiving equipment.

P241: Use explosion-proof electrical/ ventilating/ lighting/ equipment.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 175
- Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, MO, 1980
- Selected Histochemical and Histopathologic Methods, SW Thompson, RD Hunt, Editors, CC Thomas, Springfield, IL, 1966, p 453
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 4.0	2016
Rev. 5.0	2022
Rev. 6.0	2022 Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards.
Rev. 7.0	2023 Revised Intended Use statement.
Rev 8.0	2025 Added CH-Rep information.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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The Initial M, TISSUE-TROL, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

1. Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 175
2. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, MO, 1980
3. Selected Histochemical and Histopathologic Methods, SW Thompson, RD Hunt, Editors, CC Thomas, Springfield, IL, 1966, p 453
4. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
5. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
6. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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Références

- Cellular Pathology Techniques, 4th ed. C.F.A. Culling, R.T. Addison, H.T. Barr, Editors Butterworth, London, 1985, p 175
- Theory and Practice of Histotechnology, 2nd ed. D.C. Sheehan, B.B. Hrapchak, Editors, Mosby, St. Louis, MO, 1980
- Selected Histochemical and Histopathologic Methods, SW Thompson, R.D. Hunt, Editors, C.C. Thomas, Springfield, IL, 1966, p 453
- Leong A.S.-Y., Milios J.: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J. Pathol 148:183, 1986
- Brinn N.T.: Rapid metallic histologic staining using the microwave oven. J. Histotechnol 6:125, 1983
- Valle S.: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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- Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, σελ. 175
- Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, MO, 1980
- Selected Histochemical and Histopathologic Methods, SW Thompson, RD Hunt, Editors, CC Thomas, Springfield, IL, 1966, σελ. 453
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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Reference

- Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 175
- Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, MO, 1980
- Selected Histochemical and Histopathologic Methods, SW Thompson, RD Hunt, Editors, CC Thomas, Springfield, IL, 1966, p 453
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Reference

- Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 175
- Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, MO, 1980
- Selected Histochemical and Histopathologic Methods, SW Thompson, RD Hunt, Editors, CC Thomas, Springfield, IL, 1966, p 453
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986.
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983.
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986.

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Literatura

- 1. Cellular Pathology Techniques, 4th ed. CFA Culling, RT Addison, HT Barr, Editors Butterworth, London, 1985, p 175
- 2. Theory and Practice of Histotechnology, 2nd ed. DC Sheehan, BB Hrapchak, Editors, Mosby, St. Louis, MO, 1980
- 3. Selected Histochemical and Histopathologic Methods, SW Thompson, RD Hunt, Editors, CC Thomas, Springfield, IL, 1966, p 453
- 4. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- 5. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- 6. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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MDSS GmbH
MDSS CH GmbH
Chemie GmbH