## User Guide

# Centrifree® Ultrafiltration Device

4104

**IVD** For in vitro diagnostic use

### Languages available

English . . . . .	2	Lietuviškai . . . . .	57
Français . . . . .	7	Magyar . . . . .	62
Italiano . . . . .	12	Čeština . . . . .	67
Deutsch . . . . .	17	Polski . . . . .	72
Español . . . . .	22	日本語 . . . . .	77
Português . . . . .	27	中文 . . . . .	82
Ελληνικά . . . . .	32	한국어 . . . . .	87
Nederlands . . . . .	37	Norsk . . . . .	92
Dansk . . . . .	42	Slovenčina . . . . .	97
Svenska . . . . .	47	Türkçe . . . . .	102
Latviski . . . . .	52		

## Introduction

Centrifree® devices rapidly and efficiently separate free from protein-bound microsolite in small volumes (0.15–1.0 mL) of serum, plasma, and other biological samples using a method called ultrafiltration. Accurate partitioning occurs in minutes without dilution, change in physiologic pH, ion composition, or unbound microsolite concentration. These devices contain low-adsorptive hydrophilic membranes and O-rings without plasticizers to ensure excellent recovery. Hold-up volume is 10 µL or less.

In contrast to dialysis, gel filtration, or charcoal adsorption, ultrafiltration provides improved accuracy when measuring free ligand concentration, binding capacity, or affinity constants. It eliminates time-consuming methodology, the need for specialized equipment, dilution errors, and shifts in binding equilibrium.

Ultrafiltration does not change free microsolite ligand concentration. Protein becomes selectively partitioned into a fraction of the sample volume (the concentrate), while free ligand passes essentially unhindered through the membrane along with solvent.

The laws of mass action and conservation of mass for ideal protein binding predict that free ligand concentration in the ultrafiltrate will remain constant, provided molar binding capacity and affinity are independent of total protein concentration. Results for several systems showing constant free ligand concentration in successive ultrafiltrate fractions support those predictions.


Changing free ligand concentration not caused by membrane retention or adsorption is evidence of altered capacity or affinity of binding proteins due to aggregation or other non-ideal protein-protein interactions. Depending upon the application, you can analyze the resulting filtrate for free ligand concentration either quantitatively or qualitatively.

## Intended Use

Centrifree® devices are non-sterile, disposable ultrafiltration devices for in vitro diagnostic use and are intended to separate free from protein-bound microsolite in small volumes (0.15–1.0 mL) of biological samples e.g. serum, urine, cerebrospinal fluid and other body fluids prior to in vitro diagnostic analysis. Device intended to be single-use and used by a laboratory professional.

## Centrifree® Device Components

The Centrifree® ultrafiltration device provides maximum efficiency for multiple sample processing. Each unit consists of a membrane and O-ring permanently sealed between the sample reservoir and support base. A removable filtrate cup is attached to the base. The Centrifree® device is intended for single use only.

## Materials Supplied

The following components are supplied with Centrifree® ultrafiltration devices.

- 50 ultrafiltration devices
- 50 reservoir caps
- 50 filtrate cups
- 50 filtrate cup caps

**NOTE:** The red reservoir caps are provided to prevent sample evaporation and pH change due to loss of CO<sub>2</sub>.

## Equipment Required

- Centrifuge with rotor adapter or carriers that can accept 17 × 100 mm tubes and is capable of 1,000–2,000 × g.

**NOTE:** For optimum solvent flow, use a fixed-angle centrifuge rotor.

- Pasteur or fixed-volume pipette for sample delivery.

## Storage and Stability

Refer to product label for storage conditions and shelf life.

## Guidelines for Use

- Serum, plasma, or other biological fluids may be used with Centrifree® devices. For more efficient flow rates, remove fibrin by centrifuging sample before loading into reservoir.
- The recommended sample volume is 0.15–1.0 mL.
- For best results, use a fixed-angle rather than a swinging-bucket centrifuge rotor and spin at 1,000–2,000 × g.  
**CAUTION:** Do not operate at relative centrifugal force above 2,000 × g.
- Conduct preliminary ligand adsorption and protein retention studies to ensure suitability of the system for the intended application.
- Centrifree® ultrafiltration devices are designed for use with biological fluids and aqueous solutions. Do not use devices with organic solvents.
- Determine the correct temperature for your application. Maintain the desired temperature by using a heated centrifuge or by prewarming the centrifuge rotor and chamber. In a fixed-angle rotor, the contents of the device reach thermal equilibrium with the rotor in 5–10 minutes. The rapid rate of ultrafiltration achieved by Centrifree® devices minimizes the effects of temperature on ligand binding.
- Do not use for concentration and recovery of macromolecules.
- Do not autoclave component parts.
- Do not reuse Centrifree® devices.
- Do not allow the membrane in the ultrafiltration devices to dry out once wet. If you are not using the device immediately after rinsing, leave fluid on the membrane until the device is used.
- The Ultracel® membrane in the device contains trace amounts of glycerine (~2 µL). If this substance interferes with analysis, rinse device with deionized water or 0.1 N NaOH until no more interference is noted. If 0.1 N NaOH is used for glycerine removal, rinse the device thoroughly with deionized water or buffer before use. Spin for at least 15 minutes to remove as much fluid as possible. To avoid dilution error caused by entrapped fluid (~10 µL) after rinsing, discard the first aliquot of ultrafiltrate.

## How to Use Centrifree® Ultrafiltration Devices

Before using a Centrifree® device, remove the red reservoir cap.

1. Hold reservoir angled at approximately 45° with pipette tip touching reservoir wall. Add sample solution smoothly in one even flow to avoid air locking.

**NOTE:** Avoid touching membrane with pipette tip.

2. Cap sample reservoir, then place device in a centrifuge rotor with 17 × 100 mm adapters. Counterbalance centrifuge with a similar device.

**NOTE:** Red reservoir cap should extend no more than 3–4 mm down over reservoir top. Further compression can cause prefiltration, which could subsequently introduce analytical errors if sample is not at desired operating temperature.

**NOTE:** Check Centrifuge clearance before use.

3. Equilibrate device and sample to temperature required by your application.
4. Spin device at 1,000–2,000 × g for required time to obtain desired filtrate volume. For spin time guidelines, refer to the “Ultrafiltration Rate” section.
5. Carefully remove device from centrifuge rotor and disconnect filtrate cup containing the ultrafiltrate. Cover filtrate cup with supplied cap until ultrafiltrate can be analyzed.

**WARNING:** If discarding used components, be sure to follow precautions for disposal of items contaminated with potentially infectious or hazardous biological material.


**NOTE:** Filter may not function correctly if allowed to dry out after wetting.

## Performance

The following sections discuss various performance characteristics including polarization control, ultrafiltration rate, membrane performance, pH control, and nonspecific adsorption.

## Polarization Control

Use of a fixed-angle rotor provides polarization control and minimizes the potential for artifacts due to protein-protein interactions. The angle counteracts buildup of retained protein at the membrane surface because this dense layer slides outward and accumulates at the edge of the membrane. In a swinging-bucket rotor, the polarization layer is compacted over the entire membrane surface, restricting solvent flow.

## Ultrafiltration Rate


Flow rate depends on sample protein concentration, starting volume, relative centrifugal force (RCF), rotor type, and temperature. High flow rates result from maximal membrane surface area, polarization control, and optimal transmembrane pressure achieved with sample volumes of greater than 300  $\mu\text{L}$ . Optimal filtration rates are achieved at 1,000–2,000  $\times g$ . Higher centrifugal force does not significantly increase flow rate and is not recommended.

## Comparison of Fixed-angle vs. Swinging-bucket Centrifuge Rotors

Rotor	RCF ( $\times g$ )	Typical Ultrafiltrate Volume ( $\mu\text{L}$ )
Fixed-angle (33°)	1,000	140–150
Swinging-bucket	1,000	115–120


1 mL normal human serum, 10 minute spin

## Typical Ultrafiltration Rates

Conditions: 33° fixed-angle rotor, 25 °C, 1,000  $\times g$

## Effect of Relative Centrifugal Force

Conditions: 33° fixed-angle rotor, 25 °C, 10 minute spin

## Membrane Performance

The high selective permeability of the anisotropic, hydrophilic Ultracel® ultrafiltration membrane is due to narrow pore size distribution. Typically, Centrifree® devices retain 99.9% of serum protein and <5% L-thyroxine (0.1 mg/mL in 0.1 N NaOH).

Protein retention requirements for each application depend on the degree of ligand binding. For a ligand that is 90% bound, 1% serum protein leakage would result in a 9% overestimate of free ligand concentration. For measurement of free thyroxine (>99.95% bound) in undiluted serum, protein retention of >99.995% is required for <10% error. Centrifree® devices are suggested for separation of ligands that are up to 99% bound.

Occasionally, serum leakage can occur due to a defect in the membrane or a mechanical defect. If the leakage is greater than 20%, the ultrafiltrate will be yellow colored. Determination of leakage less than 20% requires measurement of the protein concentration in the ultrafiltrate with a sensitive microprotein assay. If you need greater reliability but want to avoid routine protein assay, filter the sample in duplicate devices.

## pH Control

Establish the acceptable limit of pH variation for each ligand-binding system. Experimentally, sample pH can be controlled by anaerobic specimen handling and transfer, gassing the sample with 5% CO<sub>2</sub>, 95% O<sub>2</sub>, or addition of a small volume of concentrated buffer, provided that the buffer ion does not alter ligand binding equilibria.

## Nonspecific Adsorption

Adsorptive losses depend on the concentration of ligand, its ionic and hydrophobic nature, the temperature and time of contact with component surfaces, and sample matrix. Recovery of ligand that has been added to protein-free serum ultrafiltrate is generally greater than from buffer and more predictive of behavior with whole serum. Low-molecular-weight, free fatty and amino acids in protein-free ultrafiltrate interact with and passivate component surfaces, while serum proteins generally passivate surfaces encountered by whole serum. Low recovery from either protein-free ultrafiltrate or buffer may indicate adsorptive losses and/or membrane retention of ligand. Ligand recovery in Centrifree® devices is most accurate in the presence of binding protein and competing microsolutives.

## Specifications

<b>Maximum sample volume</b>	1.0 mL
<b>Minimum sample volume</b>	0.15 mL
<b>Maximum relative centrifugal force</b>	2,000 × g
<b>Active membrane area</b>	0.92 cm <sup>2</sup>
<b>Hold-up volume (membrane and support)</b>	10 µL
<b>Dimensions</b>	
Length of capped device with filtrate cup	95 mm
Diameter	16 mm
<b>Materials of Construction</b>	
Membrane	Ultracel® regenerated cellulose
Sample reservoir	Polycarbonate
Membrane support base	Polycarbonate
Filtrate cup	Polyethylene
Filtrate cup cap	Polyethylene
O-ring	Ethylene propylene diene monomer (EPDM)

**CAUTION:** Centrifree® components are not autoclavable. Do not use with organic solvents.

## Chemical Compatibility

Centrifree® devices are intended for use with biological fluids and aqueous solutions. Before use, check the sample for chemical compatibility with the device. Go to [SigmaAldrich.com/FilterChemicalCompatibility](https://www.sigmaaldrich.com/FilterChemicalCompatibility) for more information.

## Product Ordering

Purchase products online at [SigmaAldrich.com](https://SigmaAldrich.com).

## Symbol Definitions

Symbol	Definition	Symbol	Definition
	In vitro diagnostic medical device		Do not use if package is damaged and consult instructions for use
	Catalogue number		Date of manufacture
	Batch code		Manufacturer
	Consult electronic instructions for use		Importer
	Download product documentation online		CE conformity marking
	Non-sterile		UK Conformity Assessed
	Do not reuse		Swiss authorized representative
	Use-by date		Unique device identifier
	Temperature limit		

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

## Collection and Disposal

All samples must be clearly labelled. Suitable instruments must be used for obtaining and preparing samples.

**NOTE:** Follow precautions for disposal of items contaminated with potentially infectious or hazardous biological material according to all applicable international, federal, state, and local regulations.

## Technical Assistance

Visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

Any serious incident of this device should be reported to manufacturer and competent authority of country where user is established.

## Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

## Revision History

2021-OCT	<ul style="list-style-type: none"><li>• IFU PR05782 Date of Issue OCT 2021 - Replaced PR05180.</li><li>• IFU and Packaging Damage symbols added.</li><li>• Chemical Compatibility and Ordering Information linked to website.</li><li>• Chemical Compatibility, Disposal and Complaints information added.</li><li>• UK Responsible Person and UKCA symbol information added</li></ul>
2024-OCT	<ul style="list-style-type: none"><li>• Added to Symbol Definitions table: CH-REP, Importer, UDI</li><li>• Added UKCA symbol on title page</li><li>• Corrected product name in Chemical Compatibility section</li></ul>

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an affiliate of Merck KGaA,  
Darmstadt, Germany

UK Responsible Person:

Millipore (UK) Ltd.  
Kirkton Campus  
Fleming Road,  
Livingston, UK  
EH54 7BN

Merck, Millipore, Centrifree, Ultracel, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.  
© 1998-2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck operates  
as MilliporeSigma in the U.S. and Canada.

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Merck Millipore Ltd.