3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

# **Product Information**

Luciferase from *Photinus pyralis* (firefly) recombinant, expressed in *Escherichia coli* 

Catalog Number **SRE0045** Storage Temperature –20 °C

CAS RN 61970-00-1 EC 1.13.12.7

Synonyms: Luciferin 4-monooxygenase, Firefly Luciferase

#### **Product Description**

Firefly luciferase is a 62 kDa protein that catalyzes the production of light. The enzyme requires ATP, molecular oxygen, and the heterocyclic compound, firefly luciferin, to generate light in a two-step process. The light producing reaction is initiated by luciferin activation (adenylation of its carboxylate group) and proceeds in the presence of molecular oxygen to yield a photon of yellow-green light. 1.2

Firefly luciferase is used extensively in molecular and cell biology, in particular for the efficient detection and quantitation of ATP and as a reporter for genetic function.<sup>3,4</sup>

This product is a recombinant luciferase from *Photinus pyralis* (American firefly) produced from the *luc* gene expressed in *E. coli.* It is lyophilized from a buffered solution containing HEPES, pH 7.5, NaCl, MgCl<sub>2</sub>, EDTA, DTT, and a carbohydrate stabilizer.

Specific Activity: ≥10 × 10<sup>10</sup> light units/mg protein

Unit definition: One luciferase enzyme unit will produce one Relative Light Unit (RLU) at 20–25 °C over a 10 second period, measured in 100  $\mu$ l assay mixture containing 40 pmole ATP and 15 nmole luciferin in Tris-glycine buffer, pH 7.6, using a GloMax® 20/20 Luminometer.

### **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

#### **Preparation Instructions**

To obtain maximal solubility it is important to reconstitute the enzyme at a high salt concentration, such as 1 M Tris buffer with any counter ion at pH 7–8. The enzyme can be prepared at a concentration of up to 5 mg protein/ml. Do not vortex and avoid agitation.

## Storage/Stability

Store the product at -20 °C.

After reconstitution, the enzyme solutions can kept at 4–8 °C for up to 2 days or frozen in working aliquots at –20 °C for at least one month. Repeated freezing and thawing is not recommended. Do not store in a frost-free freezer.

#### References

- 1. DeWet, J.R. et al., Firefly luciferase gene: Structure and expression in mammalian cells. Mol. Cell. Biol., **7**, 725-737 (1987).
- Stanley, P.E., A review of bioluminescent ATP techniques in rapid microbiology. J. Biolumin. Chemilumin., 4, 375-380 (1989).
- 3. Kricka, L.J., Clinical and biochemical applications of luciferases and luciferins. Analyt. Biochem., **175**, 14-21 (1988).
- 4. Chappelle, É.W. et al., Determination of bacterial content in fluids. Meth. Enzymol., **57**, 65-72 (1978).

GloMax is a registered trademark of Promega Corporation.

DS,DT,MAM 11/16-1