

Data Sheet

THP-1 NF- κ B-eGFP Human Acute Monocytic Leukemia Cell Line

Cancer Cell Line

SCC223**Pack Size: $\geq 1 \times 10^6$ viable cells/vial****Store in liquid nitrogen.**FOR RESEARCH USE ONLY**Not for use in diagnostic procedures. Not for human or animal consumption.**

Background

Microbial contamination of biological cultures is a universal concern for laboratories. Detection and monitoring of microbial contamination often require a significant commitment of time, labor and expense, especially when maintaining multiple cultures and preserving precious specimens. Methods that are at once efficient, cost-effective, and sensitive for a range of microbial contamination are of great value in facilitating detection and ensuring contamination-free cultures.

The THP-1 NF- κ B-eGFP cell line makes use of the toll-like receptor (TLR) signaling pathway as a simple and sensitive biological assay for detection of microbial contamination. The suite of TLR proteins recognize distinct molecular signatures conserved across bacterial species and activate the transcription factor NF- κ B. THP-1 NF- κ B-eGFP cells exhibit high sensitivity for mycoplasma and gram-negative bacterial LPS.¹ Incubation of THP-1 NF- κ B-eGFP cells for 24 hours with contaminated culture supernatants activates expression of the eGFP reporter which may be easily detected via flow cytometry or light microscopy.¹ Presence of amphotericin B in cultures also induces the reporter gene. This unique cell line is a simple method for detecting a range of common microbial contaminants and is suitable for both routine assays and high throughput screening of biological cultures.

Source

The THP-1 NF- κ B-eGFP cell line is derived from a single-cell clone of THP-1 cells stably transfected with an NF- κ B-eGFP reporter construct.¹ The parental THP-1 cell line was derived from peripheral blood of a 1-year-old male suffering from acute monocytic leukemia.²

Short Tandem Repeat

D3S1358: 15, 17	D16S539: 11, 12
TH01: 8, 9.3	CSF1PO: 11, 13
D21S11: 30, 31.2	Penta D:
D18S51: 13, 14	vWA: 16
Penta E: 11, 15	D8S1179: 10, 14
D5S818: 11, 12	TPOX: 8, 11
D13S317: 13	FGA: 24, 25
D7S820: 10	Amelogenin: X, Y

Cancer cell lines are inherently genetically unstable. Genetic instability may arise in the form of loss of heterozygosity of alleles at one or more genetic sites with increased passages.

Quality Control Testing

- Each vial contains $\geq 1 \times 10^6$ viable cells.
- Cells are tested negative for infectious diseases by a Human Essential CLEAR panel by Charles River Animal Diagnostic Services.
- Cells are verified to be of human origin and negative for inter-species contamination from rat, mouse, Chinese hamster, Golden Syrian hamster, and non-human primate (NHP) as assessed by a Contamination CLEAR panel by Charles River Animal Diagnostic Services.
- Cells are negative for mycoplasma contamination.
- Each lot of cells is genotyped by STR analysis to verify the unique identity of the cell line.

Storage and Handling

Store in liquid nitrogen. The cells can be cultured for at least 10 passages after initial thawing without significantly affecting the cell marker expression and functionality.

Presentation




Figure 1. (A) Days one and (B) five after thaw. (C) No LPS control. (D) Activation of GFP fluorescence by exogenously applied lipopolysaccharide (LPS, 3 $\mu\text{g}/\text{mL}$). (E) Flow analysis.

Protocols

THP-1 NF- κ B-eGFP cells grow as suspension cells and thus do not require enzymatic detachment or dissociation. Passage when the cell density reaches 1-1.5 million cells/mL. Optimal plating density should be $\sim 400,000$ cells/mL. The cells should not be grown at excessively high densities. Do not add amphotericin B to the culture medium as this will induce the reporter gene.

1. Do not thaw the cells until the recommended medium is on hand.
2. THP-1 NF- κ B-eGFP Expansion Medium: cells are thawed and expanded in RPMI-1640 (R0883) and 10% heat-inactivated FBS (ES-009-B).
3. Remove the vial of frozen THP-1 NF- κ B-eGFP cells from liquid nitrogen and incubate in a 37 $^{\circ}\text{C}$ water bath. Closely monitor until the cells are completely thawed. Maximum cell viability is dependent on the rapid and complete thawing of frozen cells.
Important: Do not vortex the cells.
4. As soon as the cells are completely thawed, disinfect the outside of the vial with 70% ethanol. Proceed immediately to the next step.
5. In a laminar flow hood, use a 1-2 mL pipette to transfer the cells to a sterile 15 mL conical tube. Be careful not to introduce any bubbles during the transfer process.
6. Using a 10 mL pipette, slowly add dropwise 9 mL of THP-1 NF- κ B-eGFP Expansion Medium (Step 1 above) to the 15 mL conical tube.
Important: Do not add the entire volume of media all at once to the cells. This may result in decreased cell viability due to osmotic shock.
7. Gently mix the cell suspension by slowly pipetting up and down twice. Be careful not to introduce any bubbles.
Important: Do not vortex the cells.
8. Centrifuge the tube at 300 $\times g$ for 2-3 minutes to pellet the cells.
9. Decant as much of the supernatant as possible. Steps 5-8 are necessary to remove residual cryopreservative (DMSO).
10. Resuspend the cells in 5 mL of THP-1 NF- κ B-eGFP Expansion Medium.
11. Transfer the cell suspension to a T25 flask.

- 
12. Incubate the cells at 37 °C in a humidified incubator with 5% CO₂. Place the flask upright for 3 days. Thawed cells may be slow to proliferate the first 3 days. From days 5-6, cell density will double. Passage cells when the cell density is at 1-1.5 million cells/mL.
 13. Replenish the media every 2-3 days. Optimal plating density should be ~400,000 cells/mL

Cryopreservation of Cells

THP-1 NF- κ B-eGFP Human Acute Monocytic Leukemia Cell Line may be frozen in the expansion medium plus 10% DMSO using a Nalgene® slow freeze Mr. Frosty® container.

References

1. PLOS One. 2017; 12(5): e0178220
2. Int J Cancer. 1980; 26(2): 171-176.

Academic Use Agreement

Subject to local law

THIS PRODUCT MAY ONLY BE USED BY INDIVIDUALS EMPLOYED BY AN ACADEMIC INSTITUTION AND IS INTENDED SOLELY TO BE USED FOR ACADEMIC RESEARCH, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU ("PURCHASER") HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER INSTITUTION, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS ACADEMIC USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.

"Product" means THP-1 NF-κB-eGFP Human Acute Monocytic Leukemia Cell Line (SCC223).

"Academic Research" means any internal in vitro research use by individuals employed by an academic institution. Academic Research specifically excludes the following uses of whatever kind or nature:

- Re-engineering or copying the Product
- Making derivatives, modifications, or functional equivalents of the Product
- Obtaining patents or other intellectual property rights claiming use of the Product
- Using the Product in the development, testing, or manufacture of a Commercial Product
- Using the Product as a component of a Commercial Product
- Reselling or licensing the Product
- Using the Product in clinical or therapeutic applications including producing materials for clinical trials
- Administering the Product to humans
- Using the Product in collaboration with a commercial or non-academic entity

"Commercial Product" means any product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to those officers, employees, and students of PURCHASER's academic institution who need access to the Product to perform Academic Research. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These use restrictions will remain in effect for as long as PURCHASER possesses the Product.

COMMERCIAL OR NON-ACADEMIC ENTITIES INTERESTED IN PURCHASING OR USING THE PRODUCT MUST CONTACT licensing@emdmillipore.com AND AGREE TO SEPARATE TERMS OF USE PRIOR TO USE OR PURCHASE.

Genetically Modified Organisms (GMO)

This product contains genetically modified organisms.
Este producto contiene organismos genéticamente modificados.
Questo prodotto contiene degli organismi geneticamente modificati.
Dieses Produkt enthält genetisch modifizierte Organismen.
Ce produit contient des organismes génétiquement modifiés.
Dit product bevat genetisch gewijzigde organismen.
Tämä tuote sisältää geneettisesti muutettuja organismeja.
Denna produkt innehåller genetiskt ändrade organismer.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.
All other trademarks are the property of their respective owners. Detailed information on
trademarks is available via publicly accessible resources.

© 2020-2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20306518 Ver 6.0

20508023 Ver 2.0, Rev 14AUG2024, LB

The Merck logo is displayed in a bold, red, sans-serif font.