User Guide

Millex® (33 mm) Sterile Filter Unit with Durapore® Membrane

SLVVM33RS (0.10 µm), SLGVM33RS (0.22 µm), SLHVM33RS (0.45 µm)

CE₀₄₅₉ Rx only

Available languages:

English	2	Svenska (Swedish). . . .	33
Français (French).	5	Latviski (Latvian).	36
Italiano (Italian)	8	Lietuviškai (Lithuanian). 39	
Deutsch (German).	11	Magyar (Hungarian). . . .	42
Español (Spanish)	14	Čeština (Czech)	45
Português (Portuguese) 17		Polski (Polish)	48
Ελληνικά (Greek).	20	日本語 (Japanese).	51
Nederlands (Dutch)	24	中文 (Chinese)	54
Dansk (Danish)	27	한국어 (Korean)	57
Norsk (Norwegian).	30		

Definition of Symbols Used

Symbol	Definition	Symbol	Definition
Rx only	Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed healthcare professional.		CE conformity marking
	Catalogue number		Date of manufacture
	Do not reuse		Manufacturer
	Use-by date		Do not use if package is damaged
	Batch code		Sterilized using Gamma irradiation

Millex® 33 mm Sterile Filter with Durapore® Membrane

- Single use only
- Sterile
- Non-pyrogenic
- Contains no natural latex rubber

Indications for Use/Purpose

The Millex®-HV/GV/VV filters are medical devices intended for use as a syringe filter to sterilize (GV/VV only) and/or clarify low volume solutions in direct patient care and pharmacy admixture applications.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician or pharmacist.

Introduction

This document provides compatibility information, operating steps, and specifications for the Durapore® hydrophilic polyvinylidene fluoride (PVDF) family of sterile Millex® filter units. The Millex® filter unit's bidirectional support of the filter membrane enables users to filter aqueous solutions in either direction; forward (from the syringe into the container) or backward (from the container into the syringe). The Millex® filter unit removes microorganisms, particles, precipitates, and undissolved powders larger than the membrane's rated pore size. These single-use filter units consist of a membrane filter sealed in an acrylic housing. They are non-pyrogenic and non-toxic.

Applications

Millex® syringe filters are suitable for use in hospital pharmacy and direct patient care applications. In the hospital pharmacy, Millex® filter units with Durapore® membrane can be used for the sterile filtration (VV/GV) and/or clarification (VV/GV/HV) of small volumes of protein pharmaceuticals, diagnostic imaging agents, chemotherapeutics, aqueous solutions, or water during admixture preparation. Direct patient care applications include sterilization (VV/GV) and/or particulate removal (VV/GV/HV) from epidural and other liquid anesthetics as well as from irrigation solutions used in ophthalmic, otic, and other surgical procedures.

Chemical Compatibility

The Millex® filter unit with Durapore® membrane is compatible with most aqueous solutions. Based on information from technical publications, materials suppliers, and laboratory tests, Merck Millipore Ltd. believes that the agents listed in the following chart are safe to use with Millex® filter units. However, because of the effects of variability in temperature, concentrations, duration of exposure, and other factors outside of our control, Merck Millipore Ltd. does not provide or imply a warranty with respect to this information. For example, the drug list refers to the compatibility between the drug and filter materials. It does not refer to the specific drug or protein binding to the filter, or the potential loss of the active drug component during filtration. You should qualify specific drugs for protein binding or loss of drug component prior to use. Agents that are not listed should be tested prior to use.

Chemicals

Acetic acid (5%)	Sodium hydroxide (10%)
Aqueous buffers	Sulfuric acid (20%)
Cell culture media	Water
Clorox® bleach (5% solution)	

Active Drug Compounds

Drug	Comments
Aminophylline	water soluble
Ampicillin	water soluble
Aspartame	water soluble
Bleomycins	water soluble
Caffeine	water soluble
Cefazolin	water soluble
Cefoxitin	water soluble
Cephalothin	water soluble
Cisplatin	water soluble
Colistin	suspension + surfactant
Cytarabine	water soluble
Dactinomycin	water soluble
Daunorubicin	water soluble
Dexamethasone	5% alcohol
Diazepam	40% alcohol
Digoxin	50% alcohol
Dobutamine	water soluble
Dopamine	water soluble
Doxorubicin	water soluble
Ergonovine	water soluble
Etoposide	30% alcohol
Factor III	water soluble
Factor IX	water soluble
Fentanyl	water soluble
Fluorouracil	water soluble
Folic Acid	water soluble
Furosemide	water soluble
Gentamicin	water soluble
Hemin	water soluble
Heparin	water soluble
Hydrocortisone 21-glycol sodium succinate	water soluble
Immunoglobulins	water soluble
Insulin	water soluble
Isoproterenol	water soluble
Lidocaine	water soluble
Mannitol	water soluble
Metronidazole	water soluble
Mitoguazone	water soluble
Mitomycins	water soluble
Mitoxantrone	water soluble
Moxalactam	water soluble
Nitroglycerin	water soluble
Norepinephrine	water soluble
Penicillin G potassium	water soluble
Phenobarbital	water soluble
Piperacillin	water soluble
Plicamycin	water soluble
Prochlorperazine	water soluble

Drug	Comments
Protamines	water soluble
Streptokinase	water soluble
Sulfamethomidine	50% alcohol
Tobramycin	water soluble
Trimethoprim	water soluble
Urokinase	water soluble
Vidarabine	water soluble
Vinblastine	water soluble
Vincristine	water soluble

How to Use the Millex® Sterile Filter

This section lists warnings and cautions and provides steps to use the Millex® sterile filter.

WARNINGS:


- To ensure sterility, do not use this product if the package is damaged.
- Do not use this product as an in-line filter for intravenous fluid administration; it was not designed for long-term continuous use.
- Do not use with syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter unit and/or personal injury.

CAUTIONS:

- Do not use the Millex® filter unit to filter fluids at temperatures above 45 °C (113 °F).
- Do not use the Millex® filter unit to filter emulsions or suspensions because it was not designed for that purpose.
- Do not use the Millex® filter unit to filter solutions containing 5 milligrams (mg) or less of active drug materials unless binding studies have been performed.
- Do not use the same Millex® filter unit to filter solutions in both directions.
- Do not re-sterilize or reuse the Millex® filter unit, as Merck Millipore Ltd. cannot assure the sterility, integrity, and performance beyond a single use.

Procedure for Using the Millex® Sterile Filter

1. Fill syringe with solution to be filtered.


2. Aseptically remove cover from package.

3. Attach syringe to filter and remove assembly from package. Attach needle to Luer-slip outlet if necessary.

4. Hold syringe with filter (and needle if attached) pointing up and top off by pushing a few drops through. Do not contaminate underside of filter with fingers.

5. Insert needle (if attached) and push plunger to deliver filtered solution.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices).

Technical Assistance

Visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://www.sigmaaldrich.com/terms).

Specifications

Materials

Membrane	Low binding hydrophilic polyvinylidene fluoride
Pore size	Millex®-VV filter: 0.10 µm Millex®-GV filter: 0.22 µm Millex®-HV filter: 0.45 µm
Housing	Acrylic

Dimensions

Inlet to outlet	27 mm (1.06 in.)
Diameter	33 mm (1.30 in.)
Filtration area	4.52 cm ² (0.70 in ²)

Temperature limit	45 °C (113 °F) maximum
--------------------------	------------------------

Housing Pressure at 25 °C	10.3 bar (150 psi) inlet maximum
----------------------------------	----------------------------------

Filtration volume	10 mL to 100 mL
--------------------------	-----------------

Hold-up volume	≤ 0.1 mL after air purge
-----------------------	--------------------------

Sterilization method	Gamma irradiation
-----------------------------	-------------------

Connections	Female Luer-Lok™ inlet; male Luer-slip outlet
--------------------	--

Flow rate at 2.1 bar (30 psi), 25 °C	Millex®-VV filter: ≥ 25 mL/min Millex®-GV filter: ≥ 40 mL/min Millex®-HV filter: ≥ 300 mL/min
---	---

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Made in Ireland
Merck Millipore Ltd.
Tullagreen, Carrigtwohill,
Co. Cork, IRL
an affiliate of Merck, KGaA, Darmstadt, Germany

Millipore, Millex, Durapore and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners.
Detailed information on trademarks is available via publicly accessible resources.
© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.