

Product Information

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103 USA

Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757

email: techservice@sial.com sigma-aldrich.com

Glutathione S-Transferase from human placenta

Catalog Number **G8642**

Storage Temperature: -20 °C

EC 2.5.1.18

CAS RN 50812-37-8

Synonyms: GST; Glutathione S-alkenetransferase;
Glutathione S-alkyltransferase;
Glutathione S-aralkyltransferase;
Glutathione S-aryltransferase;
Glutathione S-epoxidetransferase;
RX: Glutathione R-transferase

Product Description

Glutathione S-transferases (GST) are a family of proteins important in the detoxification of many xenobiotics in mammals. They catalyze the formation of thioether conjugates between reduced glutathione and hydrophobic compounds. The major biological function is believed to be a role in the elimination of reactive, electrophilic chemical species, many of which are generated by cellular oxidative reactions catalyzed by cytochrome P450 and other oxidases. This elimination defends cells against the mutagenic, carcinogenic, and toxic effects of the compounds. GST activity is present in plants, insects, yeast, bacteria, and most mammalian tissues especially in the liver, which plays a key role in detoxification.

Enzymatic activities are based on the conjugation of reduced glutathione with a second substrate. The individual proteins generally have activity with more than one class of substrate.

Molecular mass: 47 kDa (dissociates into equal subunits in the presence of SDS)

The product is supplied as lyophilized powder containing Tris buffer salts, reduced glutathione, and EDTA.

Protein: ≥25% (biuret)

Specific Activity: 25–125 units/mg-protein

Unit Definition: One unit will conjugate 1.0 μmole of 1-chloro-2,4-dinitrobenzene with reduced glutathione per minute at pH 6.5 at 25 °C.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Immediately before use, prepare a solution containing 0.075–0.15 unit/ml of Glutathione S-Transferase in 100 mM potassium phosphate buffer, pH 6.5, with 1.0 mM ethylenediaminetetraacetic acid (EDTA). Long term storage of this dilute enzyme solution is not recommended.

Storage/Stability

Store the product at -20 °C. When stored at -20 °C, the enzyme retains activity for at least 2 years.

References

1. Habig, W.H., et al., Glutathione S-transferase. The first enzymatic step in mercapturic acid formation. *J. Biol. Chem.*, **249**, 7130-7139 (1974).
2. Asaoka, K., and Takahashi, K.J., *Biochem.*, **82**, 1313 (1977).
3. Lee, C.-Y.G., *Biochem., Soc. Trans.*, **12**, 30 (1984).
4. Asoka, K.J., *Biochem.*, **95**, 685, (1984)
5. Simons, P.C., and Vander Jagt, D.L., *Anal. Biochem.*, **82**, 334,(1977).
6. Warholm, M., et al., *Biochemistry*, **25**, 4119 (1986).

KAA,MAM 01/08-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

SIGMA-ALDRICH™