


Product Information

95348 Abberior® STAR 470SX, NHS ester



Absorption & Fluorescence Spectrum



Key Features

- Designed for STED microscopy@ ~750 nm
- Long Stokes' shift (>130 nm) for 2-color applications
- Tested in the Leica TCS STED Ti:Sa 2-color system

Description


Abberior STAR 470SX is the latest development of long-Stokes-Shift dyes for STED microscopy. The dye can be excited from 450 to 480nm. It can substitute dyes like Chromeo® 494. For STED, a depletion wavelength ~750 nm is recommended. It is therefore well suited for **2-color STED imaging** as implemented in the Leica TCS STED Ti:Sa microscope.

Chemical Data : Abberior® STAR 470SX

Molecular weight:	595.57 g/mol (NHS ester)
Solubility:	water, acetonitrile, MeOH, DMSO, DMF
Polarity:	medium hydrophilic
Charge:	0 (when conjugated)
Purity:	> 90 %

Photophysical Data : Abberior® STAR 470SX

Absorption Maximum, λ_{max} :	475 nm (MeOH), 477 nm (PBS, pH 7.4)
Extinction Coefficient, $\epsilon(\lambda_{\text{max}})$:	30,400 M ⁻¹ cm ⁻¹ (MeOH), 22,700 M ⁻¹ cm ⁻¹ (PBS, pH 7.4)
Correction Factor, CF ₂₈₀ = $\epsilon_{280}/\epsilon_{\text{max}}$:	0.69
Correction Factor, CF ₂₈₀ = $\epsilon_{280}/\epsilon_{\text{max}}$:	0.47
Fluorescence Maximum, λ_{fl} :	609 nm (MeOH), 627 nm (PBS, pH 7.4)
Recommended STED Wavelength, λ_{STED} :	740 – 770 nm
Fluorescence Quantum Yield, η :	0.80 (EtOH),
Fluorescence Lifetime, τ :	3.9 (EtOH)



Applications

Abberior STAR 470SX is the dye of choice for long Stokes STED applications in the orange fluorescent regime. The dye is particularly designed and tested for 2-color STED microscopy in combination with our STAR 635 using a single STED wavelength. The dye is our recommendation for usage in the Leica TCS STED Ti:Sa 2-color system. For more information see our [2-color dye package](#) section

Trademark info: ®Abberior is a registered Trademark of Abberior GmbH, Germany

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.