

Data Sheet

SCp2 Mouse Mammary Epithelial Cell Line

Cancer Cell Line

SCC458**Pack Size: $\geq 1 \times 10^6$ viable cells/vial****Store in Liquid Nitrogen****FOR RESEARCH USE ONLY****Not for use in diagnostic procedures. Not for human or animal consumption.**

Background

Mammary epithelial cell differentiation and reorganization occurs throughout development and in the adult organ during pregnancy and menopause. During these periods in the adult life there is an increased risk for breast cancer, which is the most common cancer in women. Breast cancer can be treated when detected at early stages but is considered incurable in advance stages with metastasis. The extracellular matrix (ECM) in mammary tissue influences growth and differentiation during development, differentiation, and carcinogenesis.¹ Epithelial cells sense ECM modifications through transmembrane receptors like integrins that can influence the expression of malignant phenotypes exhibiting increased proliferation and loss of contact inhibition growth.¹

The SCp2 cell line was isolated from the heterogenous CID9 cells by limiting dilution cloning.² CID9 cells were isolated from a mid-pregnancy mouse mammary gland. SCp2 cells require exogenously added ECM and lactogenic hormones to express b-casein, under these conditions the cells arrest growth, aggregate, and form alveolar structures that express high levels of milk proteins.² SCp2 cells are useful for the study of normal and abnormal mammary epithelial cell phenotypes because their growth and differentiation can be controlled in culture. Activated mutants of M-Ras have been introduced in normal SCp2 cells to study transformation, invasion, and tumor formation *in vivo*.³

Source

SCp2 is a murine mammary epithelial cell line that originated from limiting dilution cloning of heterogenous CID9 cells.

Short Tandem Repeat

M18-3: 17, 18	M1-1: 13, 14	M11-2: 18	D8S1106: ND
M4-2: 21.3	M3-2: 14	M17-2: 16, 17	D4S2408: ND
M6-7: 12	M8-1: 13	M12-1: 16	
M19-2: 13	M2-1: 16	M5-5: 14	
M1-2: 17	M15-3: 22.3	MX-1: 24	
M7-1: 26.2	M6-4: 19	M13-1: 16.2	

Spectral Properties

Cancer cell lines are inherently genetically unstable. Instability may arise in the form of loss of heterozygosity of alleles at one or more genetic sites with increased passages.

Quality Control Testing

- Each vial contains $\geq 1 \times 10^6$ viable cells.
- Cells are tested negative for infectious diseases by a Mouse Essential CLEAR panel by Charles River Animal Diagnostic Services.
- Cells are verified to be of mouse origin and negative for interspecies contamination from human, rat, Chinese hamster, Golden Syrian hamster, and Non-human Primate (NHP) as assessed by a Contamination Clear panel by Charles River Animal Diagnostic Services.
- Cells are negative for mycoplasma contamination.

Storage and Handling

SCp2 cells should be stored in liquid nitrogen. The cells can be cultured for at least 10 passages after initial thawing without significantly affecting functionality.

Representative Data




Figure 1. Brightfield image of SCp2 cells one day after thaw in a T75 flask (A). Cell confluency (30%) was assessed throughout the culture using the Millicell® Digital Cell Imager (B, Cat. No. MDCI 10000). Cell counting was performed using Scepter™ 3.0 handheld automated cell counter using 60 mm sensors (C, PHCC360KIT). SCp2 cells express actin (D, Cat. No. 49409), estrogen receptor alpha (ERα) (E, Cat. No. 06-935) and b-catenin (F, Cat. No. ABE208). Majority of SCp2 cells (98.3%) express the epithelial marker, pan-cytokeratin (G, Cat. No. SAB4700667) and are negative for fibroblastic marker, vimentin (H, R&D Systems® IC2105G). As a control, human foreskin fibroblasts are positive for vimentin (I).

Protocols

Thawing Cells

1. Do not thaw the cells until the recommended medium is on hand. Cells can grow on normal tissue cultureware surfaces without any additional coating. Cells are thawed and expanded in SCp2 Expansion Medium comprised of DMEM/F12 medium (Cat. No. D8062) with 5 µg/mL insulin (Cat. No. I9278) and 2% FBS (Cat. No. ES-009-B).
2. Remove the vial of frozen SCp2 cells from liquid nitrogen and incubate in a 37 °C water bath. Closely monitor until the cells are completely thawed. Maximum cell viability is dependent on the rapid and complete thawing of frozen cells.
IMPORTANT: Do not vortex the cells.
3. As soon as the cells are completely thawed, disinfect the outside of the vial with 70% ethanol. Proceed immediately to the next step.
4. In a laminar flow hood, use a 1- or 2- mL pipette to transfer the cells to a sterile 15 mL conical tube. Be careful not to introduce any bubbles during the transfer process.
5. Using a 10 mL pipette, slowly add dropwise 9 mL of SCp2 Expansion Medium (Step 1 above) to the 15 mL conical tube.
IMPORTANT: Do not add the entire volume of media all at once to the cells. This may result in decreased cell viability due to osmotic shock.
6. Gently mix the cell suspension by slowly pipetting up and down twice. Be careful not to introduce any bubbles.
IMPORTANT: Do not vortex the cells.
7. Centrifuge the tube at 300 x *g* for 2-3 minutes to pellet the cells.
8. Decant as much of the supernatant as possible. Steps 5-8 are necessary to remove residual cryopreservative (DMSO).
9. Resuspend the cells in 15 mL of SCp2 Expansion Medium.
10. Transfer the cell mixture to a T75 tissue culture flask.
11. Incubate the cells at 37 °C in a humidified incubator with 5% CO₂.

Subculturing Cells

Differential detachment is required ~ every other time that the cell line is passaged to remove fibroblast-like cells from the population.

1. SCp2 cells should be passaged at ~ 80-85% confluency. Do not allow the cells to grow over 85% confluency.
2. Carefully remove the medium from the T75 tissue culture flask containing the 80% confluent layer of SCp2 cells.
3. Rinse the flask with 10 mL 1X PBS. Aspirate after the rinse.
4. Differential detachment of cells:
 - Apply 5-7 mL of Accutase® and incubate in a 37 °C incubator for ~2 minutes. Bang hard on the flask to remove cells that detach easily.
 - Aspirate the Accutase® and any dislodged fibroblastic cells.
 - Add fresh Accutase® to the flask and incubate for an additional 2-3 minutes. Cells that detach during these 2-3 minutes should be predominantly epithelial cells of the SCp2 cell line.
5. Add 10 mL of SCp2 Expansion Medium to the flask.
6. Gently rotate the flask to mix the cell suspension. Transfer the dissociated cells to a 15 mL conical tube.
7. Centrifuge the tube at 300 x *g* for 3-5 minutes to pellet the cells.
8. Discard the supernatant, then loosen the cell pellet by tapping the tip of the tube with a finger.
9. Apply 2-5 mL of SCp2 Expansion Medium to the conical tube and resuspend the cells thoroughly.
IMPORTANT: Do not vortex the cells.
10. Count the number of cells with a Scepter™ 3.0 handheld automated cell counter using 60 µm sensors.
11. Plate the cells to the desired density. Typical split ratio is 1:6.

Cryopreservation of Cells

SCp2 cells may be frozen in SCp2 Expansion Medium and 10% DMSO using a Nalgene® slow freeze Mr. Frosty™ container.

References

1. Int J Biochem Cell Biol 2007, 39(11): 1987-1994.
2. Proc Natl Acad Sci USA 1994, 91(26): 12378-12382.
3. Oncogene 2004, 23(6): 1187-1196.

Academic Use Agreement

(Subject to local law)

THIS PRODUCT MAY ONLY BE USED BY INDIVIDUALS EMPLOYED BY AN ACADEMIC INSTITUTION AND IS INTENDED SOLELY TO BE USED FOR ACADEMIC RESEARCH, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU ("PURCHASER") HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER INSTITUTION, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS ACADEMIC USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.

"Product" means SCp2 Mouse Mammary Epithelial Cell Line (SCC458).

- "Academic Research" means any internal in vitro research use by individuals employed by an academic institution. Academic Research specifically excludes the following uses of whatever kind or nature:
- Re-engineering or copying the Product
- Making derivatives, modifications, or functional equivalents of the Product
- Obtaining patents or other intellectual property rights claiming use of the Product
- Using the Product in the development, testing, or manufacture of a Commercial Product
- Using the Product as a component of a Commercial Product
- Reselling or licensing the Product
- Using the Product in clinical or therapeutic applications including producing materials for clinical trials
- Administering the Product to humans
- Using the Product in collaboration with a commercial or non-academic entity

"Commercial Product" means any product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to those officers, employees, and students of PURCHASER's academic institution who need access to the Product to perform Academic Research. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These use restrictions will remain in effect for as long as PURCHASER possesses the Product.

COMMERCIAL OR NON-ACADEMIC ENTITIES INTERESTED IN PURCHASING OR USING THE PRODUCT MUST CONTACT licensing@emdmillipore.com AND AGREE TO SEPARATE TERMS OF USE PRIOR TO USE OR PURCHASE.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrichAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany
operates as MilliporeSigma in the U.S. and Canada.

Merck Millicell, Scepter and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20306518 Ver 6.0

20711629 Ver 2.0, Rev 06OCT2023, RC/AB

