Instructions for Use

Amyloid Stain, Congo Red

Procedure No. HT60

Intended Use

Amyloid Stain, Congo Red, kit is for use in the determination of amyloid in tissue sections. Both a standard procedure and a microwave procedure are provided. Amyloid Stain reagents are for "In Vitro Diagnostic Use". For professional use only. The data obtained from this manual, qualitative procedure serves only as an aid to diagnosis of amyloidosis in human specimens and should be reviewed in conjunction with other clinical diagnostic tests or information.

Amyloid, an abnormal protein product of an immunologic disturbance, accumulates between parenchymatous cells of connective tissue resulting in amyloidosis. Amyloid may be detected in tissue with congo red, a metachromatic anionic dye. Initial methods using congo red required examination of stained sections with a polarizing microscope.¹ Puchtler's modification, on which the Sigma-Aldrich procedure is based, intensifies congo red staining by pretreatment of tissue with an alkaline solution of sodium chloride.² The section is then evaluated with a light microscope or polarizing microscope.

An amyloid stain technique for rapid staining in microwave ovens is included.³⁻⁵

Reagents

Sodium Chloride Solution, Alcoholic (Cat. No. HT601-500ML)
SD 3A alcohol, 80%, saturated with sodium chloride

Sodium Hydroxide Solution (Cat. No. HT602-12ML)
Sodium hydroxide, 1% (0.25 mol/L). Warning. Causes skin irritation. Causes serious eye irritation. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Congo Red Solution (Cat. No. HT603-500ML)
Congo red (certified), 0.2%, C.I. 22120, in SD 3A alcohol, 80% saturated with sodium chloride.

Special Materials Required but Not Provided

- Positive control slides, such as Amyloid TISSUE-TROL™ (Cat. No. TTR002) should be included in each run
- Mayer's Hematoxylin Solution or Hematoxylin Solution, Gill No. 3
- Xylene or xylene substitute
- Ethanol, Absolute

For Microwave Procedure Only

- Microwave Oven
- Staining Jars (Cat. No. S5641)

Additional Equipment

- Microscope or polarizing microscope

Storage and Stability

Amyloid Stain reagents should be stored at room temperature (18–26°C). Congo Red Solution and Sodium Chloride Solution, Alcoholic, are stable until expiration date shown on labels.

Molarity of Sodium Hydroxide Solution can be verified by titration.

Presence of precipitate in Sodium Hydroxide Solution, Alcoholic, or Congo Red Solution does not affect use.

Preparation

Sodium Hydroxide Solution is supplied ready for use.

Alkaline Sodium Chloride Solution is prepared by adding 0.4 mL Sodium Hydroxide Solution to 40 mL Sodium Chloride Solution, Alcoholic. Mix just before use. Use once and discard.

Alkaline Congo Red Solution is prepared by adding 0.4 mL Sodium Hydroxide Solution to 40 mL Congo Red Solution. Filter solution and use within 15 minutes.

Precautions

The IVDs included in this kit are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Amyloid TISSUE-TROL™ control slides are paraffin embedded human tissue containing amyloid and should be considered potentially infectious.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Paraffin embedded tissue fixed in Carnoy's fluid, absolute alcohol, Formalin Solution, 10% Neutral Buffered, or Bouin's Solution, may be used.² Cut sections 6–12 microns. Thinner sections may not show small amounts of amyloid deposits and will not exhibit birefringence.

Notes

The use of polarized light is highly recommended for the evaluation of amyloid deposits. Apple-green birefringence is considered specific for amyloid.

Prolonged fixation in formalin may diminish Congo Red staining intensity.⁶

Alcohol fixed tissue may yield more intense staining of amyloid, as compared to formalin fixed tissue.⁷

Procedure

- Deparaffinize and hydrate to deionized water.
- Stain in Mayer's Hematoxylin Solution, for **10 minutes**
OR Hematoxylin Solution, Gill No. 3, for **3 minutes**.
- Rinse in tap water for **5 minutes**.
- Place in Alkaline Sodium Chloride Solution for **20 minutes**.
- Stain in Alkaline Congo Red Solution for **20 minutes**.⁸
- Rinse in 3 changes of absolute ethanol.
- Clear in xylene and mount. Examine microscopically.

Microwave Procedure

- Deparaffinize and hydrate to deionized water.
- Place slides in **40 mL** of Hematoxylin Solution, Gill No. 3, contained in a plastic Coplin jar. Loosely cover jar with lid, or use lids with holes drilled into them.
- Microwave at **800 watts** for **5 seconds**.
- Rinse well in running water for **1–2 minutes**.
- Place slides in **40 mL** Alkaline Sodium Chloride Solution contained in a plastic Coplin jar.
- Microwave at **400 watts** for **20 seconds**.
- Place slides directly into **40 mL** of Alkaline Congo Red Solution contained in a plastic Coplin jar.
- Microwave at **400 watts** for **30 seconds**. Let incubate for **5 minutes**.
- Rinse in three changes of absolute ethanol.
- Clear in xylene and mount. Examine microscopically.

Performance Characteristics

Target Structure	Staining Result
Amyloid	Red to Pink-Red
Nuclei	Blue
Elastic fibers	Lighter Red

When examined under polarizing microscope, amyloid shows apple-green birefringence.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT601	Sodium Chloride Solution Alcoholic	Amyloid	3 of 3	3 of 3	3 of 3	3 of 3
		Elastic	3 of 3	3 of 3	3 of 3	3 of 3
HT602	Sodium Hydroxide Solution	Molarity	3 of 3	3 of 3	3 of 3	3 of 3
HT603	Congo Red Solution	Amyloid	3 of 3	3 of 3	3 of 3	3 of 3
		Elastic	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT60:

H225: Highly flammable liquid and vapor.

H290: May be corrosive to metals.

H315: Causes skin irritation.

H319: Causes serious eye irritation.

H350: May cause cancer.

P202: Do not handle until all safety precautions have been read and understood.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308 + P313: IF exposed or concerned: Get medical advice/ attention.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer

References

1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, pp 402–405
2. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, pp 177–178
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, p 311
7. Humanson’s Animal Tissue Techniques, 5th Edition, Presnell JK and Schreiberman MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, p 270

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev. 4.0	2016
Rev. 5.0	2022
Rev. 6.0	2022
Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Updated literature references. Added a revision history table. Added Warnings and Hazards.	

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271 Darmstadt,
Germany

The Initial M, TISSUE-TROL, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, pp 402-405
2. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, pp 177-178
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in microwave oven. J. Histotechnol 9:237, 1986
6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, p 311
7. Humanson's Animal Tissue Techniques, 5th Edition, Presnell JK and Schreiberman MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, p 270

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Références

1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, pp 402-405
2. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, pp 177-178
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, p 311
7. Humanson’s Animal Tissue Techniques, 5th Edition, Presnell JK and Schreibman MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, p 270

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Riferimenti

- 1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, pp 402-405.
- 2. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, pp 177-178.
- 3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986.
- 4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983.
- 5. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986.
- 6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, p 311.
- 7. Humanson's Animal Tissue Techniques, 5th Edition, Presnell JK and Schreiberman MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, p 270.

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Referencias

1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, pp 402-405
2. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, pp 177-178
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, p 311
7. Humanson’s Animal Tissue Techniques, 5th Edition, Presnell JK and Schreiber MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, p 270

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Referencer

1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, pp 402-405
2. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, pp 177-178
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, p 311
7. Humanson's Animal Tissue Techniques, 5th Edition, Presnell JK and Schreibman MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, p 270

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Referenser

1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, pp 402–405
2. Theory and Practice of Histotechnology. Redigerat av DC Sheehan och BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, s 177–178
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986
6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, p 311
7. Humanson’s Animal Tissue Techniques, 5th Edition, Presnell JK and Schreibman MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, p 270

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Referências

- 1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, pp 402–405
- 2. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, pp 177–178
- 3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- 4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- 5. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
- 6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, p 311
- 7. Humanson’s Animal Tissue Techniques, 5th Edition, Presnell JK and Schreibman MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, p 270

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1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, pp 402–405
2. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, pp 177–178
3. Leong AS-Y, Millios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, p 311
7. Humanson’s Animal Tissue Techniques, 5th Edition, Presnell JK and Schreibman MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, p 270

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Reference

1. Selected Histochemical and Histopathological Methods. Editovali SW Thompson, Charles C. Thomas, Springfield (IL), 1966, str. 402–405
2. Theory and Practice of Histotechnology. Editovali DC Sheehan a BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, str. 177–178
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986
6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, p 311
7. Humanson's Animal Tissue Techniques, 5th Edition, Presnell JK and Schreibman MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, str. 270

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Referanser

1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, s. 402–405
2. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2. utg. Mosby, St. Louis (MO), 1980, s. 177–178
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, s. 311
7. Humanson's Animal Tissue Techniques, 5th Edition, Presnell JK and Schreibleman MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, s. 270

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Referanslar

1. Selected Histochemical and Histopathological Methods. Edited by SW Thompson, Charles C. Thomas, Springfield (IL), 1966, s. 402-405
2. Theory and Practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby, St. Louis (MO), 1980, s. 177-178
3. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
4. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
5. Valle S: Special stains in microwave oven. J Histotechnol 9:237, 1986
6. Theory and Practice of Histological Techniques, 5th Edition, Bancroft JD, Gamble M, Editors Churchill Livingstone, New York (NY), 2002, s. 311
7. Humanson's Animal Tissue Techniques, 5th Edition, Presnell JK and Schreiber MP, Editors The John Hopkins University Press, Baltimore (MD), 1997, s. 270

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MDSS GmbH