

Product Information

Heat-Labile Enterotoxin, B subunit from *Escherichia coli*, recombinant expressed in *Pichia pastoris*

Catalog Number **E8656**

Storage Temperature 2–8 °C

Synonym: LTB

Product Description

Enterotoxigenic *Escherichia coli* causes diarrhea through its heat-labile enterotoxin (LT). The LT is a periplasmic protein composed of one A subunit (LTA, 27 kDa) and five non-covalently associated B subunits (LTB, 11.6 kDa each) forming a ring-like pentamer.^{1,2}

LTB has high affinity towards the toxin receptor monosialoganglioside G_{M1}, a glycosphingolipid found ubiquitously on the surface of mammalian cells. Ganglioside G_{M1} facilitates the delivery of the A subunit to the cytosol of the target cell resulting in persistent synthesis of cAMP and subsequently diarrhea.³ This characteristic makes LTB a good label for microglial cells (due to the enrichment of ganglioside GM₁ on their cell surface). In addition, studies made mostly with animal models demonstrated that recombinant LTB could stimulate strong serum and mucosal immune responses against LT. Other studies have indicated that LTB could be used as a potent mucosal adjuvant.^{1,4}

The product is lyophilized from a solution containing 0.05 M Tris buffer, pH 7.5, 0.2 M NaCl, 3 mM NaN₃, and 1 mM sodium EDTA.

Purity: >90% (SDS-PAGE)

Binding saturation: 50% is achieved with ≤0.12 µg/mL of LTB

The activity is measured by ELISA using ganglioside G_{M1}-coated plates (Catalog Number G7641). LTB at different concentrations is incubated in the ganglioside G_{M1} coated wells, followed by anti-Cholera Toxin antibody (Catalog Number C3062), and peroxidase-labeled goat anti-rabbit IgG (Catalog Number A0545) as the secondary antibody.

Note: The use of anti-Cholera Toxin for LTB detection is possible due to similar antigenic determinants.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Reconstitute with ultrapure water. Swirl the bottle gently during reconstitution to facilitate dissolution. Avoid vigorous pipetting that may lead to foaming. The solution can be filtered through a 0.2 µm filter.

Storage/Stability

Store the lyophilized powder and reconstituted solutions at 2–8 °C.

References

1. Rezaee, M.A., et al., Expression of *Escherichia coli* heat-labile enterotoxin B subunit (LTB) in *Saccharomyces cerevisiae*. *J. Microbiol.*, **43**, 354-60 (2005).
2. Hedges, P.A., and Hardy, S.J., Formation of disulfide bonded dimer of mutated heat-labile enterotoxin *in vivo*. *FEBS Letters*, **381**, 53-56 (1996).
3. de Haan, L., et. al., Role of G_{M1} binding in the mucosal immunogenicity and adjuvant activity of the *Escherichia coli* heat-labile enterotoxin and its B subunit. *Immunology*, **94**, 424-430 (1998).
4. Freytag, L.C., and Clements, J.D., Mucosal adjuvants. *Vaccine*, **23**, 1804-1813 (2005).

KAA,EB,MAM 03/08-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.