


75744 Tributyrin-Strips

Diagnostic test for the differentiation between Branhamella and Neisseria. The test principle is an enzyme hydrolysis of tributyrin. This reaction causes colour change of acidobasic indicator. The result of the reaction is read after 18-20 h.



Composition:

(1 package contains 300 test strips) contains strips saturated with tributyrin and acidobasic indicator

Storage:

Store dry at +2 to +8°C.
Expiration can be prolonged when stored at -20°C.

Attention: The strips are delivered sterilized. Observe the maintenance of sterility when the strips are used repeatedly. Prevention of moisture during storage is necessary for proper function of test. Bottle with strips should not be opened before warm up to room temperature. If it is not observed, the product gets moist and deteriorates.

Directions:

Using a sterile for ceps, throw one tributyrin strip into the suspension of tested strain in 1 ml of buffered saline (pH 7.2).
Incubate the test sample at 37°C (without CO₂).
Preliminary results can be read after several hours when the red colour changes to yellow in the case of positive result. Perform the final evaluation of result after 18-20 hours incubation.

Interpretation of results:

Negative reaction: red colour did not change to yellow (Neisseria)
Positive reaction: red colour change to yellow (Branhamella)

Quality control:

The list below illustrates control strains in routine use:

Test Organisms (ATCC)	Result
<i>Neisseria gonorrhoeae</i> (19424)	negative
<i>Branhamella catarrhalis</i> (25238)	positive

References:

- Berger U., Über die Spaltung von tributyrin durch Neisseria., Arch Hyg. Bakteriol., 146:388-391 (1962)
- Kuzmanská P., Laboratorní průkaz kokobacilli a gramnegativních koků, Avicenum Praha, 52-57, (1987)
- Janda W.M. and P. Ruther. B. CAT CONFIRM; A rapid test for confirmation of Branhamella catarrhalis. J. Clin. Microbiol., 27:1390-1391 (1989)
- August M.J., et al., Cumitech 3A; Quality Control and Quality Assurance Practices in Clinical Microbiology, Coordinating ed., A.S. Weissfeld. American Society for Microbiology, Washington D.C. (1990).



-
5. Perez J.L., et al. Butyrate esterase (tributyrim) spot test, a simple method for immediate identification of *Moraxella* (*Branhamella*) *catarrhalis*. *J. Clin. Microbiol.*, 28: 2347-2348 (1990)
 6. Murray P.R., et al. *Manual of Clinical Microbiology*, 6th ed. American Society for Microbiology, Washington D.C. (1995)
 7. Koneman E.W., et al. *Color Atlas and Textbook of Diagnostic Microbiology*, 5th ed. J.B. Lippincott Company, Philadelphia, PA, (1997)
 8. Forbes B.A., et al. *Bailey and Scott's Diagnostic Microbiology*, 10th ed. C.V. Mosby Company, St. Louis, MO, (1998)

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The vibrant M, Millipore, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.
© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

