

# Selenite Enrichment Broth acc. to LEIFSON

## Selenite-F Broth; Selenite Broth

Medium proposed by LEIFSON (1936) for the selective enrichment of *Salmonella* from faeces, urine, water, foodstuffs etc.



*In Vitro Diagnostic Medical Device –*

*For professional use only*



Version 17-10-2008

Merck KGaA, 64271 Darmstadt

*See also General Instruction for Use  
„How to use Dehydrated Culture Media“*

*For MSDS, warnings and precautions see our website:  
[www.merck-chemicals.com](http://www.merck-chemicals.com)*

### Principle

Microbiological method.

### General Information

The medium complies with the recommendations of the APHA (1992) for food examination.

### Mode of Action

Selenite inhibits the growth of enteric coliform bacteria and enterococci, mainly during the first 6-12 hours of incubation. *Salmonella*, *Proteus* and *Pseudomonas* are not suppressed.

### Typical Composition (g/litre)

Peptone from meat 5.0; lactose 4.0; sodium selenite 4.0; di-potassium hydrogen phosphate 3,5; potassium dihydrogen phosphate 6.5.

### Preparation

Suspend 23 g/litre at room temperature; if the medium does not dissolve readily, heat briefly (max. 60 °C); if the medium is to be stored for a longer period of time filter-sterilize, dispense into suitable containers.

■ **Do not autoclave.**

pH: 7.0 ± 0.2 at 25 °C.

The prepared broth is clear and yellowish.

### Storage

After a longer storage period of the dehydrated medium, the colour of the prepared broth might change to reddish/red. The microbiological performance however is not affected.

Usable up to the expiry date when stored dry and tightly closed below +15 °C. Protect from light.

After first opening of the bottle the content can be used up to the expiry date when stored dry and tightly closed below +15 °C.

■ **Storage of the dehydrated culture medium below 15 °C!**

### Specimen

e.g. Stool, urine .

Clinical specimen collection, handling and processing, see general instructions of use.

### Experimental Procedure and Evaluation

Add solid sample material to the normal-strength broth. Mix liquid samples with double-strength broth in the ratio 1:1.

Incubation: up to 24 hours at 37 °C - according to BÄNFFER (1971) and other authors, 43 °C is better.

After 6-12 hours and, if necessary, after 18-24 hours inoculate material from the resulting culture onto selective culture media.

### Literature

BÄNFFER, J.R.: Comparison of the isolation of *Salmonellae* from human faeces by enrichment at 37 °C and 43 °C. - *Zbl. Bakt. I. Orig.*, 217; 35-40 (1971).

LEIFSON, E.: New selenite enrichment media for the isolation of typhoid and paratyphoid (*Salmonella*) bacilli. - *Am. J. Hyg.*, 24; 423-432 (1936).

American Public Health Association: Compendium of methods for the microbiological examination of foods. - 3<sup>rd</sup> ed., 1992.

### Ordering Information

| Product                                   | Ordering No. | Pack size |
|---|--------------|-----------|
| Selenite Enrichment Broth acc. to LEIFSON | 1.07717.0500 | 500 g     |

### Quality control

| Test strains                             | Inoculum     | Growth after 24 hours |
|--|--------------|-----------------------|
| <i>Escherichia coli</i> ATCC 25922       | approx. 99 % | ≤ 10 %                |
| <i>Salmonella typhimurium</i> ATCC 14028 | approx. 1 %  | ≥ 90 %                |