

Technical Data Sheet

CE MUELLER-HINTON agar CLSI

Ordering number: 1.05435.0500

For the determination of antibiotic susceptibility including sulfonamides by the disc-agar diffusion method.

Mueller Hinton Agar is recommended for the diffusion of antimicrobial agents impregnated on paper disc through an agar gel as described in CLSI Approved Standard. This medium is suitable for use in routine susceptibility testing. It is a non-selective, non-differential medium and it is useful in identifying sulfonimide-resistant and responsive strains of gonococci.

IVD in vitro diagnosticum - For professional use only

Mode of Action

The medium complies with the requirement of the Clinical and Laboratory Standards Institute (CLSI) and is manufactured to contain low concentrations of thymine and thymidine as well as appropriate levels of calcium and magnesium ions.

Thymine and thymidine concentrations are determined by the disc diffusion procedure with trimethoprim and sulfamethoxazole and *Enterococcus faecalis* ATCC 29212.

Calcium and/or magnesium concentrations are controlled by obtaining the correct zone diameters with an inoglycoside antibiotics and *Pseudomonas aeruginosa* ATCC 27853.

Typical Composition

Meat infusion	2.0
Casein hydrolysate	17.5
starch	1.5
Agar-agar	17.0

Preparation

Suspend 38 g in 1 litre of demin. water and autoclave (15 min at 121°C). If required, cool to 45-50°C and add 5-10% defibrinated blood, pour plates.

pH: 7.3 ± 0.2 at 25°C

Without blood, the plates are clear to opalescent and brownish-yellow.

Specimen

e.g. Isolated bacteria from urine.

Clinical specimen collection, handling and processing. See general instructions of use.

Experimental Procedure and Evaluation

Carry out the sensitivity test acc. to CLSI-

Incubation for 24 h at 35°C aerobically.

Storage

Usable up to expiry date when stored dry and tightly closed at +15 to 25°C. Protect from light.

After first opening of the bottle the content can be used up to the expiry date when stored dry and tightly closed at +15 to +25°C.

Quality Control

Control Strains	ATCC #	Incubation	Expected Results
<i>Escherichia coli</i>	25922	16 – 18 h at 35 °C	Growth good / very good; Zone diameters within specifications
<i>Staphylococcus aureus</i>	25923	16 – 18 h at 35 °C	Growth good / very good; Zone diameters within specifications
<i>Pseudomonas aeruginosa</i>	27853	16 – 18 h at 35 °C	Growth good / very good; Zone diameters within specifications
<i>Enterococcus faecalis</i>	29212	16 – 18 h at 35 °C	Growth good / very good; Zone diameters within specifications

Please refer to the actual batch related Certificate of Analysis.

Literature

Deweese, et al. 1970. Effect of storage of Mueller Hinton Agar plates on zone sizes for antimicrobial testing. Appl. Microbiol.; 30:203.

National Committee for Clinical Laboratory Standards. Approved Standard. Performance standards for antimicrobial disc susceptibility tests, 5th ed. National Committee for Laboratory Standards, Villanova, Pa. (1993).

NCCLS Approved Standard: ASM-2, 1979, Performance Standards for Antimicrobial disc Susceptibility Tests, 2nd Ed., National Committee for Clin. Lab. Standards.

Ryan, K.J., et al. 1970. Disk sensitivity testing. Hosp. Prac.; 5:91-100.

Ordering Information

Product	Cat. No.	Pack size
MUELLER-HINTON agar acc. to CLSI	1.05435.0500	500 g

Merck KGaA, 64271 Darmstadt, Germany
Fax: +49 (0) 61 51 / 72-60 80
mibio@merckgroup.com
www.merckmillipore.com/biomonitoring

Find contact information for your country at:
www.merckmillipore.com/offices
For Technical Service, please visit:
www.merckmillipore.com/techservice



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Merck Millipore and the M logo are registered trademarks of Merck KGaA, Darmstadt, Germany. ATCC is a registered trademark of ATCC, Manassas, VA, USA. Lit. No. TN1461EN00