



Air Gap

General Description

The air gap is a small plastic device that is placed on tubing going from the water system to the drain. This creates an unobstructed vertical space, the “air gap” which will prevent contaminants in the drain from going back into the system.

Key Benefits

- Provides cost-effective protection of the water system and its accessories.
- Prevents bacteriological contamination from the drain.

Ordering Information

Description	Cat. No.
Air Gap Device	AIRGAP001



To place an order or receive technical assistance

In Europe, please call Customer Service:

France: 0825 045 645
Germany: 069 86798021
Italy: 848 845 645
Spain: 901 516 645 Option 1
Switzerland: 0848 645 645
United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: merckmillipore.com/offices

For Technical Service visit:
merckmillipore.com/techservice
merckgroup.com/life-science

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, Milli-Q, and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK_DS1252EN
2017 - 08107
1/2018