

The Future of Drug Delivery

Ready-to-use preformulated lyophilized liposomes

Liposomes are one of the most common drug delivery carriers for small molecules and nucleic acid therapeutics.


We offer a variety of preformulated liposomes including:

- Cationic liposomes
- Natural liposomes
- PEGylated liposomes

Upon hydration, the liposomes are ready to use for your drug delivery needs.

For a complete list of preformulated lyophilized liposomes visit:

SigmaAldrich.com/liposomes



Lyophilized Cationic Liposomes

Lipid Composition	Molar Ratio	Mean Size	Cat. No.
DOTAP-based			
DOTAP:DOPC	50:50	100 nm	LDP-501-5MG
			LDP-501-10MG
	40:60	100 nm	LDP-502-5MG
			LDP-502-10MG
	30:70	100 nm	LDP-503-5MG
			LDP-503-10MG
	20:80	100 nm	LDP-504-5MG
			LDP-504-10MG
	10:90	100 nm	LDP-505-5MG
			LDP-505-10MG

Lipid Composition	Molar Ratio	Mean Size	Cat. No.
DOTAP-based			
DOTAP:DOPC	5:95	100 nm	LDP-506-5MG
			LDP-506-10MG
	2:98	100 nm	LDP-507-5MG
			LDP-507-10MG
	1:99	100 nm	LDP-508-5MG
			LDP-508-10MG
	0.5:99.5	100 nm	LDP-509-5MG
			LDP-509-10MG

Lyophilized Anionic Liposomes

Lipid Composition	Molar Ratio	Mean Size	Cat. No.
Phosphatidylglycerol (PG)-based			
DOPG	100	100 nm	LPG-514-5MG
			LPG-514-10MG
DOPG:DOPC	90:10	100 nm	LPG-501-5MG
			LPG-501-10MG
	80:20	100 nm	LPG-502-5MG
			LPG-502-10MG
	70:30	100 nm	LPG-503-5MG
			LPG-503-10MG
	60:40	100 nm	LPG-504-5MG
			LPG-504-10MG
	50:50	100 nm	LPG-505-5MG
			LPG-505-10MG
	40:60	100 nm	LPG-506-5MG
			LPG-506-10MG
	30:70	100 nm	LPG-507-5MG
			LPG-507-10MG
	20:80	100 nm	LPG-508-5MG
			LPG-508-10MG
	10:90	100 nm	LPG-509-5MG
			LPG-509-10MG
	5:95	100 nm	LPG-510-5MG
			LPG-510-10MG
2:98	100 nm	LPG-511-5MG	
		LPG-511-10MG	
1:99	100 nm	LPG-512-5MG	
		LPG-512-10MG	
0.5:99.5	100 nm	LPG-513-5MG	
		LPG-513-10MG	

Lipid Composition	Molar Ratio	Mean Size	Cat. No.
Phosphatidylserine (PS)-based			
DOPS	100	100 nm	LPG-514-5MG
			LPG-514-10MG
DOPS:DOPC	90:10	100 nm	LPG-501-5MG
			LPG-501-10MG
	80:20	100 nm	LPG-502-5MG
			LPG-502-10MG
	70:30	100 nm	LPG-503-5MG
			LPG-503-10MG
	60:40	100 nm	LPG-504-5MG
			LPG-504-10MG
	50:50	100 nm	LPG-505-5MG
			LPG-505-10MG
	40:60	100 nm	LPG-506-5MG
			LPG-506-10MG
	30:70	100 nm	LPG-507-5MG
			LPG-507-10MG
	20:80	100 nm	LPG-508-5MG
			LPG-508-10MG
	10:90	100 nm	LPG-509-5MG
			LPG-509-10MG
	5:95	100 nm	LPG-510-5MG
			LPG-510-10MG
2:98	100 nm	LPG-511-5MG	
		LPG-511-10MG	
1:99	100 nm	LPG-512-5MG	
		LPG-512-10MG	
0.5:99.5	100 nm	LPG-513-5MG	
		LPG-513-10MG	

Lyophilized PEGylated Liposomes

Lipid Composition	Molar Ratio	Mean Size	Cat. No.
PEGylated Anionic-based			
DSPC:Chol:DSPG: PEG2000-DSPE	60:30:5:5	100 nm	LPCG-501-5MG LPCG-501-10MG
DOPC:Chol:DOPG: PEG2000-DOPE	65:30:5:5	100 nm	LPGG-502-5MG LPGG-502-10MG
DOPC:Chol:DSPS: PEG2000-DOPE	65:30:5:5	100 nm	LPSG-501-5MG LPSG-501-10MG
DOPC:Chol:DOPS: PEG2000-DOPE	65:30:5:5	100 nm	LPSG-502-5MG LPSG-502-10MG
PEGylated Cationic-based			
DSPC:Chol:DSTAP: PEG2000-DSPE	55:30:10:5	100 nm	LDPG-501-5MG LDPG-501-10MG
DOPC:Chol:DOTAP: PEG2000-DOPE	55:30:10:5	100 nm	LDPG-502-5MG LDPG-502-10MG

Lyophilized Liposomes from Diether Lipids

Lipid Composition	Molar Ratio	Mean Size	Cat. No.
Diether-based			
16:0 Diether PC	100	100 nm	LDE-504-5MG LDE-504-10MG
18:0 Diether PC	100	100 nm	LDE-505-5MG LDE-505-10MG
18:1 Diether PC	100	100 nm	LDE-506-5MG LDE-506-10MG
16:0 Diether PC: Cholesterol	70:30	100 nm	LDE-510-5MG LDE-510-10MG
18:0 Diether PC: Cholesterol	70:30	100 nm	LDE-511-5MG LDE-511-10MG
18:1 Diether PC: Cholesterol	70:30	100 nm	LDE-512-5MG LDE-512-10MG

Lipid Composition	Molar Ratio	Mean Size	Cat. No.
PEGylated Neutral-based			
HSPC:Chol: PEG2000-DSPE	60:30:5:5	100 nm	LPCG-501-5MG LPCG-501-10MG
DOPC:Chol: PEG2000-DOPE	65:30:5:5	100 nm	LPGG-502-5MG LPCG-502-10MG
DOPC:Chol:DSPS: PEG2000-DOPE	65:30:5:5	100 nm	LPCG-504-5MG LPCG-504-10MG

Lyophilized Liposomes from Unsaturated Lipids

Lipid Composition	Molar Ratio	Mean Size	Cat. No.
Unsaturated Phosphatidylcholine (PC)-based			
16:0-18:2 PC	100	100 nm	LUP-513-5MG LUP-513-10MG
16:0-20:4 PC	100	100 nm	LUP-514-5MG LUP-514-10MG
18:0-:18:1 PC	100	100 nm	LUP-515-5MG LUP-515-10MG
18:0-18:2 PC	100	100 nm	LUP-516-5MG LUP-516-10MG
18:0-20:4 PC	100	100 nm	LUP-517-5MG LUP-517-10MG
18:1-14:0 PC	100	100 nm	LUP-518-5MG LUP-518-10MG
18:1-16:0 PC	100	100 nm	LUP-519-5MG LUP-519-10MG
18:1-18:0 PC	100	100 nm	LUP-520-5MG LUP-520-10MG

See next page for
Lyophilized Liposomes
from Unsaturated Lipids

Lyophilized Liposomes from Unsaturated Lipids

Lipid Composition	Molar Ratio	Mean Size	Cat. No.
Natural Phosphatidylcholine (PC)-based			
Brain PC	100	100 nm	LPC-601-5MG
			LPC-601-10MG
Egg PC	100	100 nm	LPC-602-5MG
			LPC-602-10MG
Heart PC	100	100 nm	LPC-603-5MG
			LPC-603-10MG
Liver PC	100	100 nm	LPC-604-5MG
			LPC-604-10MG
Soy PC	100	100 nm	LPC-605-5MG
			LPC-605-10MG
Brain PC:Cholesterol	70:30	100 nm	LPC-606-5MG
			LPC-606-10MG
Egg PC:Cholesterol	70:30	100 nm	LPC-607-5MG
			LPC-607-10MG
Heart PC:Cholesterol	70:30	100 nm	LPC-608-5MG
			LPC-608-10MG
Liver PC:Cholesterol	70:30	100 nm	LPC-609-5MG
			LPC-609-10MG
Soy PC:Cholesterol	70:30	100 nm	LPC-610-5MG
			LPC-610-10MG



To place an order or receive technical assistance in Europe, please call Customer Service:
 France: 0825 045 645 Spain: 901 516 645 Option 1
 Germany: 069 86798021 Switzerland: 0848 645 645
 Italy: 848 845 645 United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840
 Or visit: [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices)
 For Technical Service visit: [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice)

Merck KGaA
 Frankfurter Strasse 250
 64293 Darmstadt, Germany

[SigmaAldrich.com](https://www.sigmaaldrich.com)

We have built a unique collection of life science brands with
 unrivalled experience in supporting your scientific advancements.

Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, BioReliance, Millipore, Milli-Q, SAFC, Sigma-Aldrich, and Supelco are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK_BR9872EN Ver. 0.0 44086 09/2022

