



MERCK

SAFETY



INSIDE

Daily safety, risk prevention and emergency help that meet your high expectations.



The life science
business of Merck
operates as
MilliporeSigma in
the U.S. and Canada.

Supelco[®]
Analytical Products

SAFETY INSIDE

Trust our safety expertise

We combine ingenious products and packaging solutions with individual support to keep you safe, simplify your lab work, and help you save resources. Discover daily safety, risk prevention and emergency solutions that match your high standards.

Discover more:

[SigmaAldrich.com/lab-safety](https://www.sigmaaldrich.com/lab-safety)



Safebreak bottles for acids

Our PE-coated glass bottles protect you from splinters and chemical burns in case of accidents



Integrated safety right from the start

Risk prevention

Preventing harm starts with the right packaging and chemicals. We offer 350 years of experience in handling and designing packaging for hazardous goods, so you can enjoy unrivaled safety standards and application-oriented foresight. Our specially developed packaging materials ensure you work safely even during critical applications.



Silica gel

Our orange gel contains an indicator based on iron-salt instead of the carcinogenic cobalt dichloride used in blue silica gel



HDPE bottles

Lightweight and unbreakable alternative to glass that offers easy handling



HDPE, metal & stainless steel containers

Safe handling of higher volumes



Extran® detergents

Enjoy reliable, sustainable cleaning power with our non-toxic, biodegradable detergents



Black HDPE bottle

Protecting chemicals from UV light. Each is fitted with our unique SafetyCap, which releases excess gas to decrease internal pressure, but is absolutely leakproof against liquids

Tools & products for your Daily safety



Handling hazardous goods safely is always a challenge in daily lab work. We make it easy! With our perfectly matched components – chemicals, packaging, safety tools and withdrawal systems – you'll benefit from safe, convenient and contamination-free work even during highly sensitive applications. Our specialized tools help you avoid harmful accidents, chemical burns, costly contaminations and time-consuming issue tracking.



Withdrawal systems

Ingeniously safe solutions for:

- Solvents in bottles
- Solvents, acids & bases in higher volumes



Bottle opening keys

Safe and easy handling from the moment you open a bottle



Pouring aid

Our innovative aid simplifies pouring from 2.5 l glass bottles to protect your hands from chemical burns



Safety carrier

The carrier features PE-foam bottle protection and a tightly closing screw lid to help you avoid contact with harmful chemicals, vapors and splinters



Label set

Label your own mixtures with comprehensive safety information and GHS pictograms



Swift, simple and efficient

Emergency help

There is a lot you can do to prevent safety risks. But, especially when working with glass bottles or open vessels in labs, accidents can happen anytime. To avoid serious injuries, it is important to act quickly but wisely. With the Chemizorb® range, we offer a one-stop solution for rapid, reliable absorption and neutralization of unpleasant liquid spills.

Chemizorb® absorbents

The fast, safe and easy way to clear up aggressive or unpleasant chemical spills



Supelco®

Analytical Products

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

SigmaAldrich.com

© 2020 Merck KGaA, Darmstadt, Germany and its affiliates. All Rights Reserved.
Merck, the vibrant M, Supelco, Extran and Chemizorb are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners.
Detailed information on trademarks is available via publicly accessible resources.

Lit. No. MK_BR3854EN Ver. 1.0
07/2020