



## NovaSeptum® GO Sterile Sampling System

New and improved system offers additional security and ease-of-use.

Continuing to set the sampling standard, our NovaSeptum® GO sterile sampling system is equipped with additional features for even safer sampling throughout your entire process.

Enhanced with new magazines and trigger that are compatible with NovaSeptum® stainless steel bases, along with a new locking mechanism that prevents accidental actuation, you now have the flexibility and control to sample the way you want, where you want.

Common Features	NovaSeptum® Sterile Sampling System and NovaSeptum® GO Sterile Sampling System
Connection	<ul style="list-style-type: none"> <li>• High level of customization: connector types (Tri-Clamp, InLine, Ingold)</li> <li>• Stainless-steel bases</li> <li>• Preloaded or loadable on site</li> </ul>
Sample mechanism	<ul style="list-style-type: none"> <li>• 2mm ID cannula needle</li> </ul>
Material of construction	<ul style="list-style-type: none"> <li>• Validated materials and system</li> <li>• Containers and resin</li> </ul>
Disconnection	<ul style="list-style-type: none"> <li>• NovaSeal® crimper</li> </ul>

New Features	NovaSeptum® GO Sterile Sampling System
Device security and safety: accidental actuation avoidance	<ul style="list-style-type: none"> <li>• First actuation locking tag</li> <li>• Ring actuation avoidance</li> </ul>
Port plug	<ul style="list-style-type: none"> <li>• Improved port plug, easy to spot and reduces risk of leakage</li> </ul>
Preloaded feature	<ul style="list-style-type: none"> <li>• Less plastic to recycle</li> <li>• Reusable inox base</li> <li>• Base stays on the tank: reduce NovAseptic® connect disassembly operations</li> </ul>

If you are already a NovaSeptum® system user and interested in upgrading to the new NovaSeptum® GO sterile sampling system, we can help. Contact your sales representative for more information.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2020 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, NovaSeptum, NovaSeal and NovAseptic are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK\_FL5416EN Ver. 1.0 30502 04/2020

