

## **TraceCERT<sup>®</sup> Pesticide CRMs**

For accredited testing labs, the availability of reliable certified reference materials (CRMs) is crucial for calibration and validation in compliance with ISO/IEC 17025. High-performance quantitative NMR (HP-qNMR) is applied for the certification of our *TraceCERT<sup>®</sup>* pesticide reference materials as a relative primary method to achieve traceability to the SI unit via NIST SRMs (Standard Reference Materials).

Our comprehensive portfolio of *TraceCERT<sup>®</sup>* reference materials offer the following characteristics:

- Certified content by quantitative NMR (qNMR)
- Manufactured under ISO/IEC 17025 and ISO 17034 double accreditation
- Highest level of accuracy, calculated uncertainties, lot-specific values and traceability
- Comprehensive documentation delivered with the product according to ISO Guide 31


### **Accuracy and precision made simple**

To see our full portfolio of pesticide standards visit us at:

[SigmaAldrich.com/pesticides](https://www.sigmaaldrich.com/pesticides)

For our complete offer of qNMR certified reference materials see:

[SigmaAldrich.com/organiccrm](https://www.sigmaaldrich.com/organiccrm)



Cat. No.	Product Description
08288-100mg	Alachlor
68417-50mg	Aldicarb
08573-50mg	Aldrin
90935-50mg	Atrazine
74116-50mg	Azinphos-ethyl
43854-50mg	Azoxystrobin
96996-100mg	Biphenyl
64148-50mg	Bromobutide
74047-50mg	Captan
79888-50mg	Carbendazim
94114-50mg	Chlorpyrifos
68965-50mg	Clothianidin
72765-50mg	$\lambda$ -Cyhalothrin
90024-50mg	Cymoxanil
51991-50mg	Cypermethrin
90622-50mg	Cyproconazol
76514-100mg	2,4-D
73073-50mg	2,4-DB
44012-50mg	2,4-DDD
43537-50mg	4,4'-DDE
80076-100mg	4,4'-DDT
05995-50mg	Deltamethrin
68486-50mg	Diazinon
16826-100mg	Dicamba
78336-50mg	Dichloran
74046-50mg	Dicofol
77070-50mg	Diflubenzuron
16902-50mg	Diuron
74119-50mg	$\alpha$ -Endosulfan
40828-50mg	$\beta$ -Endosulfan
79577-50mg	Ethylparaben
89929-50mg	Etofenprox
90664-50mg	Fenvalerate
16785-50mg	Fipronil
74123-50mg	Fludioxonil
49677-100mg	Glufosinate-ammonium

Cat. No.	Product Description
89432-100mg	Glyphosate
74142-50mg	$\alpha$ -HCH
90426-50mg	Hepatchlor
68694-50mg	Imidacloprid
89977-50mg	Iprodione
90053-50mg	Kresoxim-methyl
90524-50mg	Linuron
91481-50mg	Malathion
89979-50mg	Metalaxyl
89933-50mg	Methamidophos
74189-50mg	Methidathion
75296-50mg	Metribuzin
79721-50mg	Methylparaben
90521-50mg	Monocrotophos
51625-50mg	N-(6-Chloro-3-pyridylmethyl)-N-cyano-N-methylacetamide
43900-50mg	Paclobutrazol
90067-50mg	Penconazol
89931-50mg	Pendimethalin
78447-50mg	Piperonylbutoxide
64947-50mg	Prochloraz
75120-50mg	Procymidone
16897-50mg	Propylparaben
95700-50mg	Propyzamide
92597-50mg	Pyriproxyfen
43897-50mg	Tebuconazole
74124-50mg	Terbutylazine
67554-50mg	Thiabendazole
14783-50mg	Thiacloprid
68901-50mg	Thiamethoxam
43966-50mg	Thiram
90619-50mg	Tolyfluanid
68298-50mg	Triazophos
74022-50mg	Triflumuron
90413-50mg	Vinclozolin

## To place an order or receive technical assistance

Order/Customer Service: [SigmaAldrich.com/order](https://SigmaAldrich.com/order)  
 Technical Service: [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice)  
 Safety-related Information: [SigmaAldrich.com/safetycenter](https://SigmaAldrich.com/safetycenter)

[SigmaAldrich.com](https://SigmaAldrich.com)

Merck KGaA  
 Frankfurter Strasse 250  
 64293 Darmstadt, Germany

© 2019 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Supelco and TraceCERT are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

