



CHICKEN ANTI-POLYHISTIDINE TAG AFFINITY PURIFIED POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB3517
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	1.0 mg/mL
SPECIFICITY:	Polyhistidine tag.
IMMUNOGEN:	Synthetic peptide.
APPLICATIONS:	Immunoblotting: 1:1,000-1:5,000 Immunocytochemistry: 1:400-1:2,000 ELISA: 1:1,000-1:5,000 Optimal working dilutions must be determined by the end user.
FORMAT:	Affinity purified immunoglobulin.
PRESENTATION:	Liquid in PBS containing 0.02% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

Now Available - Chicken IgY Negative Control Antibody (Cat. No. PP25),
Chicken IgY Blocking Solution (Cat. No. 2180) and Chicken IgY Precipitating Agarose (Cat. No. 2185)

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.