

G3789 G.N. (Gram Negative) Spore Anaerobic Supplement

An enrichment and antibiotic supplement recommended for the selective isolation of gram-negative anaerobes.

Composition

(per vial, sufficient for 500 ml medium)

Hemin 2.5 mg

Sodium Succinate 1.25 g

Nalidixic Acid 5 mg

Vancomycin 1.25 mg

Menadione 0.25 mg

Store at 4°C.

Directions

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile filtered 0.2 N sodium hydroxide. Mix well and aseptically add it to 490 ml of sterile, molten Wilkins Chalgren Anaerobic Agar (W1761), enriched with 5% v/v sterile defibrinated blood. Mix well and pour into sterile petri plates.

G3789 G.N. (Gram Negative) Spore Anaerobic Supplement

An enrichment and antibiotic supplement recommended for the selective isolation of gram-negative anaerobes.

Composition

(per vial, sufficient for 500 ml medium)

Hemin 2.5 mg

Sodium Succinate 1.25 g

Nalidixic Acid 5 mg

Vancomycin 1.25 mg

Menadione 0.25 mg

Store at 4°C.

Directions

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile filtered 0.2 N sodium hydroxide. Mix well and aseptically add it to 490 ml of sterile, molten Wilkins Chalgren Anaerobic Agar (W1761), enriched with 5% v/v sterile defibrinated blood. Mix well and pour into sterile petri plates.