

The pHast solution for gel electrophoresis

pHast Pack™ Ready-To-Use Buffers & Ready-to-Pour Agarose Gels

The pHast Pack™ platform was developed to expedite your research and provide a convenient solution for tedious lab tasks. Our buffers and new ready-to-pour agarose gels are formulated to save time and effort in native RNA and DNA electrophoresis assays:

- **Fast, reliable buffer prep for precise accuracy** without the need for calculations, weighing, or adjusting pH
- **Accelerate gel preparation** – simply dissolve in boiling water, add your preferred nucleic acid stain, then pour the gel
- **Available in 1% agarose gel** with high gel strength and a band resolution of 300 – 8,000 bp
- **Certified performance** – suitable for electrophoresis with no trace elements, DNase, RNase, or Nickase

Simplify your work with pHast Pack™. Discover our ready-to-use buffers and ready-to-pour agarose gels that are cost effective and convenient for research use in gel electrophoresis assays.

Cat. No.	pHast pack™	Application
PPB008-20PAK	Tris Acetate EDTA buffer (TAE), pH 8.3	1X Running buffer
PPB009-20PAK	Tris Borate EDTA buffer (TBE), pH 8.3	1X Running buffer
PPB010-20PAK	10X Tris EDTA buffer (TE), pH 8.0	10X Stock buffer
PPB012-20PAK	1.0% Agarose in TAE, pH 8.3	1% Agarose gel
PPB013-20PAK	1.0% Agarose in TBE, pH 8.3	1% Agarose gel



For more information, visit
SigmaAldrich.com/pHast-Pack

The Life Science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Sigma-Aldrich, and pHast Pack are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK_FL9076EN Ver. 2.0 12/2022

