

#### Product No. A 0420 Lot 077H4863

## Anti-Human IgM (µ-chain specific) Peroxidase Conjugate

Antibody developed in Goat Affinity Isolated Antigen Specific Antibody

Antiserum is developed in goat using the human IgM as the immunogen. The antibody is isolated from goat anti-human IgM antiserum by immunospecific purification to remove essentially all goat serum proteins, including immunoglobulins, which do not specifically bind to the  $\mu$ -chain of human IgM. Goat anti-human IgM is conjugated to peroxidase by means of a two-step glutaraldehyde method. The product is purified to remove unconjugated material and supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 0.01% thimerosal as a preservative.

#### **Specificity**

Specificity of the Peroxidase Conjugated Anti-Human IgM is determined by Enzyme Linked Immunosorbent Assay (ELISA). The conjugate is specific for human IgM when tested against human IgA, IgG, IgM, Bence Jones kappa and lambda myeloma proteins.

Cross-reactivity of the antibody-conjugate is determined by ELISA. The conjugate shows no reactivity with mouse or rat IgG.

#### **Identity and Purity**

Identity and purity of the antibody is established by immunoelectrophoresis (IEP), prior to conjugation. Electrophoresis of the antibody preparation followed by diffusion versus anti-goat IgG and anti-goat whole serum results in single arcs of precipitation.

**Molar Ratio** (IgG:Peroxidase) = 1.1

# ABPT

In an agar diffusion assay the conjugate produces a precipitation arc at a dilution of 1:32 versus a 1:10 dilution of IgM enriched human serum.

Enzyme Activity: 540 purpurogallin units/ml

Enzyme activity is determined using 5% Pyrogallol (Sigma Product No. P 0381) in deionized water, pH 6.0, at 20°C. One purpurogallin unit will form 1 mg of purpurogallin from pyrogallol in 20 seconds at pH 6.0, 20°C.

#### **Titers**

#### 1. 1:50,000 (Direct ELISA)

We are now reporting lot specific information as a titer by direct ELISA rather than as a working dilution (see below). Titer is defined as the dilution of conjugate sufficient to give a change in absorbance of 1.0 at 450 nm after 30 minutes of substrate conversion at  $25\,^{\circ}\mathrm{C}$  (Voller, et al.)¹. Microtiter plates are coated with purified human IgM at a concentration of 5 µg/ml in 0.05 M carbonate-bicarbonate buffer, pH 9.6 (Carbonate-Bicarbonate Buffer Capsules are available as Sigma Product No. C 3041).

**Substrate:** *o*-Phenylenediamine Dihydrochloride (OPD, Sigma Product No. P 8287), 0.4 mg/ml in 0.05 M phosphate-citrate buffer, pH 5.0 containing 0.03% sodium perborate (Phosphate-Citrate Buffer Capsules with Sodium Perborate are available as Sigma Product No. P 4922).

### 2. Dot Blot

- A dilution of 1:6,000 was determined by direct dot blot assay using 40 ng human IgM/dot.
- A dilution of 1:100,000 was determined in a direct chemiluminescence assay using 20 ng human IgM/dot. Luminol plus enhancer was used as substrate.

## 3. Immunohistology

A dilution of 1:200 was determined in a direct assay using formalin-fixed, paraffin- embedded human tonsil sections.

## **Working Dilutions**

Working dilutions should be determined by titration assay. Due to differences in assay systems, these titers may not reflect the user's actual working dilution.

This goat antisera was maintained at pH 5.0 for 40 minutes to meet USDA requirements.

#### Storage

For continuous use, store at 2-8°C for up to one month. For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is **not** recommended. Storage in "frost-free" freezers is **not** recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

#### Reference

1. Voller, A., et al., Bulletin WHO, **53**, 55 (1976).

Sigma warrants that its products conform to the information contained in this and other Sigma products. Purchaser must determine the suitability of the products for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Issued 09/97.