

RAISING THE BAR IN
STERILITY
TESTING

With the Steritest® NEO Membrane Filtration Devices

The Life Science business
of Merck operates as
MilliporeSigma in the
U.S. and Canada.

Millipore®

Preparation, Separation,
Filtration & Monitoring Products

DESIGNED TO YOUR NEEDS

THE SUM OF EACH ELEMENT:
Steritest® NEO device

Since 2018, when we introduced the new Steritest® NEO, we have launched new elements every year to significantly improve your workflow safety, reliability and convenience.

Introducing the Velax® cutting clamp, a safe and easy way to cut the tubing. Scissors are no longer needed.

We invented the original Steritest® closed filtration system in the 1970's which changed sterility testing forever. And the Steritest® NEO devices is the 4th generation. We are continuously raising the industry standards, and driving innovation to provide you the most secure and reliable test ever.

Thanks to our strong and long experience with the Sterility Test, we have a complete offer with the Steritest® NEO devices, the Steritest® culture media and rinsing fluids, the Steritest® Symbio pump, and broad service offering with method development and consultancy, validation, training and more.

THE EVOLUTION OF SAFETY

Feel safe: tubing disconnection tool

Experience the scissor free Velax® cutting clamp and reduce the risks during the test. A safe and convenient way to cut the tubing using the Velax® cutting clamp mounted on every Steritest® NEO device.

Feel confident: colored clamps

Prevent filling errors and improve your workflow clarity, thanks to the pre-installed colored clamps and the existing blackline for accurate media filling.

Feel safer: newly designed needle guard and needle protector

Grips on the guard and ridges on the protector improve the confidence in needle manipulation.

THE EVOLUTION OF CONVENIENCE

Feel control: a brand new short needle for small sample containers

Better dexterity with the new 20 mm length needle when piercing cartridges or small soft plastic containers, without compromising the flow rate.

Feel flexible: protective caps for long needles are now in 2 parts

The protective cap in 2 parts gives access to either a short (35 mm) or a long (60 mm) needle designed to fit your sample packaging configuration. Color-coded protectors help you to differentiate the needle type once covered.

Feel free: upgraded accessory bag

Simplified opening of the accessory bag improves your workflow convenience thanks to the perforated tear strip.

Feel comfortable: new placement mark

Be sure of placing the Steritest® NEO tube in the pump head precisely by using the new placement mark.

THE EVOLUTION OF RELIABILITY

Feel secure: optimized identification and traceability

Clear packaging identification:
The selection of the appropriate box of Steritest® NEO devices is facilitated thanks to the new designed label using color coding linked to canister base color and using a needle/application drawing.

Smart Label with 2D data matrix barcode to retrieve product data, including the item code, batch number, shelf life, and global trade item number (GTIN). It allows an easy data processing in a broad range of systems (EPR, LIMS, etc...)

Learn more on
SigmaAldrich.com/smart-label

Peel-and-stick label with critical information to place in a lab notebook for accurate tracking.

Feel confident: volume graduation on the canisters

Be precise and improve your workflow accuracy through the addition of a 25 mL graduation line and volume engraved in the Steritest® NEO canisters.

Learn more at SigmaAldrich.com/steritestneo

To Place an Order or Receive Technical Assistance

Order/Customer Service: SigmaAldrich.com/order

Technical Service: SigmaAldrich.com/techservice

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

