

40568 Bolton Selective Supplement

This supplement is recommended for the selective enrichment of *Campylobacter* species from foods as recommended by ISO Committee under specifications of ISO 10272-1:2017.

Composition

(per vial, sufficient for 500 ml medium):

Cefoperazone 10 mg

Vancomycin 10 mg

Trimethoprim 10 mg

Amphotericin 5 mg

Directions

Rehydrate the contents of 1 vial aseptically with 2.5 ml of 50% ethanol. Mix well and aseptically add to 472.5 ml of sterile, molten Bolton Broth, Base (Cat. No. 67454), with 25 ml of sterile lysed defibrinated horse blood. Mix well and dispense into sterile tubes.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

