

## C3222 Campylobacter Supplement III (Skirrow)

An antibiotic supplement recommended for the selective isolation of the Campylobacter species.

### Composition

(per vial, sufficient for 500 ml medium)

Polymyxin B	1'250 IU
Vancomycin	5 mg
Trimethoprim	2.5 mg

Store at 4°C.

### Directions

Rehydrate the contents of one vial aseptically with 2 ml of sterile distilled water and mix well to dissolve. Avoid frothing of the solution. Aseptically, add the contents of the vial to 500 ml of sterile, molten Columbia Agar Base (Prod. No 27688) with 5-7% v/v lysed defibrinated horse blood. Mix gently and pour into sterile petri plates.

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

